

Data Acquisition and Analysis of 7-ethoxycoumarin and its Metabolites Using TurboDDS Software in the Agilent 500 Ion Trap LC/MS

Application Note

Clinical Research

Author

Fran Lai
Agilent Technologies, Inc.
5301 Stevens Creek Boulevard
Santa Clara, CA 95051
USA

Abstract

This note describes an MS^n data acquisition method using the 500 Ion Trap LC/MS TurboDDS software for the analysis of 7-ethoxycoumarin and two of its metabolites. TurboDDS saves time and streamlines data reporting for multiple stages of MS/MS and MS^n analysis.

Introduction

The ion trap is a powerful tool for metabolism studies, due to the MS^n capability for the identification of metabolites and structure elucidation of unknown peaks. The TurboDDS software on the Agilent 500 Ion Trap LC/MS eliminates the need to do multiple runs in order to identify metabolites. Once the expected m/z of a metabolite is entered into the method, the 500 Ion Trap will automatically look for this m/z and when found, will immediately trigger an MS^n (where n can be 2 or higher) experiment.

7-ethoxycoumarin (Figure 1A) is commonly used for metabolism studies, including hepatocyte studies [1], liver perfusion techniques [2] and cytochrome P450 studies [3]. This application note describes an MS^n data acquisition method using the Agilent TurboDDS for the analysis of 7-ethoxycoumarin and two of its metabolites: 7-hydroxycoumarin (Figure 1B) and 7-hydroxycoumarin β -D-glucuronide (Figure 1C).



Agilent Technologies

Instrumentation

The following instruments were used in this study:

- Agilent 500 Ion Trap LC/MS equipped with ESI source
- Agilent ProStar 420 AutoSampler
- Agilent 212-LC Binary Solvent Delivery Modules equipped with 150 μ L Static Mixer

Materials and Reagents

The following reagents were used and obtained from Sigma Aldrich, St. Louis, MO:

- 7-ethoxycoumarin (CAS number 31005-02-4),
- 7-hydroxycoumarin (CAS number 93-35-6)
- 7-hydroxycoumarin β -D-glucuronide sodium salt (CAS number 168286-98-4).

All other chemicals were reagent grade or HPLC grade.

Sample Preparation

To simulate a typical mixture of a drug and its metabolites from metabolic study samples, a mixture of 7-ethoxycoumarin (7-EC), 7-hydroxycoumarin (7-HC) and 7-hydroxycoumarin β -D-glucuronide sodium salt (7-HC glu) in 50:50 methanol:water was prepared, with concentrations of 0.12, 5.5 and 15 ng/ μ L respectively. These concentration levels were used in order to show comparative peak sizes in the resultant chromatograms.




Figure 1A. Structure of 7-ethoxycoumarin (7-EC).




Figure 1B. Structure of 7-hydroxycoumarin (7-HC).




Figure 1C. Structure of 7-hydroxycoumarin glucuronide (7-HC glu).

HPLC Conditions

Column: Pursuit XRs C18, 150 \times 2 mm, 5 μ m (Agilent p/n A6000150X020)

Solvent A: 0.1% formic acid in water

Solvent B: 0.1% formic acid in acetonitrile

LC program:	Time (min:sec)	%A	%B	Flow (μ L/min)
	0:00	60	40	200
	5:00	15	85	200
	5:30	15	85	200
	5:42	60	40	200
	7:00	60	40	200

Injection volume: 20 μ L

Injection solvent: 50:50 methanol:water

MS and API Parameters

API drying gas:	25 psi at 400 °C
API nebulizing gas:	25 psi
Needle:	5000 V
Capillary:	74 V
Shield:	500 V
Enhanced scan mode:	5000 Da/sec
RF loading:	72%

TurboDDS Scan Parameters

Survey scan: m/z 100–400

“Include List” for MS^n triggers:

- m/z 191 (for 7-EC)
- m/z 163 (O-deethylation [–28] to 7-HC)
- m/z 339 (glucuronidation [+176] to 7-HC glu)

Data Dependent Scans (DDS) with MS^2 and MS^3 were included in the method.

Results and Discussion

Table 1A. TurboDDS results with MS^2 :

Analyte	Precursor ion (m/z)	Product ion (m/z)
7-EC	191	163
7-HC	163	119
7-HC glu	339	163

Table 1B TurboDDS results with MS^3 :

Analyte	Precursor ion (m/z)	Product ion (m/z)
7-EC	163	119
7-HC	119	91
7-HC glu	163	119




Figure 2. Chromatogram of 7-EC, 7-HC and 7-HC glu (top). All scan descriptors, (bottom) Survey scan.

Figures 3–5 (tree reports) show MS² and MS³ data tabulated in Tables 1A and 1B.





Figure 3. Survey, MS² and MS³ spectra of 7-EC (peak #3) from TurboDDS.



Product Ion Tree Report for 163.0 m/z

File: ...x5-dds-ms2-2usc_0202-coumarin8d-dds-ms3exin.xms

Survey Scans [99.5 : 400.5] Scans: 436 - 444




MS2 [99.5 : 400.5] - 163.0 > 59.1 : 173.5

Figure 4. Survey, MS² and MS³ spectra of 7-HC (peak #2) from TurboDDS.



Product Ion Tree Report for 339.0 m/z

File: ...x5-dds-ms2-2usc_0202-coumarin8d-dds-ms3exin.xms

Survey Scans [99.5 : 400.5] Scans: 360 - 374



MS2 [99.5 : 400.5] - 339.0 > 122.2 : 349.5

Figure 5. Survey, MS² and MS³ spectra of 7-HC glu (peak #1) from TurboDDS.

These tree reports (Figures 3–5) make it very easy for an operator to view product ions and trace them back to the original precursor ions. Also, mass differences are displayed at each MS/MS step, a clear benefit for structural elucidation.

Conclusion

The Agilent 500 Ion Trap LC/MS TurboDDS software demonstrates excellent performance for data dependent scans of 7-ethoxycoumarin and two of its metabolites.

MSⁿ was performed up to MS³ in this application note. However, further MS/MS experiments could be performed in the same run as long as there is enough sample. Without TurboDDS, the same data would have required multiple runs, but with it, the time saving is tremendous and the data reporting for the multiple MS/MS stages is streamlined.

References

1. Carlile et al., Drug Metabolism and Disposition 26 (1998) 216.
2. Andersson et al., Drug Metabolism and Disposition 11 (1983) 494.
3. Yamazaki et. al., Biochemical Pharmacology 51 (1996) 313.

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2011
Printed in the USA
February 24, 2011
SI-0926



Agilent Technologies