

Detection of Recombinant Bovine Somatotropin in Milk by LC-ESI-MS/MS



Food Safety

Authors

Marie-Hélène Le Breton, Sandrine Rochereau-Roulet, Gaud Pinel,
Fabrice Monteau, and Bruno Le Bizec
Laboratoire d'Etude des Résidus et Contaminants dans les
Aliments (LABERCA)
Ecole Nationale Vétérinaire de Nantes (ENVN)
BP 50707, 44307, Nantes
France
Marie-Hélène Le Breton
2 Nestlé Research Center
Nestec Ltd., P.O. Box 44
CH-1000 Lausanne 26
Switzerland

Abstract

Recombinant bovine somatotropin (rbST), also called growth hormone, is a protein hormone used in dairy farming to enhance milk production. A method has been developed for the detection of rbST in milk by ESI(+)LC-MS/MS. This method allowed a detection limit of 20 pg of tryptic N-terminal peptide rbST in standard solution injected on-column and was successfully applied to extracts obtained from milk samples spiked with 50 ng/mL⁻¹ (2.3 pmol/mL⁻¹) rbST.

Introduction

Recombinant bovine somatotropin (rbST), also called growth hormone, is used in lactating cows to increase milk production. Different regulations

exist regarding its use, but the lack of confirmatory methods [1] for its detection makes it difficult to apply these regulations. It turns out to be an international issue in terms of animal doping and in food safety as well. Indeed, residues of rbST can be present in the milk of dairy animals treated with this hormone.

In order to detect residues of rbST in milk, the choice has been made to focus the analysis on the tryptic N-terminal peptide of the protein, specifically the difference between the endogenous and recombinant forms. The N-terminal amino acid alanine that is present in the endogenous form is replaced by a methionine in the recombinant one [2].

This application describes a method for the detection of rbST by ESI(+) LC-MS/MS. The method was successfully applied to extracts from milk samples spiked with rbST.

Experimental

Standards of Proteins and Peptides

Protein standards of rbST and recombinant equine somatotropin, reST (EquiGen-5), were obtained from the Harbor-UCLA Medical Center, National Hormone and Pituitary Program (Torrance, CA, USA) and Bresagen Limited (Thebarton, Australia), respectively.



Agilent Technologies

The peptides used as standards, with the following amino acid sequence MFPAMSLSGLFANAVLR (N-terminal tryptic rbST), MFPAMPLSSLFANAVLR (N-terminal tryptic reST), and AFPAMSLSGLFANAVLR (N-terminal tryptic bST) were synthesized from Millegen (Labege, France).

Instrumentation

The detail of the instrumentation used for the detection of the N-terminal peptides is described in the following tables.

LC

Instrument	Agilent 1200	
Column	Column Interchrom ModuloCart QS Uptisphere 3HDO 150 mm × 2 mm	
Mobile phase	A: Acetonitrile + 0.1% formic acid B: H ₂ O + 0.1% formic acid	
Flow rate	0.3 mL/min	
Injection volume	20 µL	
Gradient	Time (min)	%A
	0	10
	5	55
	10	60
	15	100
	17	10
	20	10

MS

Instrument	Agilent 6410 LC/MS Triple Quadrupole	
Ionization mode	ESI (+)	
Capillary	5000 V	
Nebulizer	55 psi	
Gas flow	13 L/min	
Gas temperature	300 °C	

Selected Reaction Monitoring (SRM) Method Parameters

In order to obtain a better specificity, the detection was performed in SRM mode. The transitions monitored are displayed in Table 1.

Table 1. SRM Method Parameters

Compound	RT	Charge	Transitions monitored	Collision energy (V)
Nterm rbST	8.33	z = 2	913.2 → 1047.7	30
			913.2 → 774.1	20
		z = 3	609.3 → 774	10
			609.3 → 643.5	20
Nterm reST	8.39	z = 2	933.2 → 1287.9	30
			933.2 → 794.1	20
Nterm bST	8.20	z = 2	883.2 → 1047.8	20
			883.2 → 774.1	20

Results and Discussion

Application of Triple Quadrupole MS-MS and Electrospray Ionization Mode Methodology

In this method, the choice has been made to use electrospray ionization in positive mode. Indeed, this ionization mode presented as a “soft” ionization technique is optimal for peptides. The ionization of the N-terminal peptide rbST leads to two main forms (z = 2 and z = 3).

This use of a triple quadrupole based methodology enabled very good sensitivity and selectivity and also a possible quantification of the monitored signals.

Separation of the Different Compounds

The detection method was developed for the detection of the tryptic N-terminal peptide of rbST and also for the tryptic N-terminal peptide of endogenous pituitary bovine somatotropin (bST) and reST as well. Due to the high homology in the amino acid sequence with rbST, reST was used as the internal standard.

The three compounds were separated chromatographically and analyzed utilizing the transitions described in Table 1. The chromatogram corresponding to the injection of 0.2 ng of N-terminal peptide bST, rbST, and reST is shown in Figure 1.




Figure 1. SRM ion chromatograms of standard solutions of tryptic N-terminal peptide of (a) reST, (b) rbST and (c) bST. The injection aliquot used corresponded to 0.2 ng on-column.

Even with an optimized gradient, due to their high homology in terms of sequence, the three compounds eluted with very similar retention times.

Linearity and Sensitivity of the Method

The method described allowed detection of the three peptides with very good sensitivity. A limit of detection of 20 pg injected on-column (~900 femtomole) was reached. Quantification was possible as shown by the good linearity of the calibration curves (Figure 2).




Figure 2. Calibration curve of tryptic N-terminal peptides rbST ($913.2 \rightarrow 774.1$), reST ($933.2 \rightarrow 794.1$) and bST ($883.2 \rightarrow 774.1$).

Results of Spiked Samples

The detection method was applied to extracts obtained from milk samples spiked with rbST. The purification procedure used is described in [3].

Figure 3 shows the chromatogram of a milk sample spiked with 50 ng.mL^{-1} of rbST, in accordance with guidelines for the identification of rbST according to 2002/657 criteria [4]. The chromatogram shows excellent peak shape, and above all, nearly null




Figure 3. SRM ion chromatograms obtained from milk samples. The different signals correspond to (a) blank milk, (b) the same milk spiked with 50 ng.mL^{-1} reST (internal standard), and (c) 50 ng.mL^{-1} rbST.

background noise, demonstrating the selectivity of the method. The intensity of the signal, although lower than the internal standard, is, however, significant, and shows a clear and distinct signal. The method clearly allows for unambiguous identification of rbST in milk.

Conclusions

The detection of rbST in milk was performed with detection by ESI(+) LC-MS/MS. The method showed very good sensitivity, specificity, and robustness. It was successfully applied to milk samples spiked with rbST at 50 ng.mL⁻¹, in accordance with criteria outlined by the 2002/657 Council Directive.

References

1. L. Bailly-Chouriberry, E. Chu-Van, G. Pinel, P. Garcia, M.-A. Popot, G. André-Fontaine, Y. Bonnaire, and B. Le Bizec, *Analyst*, 133 (2008) 270
2. G. Pinel, F. André, and B. Le Bizec (2004) *Journal of Agricultural and Food Chemistry*, 52, 407-414
3. M. H. Le Breton, S. Rochereau-Roulet, G. Pinel, L. Bailly-Chouriberry, G. Rychen, S. Jurjanz, T. Goldmann, and B. Le Bizec (2008) *Analytical and Bioanalytical Chemistry*, submitted
4. Commission Decision of 12 August 2002 implementing Council Directive 96/23/EC concerning the performance of analytical methods and the interpretation of results, Off. J. Eur. Commun. 2002/657/EC, 2002

For More Information

For more information on this application, you may also contact Agilent Technologies, Inc., Attention: Lea Bonnington, European Field Support Centre, Hewlett-Packard Strasse 8, 76337 Waldbronn, Germany.

For more information on our products and services, visit our Web site at www.agilent.com/chem.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2008

Printed in the USA
June 24, 2008
5989-8481EN



Agilent Technologies