

Determination of Polyacrylamides in Reclaimed Water

Application Note

Authors

Greg Saunders, Ben MacCreath Agilent Technologies, Inc.

Introduction

Polyacrylamides are often used as settling aids in water clarification processes. These polymers act as flocculants to help remove contaminants from the water stream. Residual polyacrylamides in the reclaimed water can affect subsequent processing operations and have an obvious environmental impact. A fast analytical method is therefore required to monitor traces of residual polyacrylamide in reclaimed water. Aqueous size exclusion chromatography (SEC) is a well documented technique for the characterization of polyacrylamide molecular weight distribution. In this particular application, the primary aim is to quantify the amount of residual polymer, whereas the molecular weight information is of little importance. A single pore size Agilent PL aquagel-OH column was selected to partially exclude the polymer peak and thereby give increased response, and to provide separation between the polymer and any lower molecular weight species, which may be present in the sample. PL aquagel-OH columns are ideal for this purpose because they combine high pore volume and high column efficiency (>35,000 plates/meter) for maximum resolution.



Conditions

Polyacrylamides can exhibit a variety of chemistries rendering the polymer nonionic, anionic or cationic in nature. In order to stabilize the polymer solution and to inhibit interaction between the polymer and the column packing material, 0.05M Na₂SO₄ at pH 3 was added to the aqueous eluent. By selection of the appropriate eluent, the full range of polymer types was evaluated.

Cal	umn	
60	umm	

Eluent:

Flow Rate:

Detection:

PL aquagel-OH 50 8 µm, 300 x 7.5 mm (p/n PL1149-6850) 0.1 mg/mL, 200 µL Sample Loading: 0.05 M Na, SO, at pH 3 1.0 mL/min

UV, 208 nm

Results and Discussion

Acrylamides exhibit strong UV absorbance at low wavelengths. However, the limiting factor in the determination of trace polyacrylamides is purely down to detector sensitivity. Figure 1 reveals the presence of nonionic, anionic and cationic polyacrylamides with good separation between the polymer and lower molecular weight contaminants.



Figure 1. Chromatogram shows the presence of nonionic, anionic and cationic polyacrylamides with good separation between the polymer and lower molecular weight contaminants

Conclusion

SEC and PL aguagel-OH columns successfully quantified polyacrylamides in a sample of reclaimed water. The 'neutral' surface and ability to operate across a wide range of eluent conditions equip PL aquagel-OH for the high performance analysis of analytes with neutral, ionic and hydrophobic moieties, singly or combined.

www.agilent.com/chem

This information is subject to change without notice. © Agilent Technologies, Inc. 2010 Published in UK, September 2, 2010 SI-01627



Agilent Technologies