# Agilent EcFWUPDATE

# **Operating Guide** 操作说明书

## Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

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#### Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

### **Technology Licenses**

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

#### **Restricted Rights Legend**

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data -Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

#### **Safety Notices**

### CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

### WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

Agilent EcFWUPDATE Operating Guide

# **EcFWUPDATE Operating**

Overview 4 Install the Software 4 Install EcFWUPDATE 4 Install the USB driver 5 Configure the USB port 7 Run the Program 9 User Interface Reference 12

Agilent Technologies

# **Overview**

The EcFWUPDATE tool allows a PC to update the firmware on an Agilent electrochemical instrument. Using the PC USB port, the software can load new firmware into the internal program memory of the instrument.

# **Install the Software**

This program is supported on the following operating systems:

- 32 bit Windows 7
- 32 bit Windows XP, Service pack 3

## Install EcFWUPDATE

- 1 Download the EcFWUPDATE.ZIP file from this Agilent website: http://www.chem.agilent.com/en-US/Support/Downloads/ Utilities/Pages/default.aspx.
- **2** Unzip the EcFWUPDATE.ZIP file to a temporary directory.
- **3** From the temporary directory location, double-click the EcFWUPDATE.EXE file and follow the instructions provided on the software setup wizard to complete the installation of the EcFWUPDATE program.

The USB driver must be installed before running this program.

### Install the USB driver

Install the provided USB driver on a PC to allow communication with the instrument. If this driver was previously installed for EcPRINT please skip this procedure.

- 1 Run **CP210xVCPInstaller.exe** which is located in the temporary directory used for unzipping the EcFWUPDATE.ZIP file.
- **2** With the instrument powered off, connect the USB cable (provided with the instrument) from the instrument port labeled USB to a USB port on the PC.
- **3** Power on the instrument.

4 Open the device manager on the PC to view the newly installed USB driver labeled Silicon Labs CP210x USB to UART Bridge (COMX), where X is the communication port number. This label indicates that the USB driver has been properly configured for the instrument.

The port number must be in the range of 1 to 16. If the port number shown is outside this range, double-click on the driver label to open the properties dialog, then click on the **Port Settings** to bring up the **Advanced Settings** dialog, and select a valid port number from the **COM Port Number** drop-down.

**5** Record this USB communication port number (12 in this example) for configuration of the EcFWUPDATE software.

## **Configure the USB port**

This procedure assigns the previously recorded USB port number to the EcFWUPDATE program.

1 Double-click the EcFWUPDATE icon on the PC desktop to start the program. Initially the PC screen displays **Connecting instrument**....

A **connecting error!** message displays because the port number is not yet configured.

2 Click Set Commport to display the Set Commport dialog.

- **3** Enter the previously recorded communication port number (12 in this example) and click **OK**.
- **4** Click **Retry** and the EcFWUPDATE software should now connect to the instrument and display the instrument type, and firmware version number.
- **5** A prompt displays if an update is required.

6 Click **Cancel** to exit or click **OK** to update the firmware now. See the next section on "Run the Program" and continue with step 2.

## **Run the Program**

Configure the USB port before running the software.

1 Double-click the EcFWUPDATE icon on the PC desktop to start the program.

When the instrument communicates with the EcFWUPDATE software, it displays the program user interface window with recent data readings in the data display area.

2 Click **Read BIN File** to select the latest firmware file.

The firmware files are named as follows where x.00 is the version number:

3200P\_ISP\_V*x.00*.Bin, for 3200P pH meter 3200C\_ISP\_V*x.00*.Bin, for 3200C Conductivity Meter 3200D\_ISP\_V*x.00*.Bin, for 3200D Dissolved Oxygen Meter

3200I\_ISP\_Vx.00.Bin, for 3200I Ion Meter

3200M\_ISP\_Vx.00.Bin, for 3200M Multi-Parameter Analyzer

#### **EcFWUPDATE Operating**

**3** Click **Update** to update the firmware.

When the firmware update is complete, the following message displays.

4 Click **Exit** to exit the program.

# **User Interface Reference**

	CP	Agilent EcFWUPDATE Control Panel	
Instrument type	3200M Multi-Para	ameter Analyzer	
File name of firmware —		BIN File Name:	
Update status		Actives.	
	Connect		
Operating menu	Read BIN File	Grifder	1
Operating steps	Vpdate	Step: Install USB driver. You any skip this step if you have installed.	
oporating oropo	Exit	Step5. Red wedste BIF file. Step4. Fress lödale butten to størt spåste.	

The software interface is composed of 5 parts:

- Instrument type
- File name of firmware (.BIN)
- Update status
- Operating menu
- Operating steps

#### Cancel

Click **Cancel** to cancel the connection.

#### Exit

Click **Exit** to exit the EcFIRMWARE program.

### Connect (2nd try)

Click **Retry** to attempt to connect the instrument.

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# Common Common Marcol Common Set Common Marcol Cancel Cancel Cancel Cancel

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	Agilent EcFWUPDATE Control Panel	
3200M Multi-Para	meter Analyzer	
:01:02	BIN File Name: G-VA2000_ISF_VI 01 Bin Status: Restyl	
	Updating is completed.	
Connect Read BIN File	Cuidas	
Vjrðate	Step: Install USB driver. You may skip this step if you have installed. Step: Em. IdfWTDATE.	
Exit	Step3. Reed update BIN file. Step4. Frees Update buttom to start update.	

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### www.agilent.com

© Agilent Technologies, Inc. 2012

Printed in China 04/2012 First Edition, April 2012

5973-1795