

Effect of Bonded Phase on Selectivity 10 Polypeptides Separated By Wide-Pore Reversed Phase

- Application Biochemical
- Biochemica

Robert Ricker



Highlights

- Improved separation is achieved by taking advantage of different bondedphase chemistries.
- Both long-chain and short-chain bonded phases are stable when utilizing StableBond sterically protected phases.

Conditions:

ZORBAX 300SB-C18, 300SB-C8, 300SB-C3, 300SB-CN (4.6 x 150 mm) (Agilent P/N: 883995-902, 883995-906, 883995-909, 883995-905) Mobile Phase: 25-70% B in 40 min, A: 0.1% TFA in H₂O, B: 0.09% TFA in 80% ACN/20% H₂O Injection 3µg each polypeptide, 1 mL/min, 60°C, Detect. UV (210 nm)



Robert Ricker is an application chemist based at Agilent Technologies, Wilmington, Delaware.

For more information on our products and services, visit our website at: www.agilent.com/chem

Copyright[©] 2002 Agilent Technologies, Inc. All Rights Reserved. Reproduction, adaptation or translation without prior written permission is prohibited, except as allowed under the copyright laws.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

Printed in the USA April 25, 2002 5988-6295EN



Agilent Technologies