


## Effect of Mobile-Phase Strength in Separations of Mono- and Disaccharides

Application
Food Analysis
Robert Ricker

Sugars are important component of many foodstuffs derived from plant and animal sources Sugars are also of great clinical importance. They have traditionally separated by normal- phase on  $\mathrm{NH}_2$  bonded phases. Effect of mobile-phase solvent strength in these separations is demonstrated for various mono- and disaccharides.



## ACN/H<sub>2</sub>O: 75/25

**Conditions:** 

ZORBAX NH<sub>2</sub> (4.6 x 250 mm) (Agilent P/N: 880952-708)

Mobile Phase: ACN: H<sub>2</sub>O, as indicated 1 mL/min, Detect. = Refractive Index

## Highlights

- Good resolution of various mono- and disaccharides is obtained using a ZORBAX NH2 column.
- Retention of sugars <u>increases</u> with increasing organic content in the mobile phase.



Robert Ricker is an application chemist based at Agilent Technologies, Wilmington, Delaware.

For more information on our products and services, visit our website at: www.agilent.com/chem

Copyright<sup>©</sup> 2002 Agilent Technologies, Inc. All Rights Reserved. Reproduction, adaptation or translation without prior written permission is prohibited, except as allowed under the copyright laws.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

Printed in the USA April 25, 2002 5988-6360EN

