



# **Agilent PLgel 10 $\mu$ m MIXED-B Columns**

## **Effective Analysis of Butyl Rubber with Hexane Solvent**

### **Technical Overview**

#### **Introduction**

Hexane is a good solvent for butyl rubber although it can be chromatographed using other solvents, such as THF. The solubility parameter of hexane is very low (7.3). However, it can be used successfully with PLgel columns, such as PLgel MIXED-B.




**Agilent Technologies**

PLgel 10  $\mu\text{m}$  MIXED-B columns are designed for high MW polymer analysis and demanding eluent conditions. The PLgel 10  $\mu\text{m}$  MIXED-B spans a wide range of molecular weights, up to 10 million, with a linear calibration curve. It is particularly useful for molecular weight distributions where slightly higher than average MWs are encountered. The 10  $\mu\text{m}$  particle size provides good resolution with relatively low pressures for enhanced lifetimes in demanding conditions.

### Conditions

Columns:	2 x PLgel 10 $\mu\text{m}$ MIXED-B, 300 x 7.5 mm (part number PL1110-6100)
Eluent:	Hexane
Flow Rate:	1.0 mL/min
Detection:	RI



*Figure 1. Analysis of butyl rubber using PLgel 10  $\mu\text{m}$  MIXED-B columns*

[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 22, 2010

SI-1783



**Agilent Technologies**