

# GC/MS

# **RTL-MSD** Pesticide Method and Database

# Efficient Screening for Pesticides and Endocrine Disrupters Using the 6890/ 5973 GC/MSD System

#### Introduction

Retention-time locking (RTL) is a new technique by means of which the GC can precisely and accurately reproduce compound retention times from one chromatographic analysis to another on the same gas chromatograph/mass spectrometer detector (GC/MSD) or on any similar GC/MSD system in any laboratory. RTL allows the chromatographer to develop and export methods that can be universally reproduced on similar systems. RTL is easy on the Agilent Technologies 6890 GC/ 5973 MSD system, and to complete the system, MSD Productivity ChemStation software provides RTL functionality as well as instrument control, data analysis, reporting, and mass spectral library matching.

Using RTL, a GC/MS method has been developed that allows any analyst to analyze for 567 of the most common pesticides and endocrine disrupters known to be in use worldwide. By using the RTL functionality with the 6890/ 5973 GC/MSD system, it has been possible to develop a new database containing both mass spectra and an index of compound retention times. When used with the compound screening tools included with the MSD Productivity ChemStation software, data files can be easily, accurately, and quickly screened for the presence of any of the 567 compounds listed in the database. In addition, software features make it possible to detect compounds that might ordinarily be masked by others.



Figure 1. Selecting the Results Screener option from the View menu in the Enhanced mode of data analysis.

# The GC/MSD RTL Method [RTLPEST.M]

The GC/MSD RTL method is a universal pesticide and endrocinedisrupter method that accommodates many different injection schemes-for example, cold on-column, pulsed split/splitless, large-volume injection with PTV and more. However, because the method is based on retention times, the column used must be an HP-5 MS column (30-m, 0.25-mm i.d., 0.250 um thickness—available as 190913-433). The HP-5MS columns are highly uniform in terms of specifications and performance. The 6890 GC oven parameters were chosen to provide the best resolution for the wide variety of chemicals that are potentially present in complex samples. The total pesticide run time for the HP-5 MS column is short, only 41.87 min. The pesticide method is locked on the compound chlorpyrifos-methyl at a retention time of 16.596 min using the automated RTL Setup feature included with the MSD Productivity ChemStation software. For more information on how to acquire RTL data, see Publication 5968-3433E.

#### Using The RTL MS Pesticide Library and Screener Software

In addition to the 6890/5973 GC/MSD system, the RTL-MSD Pesticide Library and RTL-GC/MSD acquisition screening method are required. The RTL-MSD Pesticide Library (RTLPEST.L) and acquisition method (RTLPEST.M) are ordered and installed separately (product number G1049A). Installation instructions for the RTL-MSD Pesticide Library can be found on the MS Library CD ROM. The pesticide acquisition method is loaded by using the MSD Productivity ChemStation software. For additional information regarding the use of the MSD Productivity ChemStation software. see Publication 5968-3433E and On-Line Help.

After the method is locked on chlorpyrifos-methyl and a sample acquired by using the locked RTL pesticide method, the RTL Results Screener option is selected from the View menu in the Enhanced mode of data analysis to invoke the screening program for the loaded data file. See Figure 1.

## How Does It Work?

The Results Screener software automatically scans through the sample data file, comparing the acquired MS data with those of the 567 compounds in the RTL pesticide library. The data are compared with the target ion fragment information in the library, as well as with the expected retention time of each library compound.

The upper left window in Figure 2 displays the extracted ion chromatogram (EIC) for the sample based on the ions for the indicated compound. Compound matches ("hits") that are confirmed with high confidence are indicated by an "X" in the Quick Screener compound results list (Figure 2a) Possible and questionable hits are indicated by a question mark ("?") and suggest that further review of the data is appropriate. Using the Quick Screener tool, the analyst can automatically scroll through every compound in the library by pressing Start in the Quick Screener panel or can manually scroll to individual compounds of interest for further review.

Enabling the Display Reference Spectra function allows the analyst to simultaneously view the acquired mass spectrum versus the library's reference spectrum. See Figure 2b. The window located at the lower right of Figure 2 contains the compound name, retention time,



Figure 2. Results for GC/MS data acquired by the locked method, RTLPEST.M, as presented by the Results Screener view. The top window includes the extracted ion chromatogram for an unknown compound tentatively identified as trifluralin.



deviation from the expected locked retention time, and the compound's target, actual, and expected ion ratios. See Figure 2c.

> Figure 2a. Expanded view of right-hand window in Figure 2. Compound confirmation is denoted by an "X," and uncertain identification is denoted by a question mark ("?"). The "?" notifies the operator that further evaluation of a compound may be in order.

₩indo			- 🗆 ×
Abundance	Scan 1599 (11.631 min): RTLDEMO.D		
1		30	6
1		264	
500000		1	
	43	290	
1 1	57 75 89 105 126 145 160 172 187 206 218 232	276	318 <sup>335</sup>
m(7>		<u>4 4 € , ···, 1 7 7</u> 250 300	<u> </u>
Abundance	#103: Trifluralin		
	(c)	30	6
	· v	264	
1 1		204	
5000			
1 1			
	172 186 206 <sup>24</sup>	48 290	318 335
L			┡┰┸╄╾┲┸┲╼
m/z≻	50 100 150 200 2	:50 300	

Figure 2b. Expanded view of lower left-hand window in Figure 2. Example of actual (upper) and library (lower) spectral comparison.

A[7] TIC: RTLI	DEMO.D		<u>- 0 ×</u>
(103) T	rifluralin		
11.630n	nin (-0.00	17)	
respons	se 19795	103	
lon	Exp%	Act%	
306.00	100	100	
264.00	70.20	67.41	
307.00	13.20	13.04	
290.00	12.90	12.28	
$\longleftrightarrow$			

B

Figure 2c. Expanded view of lower righthand window in Figure 2. Compound information provided includes compound name, the retention time and deviation from the expected locked retention time, and the compound's target ions and their actual and expected ratios.

As mentioned previously, positive hits are denoted by an "X" in the Quick Screener list, and possible or questionable identifications are denoted by a question mark ("?"). For example, although benfluralin co-elutes with trifluralin and contains common ions, the Screener software correctly indicates the failure of the qualifying ions. Upon closer inspection of the Quick Screener list shown in Figure 3a, the spectral comparison shown in Figure 3b, and qualifier ion table shown in Figure 3c, one can easily discern that the qualifier ions are not present. This information provides the user with high confidence that there is no detectable benfluralin in the sample. See Figure 3.



Figure 3. Evaluation of the possible presence of benfluralin.





🍋 Quio	:kScreei	ner	×
<u>S</u> tart	E <u>x</u> it	QDel	Con <u>f</u> igure
96	2	2,3,5-Trin	nethacart 🔺
97	2	2,3,5-Trin	nethylphe
98 :	x 11.27	5 Ethalf	uralin
99	0	)icrotoph	los
100	1	hiofano	к <u>—</u>
101	? 11.62	8 Bromos	kynil
102	? 11.62	3 Bendia	carb
103 :	x 11.63	O Triflura	alin
104	? 11.63	0 Benflu	ralin
105 :	x 11.62	8 Mon 🔉	rotophos
106	9	ulfotep	
107	1	ebutam	
108	0	)esbroma	o-bromobi
109	F	romecar	b
110	0	)i-allate	
111	F	Phorate	
112	1	riclopyr	methyl es
113	E	BHC alph	a isomer
114	h	lethyl-1-	naphthal
115		)i-allate	▼

Figure 3a. The Quick Screener list informs the user of a questionable ("?") hit for benfluralin.

🎇 [7] TIC: RTL	DEMO.D	_ <b>_</b> ×		
(104) E	Benfluralin			
11.630min (-0.095)				
respon	se 53151			
lon	Exp%	Act%		
292.00	100	100		
264.00	17.00	24881.36#		
276.00	12.20	250.09#		
293.00	11.80	18.68		

Figure 3c. The qualifying ion list clearly identifies the two ions that have failed to meet the expected qualifier ion ratios as indicated by the "#" sign.

#### **Co-Eluting Compounds**

Another utility of the Results Screener is discerning compounds that may be obscured by co-eluting chromatographic peaks. In Figure 4a, the Screener software suggests a hit for the compound chlorpyrifosmethyl, but there is an unexplained ion fragment at m/z = 263 and a skewing of the ion ratios for the m/z = 125 ion. By selecting a compound with a close retention time and similar ions as indicated by the "?" in the Quick Screener list, the m/z = 263 ion fragment and skewing of the m/z = 125 ion can be explained as co-elution with methyl parathion. See Figure 4b.



Figure 4a. Detection of possible compound co-elution is assisted by information presented in the Result Screener software. The Result Screener indicates a hit for the compound chlorpyrifos-methyl. However, closer examination of the spectral comparison and qualifier ion ratios indicate the possible presence of a co-eluting compound.

### **RTL-MSD** Pesticide Method and Database



Figure 4b. Selecting the compound with the closest retention time in the Results Screener window, the spectral match and qualifying ion ratios identify the co-eluting peak to be methyl parathion.

#### **Identifying Isomers**

Many pesticide formulations consist of isomer mixtures. The classic example is the hexachlorocyclohexanes (BHCs). Because these are isomers, they have identical mass spectra and can only be distinguished by retention times. See Figures 5a and 5b.

The Screening software also automatically creates a report of possible and confirmed hits. See Figure 6. The report presents the experimental retention time, deviation from the expected locked retention time, target ions, response, and qualifier ions that fall out of range.



Figure 5a (above) and 5b (next page). The Screener software correctly identifies the two isomers using their respective retention times: the BHC beta isomer at 13.203 min and lindane (gamma BHC isomer) at 13.446 min.

### **RTL-MSD** Pesticide Method and Database



Figure 5b. See caption under Figure 5a on previous page for description.



MultiVu - [C:\HPCHEM\MSDEMO\RTLDEMO.D\SCRNTEMP.TXT]					
蹭 <u>Fi</u> le <u>E</u> dit <u>Search</u> <u>Window</u>	_ 문 ×				
Screen Report	_				
	_				
Data File : C:\HPCHEM\MSDEMO\RTLDEMO.D Vial: 4					
Sample TEST MIXIA Det COMS Inc.	-				
Misc 26 pests 10 ng/ul 17 2 psi const p Multiplr: 1 00	2				
Sample Amount: 0.00					
MS Integration Params: events.e					
Screen File: screen RTLPEST.RES					
Target Qualifiers	VOD				
Compound Status ExpRI Delta m/z Resp. Out of Range A	ACR				
21 Dichlorvos x 5.831 -0.017 109 28872469 0	.88				
35 Mevinphos x 7.595 -0.013 127 28860822 1	.00				
36 Butylate ? 7.607 -0.023 57 27318 146,156,174 0	.00				
3/3, 4-Dicnioroaniline ( 7.664 +0.122 161 -47520 163 U	.09				
40 vermulate X /.000 -0.021 120 10043988 U 94 (blordimeform 2 11 201 +0.073 196 111277 191 117 152 0	.40				
95 Dibron (naled) x 11.221 -0.020 109 9078031 0	94				
98 Ethalfluralin x 11.283 -0.008 276 4696405 0	.84				
101 Bromoxynil ? 11.550 +0.078 277 42935 279,88 0	.00				
102 Bendiocarb ? 11.540 +0.083 151 157380 126,166,57 0	.00				
103 Trifluralin x 11.637 -0.007 306 19795103 0	.96				
104 Bentfuralin ? 11.725 -0.095 292 53151 264,276 0	.10				
130 Prometon v 12 990 40 011 210 13930019 0	95				
132 Sven 2 13 052 +0 115 187 19798 189 219 0					
133 Chlorbufam ? 13.063 -0.063 153 1912228 127.53 0	.08				
134 Monolinuron ? 13.098 +0.105 61 165473 153,127 0	.02				
135 Dimethipin ? 13.098 +0.067 54 103136 39,53 0	.06				
136 Atrazine x 13.159 +0.005 200 15570363 0	.91				
137 BHC beta isomer x 13.200 +0.003 219 7163009 0	.83				
141 Propazine x 13.353 -0.015 214 108048 0	./5				
179 Chorothalonil x 14 784 -0.000 266 2682153 0	66				
216 Fuberidazole ? 16.549 +0.053 155 686802 156 0	.00				
217 Methyl parathion ? 16.594 +0.010 263 17176733 125 0	. 41				
218 Chlorpyrifos Methyl x 16.593 +0.016 286 38888465 0	.52				
221 Terbucarb ? 16.686 -0.086 205 160687 220 0	.00				
225 Carbaryi 7 16.806 -0.021 144 54379 115,116,145 0	.00				
228 Ferror 2 16 801 -0.024 255 445501 257,255 U	.00				
260 Bromacil x 18.350 -0.03 205 15483082 0	95				
263 Di-n-butylphthalate ? 18,415 -0.015 149 551586 205 0	.00				
269 Oxamyl ? 18.728 +0.062 44 50069 98,72,69 0	.01				
270 Malathion x 18.800 -0.011 173 12954421 0	. 41				
273 Fenthion ? 19.120 +0.106 278 377014 125,109,169 0	.00				
275 4, 4 - Dichlorobenzophenone ? 19.201 +0.021 139 47474 111, 141, 250 0	.00				
278 Chlorpurifoe v 19 234 _0 009 197 10194202 0	50				
280 Parathion 2 19 275 -0.058 291 9017 109 97 139 0					
	•				
3	F				

Figure 6. Example of pesticide screener report.

*Harry Prest*, Ph.D., is an Application Chemist; *Phil Wylie* is an Application Chemist; *Ken Weiner* is a Software R & D Engineer; and *Doug Agnew* is a Product Manager at Agilent Technologies Company.

### Conclusions

The highly accurate and reproducible pressure and temperature control of the 6890 GC made it possible to create a mass spectral library and a retention time-locked database of 567 chemicals of concern. The MSD Productivity ChemStation software provides a Screener tool that automatically presents concise and complete information for identifying and confirming the presence of target compounds utilizing both mass spectral and retention-time information. When used together, analyst can rapidly survey for the presence of any of the 567 compounds with both mass spectral and chromatographic confidence in compound identity.

Windows NT<sup>®</sup> is a U.S. registered trademark of Microsoft Corporation.

Windows<sup>®</sup> is a U.S. registered trademark of Microsoft Corporation.

Agilent Technologies shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this material.

Information, descriptions and specifications in this publication are subject to change without notice.

Copyright © 1999 Agilent Technologies Company All rights reserved. Reproduction and adaptation is prohibited.

Printed in the U.S.A. December 1999 (23) 5968-4884E