



# Gas Clean Filters

For chromatography, ICP, and TOC

**Delivering clean gases for  
accurate analyses**

The Measure of Confidence



**Agilent Technologies**



# Enhanced gas quality for maximum productivity

The Agilent Gas Clean Filter System delivers clean gases, reducing the risks of column damage, sensitivity loss, and instrument downtime.

Contaminants in gases can significantly affect your analysis. Oxygen, hydrocarbons, and moisture can cause loss of sensitivity and accuracy of the GC and damage your column. Impurities activate glass wool in liners and accelerate septa degradation, causing high background signals and ghost peaks and leading to time-consuming troubleshooting. Oxygen in the supply gas for ICP-OES or ICP-MS can cause plasma shut-down and loss of sensitivity. Carbon dioxide in supply gas for TOC analyzers causes elevated baselines and loss of sensitivity and accuracy.

Supply gases can pick up contaminants from every part of the gas line. You therefore need a Gas Clean Filter System even if your supply gas is of the highest quality – it is not economical to buy expensive, high purity gases if their quality is downgraded by impurities in the gas line.



Examples of the Gas Clean Filter range, showing the 4-position connecting unit for ease of operation.

## Agilent introduces a range of Gas Clean Filters and related components:

GC/MS Filter – provides faster stabilization times to reduce gas consumption (see Figure 1)

Moisture Filter – delivers faster stabilization times for increased productivity in GC


Process Moisture Filter and connecting unit – designed for use with acetylene in process GC

Carbon Dioxide Filter – eliminating CO<sub>2</sub> from supply gas, to minimize loss of sensitivity and improve accuracy in TOC

High flow connection unit – for applications up to 20 L/min, including ICP-OES, ICP-MS, and LC/MS

# Fast stabilization and replacement


## Fast stabilization



**Figure 1.** Shows the fast stabilization rate (the  $N_2$  mass measured by mass spectrometry) of a GC/MS after replacement of the filter.

## Reduced impurities for improved analysis

Inserting a Gas Clean Filter System in the gas line immediately before the instrument inlet greatly reduces the level of impurities, thus improving trace analysis (see Figure 2). Contaminants entering your GC column will also be reduced, which is critical for high temperature analysis and essential for longer column lifetime.



**Figure 2.** Shows the difference in bleed levels of two GC columns due to moisture exposure with and without a filter when running a temperature program (50 °C to 350 °C, 20 °C/min). When no filter is used, an extreme rise in the bleed profile is clearly visible due to moisture in the carrier gas. By using Gas Clean Filters, a normal bleed profile is achieved with the removal of all moisture in the carrier gas.

## Fast, leak-free filter replacement for superb performance

Every Gas Clean Filter features a unique connection design that permits easy filter replacement. The connecting unit shown in Figure 3 allows the instrument to remain under pressure during filter changeover and prevents air from entering the system. Changing a Gas Clean Filter is easy, requiring no tools or gas shut-off, unlike in-line filters.

The Technischer Überwachungs-Verein (TÜV\*) certificate\*\* assures safe use up to 15 bar (217 psi). A TÜV certificate at this level is unique to our design and guarantees the highest quality of Gas Clean Filters.

\* TÜV NORD Group has expertise in technical safety, environmental protection, and the conformity assessment of management systems and products – nationally and internationally. \*\* Certificate number 973/980538 hel-CHROMP1.doc.

## Fast replacement



**Figure 3.** The unique connection design allows fast, leak-free replacement of the Gas Clean Filter.

# Highest performance

## Gas Clean Filter indicators are more sensitive

Replacing the filters when they have reached absorption capacity ensures maximum protection of your GC columns and analytical hardware. The sensitive indicators change color, alerting you that the filter needs to be replaced.


In Figure 4, the open-topped glass tube on the left contains the Gas Clean moisture indicator (green) packed below another manufacturer's indicator (yellow). On the right is the same tube after exposure to air at room temperature with average moisture content. As the moisture penetrated the tube downwards, it passed through the other manufacturer's indicator before it reached the Gas Clean indicator at the bottom. The tube on the right shows that the Gas Clean indicator changed color from green to pale brown, revealing the presence of moisture. However, the other indicator failed to change color, confirming that the Gas Clean moisture indicator is a great deal more sensitive.

With the Gas Clean Filter System, you can be completely confident that your gas supply contains less than 0.1 ppm of moisture. The system will not suffer from hydrolization damage, saving you expense and needless downtime in replacing degraded components. Our moisture filters offer unsurpassed protection against moisture contamination for your GC system.


## Economical, with immediate payback

A Gas Clean Filter System allows you to use 99.996% (4.6) pure helium instead of the more expensive 99.999% (5.0) or 99.9999% (6.0) grade, while still yielding high quality analytical results. Figure 5 compares the costs of carrier gas with the use of helium 4.6 and 6.0. The expected cost savings are 30%.

The Gas Clean Filter System, in combination with 4.6 grade helium, delivers at least the same quality gas to your GC and GC/MS as using 6.0 grade helium with respect to oxygen and water. This presents considerable savings over the use of higher quality helium.



**Figure 4.** Shows how Gas Clean filter indicators are extremely sensitive to moisture. After exposure, the Gas Clean indicator in tube 2 changes color before another manufacturer's indicator has barely begun to react.



**Figure 5.** Shows the cost savings in helium provided by the Gas Clean Filter System when using 4.6 grade rather than 6.0 grade helium.

# The right filter for every gas

## A full range of filters

Six different filters are available. Each fast-stabilizing absorbent is packed in a transparent and virtually unbreakable, heavy-walled polycarbonate housing. Our full range offers you:

- A GC/MS Filter that delivers faster stabilization times for lower gas consumption and faster set-up time. The single filter unit with a combination of filters removes oxygen, moisture, and hydrocarbons from carrier gas for MS applications.
- A Moisture Filter that also delivers faster stabilization times for increased productivity in GC. Water removal prevents needless hydrolization damage to the stationary phase, fused silica surface, glass wool, or septum.
- An Oxygen Filter that prevents oxidation of the liquid stationary phase, septum, or glass wool in GC.
- A Process Moisture Filter, process connecting units, and male connectors, free from any brass and copper, for safe use with acetylene in process GC.
- A Charcoal Filter that removes organic compounds and ensures correct performance of FID detectors in GC.
- A new Carbon Dioxide Filter that eliminates CO<sub>2</sub> from supply gas. When combined with a Moisture Filter, it is ideal for zero-air gas generators.

Improvements to the absorbance materials in the GC/MS and Moisture Filters mean that they stabilize more quickly to reduce downtime and improve productivity. We have also introduced a high flow connection unit that handles flow rates up to 20 L/min for collision gas applications, supply gas for ICP and ICP-MS, or any application where high flows are required, expanding the range of applications you can run. For operations requiring flows above 10 L/min, you can save money by using cheaper gas and by eliminating contaminants. Each Gas Clean Filter is compatible with both the high flow and standard connection units.

All GC applications benefit from the use of Gas Clean Filters, whether your GC is connected to an FID or MS or to any other type of detector, such as flame photometric, thermal conductivity, electron capture, nitrogen-phosphorous, or thermoionic.

## Modular design for easy mounting

The Gas Clean Filter System consists of two key parts: the connecting unit and the filters. The connecting unit has inlet and outlet connectors for the gas lines, and the system can be wall-mounted or fixed to a bench. Connecting units come with a 1, 2, or 4 filter capacity and are available for 1/4 in. or 1/8 in. gas lines.

## Filter selection guide

| Technique                 | Filters  | Benefit                                   |
|---------------------------|--|---|
| GC/MS                     | GC/MS Filter   | Higher data accuracy and less maintenance |
| GC column                 | Moisture Filter and Oxygen Filter                            | Longer lifetime                           |
| ECD detector (GC)         | Moisture Filter and Oxygen Filter                            | Greater sensitivity                       |
| TCD detector (GC)         | Moisture Filter and Oxygen Filter                            | Greater sensitivity and less maintenance  |
| Process GC                | Process Moisture Filter                                      | Long-term stability                       |
| FID detector (GC)         | Two Charcoal Filters (for air and hydrogen)                  | Greater sensitivity                       |
| PID detector (GC)         | Oxygen Filter and Charcoal Filter                            | Greater sensitivity                       |
| PFPD or FDP detector (GC) | Charcoal Filter, CO <sub>2</sub> Filter, and Moisture Filter | Greater sensitivity                       |
| TSD or NPD detector (GC)  | Charcoal Filter, CO <sub>2</sub> Filter, and Moisture Filter | Greater sensitivity                       |
| Total organic carbon      | CO <sub>2</sub> Filter and Moisture Filter                   | Greater sensitivity                       |
| Zero-air generator        | CO <sub>2</sub> Filter and Moisture Filter                   | Cleaner gas                               |
| ICP-OES, ICP-MS           | High Flow Connection Unit with two Oxygen Filters            | Greater sensitivity                       |

Contact your local Agilent representative for the filter sets applicable to your GC configuration.

## Gas Clean Filter technical specifications

|  | Oxygen Filter  | Moisture Filter/<br>Process Moisture<br>Filter                      | Charcoal Filter                            | GC/MS Filter  | CO <sub>2</sub> Filter  |
|--|--|---|--|---|---|
| Function   | Removes oxygen as well as traces of sulfur and chlorine compounds from carrier gas | Removes water, oil, and other foreign material from the carrier gas | Removes organic compounds from gas streams | Single combination filter; removes water, oxygen, and organic compounds | Removes CO <sub>2</sub> from gas stream; to be used in combination with Moisture Filter |
| Indicator color change                               | From green to gray   | From green to pale brown  | No indicator                               | Oxygen:<br>from green to gray<br>Moisture:<br>from green to pale brown  | From white to violet  |
| Capacity   | 150 mL oxygen  | 7.2 g water   | Approximately 7 g, depending on impurities | 100 mL oxygen, 1 g water, organics depending on impurities              | 9 g CO <sub>2</sub>   |
| Outlet concentration at operating flow of 1-10 L/min | <50 ppb  | <0.1 ppm  | <0.1 ppm                                   | Oxygen <50 ppb<br>Moisture <0.1 ppm<br>Organics <0.1 ppm                | <1 ppm  |

## Agilent CrossLab GC supplies

### Agilent CrossLab supplies are a convenient, practical, and productive choice.

Now you can simplify your ordering process and increase your overall productivity with Agilent CrossLab supplies. Agilent's new CrossLab portfolio provides a wide array of high-quality supplies designed to fit most major brands of GC systems. CrossLab lets you consolidate purchasing and reduce your administrative burden, without sacrificing quality, reliability, or reproducibility.

In choosing Agilent CrossLab supplies, you can pick up the phone, dial a single number, and be confident that you're ordering the right supplies for your instruments, regardless of manufacturer.

Each product is packaged so that it's easy to use, easy to store, and easy to find. And thanks to Agilent's global delivery infrastructure you can count on getting the parts you need when you need them.

At the same time, you have the assurance of knowing that your CrossLab supplies are manufactured with the same commitment to quality that Agilent brings to everything we make. You'll benefit from our 40 years of chromatography expertise and our ongoing innovation.

### A growing portfolio that works seamlessly with your instruments.

The growing Agilent CrossLab portfolio is manufactured to perform seamlessly with GCs from Bruker, Varian (now Bruker products), PerkinElmer, Shimadzu, and Thermo Scientific. CrossLab products include Agilent's innovative Ultra Inert inlet liners for active compounds, plasma-treated non-stick liner O-rings and inlet septa, autosampler syringes, autosampler vials and closures, capillary column ferrules and column nuts, and more.

New products critical to instrument performance and productivity are continuously being added to the CrossLab portfolio, so you are almost certain to find exactly what you need.



For more information, go to [www.agilent.com/chem/CrossLab](http://www.agilent.com/chem/CrossLab)



# Ordering information

| Gas Clean Filter Kits  | Part number |
|--|-------------|
| Agilent Gas Clean Filter kit (connecting unit for four filters, incl. four filters [1/4 in. tube*])                        | CP7995      |
| Agilent Gas Clean Filter kit (connecting unit for four filters, incl. four filters [1/8 in. tube])                         | CP736530    |
| Agilent Gas Clean GC/MS Filter kit (incl. one connecting unit 1/8 in. and two GC/MS Filters)                               | CP17976     |
| Agilent Gas Clean GC/MS Filter kit (incl. one connecting unit 1/4 in. and two GC/MS filters)                               | CP17977     |
| Agilent Gas Clean GC/MS Filter installation kit (incl. CP17976 1 m copper tubing, two nuts, and two ferrules 1/8 in.)      | CP17978     |
| Agilent Gas Clean CO <sub>2</sub> kit, 1/4 in. (incl. 2 position connecting unit and CO <sub>2</sub> and Moisture Filters) | CP17982     |
| Agilent Gas Clean CO <sub>2</sub> kit, 1/8 in. (incl. 2 position connecting unit and CO <sub>2</sub> and Moisture Filters) | CP17983     |
| TCD Filter kit (with Oxygen and Moisture Filters)  | CP738408    |

| Connecting unit for:                       | Part number |
|--|-------------|
| 1 Filter (1/4 in. tube)                    | CP7980      |
| 1 Filter (1/8 in. tube)                    | CP7988      |
| 2 Filters (1/4 in. tube)                   | CP738406    |
| 2 Filters (1/8 in. tube)                   | CP738407    |
| 4 Filters (1/4 in. tube)                   | CP7989      |
| 4 Filters (1/8 in. tube)                   | CP736520    |
| High flow connecting unit, 1/4 in.         | CP17984     |
| High flow connecting unit, 1/8 in.         | CP17985     |
| Agilent Gas Clean Process Moisture Filter: |             |
| 1 Filter (SS, 1/4 in. tube)                | CP7980P4    |
| 1 Filter (SS, 1/8 in. tube)                | CP7988P8    |
| 1 Filter (SS, 3 mm tube)                   | CP7988P3    |
| 1 Filter (SS, 6 mm tube)                   | CP7980P6    |

| Replacement Gas Clean Filters             | Part number |
|---|-------------|
| Agilent Gas Clean CO <sub>2</sub> Filter  | CP17969     |
| Agilent Gas Clean Oxygen Filter           | CP17970     |
| Agilent Gas Clean Moisture Filter         | CP17971     |
| Agilent Gas Clean Process Moisture Filter | CP17971P    |
| Agilent Gas Clean Charcoal Filter         | CP17972     |
| Agilent Gas Clean GC/MS Filter            | CP17973     |

| Accessories and Fittings  | Part number |
|---|-------------|
| Wall mounting bracket for connecting unit (for CP7980 and CP7988) | CP7981      |
| Upper part filter connecting unit                                 | CP7978      |
| Flush head for connecting unit                                    | CP7987      |
| Male connector 1/4 in. with dust filter                           | CP7986      |
| Male connector 1/8 in. with dust filter                           | CP82117     |
| Viton O-rings (two sets)  | CP7983      |
| Male connector for Gas Clean Process Moisture Filter:             |             |
| Male connector SS 1/4 in. with dust filter                        | CP7986SS    |
| Male connector SS 1/8 in. with dust filter                        | CP82117SS   |
| Male connector SS 3 mm with dust filter                           | CP82117SS3  |
| Male connector SS 6 mm with dust filter                           | CP7986SS6   |

\* For 1/8 in. tube, use reducer 1/8 in. x 1/4 in., part. no. CP4392



## A comprehensive portfolio from the leader in gas chromatography


Agilent provides the broadest selection of GC and GC/MS systems, support, and supplies in the industry. So whether you need flexible, reliable hardware and software for complex research; simple, robust systems for routine production environments; or fast, rugged portable solutions for real-time measurements in the plant or in the field, we have a GC or GC/MS to meet your analytical and business challenges.



*Agilent 7890A GC System –  
for the next level of GC  
and GC/MS performance.*

## Agilent J&W GC Columns – for utmost confidence in your results

Put 40 years of Agilent quality and innovation behind your every separation. Agilent J&W GC Columns offer the lowest bleed levels, the best inertness for acids/bases/mixed functional compounds, and the tightest column-to-column reproducibility. So when you put industry-leading Agilent J&W GC Columns to work in your lab, you can have the utmost confidence in your column, and in every separation.



### Learn more

[www.agilent.com/chem](http://www.agilent.com/chem)

### Buy online

[www.agilent.com/chem/store](http://www.agilent.com/chem/store)

### Find an Agilent office or authorized distributor

[www.agilent.com/chem/contactus](http://www.agilent.com/chem/contactus)

### U.S. and Canada

1-800-227-9770, [agilent\\_inquiries@agilent.com](mailto:agilent_inquiries@agilent.com)

### Europe

[info\\_agilent@agilent.com](mailto:info_agilent@agilent.com)

### Asia Pacific

[inquiry\\_lsca@agilent.com](mailto:inquiry_lsca@agilent.com)

### India

[india\\_lsca\\_marketing@agilent.com](mailto:india_lsca_marketing@agilent.com)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice

© Agilent Technologies, Inc., 2011  
Published in the USA, June 30, 2011  
Publication Number 5990-8243EN



**Agilent Technologies**