



Agilent Genomic Workbench Lite Edition 6.0

Data Viewing

User Guide



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2010

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G3800-90014

Edition

Revision A, May 2010

Printed in USA

Agilent Technologies, Inc.
3501 Stevens Creek Blvd.
Santa Clara, CA 95051

Trademarks

Microsoft® is a registered trademark of Microsoft Corporation in the United States and other countries.

Adobe®, Acrobat® and Reader® are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries.

Software Revision

This guide is valid for 6.0 and later revisions of the Agilent Genomic Workbench software, until superseded.

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data - Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

In This Guide...

This guide describes how to import, organize, manage, export and display data and other content (experiments, gene lists, tracks) within Agilent Genomic Workbench Lite Edition 6.0. It is targeted for users who have no DNA Analytics application license(s). If you do have a DNA Analytics license and intend to analyze your data, see the corresponding *User Guide*.

1 Getting Started

This chapter gives an overview of the capabilities you have in Agilent Genomic Workbench Lite Edition without a license, and describes the parts of the Agilent Genomic Workbench main window that you use to import, organize, manage, export and display array data and other content.

2 Importing, Managing, and Exporting Data and Other Content

This chapter describes how to import, organize, manage, and export data and other content within the user interface of Agilent Genomic Workbench.

3 Displaying Data and Other Content

This chapter shows you how to display log ratio data from imported feature extraction data files, as well as gene list and track content, in the Genomic Viewer. It also gives you instructions on how to modify the display to see the data and content the way you prefer.

4 Data Viewing Reference

This chapter describes the tab commands, shortcut menus, and dialog boxes that can appear.

Contents

1 Getting Started 11

Using Agilent Genomic Workbench Lite Edition on a Mac 13

Using Main Window Components to Display Data 14

What are the main window components? 14

What can you do with the main components to display data? 16

Switching applications 18

Using Tabs and Command Ribbons 19

Tabs 19

Commands 21

Using the Navigator to Search for Data 22

To search the Navigator 24

Using the Genomic Viewer to Display Data 26

What is the Genomic Viewer? 26

To change the size of and detach panes from the Agilent Genomic Workbench main window 28

To maximize and reattach panes to the Agilent Genomic Workbench main window 29

General Instructions for Displaying Microarray Data 30

Quick-start Instructions for Displaying Microarray Data 31

Getting Help 39

To get help within Agilent Genomic Workbench 39

To contact Agilent Technical Support 40

To learn about Agilent products and services 40

2 Importing, Managing, and Exporting Data and Other Content 41

Importing Files 42

To select a different location for data files 43

To import Agilent GEML design files	43
To import Axon design files	44
To import Agilent FE or Axon data files	45
To import a UDF file	46
To import a genome build	50
To import tracks	51
To import array attributes	51
To import an experiment file	52
To import filters	53
Working with Experiments to Organize Imported Data	54
To display the array designs and data in the program	54
To create a new experiment	55
To add arrays to an experiment	57
To change the order of arrays in an experiment	58
To change the display names for arrays in an experiment	59
To rename an array in an experiment	59
To remove arrays from an experiment	60
To display or edit the attribute values of a specific array	60
Managing Content (Data, Gene Lists, Tracks)	62
To display a list of the content stored in the program	62
To find specific content items in the Navigator	63
To display the properties of a specific design	64
To update probe annotation in design files	64
To rename an array in the Data pane	65
To remove data or design files from the program	65
To create a gene list	66
To display the genes in a gene list	67
To rename a gene list	67
To delete gene list(s)	67
To create a track (CGH only)	68
To display the details of a track	68
To rename a track	69


To delete tracks	69
Exporting and Saving Content	70
To export array attributes	70
To export experiments	71
To export filters	72
To export tracks	72
To copy what you see in the main window	73
To copy the list of array colors for an experiment	73
To save data and design information from an experiment	74
3 Displaying Data and Other Content	75
Selecting an Experiment for Displaying Data	76
To select an experiment	76
To select or deselect arrays in the experiment	77
To change the display color of an array	78
Displaying Array Data	80
To display the scatter plots	80
To show or hide data in scatter plots	80
To customize scatter plot ranges and colors	81
To change scatter plot appearance	82
To print the scatter plot	83
To create custom scales for Views	83
To locate and display data within the Views	84
To smooth and plot CGH log ratio data	86
To produce an echo example plot (CGH only)	87
To produce a moving average example plot (CGH only)	87
Displaying Content (Gene Lists/Tracks)	89
To show gene lists in Gene View	89
To change the appearance of genes in Gene View	89
To show tracks in Gene View	90
To change the appearance of tracks	91
To display tracks in UCSC Browser	92

To change the graphical display to a different genome build	93
Searching for Probe and Gene Information	95
To search Tab View for specific probe information	95
To search Agilent eArray for probe information	97
To search the Web for information on probes in Tab View	98
To create a custom Web search link	98
To update or delete a custom Web search link	99
4 Data Viewing Reference	101
Agilent Genomic Workbench Main Window	102
Command Ribbons	103
Home command ribbon	103
View Command Ribbon	106
Tool command ribbon	107
Help command ribbon	109
Switch Application Menu	111
Search pane	112
Navigator Pane	114
Data pane – icons, special text, and buttons	115
Data pane – actions and shortcut menus	116
Experiment pane – icons, special text, and buttons	117
Experiment pane – actions and shortcut menus	118
My Entity List pane – icons, buttons, and special text	122
My Entity List pane – actions and shortcut menus	122
Genomic Viewer	124
Genome View	124
Chromosome View	126
Gene View	128
The View Cursor	132
Tab View	133
Status Bar	138

Dialog Boxes	139
Agilent Feature Extraction Importer	139
CGHSmooth Parameters	141
CGHSmooth Plot	143
Chart Properties	145
Configure Coloring Ranges and Shades	150
Confirm Overwrite	151
Create Experiment	152
Create Gene List	154
Create Track	156
Customize Search Link	157
Design Properties	158
Echo Example Plot	162
Edit Array Color	165
Edit Array Order	166
Experiment Properties	167
Export	169
Export Array Attributes	170
Export Experiments	174
Export Filters	175
Export Tracks	176
Find in column	177
Go To Gene/Genomic Location	178
Import	179
Import (experiments)	181
Import (filters)	182
Import GEML design files	183
Import Genome Build	185
Import Track	186
Microarray Properties	187
MovAvg Example Parameters	190
MovAvg Example Plot	192
Scroll to Column	194

Contents

Search probes in eArray	195
Select Color	196
Select data type for experiments (UDF files – CGH or CH3)	199
Set genome build and species for Axon design files	200
Show/Hide Columns	202
Track	203
UDF Import Summary (CGH or CH3)	205
Universal Data Importer - Map Column Headers	206
User Preferences	208
View coordinates in UCSC browser	214
View Preferences	216



1

Getting Started

Using Agilent Genomic Workbench Lite Edition on a Mac	13
Using Main Window Components to Display Data	14
Switching applications	18
Using Tabs and Command Ribbons	19
Using the Navigator to Search for Data	22
Using the Genomic Viewer to Display Data	26
General Instructions for Displaying Microarray Data	30
Quick-start Instructions for Displaying Microarray Data	31
Getting Help	39

This guide describes how to use Agilent Genomic Workbench Lite Edition to display data if you do not have a CGH, ChIP, or methylation (CH3) DNA Analytics license.

This chapter gives an overview of the window components and how to use Agilent Genomic Workbench Lite Edition to view data. Without a license, you have a number of capabilities, that include the import, management, export, and display of CGH, ChIP, and CH3 data.

To display imported data, you organize the data files into logical units called *experiments*. Experiments are used to define the data you want to display using Agilent Genomic Workbench. After you create them, and add array data, you can then display the data.

For a description of the commands and dialog boxes that appear when you use the program, see [Chapter 4](#), “Data Viewing Reference”.






NOTE

Descriptions in this guide cover only the commands and options that are available for viewing data using Agilent Genomic Workbench Lite Edition without a DNA Analytics license. For information on commands and options that are available with a license, or for information on Sample Manager, Workflow, or SureSelect Target Enrichment, see the *User Guide* for the module that you want to use.

Using Agilent Genomic Workbench Lite Edition on a Mac

The content of this User Guide applies to both the Windows and Mac versions of Agilent Genomic Workbench Lite Edition. Both of these versions have the same features. However, when you use the Mac version of the program, please note the following:

Windows command	Equivalent Mac command
Right-click	<ul style="list-style-type: none">• Command-click ( -click)• On Macs with trackpads, other options are available. On certain machines, you place two fingers on the trackpad while you press the button below the trackpad. See the user guide for your specific machine.• If you have a third-party mouse that has more than one button, you may be able to use one of the buttons as a right mouse button.
Control-click	Control-click (Same as the Windows command)
Shift-click	Shift-click (Same as the Windows command)
 (Close button)	 (Close button)

Using Main Window Components to Display Data

You can use the data *viewing* capability in Agilent Genomic Workbench Lite Edition without a license. You can view data for many types of arrays, including CGH, ChIP, and Methylation (CH3). You can use the data *analysis* capability in Agilent Genomic Workbench Lite Edition only if you have a license for one or more of the DNA Analytics programs (CGH, ChIP, or Methylation).

What are the main window components?

You use four primary components of the Agilent Genomic Workbench Lite Edition main window to import, manage, export and display extracted data.

- Home tab commands – import, manage and export data
- Navigator – create and fill new experiments with array data
When you make the experiment active, the data appear in the display, called Genomic Viewer.
- Genomic Viewer – display data and content in four Views: Genomic View, Chromosome View, Gene View, and Tab View
- View tab commands – change appearance of Genomic Viewer display

[Figure 1](#) shows the main window of Agilent Genomic Workbench Lite Edition when the Genomic Viewer tab is selected, and identifies the names of its components.

To learn how to display log ratio data, content, and analyze data to show results, see the *User Guide* for which you have a DNA Analytics program license(s).




Figure 1 Agilent Genomic Workbench Lite Edition main window with major components – unlicensed CGH version

What can you do with the main components to display data?

See the table below for the parts of the main window you use to display log ratio data.

Table 1 Components of Agilent Genomic Workbench main window for display of data

To do this	Use this part of the main window
Change program to CGH, ChIP, Methylation (CH3),and SureSelect Target Enrichment	Switch Application button: Click the button and click the program you want to open. Do this to display different data types, even if you have no license. The scatter plot options are different for the different program types.
Import Agilent design files	Home tab: Click the Import button and select Design Files>GEML File to select a design file to import. See Chapter 2 , “Importing, Managing, and Exporting Data and Other Content” for more information.
Import or export data	Home tab: Click the Import or Export button to select the data you want to import or export. See Chapter 2 , “Importing, Managing, and Exporting Data and Other Content” for more information.
Select array data to display in the three graphical views or in the Tab View as a table	Experiment pane of the Navigator: Create an experiment with the imported data, select the experiment, and then select the data within the experiment to display data. See Chapter 3 , “Displaying Data and Other Content” for more information.
Display array data for only a certain portion of a chromosome	Genome View: Select a chromosome to display in Chromosome View. You cannot view log ratio data points here. Chromosome View: Select a gene region to display in Gene View. You can display log ratio data points here if you select Scatter Plot in the View Preferences dialog box. Gene View: See the log ratio data next to a selected region of a chromosome, with associated genes and track-based annotation. See Chapter 4 , “Data Viewing Reference” for details about these Views.


What can you do with the main components to display data?

Table 1 Components of Agilent Genomic Workbench main window for display of data

To do this	Use this part of the main window
Show/Hide or customize the data points for the scatter plots	<p>Gene View: Move the mouse pointer over Scatter Plot to display the options. Or, right-click and then click View Preferences.</p> <p>Chromosome View: Right-click and then click View Preferences.</p> <p>View tab: Click View Preferences.</p> <p>See Chapter 3, “Displaying Data and Other Content” for information on how to do this.</p>
Display array data next to tracks or gene lists	<p>My Entity List pane of Navigator: Add or select a track or gene list to have it appear in Gene View.</p> <p>See Chapter 3, “Displaying Data and Other Content” for information on how to do this.</p>
Change the appearance of the display	<p>View Tab: Click View Preferences. From the View Preferences dialog box, you can change the orientation, select what type of data to view, and configure scatter plot options.</p> <p>Genomic Viewer: Right-click any View except the Tab View and select View Preferences. In the View Preferences dialog box, you can select to show or hide the scatter plots and how to display them. If you have one or more DNA Analytics licenses (CGH, ChIP, or Methylation), you can show or hide the results.</p> <p>See Chapter 3, “Displaying Data and Other Content” for more information.</p>

Switching applications

You can use the Agilent Genomic Workbench to work with a variety of different data types. Because the requirements for the display of data (and calculation of results, if using a license) are different for different data types, you must switch the application for the type of data you want to display.

The Switch Applications menu, located at the upper right corner of the Agilent Genomic Workbench window, is used to change the application. The selected application is marked . The selected application is also displayed in the title bar of the Agilent Genomic Workbench main window.




Figure 2 Switch Application menu

Using Tabs and Command Ribbons

Tabs

When you click a *tab*, groups of commands or single commands appear that are specific for that tab. The tabs that are displayed change depending on what licenses you have, and what application is selected (such as CGH, ChIP, CH3). Without a license, you only use the Home and View tabs to display data.




Figure 3 Agilent Genomic Workbench Lite Edition tab menu for CGH without a license

The following table summarizes what you can do from the tabs of Agilent Genomic Workbench, with a DNA Analytics (CGH, ChIP, CH3) application selected, but without any license installed.

Table 2 Capabilities in tabs

Tabs	Capabilities
Home	Set preferences for display of tracks. Set eArray user and data locations. Set licences for analysis applications. Import array files, design files, genome builds, tracks, array attributes, and experiments. Export experiments, tracks, and array attributes. Create an experiment. Find and go to a gene or genomic location.
View	Set up preferences for display of data Copy displayed data to the Clipboard Turn on or off display of Views and Navigator Turn on or off tabular display of signal intensity and annotations Turn on or off display of Cytoband information in Gene View Turn on or off highlight of nonunique probes Turn on or off display of custom data
Tool	Set parameters for plug-ins Display plug-in examples
Help	View program information and User Guides.

Commands

The area where commands appear is called a *command ribbon*. The command ribbon that appears when you click the Home tab is shown below. The commands that appear in the command ribbon change depending on what application module is selected, and which tab in that application module is selected.



Figure 4 Command ribbon for the Home tab

For a complete description of all of the command ribbons and commands you see in Agilent Genomic Workbench, see “[Command Ribbons](#)” on page 103.

Using the Navigator to Search for Data

This section gives you instructions on how to search for design files, extracted FE data, experiments and other information in the Navigator of Agilent Genomic Workbench.




Figure 5 Navigator panes

The Navigator ([Figure 5](#)) shows the array data, experiments, and other content stored in Agilent Genomic Workbench that is available to the user for display. It contains the following panes:

Navigator Pane	Comments
Search	Lets you search within any pane of the Navigator for a specific design or content, or for items that contain a specific string of characters, when using asterisks (*) as wildcards. See “Search pane” on page 112 for more information.
Data	<p>Contains microarray data files, organized by type, then by design and genome build.</p> <p>Shows all microarray designs that are available to you, organized by folders. In general, you can:</p> <ul style="list-style-type: none"> • Expand or collapse folders to show or hide content • Right-click the name of a folder or item to open a shortcut menu that lets you take action on the item. <p>See “Data pane – icons, special text, and buttons” on page 115 and “Data pane – actions and shortcut menus” on page 116.</p>
Experiment	Contains Agilent Genomic Workbench experiments. Experiments are organizational units that contain links to microarray data and design files. See “Experiments Folder” on page 118, “Experiment pane – icons, special text, and buttons” on page 117, and “Experiment pane – actions and shortcut menus” on page 118.
My Entity List	<p>Contains gene lists and tracks:</p> <ul style="list-style-type: none"> • Gene Lists are collections of genes of interest. You can create them within the program, import and export them, and apply them to Gene View and Chromosome View. • Tracks are collections of annotation or other information that map to specific genomic locations. You can import, export, and combine tracks, and display them in Gene View with your array data and analysis results. See “My Entity List pane – icons, buttons, and special text” on page 122 and “My Entity List pane – actions and shortcut menus” on page 122.

To search the Navigator

You can search one or all of the panes of the Navigator for items that match a specific search term. Figure 6 shows the search pane of the Navigator, and identifies a couple of its elements.





Figure 6 Search pane of the Navigator

- 1 At the top of the Navigator, in the Pane list, select the pane to be searched. To search in all panes, select **All Panels**. If the pane list is not visible, click  to show it.
- 2 In the search term box, type the desired search term. The search term is not case sensitive, but it must contain the complete entry that you want to find. You can use asterisks (*) to represent one or more unspecified characters. For example, type *12345* to find any item that contains "12345".

- 3 Click .

The program searches the selected pane(s) for items that match your search term. If it finds matching items, the program expands the appropriate folders, and displays the names of the matching items in red. The first matching item is highlighted in yellow.


- 4 Do any of the following:

- To highlight the next matching item, if one is available, click



- To highlight the previous matching item, click



- 5 After you complete the search, click  to clear the results of the search, as well as your search term.

Using the Genomic Viewer to Display Data

What is the Genomic Viewer?

Genomic Viewer is the graphics and tabular display section of the Agilent Genomic Workbench main window. In the Genomic Viewer, extracted data and analysis results can be tabulated and displayed next to depictions of the genome, selected chromosome, and selected genes of the species whose array data you are analyzing.

There are four main views in the Genomic Viewer, as shown in [Figure 7](#).

- **Genome View** – A graphical representation of the entire genome for the selected species. Use this view to select the chromosome to show in the other views.
- **Chromosome View** – A graphical representation of the selected chromosome, displayed with cytobands and a plot area. Click or drag the mouse to select a region to display in the Gene View.
- **Gene View** – A more detailed view of the chromosomal region selected in the Chromosome View.
- **Tab View** – Displays design annotation and log ratio data related to the chromosome you select in Chromosome View

For more information on the Genomic Viewer and its views, see [Chapter 4](#), “Data Viewing Reference”.





Figure 7 Genomic Viewer in vertical orientation

1 Getting Started

To change the size of and detach panes from the Agilent Genomic Workbench main window

To change the size of and detach panes from the Agilent Genomic Workbench main window

- To change the size of a pane in the main window, drag one of its inside borders.
- To detach a pane from the main window and open it in a separate window, click its **Detach** button  .




Figure 8 Changing the size of and detaching panes

To maximize and reattach panes to the Agilent Genomic Workbench main window

To maximize and reattach panes to the Agilent Genomic Workbench main window

- To display a view full-screen in a separate window, click its **Maximize** button.
- To reattach a view in a separate window to the main window, click its **Close** button.




Figure 9 Maximizing and reattaching panes

General Instructions for Displaying Microarray Data

*An **experiment** is the folder that holds data from any array set you select for the experiment. The folder also holds analysis results.*

You set up experiments to display all data in the Genomic Viewer. To set up an experiment you:

- Import data
- Create a new experiment
- Add the imported data to the experiment
- Select the experiment to display data




Figure 10 Typical pathway for display of microarray data

Quick-start Instructions for Displaying Microarray Data

The instructions in [Table 3](#) show how to organize imported log ratio FE data so you can display your data next to the corresponding cytobands. Without a DNA Analytics license, only log ratio data is displayed, not results.

These instructions assume that:

- All instructions apply whether you have a license or not.
- You use only Agilent data and design files. If you choose to use the demo Agilent design and data files that come with the program, you do not need to import those files.

Table 3 How to display data in Agilent Genomic Workbench



To do this	Follow these instructions	Comments
 A blue rounded rectangle containing the text "Select the application" in white.	<ol style="list-style-type: none">1 If you are in another application, click Switch Application.2 Select the CGH, ChIP, or CH3 application type. You do not need a license to perform the following steps.	<ul style="list-style-type: none">• If you are using a licensed version of the application, make sure all the analysis options are turned off. (Clear check boxes in the Analysis command ribbon.)
 A blue rounded rectangle containing the text "Turn on scatter plot options" in white.	<ol style="list-style-type: none">1 By default, the scatter plot is turned on for Gene View. To view or change the scatter plot options, in Gene View move the pointer over the arrow next to Scatter Plot and mark one or more of the plot check boxes. See "Gene View" on page 128.2 To view or change additional options for the scatter plot, or to change the orientation of the panes in the Genomic Viewer, right-click in one of the panes and select View Preferences. See "View Preferences" on page 216.	<ul style="list-style-type: none">• The check boxes in Scatter Plot set the program to draw data points that represent log ratio or other selected values.• If you do not turn on the scatter plot functions, you will see nothing in the Genomic Viewer.• If you are using the Agilent demo files, continue to <i>Create an Experiment</i>.

Table 3 How to display data in Agilent Genomic Workbench (continued)

To do this	Follow these instructions	Comments
Import design files	<ol style="list-style-type: none"> 1 To select Agilent GEML-based microarray design files, click Home > Import > Design Files > GEML File. 2 In the dialog box that appears, select the file you want to import, then click Open. The Import GEML design dialog box appears. 3 If necessary, select the Genome Build for your files. 4 Click Start Import. Design Import Summary appears, and the design file with selected genome build appears in the Navigator in the Data pane. See “Using the Navigator to Search for Data” on page 22. 	<ul style="list-style-type: none"> • If you want to import Agilent or Axon array data files, the program requires their design file(s). If the design file(s) are not already available in the Navigator, (for example, downloaded from the eArray Web site) you must import them. • When you import a design file, the program shows the genome build(s) that can be used by the design file as nodes under the design file. • The current builds of the human genome (hg), mouse genome (mm) and rat genome (rn) are available in Agilent Genomic Workbench. Should you want to import a design file for a different genome build, you must import the genome build first.
Import data files <ul style="list-style-type: none"> • Agilent FE files • Axon files • UDFs 	<ol style="list-style-type: none"> 1 To import Agilent FE files, click Home > Import > Array Files > FE File. 2 Find and select the desired file, then click Open. To select multiple files, hold down the Ctrl key and click their names. 3 In the dialog box that appears, in Dye Flip, select either Normal or Flipped for each FE or Universal Data File (UDF). 4 Click OK. 5 In the Navigator, check the Data folder to make sure that the program imported the correct files. 	<ul style="list-style-type: none"> • In Dye-Flip, select Normal if: <ul style="list-style-type: none"> • The test samples were labeled with cyanine-5 (red). • The control samples were labeled with cyanine-3 (green). • The imported ratio (test/control) will be reported directly. • In Dye-Flip, select Flipped if: <ul style="list-style-type: none"> • The test samples were labeled with cyanine-3 (green). • The control samples were labeled with cyanine-5 (red). • The imported ratio (control/test) will be reported with the ratio inverted (test/control).

1 **Getting Started**
Quick-start Instructions for Displaying Microarray Data

Table 3 How to display data in Agilent Genomic Workbench (continued)



To do this	Follow these instructions	Comments
		<ul style="list-style-type: none">• The program automatically puts the imported data under the genome build folder that belongs to the design used with the arrays.
<div>Create an experiment</div>	<ol style="list-style-type: none">1 In the Experiment pane of the Navigator, right-click the Experiments folder, then select New Experiment. A dialog box appears.2 Type a name and an optional description for the experiment.3 Click OK.4 (optional) To add data to the experiment now, click Properties. Otherwise continue and add data, as described in the next step.	<ul style="list-style-type: none">• The new experiment appears as a node within the Experiment pane of the Navigator. The node becomes a folder once data is added to the experiment. 
<div>Add data to the new experiment</div>	<ol style="list-style-type: none">1 Fully expand the Data folder, and click the name of an array you want to add to your new experiment.2 Drag the selected arrays to the folder of the new experiment.	<ul style="list-style-type: none">• To select additional arrays within the same design, hold down the Ctrl key and click their names.• You can also right-click the name of the experiment, and select Show Properties to add arrays to an experiment.

Table 3 How to display data in Agilent Genomic Workbench (continued)


To do this	Follow these instructions	Comments
	<div><div>Select the experiment</div><div><ol style="list-style-type: none">1 In the Navigator, double-click the name of the experiment. A dialog box appears.2 Click Yes.<ul style="list-style-type: none">• In the Navigator, the name of the experiment turns blue, which indicates that it is the active experiment.• You now see the scatter plot log ratio data in the Chromosome and Gene Views.</div></div>	<ul style="list-style-type: none">• The program automatically selects the first array in the experiment for display.• A data table appears in Tab View.• At this point if the arrays within the experiment use a genome build different from the one represented in the Views, the program changes the chromosomal and gene information used in the Views.

Table 3 How to display data in Agilent Genomic Workbench (continued)


To do this	Follow these instructions	Comments
<div>Select arrays to display</div>	<p>To select an array, right-click the array name, and click Select. To clear an array selection, right-click the name of the array in the Navigator, then click Deselect.</p>  <p>You can also select or deselect several arrays at a time. Hold down the Shift key and click the contiguous arrays whose log ratio data you want to display. Hold down the Ctrl key and click the non-contiguous arrays whose log ratio data you want to display.</p>	<ul style="list-style-type: none">• In the Navigator, the icons beside the arrays become colored, when enabled for the selected experiment.• In Tab View, colored squares appear in the column heading for the arrays when selected. You can select or deselect only one array at a time in Tab View, or you can select or deselect all arrays at the same time.• The program adds the data from the array to the Chromosome and Gene views.

Table 3 How to display data in Agilent Genomic Workbench (continued)

To do this	Follow these instructions	Comments
<div>Review data</div>	<ul style="list-style-type: none"> • In Genome View, click a chromosome of interest. • In Chromosome View, drag the pointer over a region of the chromosome graph to display it with more resolution in Gene View. • In Gene View, click the + and – buttons to zoom in and out. • In Gene View, click anywhere within the scatter plot to recenter the view at that location. • To see information for the log ratio data, in Gene View, move the pointer over the arrow next to Scatter Plot to show the options. Under Configure Coloring schemes, mark the check box and select Log Ratio Values. • In Gene View, zoom in so that single data points are visible, then place the pointer over a data point. If ToolTip is enabled in View Preferences, a box appears that describes the data point. 	<ul style="list-style-type: none"> • The solid, horizontal blue lines in Chromosome and Gene views are referred to as the <i>View cursor</i>. The chromosomal location of the cursor appears in the Status bar, located on the lower left corner of the screen. • If you still cannot see the Scatter Plot data in Chromosome and Gene View or ToolTips in Gene View, do the following: <ol style="list-style-type: none"> a Right-click either View and click View Preferences. b Under Data Visibility, select All Views, then mark the Scatter Plot check box. c Under Data Visibility, select Gene View and then click Scatter Tool Tip. • When you right-click an empty area of Gene View, you can also use the shortcut menu to create a gene list or track, or to search the Agilent eArray database for probes from the selected region. See the <i>User Guide</i> for your application.

1 Getting Started

Quick-start Instructions for Displaying Microarray Data




Figure 11 Segment of Chromosome View and Gene View with scatter plot of log ratio data and ToolTip

Getting Help

To get help within Agilent Genomic Workbench

Agilent Genomic Workbench has several help resources. All help guides open in Adobe® Reader®.

Help Resource	Description/Instructions
Data Viewing User Guide	<p>This user guide, which you are now reading, supplies comprehensive help on all available Data Viewing tasks. You can access it easily from anywhere within the program.</p> <ol style="list-style-type: none">1 In any tab of Agilent Genomic Workbench, click the Help tab.2 On the Help Ribbon, click Data Viewing. The Data Viewing User Guide opens.
Other User Guides	<p>The Help tab in Agilent Genomic Workbench lets you view any of the available user guides that apply to the currently selected application type.</p> <ol style="list-style-type: none">1 Set the desired application type from the Switch Application menu.2 In the Agilent Genomic Workbench tab bar, click Help. The names of the available user guides appear in the command ribbon.3 Click the desired help guide. The selected guide opens.
Product Overview Guide	<p>An additional guide gives an overview of the capabilities within Agilent Genomic Workbench and describes how to start and find help for all of the programs.</p> <ol style="list-style-type: none">1 In any tab of Agilent Genomic Workbench, click the Open Application tab.2 At the upper right corner of the Open Application tab, click Product Overview.

To contact Agilent Technical Support

Technical support is available by phone and/or e-mail. A variety of useful information is also available on the Agilent Technical Support Web site.

Resource	To find technical support contact information
Agilent Technical Support Web site	<ol style="list-style-type: none">1 Go to http://chem.agilent.com.2 Select a country or area.3 Under Quick Links, select Technical Support.4 Select from the available links to display support information.
Contact Agilent Technical Support by telephone or e-mail (United States and Canada)	Telephone: (800-227-9770) E-mail: informatics_support@agilent.com
Contact Agilent Technical Support by telephone or e-mail (for your country)	<ol style="list-style-type: none">1 Go to http://chem.agilent.com.2 Select Contact Us.3 Under Worldwide Sales and Support Phone Assistance, click to select a country, and then click Go. Complete e-mail and telephone contact information for your country is displayed.

To learn about Agilent products and services

To view information about the Life Sciences and Chemical Analysis products and services that are available from Agilent, go to www.chem.agilent.com.



2 Importing, Managing, and Exporting Data and Other Content

Importing Files 42

Working with Experiments to Organize Imported Data 54

Managing Content (Data, Gene Lists, Tracks) 62

Exporting and Saving Content 70

This chapter describes how to import, organize, manage, and export data and other content within the user interface of Agilent Genomic Workbench. The program lets you import many different kinds of files, including array data and design files from Agilent products and other sources, and other content such as annotation tracks.

See [Chapter 4](#), “Data Viewing Reference” for a description of the Agilent Genomic Workbench main window and its contents, and descriptions of the dialog boxes that can appear.



Importing Files

You use the Home tab to import many kinds of files into Agilent Genomic Workbench. The table below summarizes the kinds of files you can import, and the topics in this section that describe how to import them.

The Data pane of the Navigator displays all of the content available for the user. See [“Navigator Pane”](#) on page 114 for more information on the Navigator panes and how to use them.

File type	Comments	See these topics
Microarray data files	<ul style="list-style-type: none">Agilent Feature Extraction (*.txt) data filesAxon (*.gpr) data filesUniversal Data Files (UDFs) (*.txt files)	“To import Agilent FE or Axon data files” on page 45 “To import a UDF file” on page 46
Microarray design files	<ul style="list-style-type: none">Agilent GEML (*.xml) design filesAxon (*.gal) design files	“To import Agilent GEML design files” on page 43 “To import Axon design files” on page 44
Genome builds	Agilent-supplied genome information for human, mouse and rat genomes	“To import a genome build” on page 50
Tracks	BED format annotation track files	“To import tracks” on page 51
Array attributes	Attribute .txt files that you have created yourself or previously exported from Agilent Genomic Workbench	“To import array attributes” on page 51
Experiments	ZIP file of experiments exported from Agilent Genomic Workbench	“To import an experiment file” on page 52

To select a different location for data files

By default, the program stores microarray and experimental data in **C:\Program Files\Agilent\Agilent Genomic Workbench Lite Edition <version>\data**. If you want, you can select a different location.

- 1 In the Home tab, click **User Preferences**.

The User Preferences dialog box appears. See “[User Preferences](#)” on page 208.

- 2 In the **Miscellaneous** tab, under **Data Location**, click **Browse**.

An Open dialog box appears.

- 3 Select a location, then click **Open**.

The selected location appears in the User Preferences dialog box, in Data Location.

- 4 Click **OK**.

CAUTION

If you change the location for data files, and there is a data folder in that location, the data will be overwritten by the current data.

To import Agilent GEML design files

The Agilent Genomic Workbench database must contain designs that match the Agilent Feature Extraction data files you want to import. Your imported GEML files contain array-specific information such as probe names, annotations, and chromosomal locations, and are associated with a specific genome build.

To import an Agilent GEML file, use the following procedure:

- 1 In the Home tab, click **Import > Design Files > GEML File**.


The Import Design Files dialog box appears. See “[Import](#)” on page 179. The dialog box shows only *.xml files.

- 2 To select a file for import, click its name. To select additional files, hold down the **Ctrl** key while you click their names.
- 3 Click **Open**.

2 Importing, Managing, and Exporting Data and Other Content

To import Axon design files

The program validates the selected file(s), and the Import GEML Design Files dialog box appears. See “[Import GEML design files](#)” on page 183.

If a design file passes validation, the Status column will show **Update** in green. If a design file fails validation, **Corrupt** appears in the Status column beside it, and the program will not import the file. To remove the corrupt array from the list, click **Remove** .

4 Click **Start Import**.

The program imports the file(s). The files appear as new design folders in the Imported External Designs folder of the Data pane of the Navigator, with the genome build as a node within the folder.

You can import two design files with the same name, but associated with different genome builds; for example, Hg17 or Hg18. If you do, the program creates a single design folder with two nodes, one for each genome build.

To import Axon design files

You can import Axon (*.gal) microarray design files into Agilent Genomic Workbench. The program requires the Axon design files that match all Axon array data files you import.


1 In the Home tab, click **Import > Design Files > Axon File**.

The Import Axon Design Files dialog box appears. See “[Import](#)” on page 179. The dialog box shows only *.gal files.

2 To select a file for import, click its name. To select additional files, hold down the **Ctrl** key while you click their names.

3 Click **Import**.

The program validates the selected file(s), and the Set genome build and species for Axon design files dialog box appears. See “[Set genome build and species for Axon design files](#)” on page 200.

If a design file passes validation, the Status column will show **Update** in green. If a design file fails validation, **Corrupt** appears in the Status column beside it, and the program will not import the file. To remove the corrupt array from the list, click **Remove** .

4 For each design file, select the appropriate **Species** and **Genome Build**.

5 Click **Start Import**.

The program imports the file(s). The files appear as new design folders in the Data pane, organized by application (CGH, ChIP, or methylation, for example).

To import Agilent FE or Axon data files

You can import several types of microarray data files into Agilent Genomic Workbench:

- Agilent Feature Extraction (FE) *.txt data files
- Axon (*.gpr) data files
- Universal Data Files (UDFs) (*.txt files) See [“To import a UDF file”](#) on page 46 for instructions on how to import this file type.

In order to import Agilent Feature Extraction files, the representative GEML array design files must be imported first. In order to import Axon data files, the representative Axon.gal design files must be imported first. See [“To import Agilent GEML design files”](#) on page 43 or [“To import Axon design files”](#) on page 44.

1 In the Home tab, do one of the following:

- To import Agilent FE data files, click **Import > Array Files > FE File**.
- To import Axon data files, click **Import > Array Files > Axon File**.

A dialog box appears. Only data files of the appropriate type appear. See [“Import”](#) on page 179.

2 To select a file for import, click its name. To select additional files, hold down the **Ctrl** key while you click their names.

3 Do one of the following:

- For Agilent FE files, click **Open**.
- For Axon files, click **Import**.

The Agilent Feature Extraction/Axon File Importer dialog box appears. See [“Agilent Feature Extraction Importer”](#) on page 139.

4 Set the following, as needed:

2 Importing, Managing, and Exporting Data and Other Content

To import a UDF file

Setting	Comments
Name	The names of imported arrays are often cryptic. You can give any array a more meaningful label. a Double-click the name of the array. b Edit the name. c Press Enter .
Dye Flip	For each array: • Select Normal if: • The test samples were labeled with cyanine-5 (red). • The control samples were labeled with cyanine-3 (green). • The imported ratio (test/control) should be reported directly. • Select Flipped if: • The test samples were labeled with cyanine-3 (green). • The control samples were labeled with cyanine-5 (red). • The imported ratio (control/test) should be reported with the ratio inverted (test/control). The program does not combine dye-flip pairs.
Overwrite arrays with duplicate names	If you mark this option, the program deletes an existing array data file if it has the same name as one you import.

5 Do one of the following:

- To import the file(s) while you wait, click **OK**.
- To import the file(s) in the background, click **Run in Background**. This lets you continue while the program imports the files.

To import a UDF file

UDF files are plain text files that contain array data in tab-delimited format. Files must contain the following six columns of information, in any order. Each column must contain the following column names, as column headers, or you must “map” the names from the file to these columns in Agilent Genomic Workbench:

- Probe name
- Chromosome name
- Start position

- Stop position
- Description
- Signal intensity data (The file can contain additional columns, each with data from an additional array.)

When you import a UDF file, the program creates a new design based on the information you enter during import, and the information in the file itself. This design contains all of the arrays represented in the file. The program also creates a new experiment that contains the arrays.

1 In the Home tab, click **Import > Array Files > UDF File**.

The UDF Files dialog box appears. See [“Import”](#) on page 179. Only *.txt files appear in the dialog box.

2 Select the UDF file, then click **Open**.

The Select data type for experiments dialog box appears. See [“Select data type for experiments \(UDF files – CGH or CH3\)”](#) on page 199.

3 For each array, set the following, as needed:

Setting	Comments
Experiment Name	By default, the program creates an experiment with the same name as the imported file. To change the name: <ul style="list-style-type: none"> a Double-click the name. b Edit the name. c Press Enter.
Data type	<ul style="list-style-type: none"> • Select the mathematical form of the signal intensity data for the array. The options are ratio, log₂ ratio, log₁₀ ratio, and ln ratio.
Design type	<ul style="list-style-type: none"> • Select cgh, expression, or CH3.

4 Click **Continue**.

2 Importing, Managing, and Exporting Data and Other Content

To import a UDF file

When you “map” a column, you assign the column heading (in an external file) to a column heading in Agilent Genomic Workbench.

The Universal Data Importer – Map column headers dialog box appears. The main table in the dialog box contains the first few rows of data from the file. Column headings derived from the first line of the file appear at the top of the table as a guide, but the program does not interpret these headings. See “[Universal Data Importer - Map Column Headers](#)” on page 206.

- 5 Below each column heading, select the label that identifies the content of the column. Use each label exactly once, except for LogRatio, which you can use many times. Alternatively, in **Select Mapping**, select a saved column map.

These options are available:

Column Label	This column contains:
ProbeName	Names of probes.
ChrName	Names of chromosomes.
Start	First chromosomal location to which each probe is designed.
Stop	Last chromosomal location to which each probe is designed.
Description	Text annotation related to the probe.
LogRatio	Array data values that correspond to each probe. You can use this label more than once.

- 6 Under **Species Info**, select the **species** and **Genome Build** appropriate to the data in the file.
- 7 If you expect to import many similar UDFs in the future, follow these steps to save the column map:
- a Under **Mapping Info**, click **Save Mapping As**.
An Input dialog box appears.
 - b Type a name for the column map, then click **OK**.
The name of the saved map appears in **Select Mapping**.

In the future, you can select this mapping and apply it to any UDF file that you import.

- 8** By default, the program creates a “Virtual Array ID” that becomes the ArrayID attribute for the array(s) in the UDF. To create your own virtual Array ID, follow these steps:

- a** Under **ArrayID Info**, clear **Use System Generated Array ID**.
- b** Double-click the number in **Virtual Array ID**, then type your own Array ID.

For more information on Array IDs, see the *Sample Manager User Guide*.

- 9** Click **Import**.

The program validates your column mapping. A dialog box appears. If you need to fix the column map, the dialog box has a list of the missing column label(s). If the column map is complete, a message asks if you want to import additional files with the same mapping.

- 10** Do one of the following:

- If you want to import additional files with the same column mapping, follow these steps to include these files in the import:
 - a** Click **Yes**.
The UDF Files dialog box appears.
 - b** Click the name of a file to select it for import. Hold down the **Ctrl** key while you click the names of additional files.
 - c** Click **Open**.
- If you do not want to include additional file(s) in the import, click **No**.

The Program imports all requested files, and the UDF Import Summary dialog box appears. This dialog box shows the imported files, the number of lines of data that were imported for each file, and the number of lines that were skipped, if any. If a file name appears in red, the program may not have imported the file. See “[UDF Import Summary \(CGH or CH3\)](#)” on page 205.

- 11** Click **OK**.

In the Data pane, in the appropriate design type folder within the Data folder, a new design folder appears. The design folder contains the imported array data.

A new experiment appears in the Experiments folder in the Experiment pane, that contains the array data. This experiment has the name of the imported UDF file, unless you changed it during import.

To import a genome build

In general, the program uses the genome build specified in the array design file, and protects it from changes. If a genome build is not available in the program, you can import one.

NOTE

Use arrays from a single genome build in an experiment.

- 1 In the Home tab, click **Import > Genome Build**.

The Import Genome Build dialog box appears. See “[Import Genome Build](#)” on page 185.

- 2 Set the following. All are required.

Setting	Instructions
Species	<ul style="list-style-type: none">Type the genome’s species of origin, as you would like it to appear within the program.
Build Name	<ul style="list-style-type: none">Type the name of the genome build you want to import, as you would like it to appear within the program.
Refseq File	<p>This file contains information on gene locations for Gene View.</p> <ul style="list-style-type: none">a Click Browse. A dialog box appears.b Select the file, then click Open.
Cyto-band File	<p>This file contains the graphic information on the cytobands for Genome and Chromosome Views.</p> <ul style="list-style-type: none">a Click Browse. A dialog box appears.b Select the file, then click Open.

- 3 Click **OK**.

To import tracks

You can import BED format track files into Agilent Genomic Workbench. Track files contain specific features correlated with chromosomal locations, and apply to a specific genome build of a given species.

- 1 In the Home tab, click **Import > Track**.

The Import Track dialog box appears. See [“Import Track”](#) on page 186.

- 2 Set the following. All are required.

Setting	Instructions
Species	<ul style="list-style-type: none"> • Select the species to which the track applies.
Build Name	<ul style="list-style-type: none"> • Select the specific genome build of the species to which the track applies.
Track Name	<ul style="list-style-type: none"> • Type a name for the track. This name identifies the track within the program, including the name that appears if you include the track in Gene View.
Track File	<ol style="list-style-type: none"> Click Browse. A dialog box appears. Select the name of the track (*.bed) file that you want to import. Click Open. The location of the file appears in Track File.

- 3 Click **OK**.

The program imports the track. To view the track in Gene View, and to manage tracks, see [“To show tracks in Gene View”](#) on page 90.

To import array attributes

An array attributes file is a tab-delimited *.txt file that contains a list of arrays by ArrayID, and values for specific array attributes. Attributes are pieces of array-specific information, such as the hybridization temperature and the name of an array set that contains the array.

2 Importing, Managing, and Exporting Data and Other Content

To import an experiment file

Although you can import array attributes with this function, the Sample Manager application lets you import and assign array attributes more easily. See the *Sample Manager User Guide* for more information.

To import an array attributes file

- 1 From the Home tab, click **Import** and then select **ArrayAttributes**.

The Import Attribute Files dialog box appears. See “[Import](#)” on page 179.

- 2 Select the microarray attributes file, then click **Import**.

The program imports the file. If the ArrayIDs in the file do not match the ArrayIDs of arrays in the program, a dialog box appears. The dialog box has a list of the ArrayIDs in the file that do not match. Click **No** to stop the import process, or click **Yes** to continue anyway.

To import an experiment file

In Agilent Genomic Workbench, an experiment is a set of links to microarray data and design files, and any associated results. An Agilent Genomic Workbench experiment file is a single ZIP file that contains the design and data files for one or more experiments. You can import

- Experiment files created in Agilent Genomic Workbench on another computer
- Agilent Genomic Workbench 5.0 and 6.0 experiment files

- 1 In the Home tab, click **Import > Experiments**.

The Import Experiments dialog box appears. See “[Import](#)” on page 179.

- 2 Select the ZIP file that contains the experiment(s) you want to import, then click **OK**.

The program imports the experiment file. Designs appear as new folders in the Data pane, in the applicable design type folder. Array data appears within the applicable design folder, organized by genome build. In addition, the experiment(s) appear in the Experiment pane, with the appropriate arrays.

NOTE

Agilent Genomic Workbench experiment files contain all of the design and array data files for an experiment, but do not include any analysis parameter settings, array selections, or analysis results. To export the data and design files from one or more experiments, see [“To export experiments”](#) on page 71.

To import filters

Filters are used in Agilent Genomic Workbench to include or exclude data from an analysis, based on filter criteria. Filters are created in the interactive CGH and ChIP applications, or in workflow setup.



- 1 In the Home tab, on the Command Ribbon, click **Import > Filters**.
The Import dialog box appears. See [“Import”](#) on page 179 for more information.
- 2 Select the file that contains the exported filter(s) for import. and then click **Import**.
- 3 In the filters Import dialog box, mark the **Import** box next to each filter you want to import, and then click **OK**.

Working with Experiments to Organize Imported Data

This section describes how to organize imported array data and designs into *experiments*. Experiments, shown in the Experiment pane of the Navigator, contain links to specific array data and design files in the Data pane. After you set up an experiment, you can then analyze selected array data within the experiment.

Because experiments only contain *links* to the actual data and design files, any number of experiments can use a given set of files. In the data analysis applications (CGH, ChIP, or methylation, for example), experiments also can contain saved experiment results.

To display the array designs and data in the program

- To display the directory of data in the program, use the Data pane (Figure 12). Double-click a folder to expand or collapse it, or click the  and  buttons.




Figure 12 Data pane of the Navigator

In the Data pane, the program organizes design files by the application (CGH, ChIP, or methylation, for example) to which they apply. It organizes array data files by genome build under the design with which they are associated.

You can right-click many elements of the Data pane to open shortcut menus. For more information, see [“Data pane – actions and shortcut menus”](#) on page 116.

Many icons can appear in the Data pane. See [“Data pane – icons, special text, and buttons”](#) on page 115 for a complete list.

The Search pane can help you find specific data files or other content. See [“To find specific content items in the Navigator”](#) on page 63.

To create a new experiment

In Agilent Genomic Workbench, *experiments* are organizational units that contain links to data and design files. To view or analyze data, you must first create an experiment and associate the data files with it. Because experiments only contain *links* to the actual data and design files, any number of experiments can use a given set of files. In data analysis applications (CGH, ChIP, or methylation, for example), experiments can also contain saved experiment results.

1 In the Home tab, click **Create Experiment**.

The Create Experiment dialog box appears. See [“Create Experiment”](#) on page 152.


2 Type a **Name** and an optional **Description** for the experiment.

3 Do one of the following:

- To create an empty experiment, and add data to it later, click **OK**. The program creates the experiment. To add arrays to the experiment later, see [“To add arrays to an experiment”](#) on page 57.

2 Importing, Managing, and Exporting Data and Other Content

To create a new experiment

- To create an experiment and add data to it now, follow these steps:
(You can add or remove data from the experiment later.)
 - a Click **Properties**.**
The Experiment Properties dialog box appears. See [“Experiment Properties”](#) on page 167.
 - b Under **Select Design**, select the design and genome build for the array data.**
The applicable arrays appear in Array List.
 - c In **Array List**, click the name of an array that you want in your experiment. Hold down the **Ctrl** key while you click the names of additional arrays.**
 - d Click .**
The program transfers the selected arrays to the Selected Array List.
The dialog box also has other options for adding arrays. See [“Experiment Properties”](#) on page 167 for more information.
 - e Click **OK**.**
The program creates the new experiment, and adds data to it from the selected arrays.
- To create an experiment and add data to it using the “drag and drop” method, follow these steps:
 - a To create an empty experiment, click **OK**.**
The program creates the experiment.
 - b From the Data pane, expand a design to see the build and array data.**
 - c Drag an array from the Data pane and drop it onto the experiment folder in the Experiment pane.**

In all cases, a folder with the name of the new experiment appears in the Experiment pane of the Navigator. For more information on the Navigator, see [“Using the Navigator to Search for Data”](#) on page 22.





Figure 13 Experiment pane of the Navigator

To add arrays to an experiment

After you create an experiment, or import one, you can add arrays to it. When you add arrays to an experiment, you create links between the experiment and the array data and design files. Because the program does not move the actual files, multiple experiments can share the same arrays.

- 1 In the **Experiment** pane, double-click the **Experiments** folder to expand it.
- 2 Right-click the name of the experiment, then click **Show Properties**.
The Experiment Properties dialog box appears. See [“Experiment Properties”](#) on page 167.
- 3 Under **Select Design**, select the design file and genome build for the arrays to add.
The arrays for the selected design file and genome build appear in Array List.
- 4 In **Array List**, select the arrays to add to the experiment. To select a single array, click its name. To select additional arrays, hold down the **Ctrl** key while you click their names.
- 5 Click .

The program transfers the selected arrays to the Selected Array List.

The dialog box also gives you other options for adding arrays. See [“Experiment Properties”](#) on page 167 for more information.

6 Click OK.

Or, to add array data to an experiment using the “drag and drop” method,

- 1** From the Data pane, expand a design to see the build and array data.
- 2** Drag an array from the Data pane and drop it onto the experiment folder in the Experiment pane.

If needed, the program adds appropriate design and genome build folders to your experiment folder in the Experiment pane. It places the arrays you selected in the appropriate genome build folder.

To change the order of arrays in an experiment

When you select an experiment, a table appears in the Tab View of Genomic Viewer that contains log ratio values and, if selected, signal intensities for arrays in the experiment. See [“Tab View”](#) on page 133. You can change the order in which the arrays appear in the table. If you display separate (stacked) scatter plots in Gene View and Chromosome View for each array, the array order also determines the order in which these plots appear. You can use this feature to organize your arrays more logically, or to make it more convenient to display certain arrays. It is especially useful if you have many arrays.



- 1** In the Experiment pane, right-click the name of the experiment, then click **Edit Array Order**.

The Edit Array Order dialog box appears. See [“Edit Array Order”](#) on page 166.

- 2** In **Design**, select the design that contains the arrays whose order you want to change.

The arrays from the selected design appear in Array Name.

- 3** Do any of the following:

- To move an array up in the list, click its name, then click .
- To move an array down in the list, click its name, then click .

- To sort the list based on a specific microarray attribute, select the attribute in **Order by**.

4 Click **OK**.

To change the display names for arrays in an experiment

You can change the name displayed for arrays in an experiment, based on array attributes. When you change the display names for arrays in an experiment, the array names are changed only for the selected experiment. The display names are unchanged in the Data pane and in the other experiments.

- 1 Expand the folders in the Experiment pane until you see the experiment you want to change.
- 2 Right-click the experiment name, and select **Show Properties**.
- 3 In the Experiment Properties dialog box, click **Display Name by** and select an attribute to use for display of array names.
- 4 Click **OK**. The names of the arrays in the experiment are changed to the selected attribute. If the attribute does not exist for an array, the Global Display Name will be displayed.

NOTE

To change the name of an array throughout Agilent Genomic Workbench, change its Global Display Name using Sample Manager. For more information, see the *Sample Manager User Guide*.

To rename an array in an experiment

When you rename an array in an experiment, you change the array's name only within the context of the selected experiment. The name of the array is unchanged in the Data pane, and in other experiments.

- 1 Expand the folders in the **Experiment** pane until you can see the array you want to rename.
- 2 Right-click the name of the array, then click **Rename**.
An Input dialog box appears.

- 3 Type the new name for the array, then click **OK**.

The name of the array in the tab view of the selected experiment is renamed. The global display name of the array is not changed.

To remove arrays from an experiment

When you remove arrays from an experiment, you only remove the links between the experiment and the data files. The files are still available in the program for use in other experiments. To completely remove files from the program, see [“To remove data or design files from the program”](#) on page 65.

- 1 In the **Experiment** pane, expand folders until you can see the experiment, and the array(s) that you want to remove from it.
- 2 In the **Arrays** or **Calibration Arrays** folder of the experiment, click the name of an array to select it for removal. Hold down the **Ctrl** key while you click the names of additional arrays.
- 3 Right-click one of the selected array names, then click **Delete**.
A Confirm dialog box appears.
- 4 Click **Yes**.

The program removes the links between the experiment and the selected array data files. If the removal of arrays leaves a design folder in the experiment empty, the program removes this folder as well.

To display or edit the attribute values of a specific array

Array attributes are pieces of information specific to an array, such as array type or hybridization temperature. In the Genomic Viewer, you can display or change attributes for each array.

NOTE

You cannot change the Array ID attribute.

- 1 Expand the folders of the Data pane or the Experiment pane until you can see the array of interest.
- 2 Right-click the name of the array, then click **Show Properties**.

The Microarray Properties dialog box appears, with a list of array attributes. See “[Microarray Properties](#)” on page 187. You can also edit the attributes of an array from this dialog box. In addition, if the array is an Agilent array, you can see header and feature information sent from the Agilent Feature Extraction program.
- 3 When you are finished, click **Close**.

NOTE

You use the Sample Manager tab to organize, create, import, and export array attributes. See the *Sample Manager User Guide*.

Managing Content (Data, Gene Lists, Tracks)

This section describes how to create, find, rename, update, combine, and/or remove content such as data, gene lists, and tracks, stored in Agilent Genomic Workbench. To display the data, gene list and track content, see [Chapter 3](#), “Displaying Data and Other Content”.

To display a list of the content stored in the program

The Data and My Entity List panes of the Navigator show the content stored in Agilent Genomic Workbench.






Figure 14 Agilent Genomic Workbench Lite Edition Navigator

Data pane – Shows all of the design and data files stored in the database. For more information, see “To display the array designs and data in the program” on page 54 and “Data pane – icons, special text, and buttons” on page 115.

My Entity List pane – Shows the gene lists and tracks stored in the program. To view the names of gene lists or tracks available in the program, double-click the names of folders to expand or collapse them, or click the  or  buttons.

To find specific content items in the Navigator

At the top of the Navigator is a search pane that can help you find specific content items. See “Search pane” on page 112.







Figure 15 Navigator search pane

- 1 Type a search term in the box at the top of the Navigator. The search term is not case-sensitive, but it must reflect the entire name of the content item that you want to find. You can use asterisks (*) as wildcards to represent a group of unspecified characters. For example, if you type *1234*, the search will find all items that contain “1234” in the name.
- 2 By default, the program searches all panes of the Navigator. To limit your search to a specific pane, click . In the list that appears, select the desired pane.
- 3 Click .

The program searches the selected pane(s). If it finds item(s) that match your search term, it expands folders so that the items are visible, and highlights them in red. You may need to scroll down to see all the search results.
- 4 To clear the results of a search, click .

To display the properties of a specific design

Design properties include general information about a design, such as its name, application type, and associated species. They also include a list of the names and chromosomal locations of probes.

- 1 Expand the folders of the Data pane until you can see the genome build folder(s) within the design folder.
- 2 Right-click the genome build folder, then click **Show Properties**.

The Design Properties dialog box appears. See “[Design Properties](#)” on page 158.

To update probe annotation in design files

Agilent regularly makes updates to probe annotations on its eArray Web portal. If you have imported Agilent array designs into Agilent Genomic Workbench, and you are a registered eArray user, you can download the updated design files from within Agilent Genomic Workbench. For more information about eArray, go to <https://earray.chem.agilent.com> and click **Help**.

- 1 In the Home tab, click **User Preferences**.
The User Preferences dialog box appears.
- 2 In the Miscellaneous tab, under **eArray User Details**, type your eArray **Username** and **Password**. See “[User Preferences](#)” on page 208.
- 3 Click **OK**.
- 4 Expand the folders of the Data pane until you can see the design that you want to update.
- 5 Right-click the design, then click **Update from eArray**. This option appears only for Agilent designs.
A confirmation dialog box appears.
- 6 Click **Yes**.

The program downloads an updated design, if one is available.

To rename an array in the Data pane

This topic describes how to rename an array in the Data pane, which changes the Global Display Name for the array. If you rename an array in this way, and subsequently add the array to an experiment, the array appears in the experiment with the new name. It also changes the array name in any experiment to which it is already linked. To rename an array only within the context of a specific experiment, see [“To rename an array in an experiment”](#) on page 59.

- 1 Expand the folders of the Data pane until you can see the array you want to rename.
- 2 Right-click the name of the array, then click **Rename**.
An Input dialog box appears.
- 3 Type a new name for the array, then click **OK**.
The program renames the array.

To remove data or design files from the program

You can delete array design and data files from the program when you are finished with them.

- 1 If an array that you want to delete is associated with an experiment, first delete it from the experiment. See [“To remove arrays from an experiment”](#) on page 60.
- 2 In the Data pane, expand folders until you can see the design folder or array that you want to delete.
- 3 Do one of the following:
 - For array data files, click the name of the first array, then hold down the **Ctrl** key while you click the names of additional arrays within the same design.
 - For array design folders, click the name of the first design folder, then hold down the **Ctrl** key while you click the names of additional ones. This selects the designs and all array data files within them for deletion.
- 4 Right-click the name of a selected design folder or array data file, then click **Delete**.

2 Importing, Managing, and Exporting Data and Other Content

To create a gene list

A confirmation dialog box appears.

5 Click Yes.

The program deletes the selected files.

CAUTION

When you delete files, you permanently remove them from Agilent Genomic Workbench. To restore deleted files, you must import them again.

To create a gene list

When you create a gene list, you create a list of the genes in a contiguous chromosomal region that you define.

1 Follow these steps to define a chromosomal region for your gene list. If you know the exact start and end locations of the chromosomal region, skip to step 2.

a In Genome View, select the chromosome.

The selected chromosome appears in Chromosome View. See [“Chromosome View”](#) on page 126,

b In Chromosome View, in the plotting area to the right of the chromosome, drag the pointer over the chromosomal region of interest.

The program draws a blue box around the region, and displays the region in greater detail in Gene View.

c In Gene View, adjust the view so only the genes of interest appear. For a description of the adjustment commands available in Gene View, see [“Gene View”](#) on page 128.

2 Right-click anywhere within the log ratio plotting area in Gene View, then click **Create Gene List**.

The Create Gene List dialog box appears. See [“Create Gene List”](#) on page 154.

3 In the dialog box set the Name, Description and Color.

4 In the dialog box select the chromosomal region for the new gene list.

5 Click OK.

The new gene list appears in the Gene List folder of My Entity List in the Navigator.

To display the genes in a gene list

You can display the genes in a gene list as a table.

- 1 Expand the folders in the **My Entity List** pane until you can see the gene list.
- 2 Right-click the gene list, then click **View In Table**.

The Gene List dialog box appears, with a table that contains the names of the genes in the gene list. You can also use this dialog box to edit the description of the gene list and its display color. See “[Gene List](#)” on page 305.

You can also create gene lists. For more information, see “[To create a gene list](#)” on page 66.

To rename a gene list

The name of a gene list identifies it within the Gene List folder of the My Entity List pane. You can rename gene lists.

- 1 Expand the folders of the **My Entity List** pane until you can see the gene list to rename.
- 2 Right-click the gene list, then click **Rename**.
- 3 Type a new name for the gene list, then click **OK**.

To delete gene list(s)

- 1 In the **My Entity List** pane of the Navigator, click to expand the **Gene List** folder.
- 2 Click the name of a gene list to delete. Hold down the **Ctrl** key while you click the names of additional gene lists.
- 3 Right-click one of the selected gene lists, then click **Delete**.
A confirmation dialog box appears.
- 4 Click **Yes**.

To create a track (CGH only)

When you create a track, you create a list of the genes in a contiguous chromosomal region that you define. To create a list of genes or other annotations, such as CNV or miRNA, in multiple regions, create additional tracks, and combine them.

- 1 Follow these steps to define a chromosomal region for your track. If you know the exact start and end locations of the chromosomal region, skip to step 2.
 - a In Genome View, select the chromosome.
The selected chromosome appears in Chromosome View.
 - b In Chromosome View, in the plot area to the right of the chromosome, drag the pointer over the approximate chromosomal region of interest.
The program draws a blue box around the region, and displays the region in greater detail in Gene View.
 - c In Gene View, adjust the view so only the genes of interest appear. For a description of the adjustment commands available in Gene View, see [“Gene View”](#) on page 128.
- 2 Right-click anywhere within the plot area in Gene View, then click **Create Track**.
The Create Track dialog box appears. See [“Create Track”](#) on page 156.
- 3 In the dialog box set the Name, Description and Color.
- 4 In the dialog box select the chromosomal region for the new track.
- 5 Click **OK**.
The new track appears in the Tracks folder of My Entity List pane in the Navigator.

To display the details of a track

You can display a table that contains the values for a list of track attributes.

- 1 In **My Entity List** pane, expand the Tracks folder to see the track.
- 2 Right-click the name of the track, then click **View Details**.

Track data appears in a Track table. See “[Track](#)” on page 203.

To rename a track

The name of a track identifies it both within the Tracks folder of the My Entity List pane, and in Gene View when you select **Show In UI** for the track. You can rename tracks.

- 1 Expand the folders of the My Entity List pane until you can see the track to rename.
- 2 Right-click the track, then click **Rename**.
- 3 Type a new name for the track, then click **OK**.

To delete tracks

- 1 In the My Entity List pane of the Navigator, expand the Tracks folder.
- 2 Click the name of a track to delete. Hold down the **Ctrl** key while you click the names of additional tracks.
- 3 Right-click one of the selected tracks, then click **Delete**.
A confirmation dialog box appears.
- 4 Click **Yes**.

Exporting and Saving Content

This section describes how to export several kinds of files from the program.

To export array attributes

You can export selected array attributes for any imported arrays. You first select the arrays and then the attributes for the selected arrays. You can export array attributes from the Home tab or from the short-cut menu for an experiment.

- 1 Click **Home > Export > Array Attributes**.

OR

In the Experiment pane of the Navigator, right-click an experiment of interest, and click **Export Attributes**.

The Export Array Attributes dialog box appears with the Array tab displayed. See [“Export Array Attributes”](#) on page 170.

If you opened this dialog box by right-clicking an experiment, only those arrays selected for the experiment appear in the Selected Array List. You can add or remove attributes from the list.

- 2 Under **Select Design**, select the design file and genome build for the arrays you want to add.

The arrays for the selected design file and genome build appear in Array List.

- 3 In **Array List**, select the arrays whose attributes you intend to export. To select a single array, click its name. To select additional arrays, hold down the **Ctrl** key while you click their names.


- 4 Click .

The program moves the selected arrays to the Selected Array List.

- 5 Click **Next** to select attributes for the selected arrays.

The Export Array Attributes dialog box appears with the Attribute tab displayed. See [“Export Array Attributes”](#) on page 170.

All of the attributes for the arrays are already located in the Selected Attribute List.

- 6 Move any attributes you don't want to export to the Available List.
 - a In the Selected Attributes List, highlight those attributes you do not want to export. To select additional attributes, hold down the **Ctrl** key while you click their names.
 - b Click .
- 7 Click **OK**.

The Export dialog box appears. See “Export” on page 169.
- 8 Select the folder in which to save the attributes, and click **Export**.

The attributes will be saved to the selected folder as a .txt file.

To export experiments

You can export experiments as a ZIP file to transfer them to another computer. Exported experiments contain the associated design and array data files, only. The program does not export information about array selections, or any analysis parameters or results.

- 1 In the Home command ribbon, click **Export > Experiments**.

The Export Experiments dialog box appears. See “Export Experiments” on page 174.
- 2 Mark the experiments that you want to export. To export all experiments, click **Select All**.
- 3 Click **OK**.

An Export dialog box appears. See “Export” on page 169.
- 4 Select a location and type a name for the exported ZIP file.
- 5 Click **Export**.

The program exports all selected experiment(s) together as a single ZIP file.

To export filters

You can export selected array, feature, design, metric, and aberration filters that are available in some data analysis applications in Agilent Genomic Workbench. The program exports all selected filters as a single *.xml file that you can import again at a later time.

- 1 In the **Home** tab, click **Export > Filters**.

The Export Filters dialog box appears. See “[Export Filters](#)” on page 175.

- 2 Under **Export**, mark the check boxes beside the filter(s) to export. To select all filters for export, click **Select All**.

- 3 Click **OK**.

An Export dialog box appears.

- 4 Select a location and type a name for the exported file, then click **Export**.

The program exports all selected filters as a single *.xml file.

To export tracks

You can export selected tracks as a BED format track file. You can then import this file into Agilent Genomic Workbench on another computer, or into a genome browser that accepts BED format files.

- 1 In the **Home** tab, click **Export > Tracks**.

The Export Tracks dialog box appears. See “[Export Tracks](#)” on page 176.

- 2 Mark the tracks to export. To select all tracks for export, click **Select All**.

- 3 Click **OK**.

An Export dialog box appears.

- 4 Select a location and type a name for the exported track file, then click **Export**.

The program exports the track(s) as a single BED format track file.

To copy what you see in the main window

You can copy panes of the main window to the Clipboard as images, and then paste them into a new document in another program (such as Microsoft® Word, or PowerPoint). The images contain only what actually appears on your screen; regions to which you must scroll are not included.

- 1 In the **View** tab, click **Copy**.
- 2 In the shortcut menu that appears, click the name of the pane that you want to copy. You can copy any view, or the Navigator. To copy all of the panes, click **All**.

The program copies the selected pane(s) to the clipboard.

- 3 Open a document in a program that accepts images. In that program, click **Edit > Paste**, or the appropriate paste command.

NOTE

To adjust how data is displayed in the panes use the View Preferences dialog box. For example, you can turn on or off the cursor. See [“View Preferences”](#) on page 216 for more information.

To copy the list of array colors for an experiment

You can copy the list of arrays in an experiment, and the colors assigned to them, to the clipboard as an image. You then paste the image into a document in another program such as Microsoft® Word or PowerPoint.

- 1 In the **Experiment** pane, expand the **Experiments** folder.
- 2 Right-click the name of the experiment, then click **Edit Array Color**.
The Edit Array Color dialog box appears. See [“Edit Array Color”](#) on page 165.

- 3 In the dialog box, click **Edit > Copy**.

The program copies the names of the arrays and their colors to the clipboard as an image.

- 4 Open a program that accepts images. Click **Edit > Paste**, or the appropriate paste command for the specific program.


2 Importing, Managing, and Exporting Data and Other Content

To save data and design information from an experiment

To save data and design information from an experiment

You can save the data and design information from a single design in an experiment as a tab-delimited text file.

- 1 In the **Experiment** pane, expand the **Experiments** folder until you see the genome build(s) for the design you want to export.
- 2 Right-click the name of the genome build, then click **Save As Text File**. A dialog box appears.
- 3 Select a location and type a name for the file, then click **Save**.



3 Displaying Data and Other Content

Selecting an Experiment for Displaying Data [76](#)

Displaying Array Data [80](#)

Displaying Content (Gene Lists/Tracks) [89](#)

Searching for Probe and Gene Information [95](#)

This chapter shows you how to display log ratio data from imported feature extraction data files, as well as gene and track content, in the Genomic Viewer. It also gives you instructions on how to customize the display of data and content to meet your needs.



Selecting an Experiment for Displaying Data

An experiment is a set of links to microarray data and design files, and any associated results. Experiments are displayed in the Experiment pane of the Navigator that appears for applicable tabs. The Experiment pane does not appear if you select the miRNA or Expression modules. See [“Using the Navigator to Search for Data”](#) on page 22.

When you select an experiment and have no CGH, ChIP, or CH3 application license, the program shows the log ratio data of selected arrays in the active experiment, if display of the data is enabled in View Preferences. See [“View Preferences”](#) on page 216 for more information.

NOTE

Without an application license (CGH, ChIP, or CH3) you cannot select an experiment that contains results.

When you select an experiment and Preprocessing and Analysis options have been turned on or set to apply, the program automatically begins the analysis of the selected array data with current settings and displays its results.

This section describes how select an experiment to make it active and select or deselect arrays for further display.

To select an experiment

When you select an experiment, the program displays log ratio data in a scatter plot, if that option is enabled.

- 1 If necessary, do one of the following to add the desired experiment to the Experiment Pane in the Navigator:
 - Create a new experiment and add data to it. See [“To create a new experiment”](#) on page 55.
 - Import a saved experiment file. See [“To import an experiment file”](#) on page 52.
- 2 In the Navigator, double-click the name of the experiment.
The Experiment Selection dialog box appears.

3 Click Yes.

In the Experiment pane of the Navigator, the name of the experiment turns blue. The name also appears in the title bar of the main window. Tables of data and design information appear in Tab View. For more information on the available tabs, see “[Tab View](#)” on page 133.

To select or deselect arrays in the experiment

To include arrays for display, you select them from the arrays available, either in an inactive experiment or the active one. When you first create an experiment, the program automatically sets the first array in the experiment for analysis. If you do not select additional arrays for display, only the first one will be shown when the experiment is selected.

To select the arrays for display before experiment selection:

- 1 Hold down the **Shift** key to highlight contiguous arrays or hold down the **Ctrl** key to highlight noncontiguous arrays.
- 2 Right-click the highlighted arrays, and click **Select**.

Even though the selected arrays do not change color, they will change color after the experiment is selected.

In the Navigator, the color of an array’s icon has the following meaning, after experiment selection:



Array is not selected.



Array is selected. The specific color matches the color of the column headings for the array in Tab View in the lower part of the window. In addition, the program displays aberration results and moving averages related to this array in this color. To configure a custom color for the array, see “[To change the display color of an array](#)” on page 78.

To select or deselect arrays in a *selected* experiment:

- 1 In the Navigator, expand the folders of the selected experiment.
- 2 Click the name of an array you want to include in the display.

3 Displaying Data and Other Content

To change the display color of an array

To include additional arrays, hold down the **Ctrl** key while you click their names. To include a contiguous block of arrays, click the name of the first array in the block, then hold down the **Shift** key while you click the name of the last one.

- 3 Right-click the name of one of the highlighted arrays, then click **Select**. After you select the arrays, the program reanalyzes the data set within the experiment and displays the data in Genome, Chromosome, and Gene Views. You can see the data for just the selected arrays in the Selected Arrays tab in Tab View.

To customize the appearance of the scatter plot in Genome, Chromosome, and Gene Views, see [“To change scatter plot appearance”](#) on page 82.

You can also use the headings of columns in Tab View that contain array data to select and deselect arrays.

- Click a column heading to select that array only.
- Hold down the **Ctrl** key while you click a column heading to select or deselect an array without changing the status of other arrays.
- Right-click a column heading to open a shortcut menu with options that let you select or deselect that array, or all arrays.

For more information on Tab View, see [“Tab View”](#) on page 133.

To change the display color of an array

The color assigned to an array sets the color of its icon when you select the array within an experiment. It also changes the colored square in the array’s column heading in Tab View.

- 1 In the Experiment pane of the Navigator, in the **Experiments** folder, expand the folder of an experiment until you can see the array of interest.

- 2 Right-click the desired array, then click **Edit Array Color**.

The Select Color dialog box appears. The dialog box gives three different ways to select the desired color. [“Select Color”](#) on page 196.

- 3 Select the desired color in one of the following ways:

Dialog box tab	Instructions
Swatches	<ul style="list-style-type: none"> Click the desired color swatch.
HSB (Hue/Saturation/Brightness)	<p>Type or adjust the values in H (Hue), S (Saturation), and B (Brightness), or alternately, follow these steps:</p> <ol style="list-style-type: none"> Select H, then drag the slider to select a hue based on the color strip to its right. Click an appropriate location in the large color box to the left of the slider to set the saturation and brightness levels of the color. Both the HSB and equivalent RGB values of the color appear in the dialog box. Note these values—they will be useful if you need to use this color in the future.
RGB (Red/Green/Blue)	<p>Do any of the following. Note the final RGB Values; they will be useful if you need to use this color in the future.</p> <ul style="list-style-type: none"> Drag the Red, Green, and Blue sliders. Type or adjust values in the boxes to the right of the sliders.

Samples of the color in different contexts appear under Preview. The upper half of the right-most color sample shows the original color for comparison.

4 Adjust the color as desired, then click **OK**.

You can also manage all of the colors for all of the arrays in an experiment. Right-click the desired experiment, then click **Edit Array Color**. For more information, see “[Edit Array Color](#)” on page 165.

Displaying Array Data

After you select an experiment, you can change how data appear within the Views or change the appearance of the Views that contain the data (or results).

To display the scatter plots

By default, display of scatter plots is turned On. If you do not see the scatter plot(s), do one of the following:

- 1 From the View tab, click **View Preferences**. See “[View Preferences](#)” on page 216 for more information.
- 2 In the View Preferences dialog box, under Data Visibility, select **All views** and then mark the box next to **Scatter Plot**.


OR

- 1 Right-click in any of the views, and select **View Preferences**. See “[View Preferences](#)” on page 216 for more information.
- 2 In the View Preferences dialog box, under Data Visibility, select **All views** and then mark the box next to **Scatter Plot**.

To show or hide data in scatter plots

- 1 In the Gene View, move the mouse pointer over the down arrow in **Scatter Plot** until the Scatter Plot box appears, and do any of the following:

To do this	Follow these steps
Show or hide data points for a selected data type	<ul style="list-style-type: none">• To show data points – Mark one or both check boxes under Configure Coloring schemes; then select the type of data from the Color by list.• To hide all data points – clear the check boxes.

- 2 Click  to close the Scatter Plot window.

To customize scatter plot ranges and colors

You can customize the display of scatter plot data. For each data type (log ratio, probe score, intensity) you can set custom ranges and colors for the display. For channels, you can set custom colors only.

NOTE

The View Preferences dialog box contents changes depending on the application type that is selected (CGH, ChIP, CH3).

Add and customize a plot

- 1 In Gene View, move the mouse pointer over **Scatter Plot** to display the options.

OR

Right-click in any of the views, and select **View Preferences**.

- 2 Mark the one or both of the check boxes under Configure Coloring schemes.

- 3 Select a data type from the list.

- 4 Click **Configure Color and Ranges**.

The Configure Coloring Ranges and Shades dialog box appears where you set ranges and colors for any of the data types. For more information, see [“Configure Coloring Ranges and Shades”](#) on page 150.

- 5 In the Configure Coloring Ranges and Shades dialog box, click one of the tabs and then select the data type to configure.
- 6 Type minimum and maximum numbers to define a range for the data type.
- 7 Click **Color** to open the Select Color dialog box. Use the tabs to select a color for the range. See [“Select Color”](#) on page 196 for more information.
- 8 Click **OK** to close the Select Color dialog box and return to the Configure Coloring Ranges and Shades dialog box.
- 9 Click **Add Range** to add the custom range to the range list.
- 10 When you are done, click **OK** to close the dialog box.

Edit or remove a range

- 1 In the Configure Coloring Ranges and Shades dialog box, click one of the tabs and then select the data type to configure.
- 2 In the range list, mark the **Edit/Delete** box to select the range. You can mark more than one range.
- 3 Click **Edit Range** to change the minimum and maximum values, or to change the color for the selected range.
- 4 Click **Delete Range** to delete the selected range.
- 5 Click **OK** to close the dialog box.

To change scatter plot appearance

You use the View Preferences dialog box to change the appearance of the scatter plots in Chromosome and Gene views.

- 1 In the Genomic Viewer, right-click in the Gene View or Chromosome View, and then click **View Preferences**.

Or, click the View tab, and then click **View Preferences**.

The View Preferences dialog box appears. See “[View Preferences](#)” on page 216.

- 2 Do any of the following:

To do this	Follow these steps
Show or hide the scatter plot	<ol style="list-style-type: none">a In the View tab under Data Visibility, in View, select All Views.b Do one of the following: To show the scatter plot, mark Scatter Plot. To hide the scatter plot, clear Scatter Plot.c Click OK.
Change the symbol that appears for data points	<p>You can select the symbol separately for each design type.</p> <ol style="list-style-type: none">a In the View tab, under Rendering Patterns, select the desired Design type.b Under Styles, select the desired symbol.c Click Apply.

To do this	Follow these steps
Show a separate scatter plot in Gene and Chromosome Views for each selected array	<p>a In the View tab, under View Alignment, under Rendering Style, select Stacked.</p> <p>b Click Apply.</p>
Show one scatter plot that contains data for selected arrays	<p>a In the View tab, under View Alignment, under Rendering Style, select Overlaid.</p> <p>b Click Apply.</p>
Enable ToolTips for the scatter plot in Gene View	<p>ToolTips show information about an individual data point when you place the pointer over it.</p> <p>a Click the View tab.</p> <p>b Under Data Visibility, in View, select Gene View.</p> <p>c Mark Scatter Tool Tip.</p> <p>d Click Apply.</p>

3 Click **OK**.

To print the scatter plot

You can print the scatter plot as it appears in Genome, Chromosome, and Gene views. Each view selected in the analysis is printed on a separate page. Chromosomes and genes appear on the printed pages, but tracks do not.

- 1 In the Home tab, click **Print**.
- 2 Set print options, as desired, then click **OK**.

To create custom scales for Views

You can customize the scale used for display in the Chromosome View and Gene View. Custom scales are applied to both views.





- 1 Click the View tab and then click **View Preferences**.
- 2 In the View Preferences dialog box, under Configure Scales, mark the box next to **Apply** for the plot for which you want to create a custom scale.

3 **Displaying Data and Other Content**
To locate and display data within the Views

In Range, enter a value to use for the range. The range you enter changes the scale for the display of the selected data.

To locate and display data within the Views

To look through the data of the selected arrays, do any of the following. In general, all views are synchronized; if you select a location or region in one view, the other views move there as well.

To do this	Follow these steps
Select a specific chromosome to display	<ul style="list-style-type: none">In Genome View, click the desired chromosome. All other views switch to the selected chromosome.
Display data in a region of the selected chromosome	<ul style="list-style-type: none">In Chromosome View, drag the pointer over the desired region. Gene View expands (or shrinks) to show only the selected region. Tab View scrolls to the new cursor location.
Zoom in and out in Gene View	<ul style="list-style-type: none">Click  to zoom in.Click  to zoom out.
Scroll through the selected chromosome	<ul style="list-style-type: none">Click  to scroll up.Click  to scroll down. <p>Note: These arrows will appear side by side for horizontal orientation.</p>
Return Gene View or Chromosome view to center	<ul style="list-style-type: none">Click anywhere in Chromosome View, or anywhere within the scatter plot in Gene View. The location you click becomes the new cursor location.
Move all Views to a specific genomic location	<ul style="list-style-type: none">a Click Home > Go To Gene/Genomic location. A dialog box appears.b Under Genomic Location, select a Chromosome, and type a Base Position.c Click Go. All views move to the selected location.

To do this	Follow these steps
Display the location of a specific gene in the center of all Views	<p>a Click Home > Go To Gene/Genomic location. A dialog box appears.</p> <p>b Under RefSeq by Symbol, either select the desired gene (if available) or type the name of the gene.</p> <p>c Click Go. All views move to the location of the selected gene.</p>
Display the data selected in Tab View in the center of Chromosome and Gene Views	<ul style="list-style-type: none"> In Tab View, click any entry in any table, except a column heading. Chromosome and Gene views: The genetic location of the selected data appears in the center of Chromosome and Gene Views.
Scroll to a specific column in Tab View	<p>a In Tab View, right-click any column heading, then click Scroll To Column. The Scroll to Column dialog box appears. See “Scroll to Column” on page 194.</p> <p>b In Select Column, select the desired column.</p> <p>c Click OK.</p>
Search for a specific column entry in Tab View, and move the cursor there	<p>a In Tab View, right-click any entry except a column heading, then click Find in column. The Find in column dialog box appears. See “Find in column” on page 177.</p> <p>b Set the desired search parameters, then click Find Next. The program searches the column using your search parameters, and highlights the row of the first entry that matches. The cursor moves to the location defined in the highlighted row. This search is only for the selected chromosome.</p>
Display the exact chromosomal location of the cursor	At the bottom of the main window, look at the first cell of the Status bar. The location appears as the chromosome followed by the base position. For more information on the status bar, see “Status Bar” on page 138.

To smooth and plot CGH log ratio data

You use a plug-in program to create separate, stacked plots of smoothed log ratio data for each of the selected CGH arrays in the current experiment. The plug-in program can perform one of several varieties of moving-average-like smoothing, and plots the data associated with the currently selected chromosome.

The Plugin Settings command lets you change the parameters when you have selected to display the plot immediately after you click Plugin.

- 1 Select an experiment.
- 2 Select the arrays that contain the log ratio data to smooth and plot.
- 3 Select the chromosome that contains the log ratio data of interest.
- 4 Click **Tool > Plugin > CGHSmooth**.

The CGHSmooth Parameters dialog box appears. See “[CGHSmooth Parameters](#)” on page 141.

- 5 Type the parameter values for the selected arrays in the active experiment and the selected chromosome.
- 6 Click **OK**.

The CGHSmooth Plot appears, showing a plot of the log ratio vs. chromosomal position for each of the selected arrays. See “[CGHSmooth Plot](#)” on page 143.

- 7 (optional) To add a title, legend and change the appearance of the plot(s), right-click the plot, and click **Properties**.

See “[Chart Properties](#)” on page 145.

- 8 (optional) To see options to save, print, zoom in or out, or change the auto range of the plot, right-click the plot.

To enter parameters from the Plugin Settings command and display the plot directly from the Plugin command

- 1 After [step 5](#) above, mark **Don't Show Again**, then click **OK**.
The plot appears.
- 2 To change the parameter values, exit the plot, and click **Plugin Settings**.
- 3 Change the values you want to change, and click **OK**.
- 4 Click **Plugin**.

The plot appears.

- 5 (optional) To show the CGHSmooth Parameters dialog box again when you click Plugin, click **Plugin Settings**, clear **Don't Show Again**, and click **OK**.

To produce an echo example plot (CGH only)

The Echo Example Plot is the output of the Echo Example plug-in. It displays the log ratio data for the selected chromosome in the active experiment. Data from all of the selected arrays in the experiment appear as a series of stacked plots, one for each array.

- 1 Select an experiment.
- 2 Select the arrays that contain the log ratio data to smooth and plot.
- 3 Select the chromosome that contains the log ratio data of interest.
- 4 Click **Tool > Plugin > EchoExample**.
- 5 Type the parameter values for the selected arrays in the active experiment and the selected chromosome.
- 6 Click **OK**.

The Echo Example Plot appears, showing a plot of the log ratio vs. chromosomal position for each of the selected arrays. See “[Echo Example Plot](#)” on page 162.

- 7 (optional) To add a title, legend and change the appearance of the plot(s), right-click the plot, and click **Properties**.

See “[Chart Properties](#)” on page 145. To see options to save, print, zoom in or out, or change the auto range of the plot, right-click the plot.

To produce a moving average example plot (CGH only)

The MovAvgExample plug-in program calculates a 10-point moving average of the data from the selected chromosome in the selected arrays of the active experiment. It displays stacked plots of moving averages for all the arrays, one plot per array.

3 Displaying Data and Other Content

To produce a moving average example plot (CGH only)

The plug-in program itself (**MovAvg Example.pl**, in the Plugins folder of the Agilent Genomic Workbench installation folder on your computer) is a short Perl program. It is a good example of how computed columns are processed. You must have Perl installed on your computer to use this plug-in.

- 1 Select an experiment.
- 2 Select the arrays that contain the log ratio data to smooth and plot.
- 3 Select the chromosome that contains the log ratio data of interest.
- 4 Click **Tool > Plugin > MovAvg Example**.
- 5 Type the parameter values for the selected arrays in the active experiment and the selected chromosome.
- 6 Click **OK**.

The MovAvg Example Plot appears, showing a plot of the log ratio vs. chromosomal position for each of the selected arrays. See “[MovAvg Example Plot](#)” on page 192.

- 7 (optional) To add a title, legend and change the appearance of the plot(s), right-click the plot, and click **Properties**.

See “[Chart Properties](#)” on page 145. To see options to save, print, zoom in or out, or change the auto range of the plot, right-click the plot.

To enter parameters from the Plugin Settings command and display the plot directly from the Plugin command

- 1 After [step 5](#) above, mark **Don't Show Again**, then click **OK**.

The plot appears.

- 2 To change the parameter values, exit the plot, and click **Plugin Settings**.
- 3 Change the values you want to change, and click **OK**.
- 4 Click **Plugin**.

The plot appears.

To show the MovAvg Example Parameters dialog box again when you click Plugin, click **Plugin Settings**, clear **Don't Show Again**, and click **OK**.

Displaying Content (Gene Lists/Tracks)

To show gene lists in Gene View

A gene list defines a set of genes of interest.

You cannot show gene lists without a license. With a license you can highlight the genes in the gene list in Gene View, or limit the display of data, genes, and tracks to the regions defined by a gene list.

You also cannot import or export a gene list without a license, but you can create a gene list in the program. See [“To create a gene list”](#) on page 66.

To change the appearance of genes in Gene View

- You use the User Preferences dialog box to change the appearance of the genes in Chromosome and Gene views.
- 1
Right-click any part of the Gene View, then click **User Preferences**.
The User Preferences dialog box appears.
- 2
Click **Tracks**.
See [“User Preferences”](#) on page 208.
- 3
Do any of the following:

To do this	Follow these steps
Show or hide genes in Gene View	<div> <div>a</div> Under Visualization Parameters: To show genes – Under Genes, mark Show Gene Symbols. To hide genes – Under Genes, clear Show Gene Symbols. </div> <div> b Click Apply. </div>
Change the display font for genes (and track annotations) in Gene View	<div> <div>a</div> In the Gene Symbols tab, under Font, select a new Font, Font Style, and Font Size. </div> <div> b Click Apply </div>

3 Displaying Data and Other Content

To show tracks in Gene View

To do this	Follow these steps
Change the display angle for genes (and track annotations) in Gene View	<p>a Under Visualization Parameters, under Genes, in Orientation (Degrees), type a new orientation in degrees. 0° is horizontal.</p> <p>b Click Apply.</p>

4 Click **OK**.

To show tracks in Gene View

Tracks contain information for specific genomic locations. A multitude of tracks from diverse sources is available for many species. You can display tracks next to genes and microarray data in Gene View.

- 1** Select and show microarray data. See [“To select an experiment”](#) on page 76.
- 2** In the My Entity List pane, open the Tracks folder.
- 3** Right-click the track you want to display, and click **Show In UI**.

Or, you can do this:

- 1** In Gene View, right-click anywhere within the scatter plot, then click **User Preferences**.
The User Preferences dialog box appears. See [“User Preferences”](#) on page 208.
- 2** Click **Tracks**.
- 3** Mark the **Show In UI** check box of each desired track.
- 4** Click **OK**.

The program displays the selected tracks in Gene View.

To change the appearance of tracks

Within the Tracks tab of the User Preferences dialog box, you can change the appearance of tracks, as described in the table below.

To do this	Follow these steps
Include track information in reports	<p>a In the list of tracks, in the Show in Report column, mark the check boxes of the desired tracks.</p> <p>b Click Apply.</p> <p>Doing this adds a column with the hits from the track file. For each aberrant interval, it reports the entries from the track file for that interval in that separate column.</p>
Show or hide annotations in all tracks	<ul style="list-style-type: none"> To show annotations in all tracks: under Tracks, mark Show Annotations. To hide annotations in all tracks: under Tracks, clear Show Annotations.
Display all selected tracks as a single track	<ul style="list-style-type: none"> Under Tracks, mark Show Overlaid. The program combines the annotations of all selected tracks into a single track named Overlaid Track. To show tracks individually again, clear Show Overlaid.
Display the parameters and the list of annotations of a track	<ul style="list-style-type: none"> In the list of tracks, for the desired track, click Details.
Change the display font for track annotations (and genes)	<p>a Under Font, select a new Font, Font Style, and Font Size for track annotations.</p> <p>b Click Apply.</p> <p>The program changes the display font of track annotations and genes in Gene View.</p>

3 Displaying Data and Other Content
To display tracks in UCSC Browser

To do this	Follow these steps
Change the order in which tracks appear in Gene View.	<p>The order of tracks in the Gene Symbols tab controls the left-to-right order of tracks in Gene View.</p> <p>a Click the name of the track you want to move.</p> <p>b Do one of the following:</p> <ul style="list-style-type: none">• To move the track up in the list of tracks (and farther left in Gene View), click its name, then click Up.• To move the track down in the list of tracks (and farther right in Gene View), click its name, then click Down. <p>c Click Apply.</p>
Change the display angle of track annotations (and genes)	<ul style="list-style-type: none">• Under Genes, in Orientation, type a new orientation (in degrees). 0° is horizontal. <p>The program changes the display angle of track annotations and genes in Gene View.</p>

To display tracks in UCSC Browser

- 1 Right-click Gene View, and click **Show in UCSC**.
The View coordinates in UCSC browser dialog box appears. See [“View coordinates in UCSC browser”](#) on page 214.
- 2 Complete the dialog box with the track parameters, and click **OK**.
The UCSC Browser appears, if you are connected to the Internet.

To change the graphical display to a different genome build




Figure 16 Track displayed in UCSC browser

3 Follow the instructions on the Web site for what you want to do.

To change the graphical display to a different genome build

The default graphical display for Genome, Chromosome and Gene Views represents human genome build 18.

- To change the graphical display to a different genome build, select an experiment whose data are based on a design file of a different genome build.

The display automatically changes when you select an experiment that contains a design file with a different genome build, such as human genome build 17, or a mouse or rat genome build.

If a genome build is not available for the design file you import, you must import the genome build first. See “To import a genome build” on page 50.

3 Displaying Data and Other Content

To change the graphical display to a different genome build

Searching for Probe and Gene Information

To search Tab View for specific probe information

You can find a specific entry in a column of a data table in Tab View. For more information on Tab View, see “[Tab View](#)” on page 133.

- 1 In Tab View, right-click anywhere in the column you want to search, then click **Find in column**. See “[Find in column](#)” on page 177.

The Find in column dialog box appears. The column to be searched also appears in the title bar of the dialog box.

NOTE

The Find in column function works within the selected chromosome.

- 2 Set the search parameters, as described below.

Parameter	Comments/Instructions
Find in column	<ul style="list-style-type: none">Type the text you want to find (the <i>search term</i>). This can be an entire entry, or part of one.
Direction	<ul style="list-style-type: none">Select one of these options:<ul style="list-style-type: none">Up – Search the column upwards from the current cursor location (the highlighted row of the table).Down – Search the column downwards from the current cursor location (the highlighted row of the table). <p>Tip: Click a row in Tab View to highlight it.</p>
Conditions	<ul style="list-style-type: none">Mark any of these, as desired:<ul style="list-style-type: none">Match Case – Find entries that match upper and lower case characters in the search term.Match whole word – Find an entry only if the entire entry matches the search term.

- 3 Click **Find Next**.

If the program finds a match, it highlights the row that contains the matching entry, and resets the cursor to the corresponding position. You can click **Find Next** as many times as you like, and the program

3 Displaying Data and Other Content

To search Tab View for specific probe information

continues to search for additional matching entries in the column. If it finds no match, the message: **String not found** appears in black in the lower part of the dialog box.

- 4 When you complete your search, click **Cancel**.

To search Agilent eArray for probe information

You can use the chromosomal region that appears in Gene View, or another chromosomal region as the basis for a probe search on the Agilent eArray Web site. eArray is a powerful microarray design system for CGH, ChIP and gene expression applications. It contains a massive database of validated, annotated probes, and a full complement of tools for custom microarray design.

Before you can search for probes in eArray, you must be a registered eArray user. For more information, go to eArray.chem.agilent.com. You must also provide your eArray user name and password in the Miscellaneous tab of the User Preferences dialog box. See “[User Preferences](#)” on page 208.

- 1 In Gene View, right-click anywhere in the plotting area, then click **Search probes in eArray**.

The Search probes in eArray dialog box appears. See “[Search probes in eArray](#)” on page 195.

- 2 Do one of the following to define the chromosomal region for your search:
 - To set the region to the one that currently appears in Gene View, select **For complete gene view**.
 - To set the region numerically, select **User Defined**, then select a **Chromosome** and type **Start** and **Stop** locations for the desired region.

- 3 Click **OK**.

The eArray Web portal opens in your internet browser.

To search the Web for information on probes in Tab View

You can use any entry in a table in Tab View as the basis for a Web search.

- 1 In Tab View, right-click any data table entry other than a column heading.
- 2 Click one of the available sites.

If the site of interest does not appear in the shortcut menu, you can create a custom search link. See [“To create a custom Web search link”](#) below.

The selected site opens in your Internet browser. The program sends the table entry to the site as a search string.

To create a custom Web search link

If you need to search a different database or site based on data table entries, you can create your own custom search link. When you right-click a table entry in Tab View, a shortcut menu opens, and your custom link appears in it. If you select this link, Agilent Genomic Workbench opens the site in your Web browser and sends the table entry to the site as a search string.

- 1 Right-click any data table entry in Tab View, except a column heading, then click **Customize Link**.

The Customize Search link dialog box appears. [“Customize Search Link”](#) on page 157.

- 2 Click **New**.
- 3 In the Input dialog box, in **URL name**, type a name for the link.
This name will appear in the shortcut menu that opens when you right-click a data table entry.
- 4 Click **OK**.
- 5 In **URL**, type the complete URL needed to send a search string to the site. Use <target> as the query string value.
For example, this URL sends selected table entries to Google.com:
`http://www.google.com/search?hl=eng&q=<target>`
- 6 Click **Update**, then click **Yes**.

To update or delete a custom Web search link

- 1 Right-click any data table entry in Tab View other than a column heading, then click **Customize Link**.
The Customize Search link dialog box appears.
- 2 In **URL Name**, select the custom search link to update or delete.
- 3 Do one of the following:

To do this	Follow these steps
Update a Web search link	<div><div>a Edit the URL name and the URL as needed.</div><div>b Click Update. A Confirm dialog box appears.</div><div>c Click Yes.</div></div>
Delete a Web search link	<div><div>• Click Delete.</div></div>

- 4 Click **Close**.

3 Displaying Data and Other Content

To update or delete a custom Web search link



4 Data Viewing Reference

Agilent Genomic Workbench Main Window	102
Command Ribbons	103
Switch Application Menu	111
Search pane	112
Navigator Pane	114
Genomic Viewer	124
Status Bar	138
Dialog Boxes	139

This chapter describes the command ribbons, Navigator panes, and dialog boxes that can appear when you are using Agilent Genomic Workbench without analysis licenses.



Agilent Genomic Workbench Main Window

The sections that follow describe the main components of the Agilent Genomic Workbench main window – Switch Application Menu, the command ribbons, the Navigator and the Views. You use these to import, organize, manage, export and display data and other content. For descriptions of the dialog boxes for these elements, see “Dialog Boxes” on page 139. Figure 17 shows the main window of Agilent Genomic Workbench, and identifies its main parts.




Figure 17 Agilent Genomic Workbench Lite Edition - major components

Command Ribbons

When you click a tab at the top of the Agilent Genomic Workbench main window, groups of commands appear below the tab bar. This group of commands is called a command ribbon, and the commands that appear are available only for the selected tab. The tabs that are displayed change depending on what application is selected (such as CGH, ChIP, CH3). This section describes the ribbon commands used to import, manage, export and display data in Agilent Genomic Workbench. For command ribbons that appear in the Sample Manager and Workflow tabs, see the User Guides for those applications.




Figure 18 Tab bar and command ribbon for unlicensed CGH application

Home command ribbon



Figure 19 Command ribbon in the Home tab of Agilent Genomic Workbench

User Preferences Opens the User Preferences dialog box with the following tabs:

Tab	Description
Tracks	Opens a dialog box that lets you manage which tracks to display in Genomic Viewer and how they appear. See “Tracks tab” on page 209.

Tab	Description
Miscellaneous	Opens a dialog box where you can select a new location for your data files and set up access to the eArray web site. See “Miscellaneous tab” on page 211.
License	Opens a dialog box where you can add a CGH, ChIP, or CH3 application license, if you want to purchase one after using the unlicensed version. See “License tab” on page 212.

Import Opens a menu of file types that you can import:

Option	Description
Array Files	<p>Opens a menu with these options:</p> <ul style="list-style-type: none"> • FE File – Opens the Import FE Files dialog box, where you can select an Agilent Feature Extraction array data file to import. See “Import” on page 179 and “To import Agilent FE or Axon data files” on page 45. • Axon File – Opens the Import Axon Files dialog box, where you can select Axon (*.gpr) files for import. See “Import” on page 179 and “To import Agilent FE or Axon data files” on page 45. • UDF File – Opens the UDF Files dialog box, where you can select a Universal Data File (UDF) to import. See “Import” on page 179 and “To import a UDF file” on page 46.
Design Files	<p>Opens a menu with these options:</p> <ul style="list-style-type: none"> • GEML File – Opens the Import Design Files dialog box, where you can select Agilent GEML-based (*.xml) array design files for import. See “Import” on page 179 and “To import Agilent GEML design files” on page 43. • Axon Design File – Opens the Import Axon Design Files dialog box, where you can select Axon (*.gal) array design files for import. See “Import” on page 179 and “To import Axon design files” on page 44.
Genome Build	Opens the Import Genome Build dialog box, where you can import Agilent-supplied genome build files. See “Import Genome Build” on page 185 and “To import a genome build” on page 50.
Track	Opens the Import Track dialog box, where you can select a BED format track file for import, and create a display name for the track. See “Import Track” on page 186 and “To import tracks” on page 51.

Option	Description
ArrayAttributes	Opens the Import microarray attributes dialog box, where you select a sample attributes file to import. See “Import” on page 179 and “To import array attributes” on page 51 for more information.
Experiments	Opens the Import Experiments dialog box, where you select an exported experiment .zip file, from which you can select experiments to import. See “Import” on page 179 and “To import an experiment file” on page 52 for more information.
Filters	Opens the Import dialog box, where you select a filter file to import. For more information, see “Import” on page 179 and “To import filters” on page 53.

Export Opens a menu that lets you export several kinds of files.

Option	Description
Experiments	Opens the Export Experiments dialog box, where you can select one or more experiments for export as a single ZIP file. See “Export Experiments” on page 174 and “To export experiments” on page 71.
Filters	Opens the Export Filters dialog box, where you can select one or more filters for export as a single *.xml file. See “Export Filters” on page 175 and “To export filters” on page 72.
Tracks	Opens the Export Tracks dialog box, where you can select one or more tracks to export as a single BED format file. See “Export Tracks” on page 176 and “To export tracks” on page 72.
Array Attributes	Opens the Export Array Attributes dialog box, where you can select arrays and their attributes for export. See “Export Array Attributes” on page 170.

Create Experiment Opens the Create Experiment dialog box, where you can create a new, empty experiment and add data to it. See [“Create Experiment”](#) on page 152 and [“To create a new experiment”](#) on page 55.

Save Experiment Result (Not available if you do not have a CGH, ChIP, or CH3 application license)

- Go to Gene/Genomic Location** Moves the cursor to the location in Chromosome and Gene Views that you select. See [“Go To Gene/Genomic Location”](#) on page 178.
- Print** Opens the Print window to print the display.
- Exit** Closes the program.

View Command Ribbon




Figure 20 View command ribbon for CGH application

- View Preferences** Opens the View Preferences dialog box where you can customize the display of data and results in the Genomic Viewer. For more information, see [“View Preferences”](#) on page 216.
- Copy** This command opens a menu with the options listed below. In general, the Copy command copies pane(s) of the main window to the Clipboard as an image. You can then paste the image into a document in another program. See [“To copy what you see in the main window”](#) on page 73.

Option	Description
All	Copies all panes of the main window to the Clipboard as an image.
Navigator	Copies only the Navigator to the Clipboard as an image.
Tab View	Copies only the Tab View to the Clipboard as an image.
SampleBySample view	(Available only in data analysis modules, when selected) Copies only the Sample-by-sample View to the Clipboard as an image.
Genome view	Copies only the Genome View to the Clipboard as an image.
Chromosome view	Copies only the Chromosome View to the Clipboard as an image.
Gene view	Copies only the Gene View to the Clipboard as an image.

Show Opens a menu with all available elements of the main window. Mark the check box for the one or ones you want to display.

View In Table

Signal Intensity Mark the check box to see the red and green raw signal intensities of the log ratio data in the Tab View.

Annotation Mark the check box to show annotations in the Tab View.

Cyto band info

View In Gene View Mark the check box to display cytobands in the Gene View.

NonUnique Probes

Highlight Nonunique probes in a microarray design have more than one mapping in the genome that is a perfect match. Because the probes represent the same sequence, the probe log ratio reflects a combination of log ratios from the redundant locations. Mark the check box to display nonunique probes in a different color.

Custom Data

Show Mark the check box to display custom data in the Genomic Viewer.

Tool command ribbon




Figure 21 Tool command ribbon

Plugin

Plugins are ancillary programs that operate on the selected array data in the active experiment in specific ways.

Opens a menu with the options described below. Custom plugins also appear in this menu.

CGHSmooth Opens the CGHSmooth Parameters dialog box. See [“CGHSmooth Parameters”](#) on page 141. You can set the parameters of the CGHSmooth plug-in, and create separate, stacked plots of smoothed log ratio data for each of the selected arrays in the current experiment. The plug-in plots the data associated with the currently selected chromosome.

When you open this dialog box, you see the default parameters you enter under Plugin Settings.

Echo Example Creates separate, stacked plots of log ratio data for each of the selected arrays in the current experiment. The plug-in plots the data associated with the currently selected chromosome. The plot appears in a new window. Although simple, this plug-in gives you a convenient way to view the log ratio data for selected arrays as separate plots. See [“Echo Example Plot”](#) on page 162.

MovAvg Example Opens the MovAvg Example Parameters dialog box. See [“MovAvg Example Parameters”](#) on page 190. You can set the parameters of the MovAvg Example plug-in. The plug-in calculates a 10-point moving average of each column of selected microarray data, and produces stacked plots of all of the input data and moving averages. To use this plug-in, you must have Perl installed on your computer.

When you open this dialog box, you see the default parameters you enter under Plugin Settings.

Plugin Settings

Opens another menu with these options:

CGHSmooth Opens the CGHSmooth Parameters dialog box, where you can set the parameters for the plug-in when you have selected to not show the parameters dialog box when you click Plugin. See [“CGHSmooth Parameters”](#) on page 141.

MovAvg Example Opens the MovAvg Example Parameters dialog box, where you can set the parameters for the plug-in when you have selected to not show the parameters dialog box when you click Plugin. See “[MovAvg Example Parameters](#)” on page 190.

Help command ribbon

The Help command ribbon lets you display the available Agilent Genomic Workbench help guides, and get information about software version. Help guides are opened in Adobe® Reader®.




Figure 22 Help command ribbon for unlicensed CGH application

Table 4 Table of Help for unlicensed version data viewing

Help Command	Action
Application Guide	Opens the Agilent Genomic Workbench application user guide for the selected application.
Sample Manager	Opens the <i>Sample Manager User Guide</i> , that shows how to use the Sample Manager module of Agilent Genomic Workbench to organize microarrays and edit their attributes. Sample Manager features are available if you have one or more DNA Analytics licenses.
Workflow	Opens the <i>Workflow User Guide</i> , that describes how to use the Workflow module of Agilent Genomic Workbench to extract image files with Agilent Feature Extraction software and/or analyze data using CGH and ChIP analysis software. Workflow features are available if you have a CGH and/or ChIP license.

Table 4 Table of Help for unlicensed version data viewing (continued)

Help Command	Action
Data Viewing	Opens the <i>Data Viewing User Guide</i> that describes how to import, organize, manage, export and display data and other content (experiments, gene lists, tracks) within Agilent Genomic Workbench. It is targeted for users who have no DNA Analytics application license(s).
About	Opens a message with information about the version number and copyright of the program.

An additional guide is available in the Open Application tab of the program. The *Agilent Genomic Workbench Product Overview Guide* gives an overview of the capabilities of Agilent Genomic Workbench. It also describes how to start each of the component programs and find help, and how to enter your license information. To open this guide, click the **Open Application** tab, then click **Product Overview**.

Switch Application Menu




Figure 23 Switch Application menu

The Switch Application menu lets you change to the other data display and analysis application types in Agilent Genomic Workbench. Select the desired application type.

CGH (Separate license required) Import, display, and analyze array-based comparative genomics hybridization (aCGH) data in both an interactive “analyze as you go” mode, and an automated workflow mode.

ChIP (Separate license required) Import, display, and analyze ChIP-on-Chip microarray data in both an interactive “analyze as you go” mode, and an automated workflow mode.

CH3 (Separate license required) Import and display data from microarray-based studies of genomic methylation patterns.

SureSelect Target Enrichment Use the Quality Analyzer function for SureSelect Target Enrichment. See the *Target Enrichment User Guide* for more information.

Search pane

The Search pane lets you find all occurrences of a specific search term in the Data, Experiment, and/or My Entity List panes. See [“To find specific content items in the Navigator”](#) on page 63. It also contains several buttons that you can use to move, hide, show or resize the Navigator.




Figure 24 Navigator – Search pane

Detach button Click to move the Navigator from the main window of the program and open it in a new, separate window.

Resize buttons Click to hide, show, or expand the Navigator.

Search term box The place where you type your desired search term. Search terms are not case-sensitive, but they must reflect the entire name of an array or other content item that you want to find. You can use asterisks (*) as wildcards to represent groups of unspecified characters. For example, a search term *25887* searches for any content that contains the string “25887”.

Pane list Lets you limit a search to a specific pane. Select the name of the desired pane from the list. To select all panes, select **All Panels**. By default, the program searches all panes.




Figure 25 Search Pane list



(Show Pane List button, available only if the Pane list is not visible) Makes the Pane list visible.



(Hide Pane List button, available only if the Pane list is visible) Hides the Pane list.



(Search button) Searches the pane(s) selected in the Pane list for all occurrences of the term you typed in the Search term box. If the program finds a matching item, it expands the folder structure to make the matching item(s) visible, makes the lettering of each item red and highlights the item in yellow. Note: The search term is not case-sensitive, but it must reflect the entire name of the desired items. You can use asterisks (*) as wildcards to represent groups of unspecified characters.

Scroll buttons

(Available only after a search) Lets you scroll up and down the lists of highlighted search items after a search.



(Clear button, available only after a search) Clears the search term from the Search term box, and resets the color of any matching item to its original color.

Navigator Pane

The Navigator contains several panes where you can look at program designs, experiments, data, or the status of tasks. Within each pane, you will see icons that tell you the status of the content. In addition, shortcut menus are available to let you perform tasks within the pane. These icons and shortcut menus are described in this section.


















Figure 26 Navigator panes

Data pane – icons, special text, and buttons

Item	Comments
	An unexpanded folder (domain) that contains subfolders or other items.
	An expanded folder. The items that it contains are visible in the Navigator.
	Expands a folder to show its contents.
	Collapses a folder to hide its contents.
	A methylation array design. This folder contains array data associated with the design, organized by genome build.
	A CGH array design. This folder contains array data associated with the design, organized by genome build.
	A gene expression array design. This folder contains array data associated with the design, organized by genome build.
	A ChIP array design. This folder contains array data associated with the design, organized by genome build.
	A genome build folder within a specific design folder. This folder contains arrays associated with the specific genome build and design.
	A single array data file.
	Data created from a multi-pack array.
text	An item that matches the search term in a search.
	(Dock out button) Moves the Data pane from the Navigator, and opens it in a, separate window.
	(Collapse button, available only if the Data pane is not collapsed) Collapses the Data pane, and shows its title bar at the bottom of the Navigator.
	(Expand button, available only if the Data pane is collapsed) Expands the Data pane.

Data pane – actions and shortcut menus

The Data pane of the Navigator shows available content items that are stored on your server for the selected application type, and any external content that you imported.

- Double-click any folder to expand or collapse it.

Data Folder

- Double-click any folder to expand or collapse it.
- Double-click a designs folder (ChIP, Expression, CGH, CH3) to display the imported designs for that data type.
- Double-click the name of a genome build folder to display imported arrays for that build.

Genome Build Folder

- Right-click the name of a genome build folder to display the following options:

Option	Description
Show Properties	Opens the Design Properties dialog box. See “Design Properties” on page 158.
Delete	Opens a Confirm dialog box. If you click Yes , the program permanently deletes all of the arrays in this genome build folder. (Not available for read-only builds.)

Specific Arrays










- Right-click the name of an array to display the following options:







Option	Description
Show Properties	Opens the Microarray Properties dialog box. See “Microarray Properties” on page 187 and “To display or edit the attribute values of a specific array” on page 60.

Option	Description
Rename	Opens an Input dialog box, where you can type a new name for the array. Click OK to rename the array. (Not available for read-only builds.)
Delete	Opens a Confirm dialog box. If you click Yes , the program permanently deletes the array. (Not available for read-only builds.)

- Drag an array from the Data pane to an experiment folder in the Experiment pane to associate it with an experiment. You can drag multiple arrays at once from one genome build in a design. Hold down the **Ctrl** key while you click the additional arrays to select them. You can also select a contiguous block of arrays; click the first array in the block, then hold down the **Shift** key and click the last one.

Experiment pane – icons, special text, and buttons

Item	Comments
	Click to expand a folder and display its contents.
	Click to collapse a folder and hide its contents.
	A folder that contains files or other folders.
	A methylation array design. This folder contains array data associated with the design, organized by genome build.
	A CGH array design. This folder contains array data associated with the design, organized by genome build.
	A gene expression array design. This folder contains array data associated with the design, organized by genome build.
	A ChIP array design. This folder contains array data associated with the design, organized by genome build.
	A genome build folder within a specific design folder. This folder contains arrays associated with the specific genome build and design.
	An array that is not selected for view

Item	Comments
	An array that is selected for view and analysis. The specific color of this icon can vary.
	An empty folder.
	Data created from a multi-pack array.
blue text	The currently active experiment. All data that appear in Chromosome, Gene, and Tab Views come from this experiment.
red text	An item that matches the search term in a search.
	(Dock out button) Moves the Experiment pane from the main window, and opens it in a separate window.
	(Collapse button, available only if the Experiment pane is not collapsed) Collapses the Experiment pane, and shows its title bar at the bottom of the Navigator.
	(Expand button, available only if the Experiment pane is collapsed) Expands the Experiment pane.

Experiment pane – actions and shortcut menus

Only the options that are not grayed out are discussed in this section. To enable the grayed-out options, you must have the appropriate license for the CGH, ChIP, or CH3 application you are using. These inactive options are explained in the *User Guide* for the application.

- In general, double-click the Experiments folder, and the folders within it, to expand and collapse them. Exception: double-click the name of an unselected experiment to select it for display. Without a license, you cannot select an experiment that contains results.

NOTE

The displayed options change depending on the user and status of the designs, builds, and arrays. You may not see all of the options that are described below.

Experiments Folder

- Right-click the **Experiments** folder to display the following options:

Option	Description
New Experiment	Opens the Create Experiment dialog box, where you can name the new experiment, and open another dialog box that lets you add microarray data to the experiment. See “Create Experiment” on page 152.
Export	Opens the Export Experiments dialog box, where you can export one or more experiments as a single ZIP file. See “Export Experiments” on page 174 and “To export experiments” on page 71.

Specific Experiment Folder

- Right-click the name of an experiment to display the following options:

Option	Description
Select Experiment	(Appears only if the experiment is not selected.) Opens the Experiment Selection dialog box, which asks if you want to select the experiment. Click Yes to select the experiment for display and analysis. Or In the Experiments folder, double-click the name of an experiment that is not selected to open the Experiment Selection dialog box. To select the experiment for analysis, click Yes .
Deselect Experiment	(Appears only if the experiment is selected.) Removes the experiment data from display.
Show Properties	Opens the Experiment Properties dialog box. Use this dialog box to see the names of the arrays in the experiment, and also to add or remove arrays from the experiment. See “Experiment Properties” on page 167.
Export	Opens the Export Experiments dialog box, where you can export this and other experiments as a single ZIP file. See “Export Experiments” on page 174 and “To export experiments” on page 71.
Export Attributes	Opens the Export Attributes dialog boxes, one for selecting arrays for which you want attributes exported and one for selecting the attributes you want to export with the selected arrays. See “Export Array Attributes” on page 170.
Edit Array Color	Opens the Edit Array Color dialog box, where you can select a display color for each of the arrays in the experiment. For more information see “Edit Array Color” on page 165.

4 Data Viewing Reference

Experiment pane – actions and shortcut menus

Option	Description
Edit Array Order	Opens the Array Order dialog box, where you can change the order of the arrays in the Experiment and in Tab View. See “Edit Array Order” on page 166.
Rename	Opens an Input dialog box, where you can type a new name for the experiment. Click OK to rename the experiment.
Delete	Opens a Confirm dialog box that asks if you want to delete the Experiment. Click Yes to delete it. Note: You can delete any experiment except the selected one.
Expand Node	Expands the selected node to display all folders and their contents.
Collapse Node	Closes all folders for the selected node.

Design Folder

- Right-click the name of a design to open a shortcut menu with a Delete option. If you select this option, a Confirm dialog box opens. If you click **Yes**, the program removes the links to all of the arrays under the design from the experiment.

Genome Build Folder

- Right-click the name of a genome build within a design to display the following options:








Option	Description
Set for Calibration	Agilent does not recommend using another array to calculate noise for the sample array.
Delete	Opens a Confirm dialog box that asks if you want to disassociate all arrays under the design from the experiment. Click Yes to remove the links between the arrays and the experiment. <ul style="list-style-type: none">If you delete a design from an experiment, the program removes the links between the experiment and the design and its arrays. The actual design and array data stay in the Data folder.

Individual Arrays

- Within the folder of a specific experiment, in the **Arrays** folder of a design, right-click the name of an individual array display the following options:

Option	Description
Select	(Available if the array is not selected) Selects the array for display.
Deselect	(Available if the array is selected) Removes the array data from Genome, Chromosome, and Gene views. Also removes the array from the Selected Arrays tab in Tab View.
Rename	Opens an Input dialog box, where you can type a new name for the array. Click OK to accept the new name for the array.
Delete	<p>Opens a Confirm dialog box that asks if you want to disassociate the array from the experiment. Click Yes to remove the link between the array and the experiment. See “To remove arrays from an experiment” on page 60.</p> <ul style="list-style-type: none"> • If you delete an array from an experiment, the program removes the link between the experiment and the array. The actual array data stays in the Data folder.
Show Properties	<p>Opens the Microarray Properties dialog box, where you can display and edit microarray attributes.</p> <p>For array files from the Agilent Feature Extraction program, you can also display the headers and feature data from the file.</p> <p>See “Microarray Properties” on page 187 and “To display or edit the attribute values of a specific array” on page 60.</p>
Edit Array Color	Opens the Edit Array Color dialog box, where you can select a display color for the array. See “Edit Array Color” on page 165 and “To change the display color of an array” on page 78.
Edit Array Order	Opens the Array Order dialog box, where you can change the order of the arrays in the Experiment and in Tab View. See “Edit Array Order” on page 166 and “To change the order of arrays in an experiment” on page 58.

My Entity List pane – icons, buttons, and special text

Item	Comments
	Click to expand a folder and display its contents.
	Click to collapse a folder and hide its contents.
	A folder that contains files or other folders.
	An individual gene list or track.
red regular text	An item that is an exact match with the search term in a search, or a gene list that has not been applied and has red assigned as its custom color.
colored italics	A gene list that has been applied.
red bold italics	A track that is selected for display in Gene View.
black bold italics	A “combined” track that is selected for display in Gene View. A combined track contains information from two or more individual tracks associated by logical criteria.
	(Dock out button) Moves the My Entity List pane from the main window, and opens it in a, separate window.
	(Collapse button, available only if the My Entity List pane is not collapsed) Collapses the My Entity List pane, and shows its title bar at the bottom of the Navigator.
	(Expand button, available only if the My Entity List pane is collapsed) Expands the My Entity List pane.

My Entity List pane – actions and shortcut menus

Only the options that are not grayed out are discussed in this section. To enable the grayed-out options, you must have the appropriate license for the CGH, ChIP or CH3 application you are using. These options are explained in the *User Guide* for the application.

- Double-click the **Gene List** folder to show or hide its gene lists.

Gene List Folder

- In the **Gene List** folder, right-click the name of a gene list to display the following options:

Option	Description
Rename	Opens an Input dialog box, where you can type a new name for the gene list. Click OK to accept the new name.
Delete	Opens a confirm dialog box that asks if you are sure you want to delete the gene list. Click Yes to confirm.

Tracks Folder

- Right-click the name of a track to display the following options:

Option	Comments
Show in UI	Mark this option to display the track in Gene View next to the data and results of the selected experiment. See “To show tracks in Gene View” on page 90 and “User Preferences” on page 208.
Show in Report	Mark the check box to show the track information in all the reports.
Genomic Boundaries	Click to use the genome track to define only the regions that aberration detection algorithms will run. You can select this for only one track.
Show in UCSC	Opens the UCSC (University of California at Santa Cruz) Genome Browser in your Web browser and uploads the track. You can then see information for the track.
View Details	Opens a table that shows all the chromosome locations defined in the track.
Rename	Opens an Input dialog box, where you can type a new name for the track. Click OK to rename the track.
Delete	Opens a Delete Track dialog box that asks if you are sure you want to delete the track. Click Yes to delete the track.

Genomic Viewer

This section describes the display areas that appear when you click the Genomic Viewer tab. The orientation of these views (vertical or horizontal) can be changed from View Preferences located in the View tab. See [“View Preferences”](#) on page 216 for more information.

Genome View








Figure 27 Genome View, vertical orientation, with human chromosomes. The X chromosome is selected.

Genome View shows pictures of each of the distinct types of chromosomes in the selected genome. A blue box is drawn around the currently selected chromosome, and the cursor appears as a blue line across the chromosome.

Genome View actions and shortcut menus

- Click a chromosome to select it. When you select a chromosome, Chromosome, Gene, and Tab Views show only genomic regions, genes, and data associated with it. The specific location where you click the chromosome sets the position of the cursor. See [“The View Cursor”](#) on page 132.
- On the selected chromosome, click anywhere to move the cursor. See [“The View Cursor”](#) on page 132. This also moves the cursor in Chromosome, Gene, and Tab Views.
- Right-click anywhere within Genome View to display a menu. If you click **View Preferences**, the View Preferences dialog box opens, where you can set preferences for the display. See [“View Preferences”](#) on page 216.
- Click the **Detach** button  (located at the top center of the pane) to remove Genome View from the main window and open in a separate window. To reattach the view, click its **Close** button . Drag the side or bottom borders of the pane to resize them.
- Drag the side or bottom borders of the pane to resize it.
- On a border of the pane, click a resize button (for example,  or ) that points away from the pane to move that border all the way to the edge of the main window. To move the border back to its previous location, click the other resize button.

Chromosome View








Figure 28 Chromosome View, human X chromosome shown

Chromosome View shows a more detailed diagram of the chromosome you select in Genome View.

- Cytobands and a plotting area appear next to the chromosome.
- When you select arrays for display, their data appear in the plotting area.
- The View cursor appears as a solid blue line across the chromosome and the plotting area.
- The selected region of the chromosome (if any) appears as a dotted blue box in the plotting area.

Chromosome View actions and shortcut menus

- Click a cytoband, any part of the chromosome, or anywhere in the plotting area to move the View cursor to that location. See [“The View Cursor”](#) on page 132.
- Drag the pointer over any part of the plotting area to select a chromosomal region for display in Gene View. Drag parallel to the chromosome. This also moves the cursor to the center of the selected region. See [“The View Cursor”](#) on page 132.
- Right-click anywhere within Chromosome View to display a menu. If you click **View Preferences**, the View Preferences dialog box opens, where you can set preferences for the display. See [“View Preferences”](#) on page 216.
- Click the **Detach** button  (located at the top center of the pane) to remove Chromosome View from the main window and open in a separate window. To reattach the view, click its **Close** button . Drag an inside border of Chromosome View to resize the view.
- Drag the side or bottom borders of the pane to resize it.
- On a border of the pane, click a resize button (for example,  or ) that points away from the pane to move that border all the way to the edge of the main window. To move the border back to its previous location, click the other resize button.

Gene View




Figure 29 Gene View, with log ratio data from an experiment and CNV tracks.

Gene View shows a more detailed view of the chromosomal region you select in Chromosome View. See [“Chromosome View”](#) on page 126.

- Regions that contain genes appear as small blue boxes. Gene names appear nearby. You can customize the appearance of gene names. Also, you can use a gene list to highlight genes of interest, or to display only the genes in the list. See [“To change the appearance of genes in Gene View”](#) on page 89, and [“To show gene lists in Gene View”](#) on page 89.

- Log ratio data from selected arrays in the active experiment appear as a scatter plot. You can also customize the scatter plot. See [“To customize scatter plot ranges and colors”](#) on page 81.
- The location of the cursor matches the location of the cursors in other views. See [“The View Cursor”](#) on page 132.
- The name of the chromosome, and the coordinates and size of the displayed chromosomal region appear at the top of the pane.
- Imported tracks can also appear in Gene View. See [“To show tracks in Gene View”](#) on page 90.

Scatter Plot




Figure 30 Scatter Plot command group in CGH Gene View

The scatter plot command group is available in Gene View or the View tab. The commands differ depending on the DNA Analytics application you are using. All the scatter plot command groups contain the commands for log ratio data in [Figure 30](#). Scatter plots appear in the Chromosome and Gene Views but only if they have been selected in the View Preferences dialog box.

The drop down lists let you select the type of data to display in the plot. For more information, see [“To show or hide data in scatter plots”](#) on page 80 and [“To customize scatter plot ranges and colors”](#) on page 81.

Gene View buttons



Zooms in to see a smaller region in more detail.



Zooms out to see a larger region in less detail.



In vertical orientation, scrolls up through the genes and data to lower-numbered chromosomal coordinates.



In vertical orientation, scrolls down through the genes and data to higher-numbered chromosomal coordinates.



In horizontal orientation, scrolls left through the genes and data to lower-numbered chromosomal coordinates.



In horizontal orientation, scrolls right through the genes and data to higher-numbered chromosomal coordinates.



(**Resize** buttons) The button that points away from Gene View expands the view. The other button restores the view to its original size.



(**Detach** button) Removes Gene View from the main window, and opens it in a separate window.

Gene View shortcut menu and other actions

- Click anywhere in the plotting area of Gene View to move the cursor to that location. See [“The View Cursor”](#) on page 132.
- Drag an inside border of Gene View to resize the View. Right-click anywhere in the plotting area of Gene View to display the following options:


Option	Description
Create Gene List	Opens the Create Gene List dialog box, where you can create a new gene list based on the currently selected (or another) chromosomal region. See “Create Gene List” on page 154 and “To show gene lists in Gene View” on page 89.
Create Track	Opens the Create Track dialog box, where you set the chromosome locations for the track. See “To create a track (CGH only)” on page 68 and “Create Track” on page 156.

Option	Description
Show in UCSC	Opens the View Coordinates in UCSC Browser dialog box where you select track information for display in the UCSC (University of California at Santa Cruz) Genome Browser. You can then view the track.
User Preferences	Opens the User Preferences dialog box, where you can set user preferences on three separate tabs. See “User Preferences” on page 208 and the related pages that follow.
View Preferences	Opens the View Preferences dialog box, where you can set the preferences for viewing data in the Genomic Viewer. See “View Preferences” on page 216.

The View Cursor

The View cursor reflects the center of the current chromosomal location of interest. It appears in several views:

- In Genome View, it appears as a blue bar across the selected chromosome.
- In Chromosome View, it is a blue bar that appears across the chromosome and across the plotting area of the view.
- In Gene View, it is a blue bar that appears across the plotting area and tracks of the view.



The position of the cursor in one View is also the position of the cursor in all Views. The exact chromosomal location of the cursor appears in the first cell of the Status bar. Several actions change the position of the View cursor:

- In Genome View, click anywhere on a chromosome to move the cursor to that location.
- In Chromosome View, click a cytoband name, part of the chromosome, or anywhere in the plotting area to move the cursor to that location.
- In Gene View, click anywhere in the plotting area to move the cursor to that location.

The cursor used in Gene View is the same cursor used for the tracks.

- In Tab View, click a row of a data table to move the cursor to the chromosomal location associated with that row.

Tab View




Figure 31 Tab View


Tab View displays design annotation and log ratio data related to the chromosome you select in Chromosome View.

- The exact column content of the tables depends on the specific tab and design, but it always includes chromosomal locations of probes
- The selected row of data appears highlighted in blue. This row represents data that corresponds approximately with the location of the cursor.
- Columns of log ratio data appear below the names of the specific arrays to which they correspond. If an array is selected for display in Chromosome and Gene views, a colored square appears next to its name.

Tab View tabs and buttons


You can see the following tabs and buttons in Tab View. See [Figure 31](#) for a diagram that identifies some of these elements.

Design tabs A separate tab appears for each microarray design included in the active experiment. The name of the design appears on each tab, along with an icon:

 **CHB** – A methylation array design

 **CGH** – An aCGH array design.

 **EXP** – A gene expression array design.

 **CHIP** – A ChIP-on-Chip array design.

When you click a design tab, the data and annotation for the arrays in the design appear in Tab View. The program separates the arrays of the design into the Arrays tab and the Calibration Arrays tab (see below).

Arrays tab (Available when you click a specific design tab.) Contains a table of data and annotation for all arrays in a design that contain biological data.

Selected Arrays tab Contains a table of data and annotation for the selected arrays from all designs in the active experiment.



(Resize buttons) The button that points away from Tab View expands the view. The other button restores the view to its original size.



(Detach button) Removes Tab View from the main window, and opens it in a separate window.

Tab View actions and shortcut menus

- Click the name of an *array in a column heading* to select the array data for display.
- Right-click the name of an *array in a column heading* to open a display the following options:

Option	Description
Rename Array	Opens an Input dialog box, where you can type a new name for the array. This only changes the name of the array within the active experiment.
Remove Array From Experiment	Opens a confirmation dialog box. Click Yes to remove the link between the array and the active experiment. This command does not delete the data file from the program. To do this, see “To remove data or design files from the program” on page 65.
Select Array	(Available if the array is not selected.) Selects the array for display. A colored square appears next to the name of the array.
Deselect Array	(Available if the array is selected.) Removes the array data from scatter plots, and removes the column of the array from the Selected Arrays tab.
Edit Array Color	Opens the Select Color dialog box, where you can change the display color of the array. See “Edit Array Color” on page 165 and “To change the display color of an array” on page 78.
Edit Array Order	Opens the Edit Array Order dialog box, where you can change the order in which the names of the arrays in a given design of the active experiment appear in Tab View and in the Data Navigator. In Gene View, when you display separate scatter plots for each array, the plots also appear in this order. See “Edit Array Order” on page 166 and “To change the order of arrays in an experiment” on page 58.
Select All Arrays	Selects all arrays in all designs in the active experiment for display. All arrays appear in the Selected Arrays tab.
Deselect All Arrays	Removes all arrays from display, and from the Selected Arrays tab.
Scroll to Column	Opens the Scroll to Column dialog box, where you can select a column in the current tab. The program then scrolls the data table in the tab so you can see the selected column.

- Right-click a *heading of a column other than an array data column* to open a shortcut menu with a Scroll To Column option. If you click this option, the Scroll to Column dialog box appears, where you can select a column in the current tab. The program then scrolls the data table in the tab so you can see the column. See “[Scroll to Column](#)” on page 194.
- Click a *data table entry* to select the row in which it appears. This also moves the cursor to the location of the data point corresponding to the selected row.
- Right-click a *data table entry* to display the following options:

Option	Description
Find in Column	Opens the Find in column dialog box, where you can search for a specific text string within the column you clicked. See “ Find in column ” on page 177.
Google LocusLink PubMed UCSC HG15(April ‘03) UCSC HG16(July’03) UCSC HG17(May’04) UCSC HG18(March’06) UCSC mm8(Feb’06) UCSC mm9(July’07) DGV(hg18) GO KEGG(HUMAN)	Opens your Web browser, and sends the column entry you clicked as a search string to the selected site. The UCSC links search the indicated University of California, Santa Cruz database related to the indicated genome build. See “ To search the Web for information on probes in Tab View ” on page 98.
Customize Link	Opens the Customize Search link dialog box, where you can create or edit a custom Web link that appears in this shortcut menu. When you click a custom link, the program opens your Web browser, and sends the column entry you clicked as a search string to the site. See “ Customize Search Link ” on page 157 and “ To update or delete a custom Web search link ” on page 99.

Option	Description
(other options)	If other options appear in this shortcut menu, they are custom Web search links. Click them to open your Web browser, and send the column entry you clicked as a search string to the site.

Status Bar




Figure 32 Status bar

The Status Bar displays information related to the currently displayed data. There are other items on the status bar that only become active if you have a DNA Analytics application license.

Cursor position	The chromosomal location of the cursor. See “The View Cursor” on page 132.
Genome build	The genome build associated with the currently displayed data.
Ratio type	The mathematical type of the array data. The possible types are: <ul style="list-style-type: none">• ratio• log₂ ratio• log₁₀ ratio• ln (natural log) ratio
Selected Row	The row in the currently displayed data table that is selected. The location of the cursor is approximately the chromosomal location associated with this row.
Table size	The number of rows and columns in the currently displayed tab. The size appears as <# of rows> x <# of columns>.

Dialog Boxes

This section describes the dialog boxes that can appear when you import, organize, manage, export and display array data and other content in Agilent Genomic Workbench. The dialog boxes appear in alphabetical order by name.

Agilent Feature Extraction Importer




Figure 33 Agilent Feature Extraction Importer dialog box

Purpose: Lets you edit the name of the FE data file you will import and to indicate whether you want to flip the red/green ratio for the data.

To open: In the Home tab, click **Import > Array Files > FE File**, select the desired FE data file(s), then click **Open**.

Name Lets you edit the names of the FE files. You can change the names of the files to names that are easier to recognize or remember.

Dye Flip For each array:

Select **Normal** if:

- The test samples were labeled with cyanine-5 (red).
- The control samples were labeled with cyanine-3 (green).
- The imported ratio (test/control) should be reported as-is.

Select **Flipped** if:

- The test samples were labeled with cyanine-3 (green).
- The control samples were labeled with cyanine-5 (red).
- The imported ratio (control/test) should be reported with the ratio inverted (test/control).

The program does not combine dye-flip pairs.

Overwrite arrays with duplicate names Mark this option to replace existing file(s) in the program with the imported one(s), if they have the same name(s).

Run in Background Imports the files, and lets you use your computer for other purposes while the import occurs. This is especially useful if you have many files to import.

OK Imports the files in the foreground. You cannot use your computer for other purposes while the import occurs.

Cancel Cancels the entire import process without importing anything.

CGHSmooth Parameters




Figure 34 CGHSmooth Parameters dialog box

Purpose: The CGHSmooth Parameters dialog box lets you configure the CGHSmooth plug-in. It can perform one of several varieties of moving-average-like smoothing, and plots the data associated with the currently selected chromosome.

To open: Click **Tool > Plugin > CGHSmooth**.

Parameters Set any of these parameters:

Parameter	Description
Smoothing Function	<p>A number from 0 to 5. The number sets one of the following options as the weighting function used by the moving average algorithm. In general, the options weight measurements closer to the center position more heavily than those more distant from it.</p> <p>0 – None. The plug-in applies no smoothing, and returns the original data. In some cases, the plug-in averages data points with identical positions. This sets, in effect, a window size of 0.</p> <p>1 – Rectangular. The plug-in performs a standard moving average. All points within the rectangle (the window) receive the same weight.</p> <p>2 – Gaussian. Applies a Gaussian weighting function.</p> <p>3 – Triangular. Applies a triangular weighting function.</p> <p>4 – Lorentzian. Applies a Lorentzian weighting function.</p> <p>5 – Biexponential. Applies a biexponential weighting function.</p>
Output Options	<p>A number from 0 to 2. The number sets one of the following options:</p> <p>0 – Overlays the unsmoothed plot of each array on the smoothed plot.</p> <p>1 – Displays smoothed and unsmoothed plots for each array.</p> <p>2 – Displays smoothed, unsmoothed, and error plots for each array.</p>
X-axis Label	The text that appears under the X-axis of the plot as a label.
Y-axis Label	The text that appears next to the Y-axis of the plot as a label.
Y-axis Range (min)	The minimum value on the Y-axis.
Y-axis Range (max)	The maximum value on the Y-axis.
X-axis Range (min)	The minimum value on the X-axis.
X-axis Range (max)	The maximum value on the X-axis.

Don't show again Mark this option to prevent the appearance of this dialog box in the future when you click **Tool > Plugins > CGHSmooth**. To restore the dialog box so it appears again, click **Tool > Plugin Settings > CGHSmooth**, then clear **Don't show again**.

OK Accepts the parameters and prepares the plot. You can make further changes to the appearance of the plot once the plug-in displays it.

Cancel Ignores any changes you made, and closes the dialog box.

CGHSmooth Plot




Figure 35 CGHSmooth Plot

Purpose: The CGHSmooth Plot is the output of the CGHSmooth plug-in. It contains separate, stacked plots of smoothed log ratio data for each of the selected arrays in the current experiment.

To open: Click **OK** in the CGHSmooth Parameters dialog box. See [“CGHSmooth Parameters”](#) on page 141.

Plot(s) Depending on the selected output option, the main plotting area shows up to three plots for each array in the active experiment. The plots can include unsmoothed and smoothed log ratio plots, and an error plot.

Edit Opens a menu with a Copy Plot(s) to Clipboard option. If you select this option, the program copies the plots to the clipboard as an image. You can then paste the image into a document in another program.

When you right-click anywhere within the plotting area, the following options are displayed:

Option	Description
Properties	Opens the Chart Properties dialog box, where you can change the appearance of many aspects of the plot. See “ Chart Properties ” on page 145.
Copy	Copies the chart to the Clipboard, where you can paste it into a word processing or other program.
Save as	Opens a Save dialog box, where you can select a location for the *.png image file of the plots.
Print	Opens the Windows Page Setup dialog box, where you can set page size, margins, and orientation, and select a printer. From there, you click OK to open the Windows Print dialog box, where you can set print options and print the plot.
Zoom In	<p>Opens another menu that lets you zoom in the plot. You can zoom in several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms in the Domain (X) axis for all stacked plots, and zooms in the range Range (Y) axis for the specific plot in which you clicked.• Domain Axis – Zooms in the Domain (X) axis for all stacked plots.• Range Axis – Zooms in the Range (Y) axis for the specific plot in which you clicked. <p>You can also drag across an area of one of the plots to select an area to expand.</p>
Zoom Out	<p>Opens another menu that lets you zoom out the plot. You can zoom out several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms out the Domain (X) axis for all stacked plots, and zooms out the range Range (Y) axis for the specific plot in which you clicked.• Domain Axis – Zooms out the Domain (X) axis for all stacked plots.• Range Axis – Zooms out the Range (Y) axis for the specific plot in which you clicked.

Option	Description
Auto Range	<p>Opens another menu that lets you zoom the plot to show the full range of the data. You can zoom in several ways:</p> <ul style="list-style-type: none"> • Both Axes – Appropriately zooms both axes of the specific plot to show the full set of data. • Domain Axis – Zooms in the Domain (X) axis for all stacked plots so this axis displays the full range of X values of the data. • Range Axis – Zooms in the Range (Y) axis for the specific plot in which you clicked to display the full range of Y values of the data.

Chart Properties

Purpose: The Chart Properties dialog box lets you create titles and legends, as well as change the appearance, for the CGHSmooth, Echo Example, and MovAvg Example plots.

To open: Use the CGHSmooth, Echo Example, or MovAvg Example plug-in to draw a plot. Right-click within the plotting area, then click **Properties** in the shortcut menu.

This dialog box has four tabs. At any point, click **OK** to accept the settings in all four tabs, or click **Cancel** to close the dialog box without making any changes to the settings.

Title Tab




Figure 36 Chart Properties dialog box – Title tab

- **Show Title** – Mark this option to display a title across the top of the chart.
- **Text** – Type a title for the chart.
- **Font** – (Available if you mark **Show Title**) Click **Select** to open the Font Selection dialog box. Select the desired font attributes, then click **OK**.
- **Color** – (Available if you mark **Show Title**) Click **Select** to open the Title Color dialog box. Select or configure a color for the title, then click **OK**. This dialog box is identical to the Select Color dialog box. See [“Select Color”](#) on page 196.

Plot Tab




Figure 37 Chart Properties dialog box – Plot tab

- Within the Plot tab, you can set these properties in the Domain Axis tab (“X” axis):

Property	Description
General	
Label	A custom label for the Domain (X) axis of the chart. Type the desired label.
Font	The font for the custom label on the Domain (X) axis. Click Select to open the Font Selection dialog box. Select the desired font attributes, then click OK .
Paint	The color of the custom label on the Domain (X) axis. Click Select to open the Label Color dialog box. Select the desired color, then click OK . This dialog box is identical to the Select Color dialog box. See “ Select Color ” on page 196.

Property	Description
Other – Ticks tab	
Show tick labels	Mark this option to show, or clear it to hide, the numerical values on the domain axis.
Tick label font	The font for the numerical values on the Domain (X) axis. Click Select to open the Font Selection dialog box. Select the desired font attributes, then click OK .
Show tick marks	Mark this option to show, or clear it to hide, tick marks on the Domain (X) axis.
Other – Range tab	
Auto-adjust range	Mark this option to automatically set the range of values on the X-axis to include all data.
Minimum range value	(Available if you do not mark Auto-adjust range) The lowest value represented on the X-axis.
Maximum range value	(Available if you do not mark Auto-adjust range) The highest value represented on the X-axis. The program automatically converts large numbers to scientific “E” notation – for example, 1.22E8 .

- Within the Plot tab, you can set the following properties in the Appearance tab:

Property	Description
Outline stroke	The thickness of the lines around each plot. Click Select to open the Stroke Selection dialog box. Select the desired line thickness, then click OK .
Outline paint	The color of the lines around each plot. Click Select to open the Outline Color dialog box. Select the desired color, then click OK . This dialog box is identical to the Select Color dialog box. See “ Select Color ” on page 196.
Background paint	The color of the background within each plotting area. Click Select to open the Background Color dialog box. Select the desired color, then click OK . This dialog box is identical to the Select Color dialog box. See “ Select Color ” on page 196.
Orientation	Select either Vertical (X-axis on the bottom of the chart) or Horizontal (X-axis on the left side of the chart).

Other tab




Figure 38 Chart Properties dialog box – Other tab

The Other tab has these options:

- **Draw anti-aliased** – Mark this option to minimize distortion and visual artifacts in the plot image. This will create a smoother image, but it can be less sharp than the original one.
- **Background paint** – The color of the chart outside of the plotting area and legend. Click **Select** to open the Background Color dialog box. Select the desired color, then click **OK**. This dialog box is identical to the Select Color dialog box. See [“Select Color”](#) on page 196.

The other options are for future expansion, and are not available in the current release of Agilent Genomic Workbench.

Configure Coloring Ranges and Shades




Figure 39 Configure Coloring Ranges and Shades dialog box for CGH

Purpose: This dialog box is used to enter ranges and select colors for scatter plot options. Tabs show scatter plot selections for the selected application type (CGH, ChIP, or CH3).

To open: In Gene View, move the mouse pointer over **Scatter Plot** to display the scatter plot options and then click **Configure Color and Ranges**. Or, click the **View** tab and click **View Preferences**. Then, under Configure Coloring schemes, click **Configure Colors and Ranges**.

For information on the contents of the various tabs, see the *User Guide* for the selected application.

Confirm Overwrite




Figure 40 Confirm overwrite dialog box

Purpose: When you import an experiment, it can contain designs and/or arrays that have the same names as those already available in Agilent Genomic Workbench. This dialog box lets you select which designs and/or arrays to overwrite.

To open: This dialog box appears when you import a ZIP format experiment file, and it contains designs and/or arrays that are already available in Agilent Genomic Workbench. See [“To import an experiment file”](#) on page 52.

Select the designs to overwrite

- Design** The names of the designs in the imported file that have the same names as designs that are already available in Agilent Genomic Workbench.
- Overwrite** Mark the check box for each existing design that you want to overwrite.
- Select All** Marks all of the check boxes under Overwrite.
- Deselect All** Clears all of the check boxes under Overwrite.

Select the microarrays to overwrite

- Array** Identification number or barcode of the array
- Name** The name of the array in the imported file that has the same name as array that is already available in Agilent Genomic Workbench.
- Overwrite** Mark the check box for each existing array that you want to overwrite.
- Select All** Marks all of the check boxes under Overwrite.
- Deselect All** Clears all of the check boxes under Overwrite.
- OK** Overwrites the selected files (both designs and arrays) and closes the dialog box.
- Cancel** Closes the dialog box, and returns you to the Import (experiments) dialog box. See [“Import \(experiments\)”](#) on page 181.

Create Experiment




Figure 41 Create Experiment dialog box

Purpose: Creates an organizational unit (an *experiment*) that lets you display and analyze array data in Agilent Genomic Workbench. You add data to the experiment with links to array data files that are available in the program, a process that you can start from this dialog box. See [“To create a new experiment”](#) on page 55.

To open: In the Home tab of Agilent Genomic Workbench, click **Create Experiment**.

- Name** (Required) The name of the new experiment. This name identifies the experiment within the program and in exported reports and files.
- Description** (Optional) Brief information that will later help to identify the experiment.
- Properties** Opens the Experiment Properties dialog box, where you can select array data files to add to the new experiment. See [“Experiment Properties”](#) on page 167.
- OK** Closes the dialog box and creates the new experiment.
- Cancel** Closes the dialog box without creating an experiment.

NOTE

Click **Properties** to open the Experiment Properties dialog box to add array data to your new experiment. Otherwise, the program creates an empty experiment. You can also add arrays to the experiment later. See [“To add arrays to an experiment”](#) on page 57.

Create Gene List


The image shows a 'Create Gene List' dialog box with a blue title bar. It contains several input fields and options. At the top, there is a 'Name' field and a 'Build' dropdown menu currently set to 'hg18'. Below these is a large 'Description' text area. A section titled 'Set Chromosome Start-Stop' contains a 'Chromosome' dropdown set to 'chr8', a 'Start' text field with '0', and a 'Stop' text field with '549999'. Underneath this section are three radio button options: 'User Defined', 'For complete gene view' (which is selected), and 'For aberrant region below cursor'. At the bottom of the dialog is a 'Color' section with a red square icon and a 'Change...' button. Finally, there are 'OK' and 'Cancel' buttons at the very bottom.

Figure 42 Create Gene List

Purpose: To limit the genes presented in Gene View to a preselected number valuable for interpreting data

To open: Right-click Gene View, and click **Create Gene List**.

Name Type in name of gene list.

Build Select the genome build for the genes to be selected for list.

Description Describe the type or nature of the genes in the list.

Set Chromosome Start-Stop

Select a chromosome and a region in Chromosome View for selecting the genes in the list before you open the Create Gene List dialog box.

- User Defined** Lets you select a region from which the genes in Gene View will be selected. The chromosome selection list and the Start and Stop positions on the Y axis are enabled when this option is selected. With this option you can override the selections you made before opening Create Gene List.
- For complete gene view** Select all the genes in Gene View.
- For aberrant region below cursor** Select those genes that appear in the aberrant region just below where the cursor sits in Gene View. Not operational in Genomic Viewer; depends on analysis.
- Chromosome** If you select User Defined, you can select a different chromosome than had been selected before opening the Create Gene List dialog box.
- Start** If you select User Defined, you can type in a Start position for defining the region contained the genes to be in the list.
- Stop** If you select User Defined, you can type in a Stop position for defining the region contained the genes to be in the list.
- Color**
- Change** Click to change the color of the gene list name in Data Navigator. See [“Select Color”](#) on page 196.

Create Track




Figure 43 Create Track dialog box

Purpose: The Create Track dialog box lets you create a track for a chromosomal region based on an assigned chromosomal region. You can display one or more tracks next to the genes and data in Gene View. See [“To show tracks in Gene View”](#) on page 90.

To open: Right-click in the plotting area of Gene View for the CGH or CH3 application, then click **Create Track** in the shortcut menu.

Name Type a name for the track. This name identifies the track when it appears in views and lists.

Build (Available if you select **User Defined** in **Set Chromosome Start-Stop**.)
Select the genome build for the track.

Description Type descriptive text to attach to the track for reference.

- Set Chromosome** Defines the region of the chromosome for which the track will be defined.
- Start-Stop** Select one of these options:
- User Defined** – Lets you define an arbitrary region of any chromosome. If you select this option, select the desired chromosome in **Chromosome**, then type the beginning (**Start**) and end (**Stop**) locations for the interval.
 - **For complete gene view** – The chromosomal region that currently appears in Gene View.
- OK** Creates the track. To display the track in Gene View, use the **Tracks** tab of the User Preferences dialog box to enable it. See “[User Preferences](#)” on page 208. To export the track, see “[To export tracks](#)” on page 72.
- Cancel** Closes the dialog box without creating a track.

Customize Search Link




Figure 44 Customize Search Link dialog box

Purpose: This dialog box lets you create a custom Web search link in the shortcut menu that appears when you right-click an entry in the Tab View. The link opens the URL of your choice, and sends the selected entry to it as a search string. See “[To create a custom Web search link](#)” on page 98.

To open: Right-click any entry in a table in Tab View, other than a column heading, then click **Customize Link**.

URL Name The name of the custom Web search link that appears in the shortcut menu (see above). To edit an existing custom Web search link, select it from the list.

URL The full uniform resource locator (URL) of the desired search page. For the query string value, type <target>

For example, this URL sends the selected Tab View entry to google.com:

`http://www.google.com/search?hl=eng&q=<target>`

New Opens an Input dialog box, where you can type a name for a new custom Web search link. Click **OK** to accept the name and add it to the URL name list.

Update Saves the settings in the dialog box.

Delete Deletes the currently selected custom Web search link.

Close Closes the dialog box.

Design Properties

Purpose: Gives general and detailed information about a given microarray design. See “[To display the properties of a specific design](#)” on page 64.

To open: In the **Data** pane of the Navigator, right-click the name of a genome build within a design folder, then click **Show Properties**. Several tabs are available.

Attribute tab

Displays general identifying attributes of the array design, and statistics such as the total number of features in the design, or the date the design was last modified.





Figure 45 Design Properties dialog box – Attribute tab

Non Unique Probes tab

Shows the nonunique probes in the design. Nonunique probes have more than one mapping in the genome that is a perfect match.



The image shows a software dialog box titled "Design Properties" with a close button in the top right corner. It has three tabs: "Attribute", "Non Unique Probes" (which is selected), and "Data". The "Non Unique Probes" tab contains a table with three columns: "S.No", "Probe", and "Value". The table lists 16 rows of data, each representing a probe and its corresponding chromosomal locations. A "Close" button is located at the bottom right of the dialog box.


S.No	Probe	Value
1	A_18_P10316231	chr1:142889397-142889343 chr...
2	A_18_P13304487	chr2:87290171-87290230 chr2:...
3	A_18_P16129688	chr7:72252622-72252564 chr7:...
4	A_18_P14749132	chr4:75511720-75511778 chr4:...
5	A_18_P12210664	chr7:57795297-57795241 chr15:...
6	A_18_P16127424	chr7:72126068-72126023 chr7:...
7	A_18_P16194541	chr7:5870637-5870589 chr7:97...
8	A_18_P10521542	chr1:220712195-220712239 chr...
9	A_18_P16198767	chr7:101111563-101111622 chr...
10	A_18_P16194545	chr7:73971446-73971387 chr7:...
11	A_18_P10305471	chr1:142503907-142503856 chr...
12	A_18_P16839817	chr4:132815293-132815234 chr...
13	A_18_P10309255	chr1:142494051-142494110 chr...
14	A_18_P15312171	chr5:98889559-98889603 chr5:...
15	A_18_P13397924	chr2:130944748-130944689 chr...
16	A_18_P16841050	chr2:94864734-94864793 chr9:...

Figure 46 Design Properties dialog box – Non-Unique Probes tab

- S. No** The sequence order of the probes within the table.
- Probe** The name of the each nonunique probe.
- Value** The chromosomal locations to which each of the probes binds. Because these are nonunique probes, multiple locations appear for each probe.

Data tab

Displays the names of the probes in the design and their target genomic locations. The tab displays the probes for one chromosome at a time.



Probe	Chromosome	Start	Stop
A_18_P10000009	chr1	3179	3223
A_18_P10000021	chr1	4753	4804
A_18_P10000034	chr1	8673	8717
A_18_P10000124	chr1	41899	41955
A_18_P10000236	chr1	67372	67431
A_18_P10000241	chr1	77534	77578
A_18_P10000243	chr1	77808	77864
A_18_P10000247	chr1	79598	79657
A_18_P10000248	chr1	79953	79997
A_18_P10000270	chr1	87402	87461
A_18_P10000277	chr1	88790	88842
A_18_P10000285	chr1	92164	92218
A_18_P10000290	chr1	95834	95886
A_18_P10000306	chr1	101360	101419

Figure 47 Design Properties dialog box – Data tab

- Select Chromosome** The chromosome whose probes appear in the list. To view the probes for another chromosome, select one from this list.
- Probe** The name (Probe ID) of each probe.
- Chromosome** The name of the probe chromosome.
- Start** The location on the selected chromosome of the first base pair for the probe.
- Stop** The location on the selected chromosome of the last base pair for the probe.

Echo Example Plot




Figure 48 Echo Example Plot

Purpose: The Echo Example Plot is the output of the Echo Example plug-in. It displays the log ratio data for the selected chromosome in the active experiment. Data from all of the selected arrays in the experiment appear as a series of stacked plots, one for each array.

To open: Select the desired experiment, select the desired chromosome in Genome View, then click **Tool > Plugin > Echo Example**.

Edit Opens a menu with a **Copy plots to clipboard** command. This command copies all of the plots to the clipboard as an image. You can then paste the image into a document in another program.

Plots Each plot displays the log ratio data for the selected chromosome from an individual array in the experiment.

You can right-click anywhere within each plot to display the following options:

Option	Description
Properties	Opens the Chart Properties dialog box, where you can change the appearance of many aspects of the plot. See “ Chart Properties ” on page 145.
Copy	Copies the chart to the Clipboard, where you can paste it into a word processing or other program.
Save as	Opens a Save dialog box, where you can select a location for a saved *.png image file of the plots.
Print	Opens the Windows Page Setup dialog box, where you can set page size, margins, and orientation, and select a printer. From there, you click OK to open the Windows Print dialog box, where you can set print options and print the plot.
Zoom In	<p>Opens another menu that lets you zoom in the plot. You can zoom in several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms in the Domain (X) axis for all stacked plots, and zooms in the range Range (Y) axis for the specific plot in which you clicked.• Domain Axis – Zooms in the Domain (X) axis for all stacked plots.• Range Axis – Zooms in the Range (Y) axis for the specific plot in which you clicked. <p>You can also drag across an area of one of the plots to select an area to expand.</p>

Option	Description
Zoom Out	<p>Opens another menu that lets you zoom out the plot. You can zoom out several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms out the Domain (X) axis for all stacked plots, and zooms out the range Range (Y) axis for the specific plot in which you clicked.• Domain Axis – Zooms out the Domain (X) axis for all stacked plots.• Range Axis – Zooms out the Range (Y) axis for the specific plot in which you clicked.
Auto Range	<p>Opens another menu that lets you zoom the plot to display the full range of the data. You can zoom in several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms both axes of the specific plot to show the full set of data.• Domain Axis – Zooms in the Domain (X) axis for all stacked plots so this axis displays the full range of X values of the data.• Range Axis – Zooms in the Range (Y) axis for the specific plot in which you clicked to display the full range of Y values of the data.

Edit Array Color




Figure 49 Edit Array Color dialog box

Purpose: Lets you show, change, and/or export the color(s) assigned to the arrays in an experiment.

To open: In the **Experiment** pane, right-click the name of an experiment, then click **Edit Array Color**.

Edit Opens a menu with a Copy command. If you click **Copy**, the program copies the list of arrays and their assigned colors to the Clipboard. You can then paste the list into a document in another program such as Word or PowerPoint.

Select Array Mark the check box for the array(s) whose color you want to change.

Color Opens the Select Color dialog box, where you can select a new color for the selected array(s). If more than one array is selected, all of the selected arrays assume the new color. For more information about selecting array colors, see [“To change the display color of an array”](#) on page 78.

4 Data Viewing Reference

Edit Array Order

- Select All** Marks all of the check boxes.
- Deselect All** Clears all of the check boxes.
- Edit Color** Opens the Select Color dialog box, where you can select a new color for the selected array(s). (Same function as the buttons under Color)
- Restore default** Restores the system default colors to all arrays.
- OK** Saves all assigned array colors and closes the dialog box.
- Cancel** Closes the dialog box without saving any changes.

Edit Array Order






Figure 50 Edit Array Order dialog box

Purpose: Changes the display order of the arrays in an experiment. This can change the order in which array data appear in Gene View and Tab View.

To open: In the Experiment pane, right-click the name of an experiment, then click **Edit Array Order**.

- Array Name** The arrays in the selected design, shown in the order that they currently appear in the Experiment.
- Design** Select a design from the list. The arrays from the selected design appear under Array Name.
- Order by** (Optional) Select an array attribute. The program can set the order of arrays based on their respective values for the selected attribute.
-  Moves a selected array up in the list. To select an array, click its name.
-  Moves a selected array down in the list. To select an array, click its name.
- OK** Sets the new order of the arrays and closes the dialog box.
- Cancel** Closes the dialog box without changing the order of any arrays.

Experiment Properties




Figure 51 Experiment Properties dialog box

Purpose: Lets you select array designs and data to link to an experiment. See “[To add arrays to an experiment](#)” on page 57.

To open: In the Create Experiment dialog box, click **Properties**, or in the **Experiment** pane of the Navigator, right-click the name of an experiment, then click **Show Properties**.

Experiment Name (Read-only) The name of the selected experiment.

Description Description that was typed when the experiment was created.

Select Design

Designs Shows all of the designs available in the program. Select the design associated with arrays that you want to add to the experiment.

Genome Builds Shows the genome build(s) that are associated with the design. Select the desired genome build to display the arrays that are associated with a single genome build.

Arrays

Array List Shows the arrays in the selected design that are available for this experiment.

- To select an array to move to the Selected Array List, click its name.
- To select additional arrays, hold down the **Ctrl** key while you click their names.
- To select a contiguous block of arrays, click the name of the first array, then hold down the **Shift** key and click the name of the last one.

Selected Array List Shows the arrays that you have selected for this experiment.



Moves the selected arrays in Array List to the Selected Array List. You can move arrays from as many designs as you like, if they are all associated with the same genome build.



Moves all of the arrays in Array List to the Selected Array List.



Removes an array from the Selected Array List. To select an array for removal, click its name. If desired, you can re-add an array.



Clears the Selected Array List.

- Display name by** Click to select an attribute to be used for display of the names of arrays in the experiment. The Global Display name is the name assigned in Sample Manager for the array. See the *Sample Manager User Guide* for more information.
- OK** Adds the arrays in the Selected Array list to the experiment and closes the dialog box.
- Cancel** Closes the dialog box without adding any arrays to the experiment.

Export











Figure 52 Export dialog box – Several types of file exports use this dialog box. This specific example exports selected experiment(s) as a ZIP format file.

Purpose: Lets you select a location for an exported file.

To open: This dialog box appears after you select specific experiment(s), track(s), filter(s) or array attributes to export. See [“To export experiments”](#) on page 71, [“To export tracks”](#) on page 72, [“To export filters”](#) on page 72 or [“To export array attributes”](#) on page 70.

- Look in** Displays the folder or other location whose contents appear in the main pane of the dialog box. To select another folder or other location, click .
-  Moves to the next higher folder level.
-  Opens the Desktop.
-  Creates a new folder in the selected location in *Look in*.
-  Displays the names, only, of folders, files, and other locations in the main pane of the dialog box.
-  Displays both the names and more information about folders, files, and other locations in the main pane of the dialog box.
- Main pane** Displays the folders, files, and other locations in the selected location in *Look in*. Only files of the selected file type are displayed. To select file, click its name. To open a folder or other location, double-click its name.
- File name** Displays the name of the file to which the exported content will be saved. To change the name, you can either select a file in the main pane of the dialog box, or type a new name.
- Files of type** Sets the type of files that are displayed. To show all files, click , then select **All Files**.
- Export** Saves the selected content to the location given in the dialog box.
- Cancel** Cancels your selections and closes the dialog box.

Export Array Attributes

Purpose: This dialog box lets you select arrays whose attributes you want to export. It contains two tabs: an Array tab where you select the arrays, and an Attribute tab where you select the attributes of the selected arrays to export. See [“To export array attributes”](#) on page 70.

To open: In the Home command ribbon, click **Export > Array Attributes**, or in the **Experiment** pane of the Navigator, right-click the name of an experiment, then click **Export Attributes**.

Array tab




Figure 53 Export Array Attributes – Array tab dialog box

Select Design

Designs Shows all of the designs available in the program. Select the design associated with arrays whose attributes you want to export.





Genome Builds Shows the genome build(s) associated with the design. Select the desired genome build to display the arrays.

Arrays

Array List Shows the arrays in the selected design.

- To select an array to move to the Selected Array List, click its name.
- To select additional arrays, hold down the **Ctrl** key while you click their names.
- To select a contiguous block of arrays, click the name of the first array, then hold down the **Shift** key while you click the name of the last one.

Selected Array List Shows the arrays that you have selected for this experiment.

-  Moves the selected arrays in Array List to the Selected Array List. You can move arrays from as many designs as you like, if they are all associated with the same genome build.
-  Moves all of the arrays in Array List to the Selected Array List.
-  Removes an array from the Selected Array List. To select an array for removal, click its name. If desired, you can re-add an array.
-  Clears the Selected Array List.
- Next** Moves to Attribute tab for attribute removal.
- Cancel** Closes the dialog box without selecting any arrays and their attributes to be exported.

Attribute Tab




Figure 54 Export Array Attributes – Attribute tab dialog box

Attributes

- Selected Attribute List** Shows the attributes for the selected arrays to be exported.
- To select an attribute to move to the Attribute List, click its name.

- To select additional attributes, hold down the **Ctrl** key while you click their names.
- To select a contiguous block of attributes, click the name of the first attribute, then hold down the **Shift** key while you click the name of the last one.

NOTE

Because certain attributes are required for importing an attributes file, it is important that you select these attributes when you export an attributes file. Required attributes are: Array ID, Global Display Name, Green Sample, Red Sample, and Polarity.

Attribute List Shows the attributes that will not be exported for the selected arrays.



Removes an attribute from the Selected Attribute List. To select an attribute for removal, click its name. You can add the attribute to the Selected Attribute List at a later time.



Clears the Selected Attribute List.



Moves the selected attributes in the Attribute List to the Selected Attribute List.



Moves all of the attributes in the Attribute List to the Selected Attribute List.

Back Moves back to the Array tab for array selection or removal.

OK Opens the Export dialog box. See [“Export”](#) on page 169.

Cancel Closes the dialog box without exporting any attributes.

Export Experiments




Figure 55 Export Experiments dialog box

Purpose: Lets you select experiments for export. The program exports all array designs and data associated with the experiments as a single ZIP file. This file does not include any parameter settings, array selections, or results. See [“To export experiments”](#) on page 71.

To open: In the Home tab, click **Export > Experiments**.

Select experiments to export Shows all experiments available for export. Mark each experiment you want to export.

Select All Selects all experiments for export.

Deselect All Clears all check boxes under Select experiments to export.

OK Opens an Export dialog box. See [“Export”](#) on page 169.

Cancel Cancels the export and closes the dialog box.

Export Filters




Figure 56 Export Filters dialog box

Purpose: Lets you select feature-level, array-level, design, and/or aberration filters, to export as a single *.xml file. You can create and use filters only if you have a DNA Analytics application license. See [“To export filters”](#) on page 72.

To open: In the **Home** tab, click **Export > Filters**.

Select filters to export Displays all of the filters available in the program. The table has these columns:

- **Export** – Mark the check box for each filter to export.
- **Filter** – The name of each filter.
- **Type** – The type of content to which the program applies each filter.

Select All Selects all available filters for export.

Deselect All Clears all of the check boxes under Select filters to export.

- OK** Opens the Export dialog box, where you can select a location for the exported *.xml file of filters. See [“Export”](#) on page 169.
- Cancel** Cancels the export and closes the dialog box.

Export Tracks




Figure 57 Export Tracks dialog box

Purpose: Lets you select tracks to export as a single BED format file. See [“To export tracks”](#) on page 72.

To open: In the **Home** tab, click **Export > Tracks**.

Select tracks Shows all of the tracks available in the program. Mark the check box for each track to export.

For more information about tracks, see [“To create a track \(CGH only\)”](#) on page 68 and [“To show tracks in Gene View”](#) on page 90.

Select All Selects all available tracks for export.

- Deselect All** Clears all of the check boxes under Select Tracks.
- OK** Opens the Export dialog box, where you can select a location for the exported BED format file. See “Export” on page 169.
- Cancel** Cancels the export and closes the dialog box.

Find in column




Figure 58 Find in column dialog box

Purpose: This dialog box lets you set search parameters for a specific column entry for the selected chromosome. Based on these parameters, the program can highlight the row of the first entry that matches. The cursor then moves to the location defined in the row.

To open: Right-click any entry in a tab in Tab View other than a column heading, then click **Find in column** in the shortcut menu.

Find in column Type all or part of the entry you want to find.

Direction Select a search direction:

- **Up** – Sets the search to move up in the selected column from the currently highlighted row.
- **Down** – Sets the search to move down in the selected column from the currently highlighted row.

Conditions Mark any of these search options:

- **Match Case** – Mark this option to take case into account. For example, if you mark **Match Case**, and you type aa351 in Find in column, the search finds the next entry in the column that contains **aa351**. It does *not* find entries that contain **AA351** or **Aa351**.
- **Match whole word** – Mark this option to only find entries in which the complete entry matches what you type in Find in column. For example, if you type AA351 in Find in column, and mark **Match whole word**, the program finds the next **AA351** entry. It does not find entries such as **AA3512** or **AA351992**.

Find Next Finds the next matching entry in the selected column, and moves the cursor to the location defined in the row that contains the entry. The search is performed only for the chromosome selected in the Genome View.

Cancel Closes the dialog box.

Go To Gene/Genomic Location




Figure 59 Go To Gene/Genomic location dialog box

Purpose: To find a specific gene location in Gene View by either selecting the RefSeq by Symbol or by selecting the Genomic Location.

To open: Click **Home** > **Go to Gene/Genomic location**.

RefSeq by Symbol Select the Reference Sequence accession symbol from NCBI, and click **Go**.

Genomic Location

- Chromosome – The chromosome number.
- Base Position – The position on the chromosome.

Click **Go** after selecting the chromosome number and the position of the gene on the chromosome.

Cancel Closes the dialog box.

Import





Figure 60 Import dialog box

Purpose: Lets you select files and import them into Agilent Genomic Workbench. The title of this dialog box changes depending on the type of file to import.

To open: In the **Home** tab, click **Import**, then select any kind of import except Genome Build or Track. The type of file to be imported appears in the title of the dialog box.

Use the standard Windows® Explorer commands in the dialog box to select a file for import.

For some imports, you can select multiple files. Click the name of the first file, then hold down the **Ctrl** key while you click the names of additional files. To select a contiguous block of files, click the name of the first file in the block, then hold down the **Shift** key while you click the name of the last one.

- File name**
- Displays the name of a file you select for import.
- Files of type**
- Lets you select the types of files to display from the types shown in the table below. To display all files, click , then select **All Files**.

File type	Extension
FE array File	*.txt
Axon array file	*.gpr
UDF file	*.txt
Design file (GEML)	*.xml
Axon design file	*.gal
Array attributes	*.txt
Experiments	*.zip
Filters	*.xml
Gene list	*.txt

- Import or Open**
- Imports the file into the program. In some cases, the name of this button is *Open*, rather than *Import*. Also, when you click **Import**, in many cases one or a series of additional dialog box(es) lets you further define the content for import. See the instructions for each specific type of import in [Chapter 2](#).
- Cancel**
- Cancels the import and closes the dialog box.

Import (experiments)




Figure 61 Import dialog box (for experiments)

Purpose: Lets you select the specific experiments within a .zip experiment file to import into the program. See [“To import an experiment file”](#) on page 52.

To open: In the **Home** tab, click **Import > Experiments**. In the dialog box that appears, select the desired .zip experiment file, then click **Import**.

Select experiments to import

These columns appear:

- **Import** – Mark the check box for the experiment(s) to import.
- **Experiment** – The names of the experiments available for import in the ZIP format experiment file.
-

Select All

Selects all of the experiments in the .zip file for import.

Deselect All

Clears all of the check boxes under Import.

OK

Imports the selected experiments into the program. If the name of an imported array design or data file matches one that is already available in the program, the Confirm overwrite dialog box appears, where you can select the data and/or design files that you want to overwrite. See [“Confirm Overwrite”](#) on page 151.

Cancel Cancels the import and closes the dialog box.

Import (filters)




Figure 62 Import (for filters) dialog box

Purpose: Lets you select the specific filters within a .zip exported filter file to import into the program. See “[To import filters](#)” on page 53.

To open: In the **Home** tab, click **Import > Filters**. In the dialog box that appears, select the desired ZIP exported filter file, then click **Import**.

Select experiments to import

These columns appear:

- **Import** – Mark the check box for the experiment(s) to import.
- **Filter** – The names of the filters available for import in the .zip filter file.
- **Type** – The type of filter

Select All Selects all of the filters in the .zip file for import.

Deselect All Clears all of the check boxes under Import.

- OK** Imports the selected filters into the program. If the name of a filter matches one that is already available in the program, the Confirm overwrite dialog box appears, where you can select the filters that you want to overwrite. See “Confirm Overwrite” on page 151.
- Cancel** Cancels the import and closes the dialog box.

Import GEML design files





Figure 63 Import GEML design files dialog box

Purpose: To display information in the design file and to remove any files that you don't want to import.

To open: In the Home tab, click **Import > Design Files > GEML File**. Select the desired *.xml design files, then click **Open**.

File Name The name(s) of the design file(s) to be imported.

ID	The Agilent ID number for the design file
Type	The application type, which can be CGH, ChIP, miRNA, or gene expression.
Species	The species for the genome build. This appears automatically when the Genome Build is selected.
Genome Build	The genome build for the design. If the genome build is not read automatically, a “?” appears. Click Genome Build and select the correct value from the list.
Status	<ul style="list-style-type: none">• Not Set – Appears if Genome Build and Species information is not shown.• Not Allowed – Appears if a Genome Build is selected that does not match the design, or if the design is a catalog or read-only design.• Overwrite – Appears when the design file has been updated and will overwrite any existing one of the same name.• Valid – Appears when the file is new.• Corrupt – Appears when the file is corrupt.
Remove	Click  to remove a specific design file from the list.
Start Import	Starts the import of the design files in the list.
Cancel	Cancels the upload and closes the dialog box.

Import Genome Build




Figure 64 Import Genome Build dialog box

Purpose: To import a new set of genome build files into Agilent Genomic Workbench. See [“To import a genome build”](#) on page 50.

To open: In the Home tab, click **Import > Genome Build**.

Species The genome’s species of origin.

Build Name The name of the build to be imported.

Refseq File The location of the RefSeq database file. This file contains chromosomal locations of genes. To select a Refseq file, click **Browse**.

CytoBand File The location of the applicable cytoband file. This file contains graphical cytoband information for Gene View and Chromosome View. To select a cytoband file, click **Browse**.

OK Imports the genome build and closes the dialog box.

Cancel Cancels the import and closes the dialog box.

CAUTION

Import only Agilent-provided genome build files.

Import Track



Figure 65 Import Track dialog box

Purpose: Lets you import a BED format track file. See [“To import tracks”](#) on page 51. Track information can appear in Gene View. See [“User Preferences”](#) on page 208.

To open: In the **Home** tab, click **Import > Track**.

- | | |
|-------------------|--|
| Species | Select the species to which the track relates. |
| Build Name | This list contains the available genome builds for the selected species. Select the desired genome build. |
| Color | Shows the currently assigned display color for the track. To change this color, click Change . For more information, see “Select Color” on page 196. You select track colors in the same way as gene list colors. |
| Track Name | Type a name to identify the imported track. |
| Track File | Type the location of the BED track file to import, or click Browse to select a file. |
| Browse | Opens an Open dialog box, where you can select the BED track file to import. |
| OK | Imports the track into the program. |
| Cancel | Cancels the import and closes the dialog box. |

Microarray Properties

Purpose: Displays the properties associated with an array. You can also edit the values of specific attributes. To add attributes to the list, see the *Sample Manager User Guide*.

To open: For any array in the **Data** folder or **Experiments** folder, right-click the array name, then click **Show Properties**. For non-Agilent arrays, only the Attribute tab appears.

Attribute tab






Figure 66 Microarray Properties dialog box with list of Attributes and their values

- **Attribute** – Displays the attributes in the array by name.
- **Value** – Indicates the values, if any, for each array. To edit the value of an attribute, select a new value for it under Value. Alternatively, click , then type or edit the value.

Close Closes the dialog box.

FE Headers tab



Attribute	FE Headers	FE Features
Index	Name	Value
1	Metric_ReproducibilityGreen_BG...	1
2	FeatureExtractor_ScanFileGUID	b4136cfe-2693-4b6c-be06-06e2...
3	OutlierFlagger_IQRatio	1.42
4	rMultDetrendSurfaceAverage	249.128
5	gOutlierFlagger_Auto_FeatB_Term	352.917
6	rAveNumPixOLLo	0.319777
7	QCMetrics_UseSpikeIns	0
8	gNegCtrlNumInliers	1467
9	AnyColorPrntSat	0.00903476
10	gDarkOffsetAverage	24.303
11	SpotAnalysis_kmeans_moi_rejec...	2.5
12	FeatureExtractor_SingleTextFile...	1
13	AnyColorPrntBGNonUnifOL	0.0217656
14	DyeNorm_RankTolerance	0.05
15	BGSubtractor_AdditiveDetrendF...	1
16	Grid_ConSnarinn	63.5

Figure 67 Microarray Properties dialog box with list of FE Headers their values

- Index** Displays a sequential index to help identify FE properties.
- Name** Displays feature parameters, statistics, and constants for the whole array.
- Value** Displays the value for each parameter, statistic, and constant.
- Close** Closes the dialog box.

FE Features tab




Figure 68 Microarray Properties dialog box with list of FE Features and associated data

Selection List Select the chromosome whose feature information you want to display.

List Box Displays FE features and the associated data. The columns are:

Index	FeatureNum	ProbeName
gIsPosAndSignif	LogRatioError	PValueLogRatio
gProcessedSignal	rProcessedSignal	gMedianSignal
rMedianSignal	gBGSubSignal	rBGSubSignal
gIsSaturated	rIsSaturated	gIsFeatNonUnifOL
rIsFeatNonUnifOL	gIsBGNonUnifOL	rIsBGNonUnifOL
rIsPosAndSignif	gIsWellAboveBG	rIsWellAboveBG

MovAvg Example Parameters




Figure 69 MovAvg Example Parameters dialog box

Purpose: This dialog box lets you set display parameters for the MovAvg Example plug-in. The plug-in calculates a 10-point moving average of the data from the selected chromosome in the selected arrays of the active experiment. It displays stacked plots of moving averages for all the arrays, one plot per array. You must have Perl installed on your computer to use this plug-in.

To open: Click **Tool > Plugin Settings > MovAvg Example**. This dialog box also opens when you click **Tool > Plugin > MovAvg Example**, if **Don't show again** is cleared.

Parameters Set any of these parameters:

Parameter	Description
X-axis Label	The text that appears under the X-axis of the plot as a label.
Y-axis Label	The text that appears next to the Y-axis of the plot as a label.
Y-axis Range (min)	The minimum value on the Y-axis.
Y-axis Range (max)	The maximum value on the Y-axis.
X-axis Range (min)	The minimum value on the X-axis.
X-axis Range (max)	The maximum value on the X-axis.

- Don't show again** Mark this option to keep this dialog box from being displayed in the future when you click Tool > Plugin > MovAvg Example. To restore the dialog box so it appears again, click **Tool > Plugin Settings > MovAvg Example**, then clear **Don't show again**.
- OK** Click to accept the parameters and prepare the plot. You can further make additional changes to the appearance of the plot once the plug-in displays it.
- Cancel** Ignores any changes you made, and closes the dialog box.

How to modify the plugin

The plug-in program (**MovAvg Example.pl**, located in the Plugins folder of the Agilent Genomic Workbench installation folder on your computer) is a short Perl program. It is a good example of how calculated columns are processed.

The plotting is very simple, but the simple plug-in architecture of MovAvg Example.pl lets you write your own computational methods to analyze data from selected arrays in the CGH application.

- Within the code of the plug-in, you can add text strings to column headers to set the format.
- To create a line graph instead of a scatter plot, you append `-plotline` to a column header.
- To prevent the plug-in from plotting a specific column, you append `-noplot` to the column heading. Note that the plug-in removes this extra text from the header before it displays it on the plot. The extra text does not appear in figures, and is only used to set the format of the plot.

MovAvg.pl shows how column-naming can be used. As you read the first line (which contains the header text), you can add text to the existing headers or add text to the headers for your generated columns, as well, to give you a small amount of formatting control.

MovAvg Example Plot




Figure 70 MovAvg Example Plot

Purpose: This plot displays the output of the MovAvg Example plug-in. The plug-in calculates a 10-point moving average of the data from the selected chromosome in the selected arrays of the active experiment.

To open: Click **OK** in the MovAvg Example Parameters dialog box. See [“MovAvg Example Parameters”](#) on page 190.

Plot(s) The main plot area shows moving average line plots for the selected chromosome. A separate plot appears for each array.

Edit Opens a menu with a Copy Plot(s) to Clipboard option. If you select this option, the program copies the plots to the clipboard as an image. You can then paste the image into a document in another program.

When you right-click anywhere within the plot area, the following options are displayed:

Option	Description
Properties	Opens the Chart Properties dialog box, where you can change the appearance of many aspects of the plot. See “Chart Properties” on page 145.
Save as	Opens a Save dialog box, where you can select a location for the *.png image file of the plots.
Print	Opens a Page Setup dialog box, where you can set page size, margins, and orientation, and select a printer. From there, you click OK to open the Print dialog box, where you can set print options and print the plot.
Zoom In	<p>Opens another menu that lets you zoom in the plot. You can zoom in several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms in the Domain (X) axis for all stacked plots, and zooms in the range Range (Y) axis for the specific plot in which you clicked.• Domain Axis – Zooms in the Domain (X) axis for all stacked plots.• Range Axis – Zooms in the Range (Y) axis for the specific plot in which you clicked. <p>You can also drag across an area of one of the plots to select an area to expand.</p>
Zoom Out	<p>Opens another menu that lets you zoom out the plot. You can zoom out several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms out the Domain (X) axis for all stacked plots, and zooms out the range Range (Y) axis for the specific plot in which you clicked.• Domain Axis – Zooms out the Domain (X) axis for all stacked plots.• Range Axis – Zooms out the Range (Y) axis for the specific plot in which you clicked.
Auto Range	<p>Opens another menu that lets you zoom the plot to display the full range of the data. You can zoom in several ways:</p> <ul style="list-style-type: none">• Both Axes – Zooms both axes of the specific plot to show the full set of data.• Domain Axis – Zooms in the Domain (X) axis for all stacked plots so this axis displays the full range of X values of the data.• Range Axis – Zooms in the Range (Y) axis for the specific plot in which you clicked to display the full range of Y values of the data.

Scroll to Column



Figure 71 Scroll to Column dialog box

Purpose: This dialog box lets you select a column. The program then scrolls the tab so that you can see the selected column.

To open: Right-click a column heading in Tab View, then click Scroll To Column in the shortcut menu.

Select column Lists the columns available in the currently selected tab. Select the one you want to view.

OK Scrolls the current tab so that you can see the selected column.

Cancel Closes the dialog box.

Search probes in eArray




Figure 72 Search probes in eArray

Purpose: To select the probes you want to update in eArray

To open: Right-click Gene View, and click **Search probes in eArray**.

Select a chromosome and a region in Chromosome View for selecting the probes related to the genes in this region.

- User Defined** Select to choose the region from which the probes to be searched in eArray will be selected. The chromosome selection list and the Start and Stop positions on the Y axis are activated when this option is selected.
- For complete gene view** All the probes related to the genes in Gene View will be searched.
- For aberrant region below cursor** Selects those probes for the genes that appear just below where the cursor sits in Gene View. Not operational without a license.
- Chromosome** If you select User Defined, you can select a different chromosome than had been selected before opening this dialog box.
- Start/Stop** If you select User Defined, you can type in Start and Stop positions for defining the region contained the genes to be in the list.

Select Color

Purpose: To select a color. Three tabs are available for selecting colors:

- Swatches tab - select colors based on samples (swatches)
- HSB tab - select colors based on an HSB schema (Hue, Saturation, and Brightness)
- RGB tab - select colors based on an RGB schema (Red-Green-Blue)

To open: This dialog box opens when a function allows you to change a color. For example, right-click on an array in an experiment, click **Edit Array Color** and click the **Swatches**, **HSB**, or **RGB** tab.

Swatches tab




Figure 73 Select Color - Swatches Tab

This tab is used to select a color based on color samples (swatches).

Preview The Preview area shows how the selected color appears. When you change the color, the original color appears at the top of the color box on the right.

Recent: Choose a recent color selection.

OK Click to select the color and close the dialog box.

Cancel Click to close the dialog box without changing the color.

Reset Click to change swatches, HSB, and RGB colors back to the default colors.

HSB Tab




Figure 74 Select Color - HSB Tab

In this tab, you can select a color based on an HSB schema (Hue, Saturation, and Brightness).

Hue Click the **H** button, and move the slider up and down, or go up and down the list of numbers, to select the hue or color of the array.

Saturation Click the **S** button, and move the slider up and down, or go up and down the list of numbers, to select the saturation level for the color.

Brightness Click the **B** button and move the slider up and down, or go up and down the list of numbers, to select the brightness level for the color.

4 Data Viewing Reference

Select Color

RGB Numbers Reflect the amount of red, green and blue in the resulting color.

Preview The Preview area shows how the selected color appears. When you change the color, the original color appears at the top of the color box on the right.

OK Click to select the color and close the dialog box.

Cancel Click to close the dialog box without changing the color.

Reset Click to change the swatches, HSB, and RGB colors back to default values.

RGB Tab




Figure 75 Select Color - RGB Tab

This tab is used to select a color based on an RGB (Red, Green Blue) schema.

Red Move the slider to change the amount of red in the color. Or, click the up or down arrow to select a number.

Green Move the slider to change the amount of green in the color. Or, click the up or down arrow to select a number.

Select data type for experiments (UDF files – CGH or CH3)

- Blue** Move the slider to change the amount of blue in the color. Or, click the up or down arrow to select a number.
- Preview** The Preview area shows how the selected color appears. When you change the color, the original color appears at the top of the color box on the right.
- OK** Click to select the color and close the dialog box.
- Cancel** Click to close the dialog box without changing the color.
- Reset** Click to return the swatches, HSB, and RGB colors back to default values.

Select data type for experiments (UDF files – CGH or CH3)




Figure 76 Select data type for experiments dialog box

Purpose: Lets you specify the mathematical form of the data in an imported UDF file, and its associated application type. See [“To import a UDF file”](#) on page 46.

To open: In the **Home** tab, click **Import > Array Files > UDF File**. In the dialog box that appears, select the desired UDF file, then click **Open**.

- Experiment Name** By default, the experiment name is the name of the imported UDF file. To change the name, double-click it, then edit it as desired.
- Data Type** Select the mathematical form of the array data in the UDF file. The options are:
- **ratio**

Data Viewing Reference
Set genome build and species for Axon design files

- **log₂ ratio**
- **log₁₀ ratio**
- **ln ratio** (base e)

Design type Select the application type (CGH or CH3, for example) associated with the array data in the UDF file.

Continue Accepts your selections, and goes to the next step in the UDF import process.

Cancel Cancels the UDF import.

Set genome build and species for Axon design files





Figure 77 Set genome build and species for Axon design files dialog box

Purpose: Lets you set the species and genome builds associated with imported Axon design file(s), and to remove specific designs files from the import, if necessary. See [“To import Axon design files”](#) on page 44.

To open: In the **Home** tab, click **Import > Design Files > Axon File**. In the dialog box that appears, select at least one Axon design file, then click **Import**.

No.	An index number within the dialog box for each Axon file.
File Name	The names of each Axon design file selected for import.
Species	The species associated with each design file. If a species is incorrect, select the correct one from the appropriate list.
Genome Build	The genome build associated with each of the design files. If a genome build is incorrect, select the correct one from the appropriate list.
Status	<p>The status of the file is one of the following:</p> <ul style="list-style-type: none"> • Valid – The file is a new file that can be imported. • Overwrite – The file is a valid design file, but when you import it, it will replace an existing design that has the same name. • Corrupt – The file failed validation. When you start the import process, the program ignores the file.
Remove	Click  to remove a specific design file from the list. This can be useful if you select a design file in error, or if you do not want to overwrite an existing one.
Start Import	Imports the file(s) and closes the dialog box.
Cancel	Cancels the import and closes the dialog box.

Show/Hide Columns




Figure 78 Show/Hide Columns dialog box

Purpose: Used to select the attributes to be displayed in the Experiment Attributes dialog box and the Sample Utility tab. The Sample Utility tab is available when you go to Sample Manager. See the *Sample Manager User Guide* for information about Sample Manager.

To open: This dialog box appears when you click **Show/Hide Attributes** at the bottom of the Experiment Attributes dialog box.

All available attributes are shown in the Attributes column. Attributes with a check-mark next to them will be displayed in the Experiment Attributes and Sample Utilities tab for each sample. To select an attribute for display, mark the **Show in Table** box next to it. To deselect an attribute, clear the **Show in Table** box again.

Save Saves the current list of selected attributes and updates the Sample Utilities table based on the selections.

Select All Selects all the attributes in the list.

Deselect All Clears all check marks from attributes in the list.

Close Closes the dialog box. If changes have been made, the program asks if you want to save your changes before closing.

Track




Figure 79 Track details

Purpose: This dialog box lets you view the chromosome locations in the track.

To open: Click the **Details** link for the desired track in the **Tracks** tab of the Preferences dialog box. See “[User Preferences](#)” on page 208.

Track Parameters These parameters appear:

Parameter	Description
Name	The name of the track.
Species	The species to which the track applies.
Format	The format of the track data. Agilent Genomic Workbench supports the BED format.
Genome Build	The specific genome build of the species to which the track applies.
Description	Descriptive text saved with the track.

Data Tracks must contain entries for at least these four columns in the table:

Column	Description
Chromosome	The name of the chromosome
Start	The first base pair of the particular feature in the chromosome.
Stop	The last base pair of the particular feature in the chromosome.
Name	The name of the feature. This name appears next to the defined region for the feature.

The other columns are additional BED track file columns that can appear for some tracks. Agilent Genomic Workbench does not display these.

Close Closes the Track dialog box.

UDF Import Summary (CGH or CH3)




Figure 80 UDF Import Summary dialog box

Purpose: Reports how many lines of data were successfully imported from a UDF file, and how many lines were skipped. Skipped lines can be caused by missing chromosome mapping information, or improper formatting of the UDF file.

To open: Import a UDF file (see [“To import a UDF file”](#) on page 46). This dialog box appears after you map the columns of the UDF file.

Table Displays the file name of the imported UDF file, the number of lines that were successfully imported, and the number of lines, if any, that were skipped during import. If many lines were skipped, review the data for improper formatting or missing chromosome mapping information.

OK Closes the dialog box.

Universal Data Importer - Map Column Headers

chr1	727595	754477	LOC643837	1000	+	754477	754477	0
Select	Select	Select	Select	Select	Select	Select	Select	Select
chr1	835123	855339	SAMD11	1000	+	835324	854913	0
chr1	854965	870958	NOC2L	1000	-	855451	870899	0
chr1	854969	870742	LOC401010	1000	-	870742	870742	0
chr1	872228	877875	KLHL17	1000	+	872334	877350	0
chr1	878658	885682	PLEKHN1	1000	+	878693	885159	0
chr1	904335	905548	HE54	1000	-	904431	905354	0
chr1	904337	905548	HE54	1000	-	904431	905354	0
chr1	922073	923139	ISG15	1000	+	922582	923078	0
chr1	929321	961320	AGRN	1000	+	929321	960189	0
chr1	999846	1001086	CTSL	1000	+	999930	1000534	0

Figure 81 Universal Data Importer - Map Column Headers dialog box

Purpose: Lets you set up a universal data file (UDF) for import. You define several properties associated with the UDF, and identify the contents of each column of data in the file. You can also save column mappings for re-use.

To open: As you go through the UDF import process (see “To import a UDF file” on page 46), in the Select data type for experiments dialog box, click **Continue**. See “Select data type for experiments (UDF files – CGH or CH3)” on page 199.

Species Info

Select Species Select the species associated with the array data in the UDF. The program supports these species:

Select Genome Build Sets the species-specific build to use.

Mapping Info

Select Mapping Applies an existing column map to the current UDF. A column map identifies the contents of each column of data. To create a new column map for the current UDF, select **CUSTOM**.

Save Mapping As Saves the column map under a new name. Opens an Input dialog box, where you can type a name for the new map.

ArrayID Info

Virtual Array ID A number that uniquely identifies the data in the UDF. Typically, an Agilent microarray slide has a physical Array ID that enables Agilent Genomic Workbench to track the data from the slide as it goes through the steps of an analysis workflow. A “virtual” Array ID is, by default, a system-generated ID that serves the same purpose for data from UDFs. You can also create your own virtual Array ID.

Use System Generated Barcode By default, the virtual Array ID assigned to the array data in a UDF is a number that is created by the program. To create your own Array ID, clear **Use System Generated Array ID**, then type a new number in **Virtual Array ID**.

Table

This table lets you identify the contents of the columns of data in the UDF. The first row of the table gives the column heading information from the UDF. The second row contains lists of labels that you apply to each column, and the rest of the table displays lines of data from the UDF. If the UDF contains data from Agilent CGH arrays, the column headings will exactly match the labels in the lists.

In the list below each column heading, select the applicable label. You must use each of the labels exactly once, except LogRatio, which you can use more than once. These labels are available:

Column Label	This column contains:
ProbeName	Names of probes.
ChrName	Names of chromosomes.

Column Label	This column contains:
Start	First chromosomal location for each probe.
Stop	Last chromosomal location for each probe.
Description	Text annotation for the probe.
LogRatio	Array data values that correspond to each probe. You can use this label more than once.

NOTE

If you select a saved column mapping, then change or reset the column labels in the table, the program changes or resets the saved column map as well.

- Reset** Clears all the column labels in the second row of the table. If you have selected a saved column mapping, this command also clears the labels in the saved map.
- Import** Imports the UDF file with the specified parameters, and opens the UDF Import Summary dialog box (see [“UDF Import Summary \(CGH or CH3\)”](#) on page 205).
- Cancel** Cancels the import and closes the dialog box.

User Preferences

Purpose: This dialog box is used to set up preferences for display of tracks, data storage locations, and licenses.

To open: From the Home tab, click **User Preferences**. Or, right-click in the Gene View, Chromosome View, or Genome View, and click **User Preferences**.

Tracks tab




Figure 82 User Preferences dialog box - Tracks tab

Purpose: To import and set up the appearance of tracks next to the Gene View. Tracks are additional graphic displays of genomic information loaded from an external file. They align with genomic coordinates in Gene View.

To open: In the User Preferences dialog box, click the **Tracks** tab.

Font Options

Select the font type, style and size for the gene annotations that appear in the selected tracks.

Tracks List

Track Name	Name of the track already loaded or imported
Show in UI	Mark the check box to display the track next to Gene View.
Show in Report	Mark the check box to display the track information in all the reports.
Genomic Boundaries	Click to use the track to define only the regions that aberration detection algorithms will run. You can choose to do this for only one track.
Delete	Mark the check box to delete the track from the list. Then, click Delete to delete the track from the list.
Details	Click to display all the chromosome locations defined in the track.
Import	Click to import new tracks.
Delete	Click to delete the tracks selected in the Delete column.
Up	Click to move a track up the list.
Down	Click to move a track down the list.

Visualization Parameters

Genes	<p>These options affect the appearance of the Track and Gene View.</p> <ul style="list-style-type: none">• Orientation – Type a number to set the angle at which the Gene Symbols will appear in Gene View and the Track Annotations appear in the tracks.• Show Gene Symbols – Mark to show gene symbols in Gene View, and clear the check box to hide them.
Genomic Boundaries	<p>These options let you include or exclude the Genomic Boundaries from the analysis.</p>
Tracks	<p>These options affect the appearance of the Track Views.</p> <ul style="list-style-type: none">• Show Annotations – Mark to show the names of the gene regions for the tracks, and clear to hide them.• Show Overlaid – Mark to overlay all the tracks that appear next to Gene View, and clear the check box to display the information in separate tracks.

Miscellaneous tab




Figure 83 User Preferences dialog box – Miscellaneous tab

Purpose: For data/content set-up, this dialog box allows you to set up eArray access and to change the location for data.

To open: In the User Preferences dialog box, click the **Miscellaneous** tab.

eArray User Details

Sets login details for the Agilent eArray Web site.

- **URL** – At present, <https://earray.chem.agilent.com>
- **Username** – The name registered on the eArray site.
- **Password** – The password registered on the eArray site.

Error Model

The DLRErrorModel (Derivative Log Ratio) is the only selection. This measures noise in the data for CGH analyses.

Data Location

The folder where the program stores array data and design files. To select a location, click **Browse**.

4 Data Viewing Reference

User Preferences

Apply Applies any changes to the preferences.

OK Accepts any changes and closes the dialog box.

Cancel Cancels all changes and closes the dialog box.

License tab




Figure 84 User Preferences dialog box – License tab

Purpose: The License tab allows you to display and update your DNA Analytics application license(s). The license enables the analysis application, and allows you to use it to analyze array data.

To open: In the User Preferences dialog box, click the **License** tab.

Host Name Displays the host computer name automatically.

- Select Analysis Application** Select the Agilent Genomic Workbench application for which you have a license.
- Server Location** Select this option if you have a concurrent user license. To edit this name, select **Server Location**, then type the path for the folder where your licenses are located. If you select this option, the Text License option is unavailable.
- Text License** Select this option if you have an application license (CGH, ChIP, CH3). To change the license, delete the old license text, and paste the new license text in the box.
- OK** Accepts any changes you have made, and closes the dialog box.
- Cancel** Closes the dialog box without changing any license information.
- Apply** Accepts any changes you have made, but does not close the dialog box.

View coordinates in UCSC browser




Figure 85 View coordinates in UCSC browser

Purpose: Defines a track to upload to the UCSC Web site so that you can see the information in the UCSC Genome Browser.

To open: Right-click in the Gene View, and select **Show in UCSC**.

Name Type a name for the track. This name identifies the track when it appears in lists and displays.

Build (Available if you select **User Defined** in **Set Chromosome Start-Stop**.)
Select the genome build with which to associate the track.

Description Type descriptive text to attach to the track for reference.

Set Chromosome Start-Stop This parameter defines the region of the chromosome for which the track will be defined. Select one of these options:

- **User Defined** – Lets you define an arbitrary region of any chromosome. If you select this option, select the desired chromosome in **Chromosome**, then type the beginning (**Start**) and end (**Stop**) locations of the desired interval.
- **For complete gene view** – The chromosomal region that appears in Gene View.

Save as Track in Genomic Workbench	Mark the check box to save this track in the Tracks folder in the My Entity List pane of the Navigator.
Change	Click to open the Choose Track Color dialog box to select the color to use for display of the track in the Tracks folder. See “ Select Color ” on page 196.
OK	Creates the track and opens the UCSC Web site, where you can display the track and associated information. For information on using the UCSC Web site, see the help and information provided there.
Cancel	Closes the dialog box without creating a track.

View Preferences




Figure 86 View Preferences dialog box for CGH

Purpose: This dialog box allows you to configure how data and results appear in Genome, Chromosome, and Gene views.

To open: In the **View** tab, click **View Preferences**.

NOTE

The View Preferences dialog box contents changes depending on what application is selected. For information on View Preferences for ChIP and CH3 applications, see the User Guide for the application.

View Alignment Selects the orientation and rendering style (described below).

Option	Description
Orientation	
Horizontal	Stacks Genome, Chromosome, and Gene views horizontally in the main program window. Genomic locations appear across the bottom of each view.
Vertical	Displays Genome, Chromosome, and Gene views from left to right as side-by-side panes in the main program window.
Rendering Style	
Overlaid	In Chromosome View and in Gene View, displays data and results as a single, combined pane for all arrays. (Default)
Stacked	In Chromosome View and in Gene View, displays a separate pane for each array.

Data Visibility For each view, or all views, selects the kind(s) of data and results to display.

In **View**, select the view you want to configure. To set availability of display items for all views, select **All views**. Some display items are only available for certain views. When you select a display item, it enables the item for display – for some items, you must take additional steps to display them. For example, you may need to configure a specific algorithm in the toolbar.

Mark any of the following options, as available:

Display item	Description/Comments
Scatter Plot	The plot(s) of individual log ratio data points.
Scatter Tool Tip	The ToolTips that appear when you place the pointer over specific data points on the scatter plot(s) in Gene View. The tool tip shows the array of origin and the numerical log ratio value for the data point.
Moving Average	The result of the Moving Average algorithm. See the <i>CGH Interactive Analysis User Guide</i> for more information.
Aberration	The result of the selected aberration detection algorithm. See the <i>CGH Interactive Analysis User Guide</i> for more information.

Display item	Description/ Comments
CNVR	Detected copy number variant regions. See the <i>CGH Interactive Analysis User Guide</i> for more information.
Log Ratio Error Envelope	The log ratio error envelope is a visual representation of the log ratio error calculated by Feature Extraction.
Penetrance plot	The probe penetrance plot for the active experiment. If you select this option, all other display items are unavailable. In addition, because the probe penetrance plot takes into account all arrays, this option overrides the <i>stacked</i> rendering style.
Common Aberration	The results of a common aberration analysis. To display this, you must first perform a common aberration analysis. See the <i>CGH Interactive Analysis User Guide</i> for more information.
Green Intensity	Mark the check box to display green raw signal intensity.
Red Intensity	Mark the check box to display red raw signal intensity.

**Rendering
Patterns**

These options control the specific appearance of data and results in Genome, Chromosome, and Gene views. You configure these options separately for each type of application design.

- **Design Type** – Select the application design type for which you want to define rendering patterns.
- **Styles** – Select the display style for each of these elements:

Display element	Details
Scatter Plot	Select the symbol used for log ratio data points in the scatter plots in Chromosome and Gene views.

Display element	Details
Moving Average	Select the line style for the moving average display. Lines appear in the display color defined for each array. See the <i>CGH Interactive Analysis User Guide</i> for more information. <ul style="list-style-type: none">• Continuous – A solid line.• Dashed – A dashed line.• Dotted – A dotted line.• Do not show area – No line.
Aberration	Select the rendering style for detected aberrations. <ul style="list-style-type: none">• Semi transparent filled – Solid, colored regions (in the display colors defined for each array, if applicable).• Hatched – Cross-hatched colored lines (in the display colors defined for each array, if applicable).• Do not show area – Aberrations do not appear.

Scatter Plot (Chr View) Point Size Select a point size to use for display of scatter plot data points in the Chromosome View.

NOTE Rendering scatter plots for more than 10 high density arrays in the Chromosome View may take significant time. Selecting filled circles as the rendering style for CGH scatter plots can also decrease performance. For faster performance, change the rendering style for CGH data from the filled circle to the plus (+) or cross hair sign.

Configure Scales For Log Ratios or Signal Intensities plots, mark **Apply** to enable the custom scale. In Range, type the value to use as the range for the scatter plot.

Configure Coloring schemes Use these options to change the display of the scatter plot in the Gene View. These options are the same as those displayed in the Scatter Plot ToolTip in the Gene View.

Show Memory Monitor in Status Bar Displays a memory usage monitor in the eighth cell of the status bar. For information about the Status Bar, see “[Status Bar](#)” on page 138.

OK Applies the changes you made to all preferences and closes the dialog box.

Cancel Closes the dialog box without applying changes.

4 Data Viewing Reference

View Preferences

Apply Applies changes without closing the dialog box.

www.agilent.com

In this book

This guide describes how to import, organize, manage, export and display data and other content within Agilent Genomic Workbench if you don't have any DNA Analytics application license(s).

© Agilent Technologies, Inc. 2010

Revision A, May 2010



G3800-90014



Agilent Technologies