




High-Pressure Injection Device for the Agilent 7890A and 6890 Series Gas Chromatographs

Accessory G3505A

Introduction

Gas chromatography sampling and representative analysis of highly volatile liquefied hydrocarbons with high precision and accuracy can be challenging. In the solution described here, a unique sample injection device based on a needle interface and liquid rotary valve, has been designed for sampling light petroleum matrices with broad boiling point distributions. The 7890A GC-based system consists of a 4-port liquid valve, a deactivated removable needle, and an auxiliary flow. The needle is directly installed on one port of the valve. This compact device is installed directly over the top of a split/splitless inlet. The unit is operated automatically just like a typical liquid autosampler; however, the needle is not withdrawn. Various pressurized liquid samples have been run on this device, such as liquefied natural gas (calibration standard), ethylene, propylene, and butadiene. Excellent repeatability is obtained with RSDs typically below 1% in quantitative analyses.

Injection Device

The high-pressure injection device (HPLI) consists of components as shown in Figure 1.

- **Valve:** Internal sample valve from Valco Instruments Co. Inc. 4-port equipped with a sample volume of 0.06 μ L. Other rotor sizes are available from Valco Instruments Company.
- **EPC:** An auxiliary flow from a 7890A Aux module is connected to port P. In sample analysis, the flow can be set at 50 mL/min to 200 mL/min. The higher auxiliary flow gives better peak shape.

Ordering Information

Order accessory G3505A. The accessory is compatible with both the 7890A and 6890 series GCs.

The following components are recommended. These are not supplied in the accessory kit.

- **Filter:** To remove particles from samples.
- **Restrictor:** To maintain sample pressure, a metering valve (Agilent PN 101-0355) is connected to the end of the sample exit line tubing. Restrictor is not included in accessory kit.

Guideline for choosing Aux flow source

7890AGC

G3471A Pneumatic Control Module (PCM) or
G3470A Aux EPC module

6890GC

G1570A Aux EPC or
G2317A PCM module

The PCM is the preferred source for both GCs.



Agilent Technologies




Figure 1. Flow diagram of the high-pressure injection device (HPLI).

Typical Instrumental Conditions

Gas chromatograph	Agilent 7890A
Injection source	High-pressure injection device (HPLI) at near ambient temperature
Injection port	Split/splitless, 250 °C (350 °C for C5–C40)
Sample size	0.06 µL
Carrier gas	Helium
Aux or PCM	150 mL/min (Helium)
FID	250 °C (350 °C for C5–C40) H ₂ , 35 mL/min Air, 400 mL/min



Agilent pneumatic air actuator/valve assembly installed on the 7890A.

Sample Chromatograms

Pressurized Propylene

This sample is analyzed by the same conditions as in ASTM D6159. A typical chromatogram is shown in Figure 2.




Figure 2. Chromatogram of pressurized propylene.

Pressurized 1,3-Butadiene

Figure 3 is an example of C4 hydrocarbons analysis showing 1,3 butadiene purity.




Figure 3. Chromatogram of pressurized 1,3-butadiene.

Summary

A unique sample injection device for the Agilent 7890A GC based on a unique deactivated interface and liquid rotary valve has been designed for sampling light petroleum matrices with broad boiling point distributions from methane to as high as C40. It is installed directly over a split/splitless GC split/splitless inlet in a few minutes. The maximum sample pressure is 3,000 psig, although typical samples will have pressures under 1,500 psig. Various pressurized liquid samples have been tested on this device with high accuracy and precision. The sampler is quick to install and easy to operate. As with all high-pressure sampling systems, appropriate safety precautions must be followed.

Competitive Advantages

The HPLI can be used with a wide variety of sample streams or pressurized vessels. Because the sampling valve is interfaced directly to the inlet with an inert needle, loss or adsorption of

analytes is minimized compared to conventional liquid sample valve systems. Compared to other gas chromatographic vaporizers for handling pressurized or nonpressurized samples, the Agilent HPLI has the following advantages:

- Better results with polar analytes such as glycols
- Superior inertness
- Low discrimination (no discrimination up to C₁₆)
- Flexibility: Install or uninstall in less than 10 minutes
- Good for trace impurity analysis with 0.5 µL rotor
- Excellent repeatability, typically RSDs below 1 %

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2008

Printed in the USA
February 25, 2008
5989-8037EN



Agilent Technologies