Agilent I9131
Sodium
Combination ISE
钠离子复合电极

Operating Guide
用户手册

Agilent Technologies

The refillable Sodium Combination ISE can measure sodium ion (Na^+) concentration in laboratory. This electrode is used with the 3200I Ion Meter or similar meters.

WARNING

Use this electrode according to the operating manual to avoid personal injury.

WARNING

The electrode solution can cause chemical burns or illness if it is taken orally or contacted by human skin. Use protective clothing or gloves to avoid contact. In case of contact, rinse contacted area with tap water or DI water thoroughly.

CAUTION

The electrode body material is glass. Handle with care to avoid damage to the instrument.

Table 1 Sodium electrode specifications

Specification	Value
Concentration range	Saturated, 10^{-6} mol/L
Temperature range	0 to 60 °C
Electrode impedance	≤ 250 M Ω (25 °C)
Reference system	Ag/AgCl
Liquid junction material	Ceramic
Electrode diameter	12 mm
Electrode length	120 mm
Cable connector	BNC
Cable length	1000 mm


Figure 1 Electrode assembly

Operation

Preparing the electrode

The electrode should be stored wet in the storage solution. It is in working condition when you receive it. Observe visually if there is any mechanical damage. If the storage bottle for this electrode is dry, soak the electrode in sodium ISE reference solution for at least 24 hours before use.

- 1** Add diisopropylamine (AR grade) to the calibration solution until its pH value is over 10.
- 2** Take the storage bottle off and store it upright for future use.
- 3** Rinse the measuring tip with distilled or deionized water.
- 4** Unplug the filling hole.
- 5** Hold the electrode measuring tip downwards and swing it for several times to remove air bubbles near the sensitive glass bulb.
- 6** Connect the electrode to the meter. Fix the electrode vertically at the electrode holder.
- 7** Revitalize the electrode. Soak the measuring tip in distilled or deionized water. Turn on the magnetic stirrer to stir the deionized water with a steady speed. Soak for 8 hours.
- 8** Check that the mV reading is higher than 310 mV. If not, repeat [step 7](#) with new deionized water until the mV reading is higher than 310 mV.
- 9** Remove the electrode from the deionized water solution and absorb away the water from the measuring tip. Do not rub.
- 10** Refer to relative professional standards to calibrate the instrument.

Calibration

If the meter is an ion meter:

- 1 Add diisopropylamine (AR grade) into the calibration solution until its pH value is over 10.
- 2 Calibrate the electrode according to the operation manual of the meter. Choose more than two NaCl calibration solutions.
- 3 Calibrate the electrode in order from the low to the high concentration.
- 4 Exchange the calibration solution during calibrating, rinse the measuring tip thoroughly and adsorb away the water on the measuring tip.

If the meter is a mV meter (such as pH meter):

- 1 Add diisopropylamine (AR grade) into the calibration solution until its pH value is over 10.
- 2 Calibrate the electrode in more than two NaCl calibration solutions in order from the low to the high concentration.
- 3 Record mV readings in each calibration solution.
- 4 Plot 'mV-pNa' linear diagram or calculate the linear equation for the electrode (pNa is negative logarithm of sodium concentration).
- 5 Exchange the calibration solution during calibrating, rinse the measuring tip thoroughly and adsorb away the water on the measuring tip.

After electrode calibration, repeat step 2 in [“Preparing the electrode”](#) on page 6.

Measurement

- 1 Soak the measuring tip in a sample solution and add diisopropylamine (AR grade) into the sample solution until its pH value is over 10.
- 2 Record the readings of the meter when the readings become stable.
 - a If the meter is an ion meter, read the concentration value of sample solution directly.
 - b If the meter is a mV meter (such as pH meter), input the mV reading into the 'mV-pNa' linear diagram or the linear equation for the electrode to calculate the pNa value of sample solution.
- 3 Rinse the measuring tip with DI water. Adsorb away the water on the measuring tip (do not rub hard).
- 4 Plug the filling hole.

Storage

Short-term

Soak the measuring tip in DI water.

Long-term

- 1 Place the electrode inside the storage bottle.
- 2 To protect the measuring tip, keep 5-10 mm between the bottom of storage bottle and the measuring tip of ISE when placing the ISE electrode.
- 3 Screw the cap onto the electrode storage bottle firmly and place it in the package box.
- 4 Store the electrode in a dry area.

Operating Hints

- The main material of the measuring tip is glass. Ensure the sample solution will not damage the measuring tip before measurement.

- During measurement, the level of filling solution in the sodium combination ISE must be higher than the level of sample solution.
- During calibration or measurement, rinse the measuring tip with DI water and adsorb away the water on the measuring tip (do not rub hard) before measuring a new sample solution.
- Control the calibration solution and sample solution at a uniform flow rate and temperature by swinging the electrode gently or by turning on the magnetic stirrer. If not the measuring accuracy will be affected.
- Do not soak the electrode in sample solution for a long time. Rinse the electrode with DI water thoroughly as soon as the measurement is completed.
- The sodium concentration range of the calibration solutions should cover the sodium concentration range of sample solution. If not, the measurement errors will occur.
- Do not apply force onto the electrode cap, cable or cable interface. This could cause damage.
- Keep the cable interface dry.

Maintenance

Reference filling solution

- 1 Drain the electrode filling solution away.
- 2 Add fresh reference filling solution (p/n 5190-0545) through the filling hole until the level is 5 mm lower than the filling hole.

Repeat steps 1 and 2 to ensure a fresh solution.

Cleaning of inorganics

Soak the measuring tip in 0.1 mol/L HCl or EDTA solution for 15 minutes.

Cleaning of organics

Soak the measuring tip in absolute ethyl alcohol or another solvent that can dissolve organics for 15 minutes.

Cleaning of grease

Soak the measuring tip in warm weakly alkaline detergent for 15 minutes.

Cleaning of protein precipitation

Soak the measuring tip in 0.1 mol/L HCl solution that contains 1% pepsin for 15 minutes.

Electrode restoration

Rub the flat surface of the measuring tip gently with a fine gauze soaked with absolute ethyl alcohol.

After any of the above cleaning procedures, replace the filling solution. Soak the measuring tip in deionized water for 2 hours.

Troubleshooting

Meters

Refer to the meter operating manual. Check all relevant parts, such as electrode, calibration solution, and samples.

Electrode

- 1 Check the mV value of the calibration solution.
- 2 Adjust the concentration of the calibration solution to ten times the original value, and measure mV again.

The difference in mV should be more than 52 mV. If not, perform maintenance as described above. Make sure the calibration solution is effective. Use distilled or deionized water that meets requirements to make sure the calibration solution is not contaminated or beyond its shelf life.

Calibration solution

Check the calibration solution. The distilled or DI water used in the solution should meet applicable requirements. Check if the solution is expired or contaminated. Prepare fresh calibration solution if needed and re-test.

Method

Refer to relative documents to ensure that the analytical method is suitable for the sample solution.

If any other problems are observed during electrode use, contact our customer service department.

To purchase a new electrode, contact distributors of Agilent Technologies or log on the Agilent official website.

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© Agilent Technologies, Inc. 2012