

Agilent ICP-MS ChemStation – Complying with 21 CFR Part 11

Application Note

Overview

Part 11 in Title 21 of the Code of Federal Regulations includes the US Federal guidelines for storing and protecting electronic records and applying electronic signatures. The intent of these guidelines is to ensure that electronic records subject to these guidelines are reliable, authentic and maintained with high integrity. This document examines each section of 21 CFR Part 11 and provides a recommended remediation approach using the Agilent ICP-MS ChemStation with User Access Control Pack enabled, in combination with the Agilent OpenLAB ECM solution for electronic records management. The User Access Control Pack provides the necessary controls for managing system access and audit trail functionality while Agilent OpenLAB ECM ensures secure record keeping and data archival. Agilent OpenLAB ECM is proven and has been deployed at many leading life sciences companies to satisfy compliance mandates for 21 CFR Part 11.

	21 (CFR 11	Sections	s (🖌 = a	pplical	ble / N	I/A = I	not appli	cable)	
Possible scenarios with ChemStation operated in a closed system	1.1 , 11.2 , 11.3	11.10	11.30	11.50	11.70	11.100	11.200 (a)	11.200 (b)	11.300 (a) , (b) , (d)	11.300 (c) , (e)
Electronic Record only (without signature)	V	~	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Electronic Signature Based upon ID Code & Password	~	~	N/A	~	~	~	~	N/A	~	N/A

Figure 1

Applicable sections of 21 CFR Part 11



Section	Result	Question	Answer
11.10a	Yes	Has the system been validated in order to ensure accuracy, reliability, consistent intended performance, and the ability to discern invalid or altered records?	Agilent develops its products according to the well estab- lished "product lifecycle" concept, which is a phase review process for soft- and hardware development, in order to ensure consistent product quality. As a result a fully quali- fied data handling system is delivered together with all nec- essary services, which are needed to implement such a sys- tem to meet the requirements of the FDA regarding 21CFR Part11. As part of this service a document is compiled that describes the installed configuration and documents the results of the executed IQ and OQ procedures.
			Electronic records generated by the ICP-MS ChemStation are securely stored in OpenLAB ECM. OpenLAB ECM's per- formance has been extensively validated with tests written to specifically evaluate accuracy, reliability and consistent performance. OpenLAB ECM incorporates the use of byte order dependant check sums at each file transfer operation to ensure that records are valid and unaltered.
11.10b	Yes	Is the system capable of generat- ing accurate and complete copies of all required records in both human readable and elec- tronic form suitable for inspec- tion, review and copying by the FDA?	Electronic records are created in electronic as well as human-readable form. OpenLAB ECM stores all data types, from raw machine data to printable reports. All files are stored complete and unaltered in the original format. To read the electronic format the Agilent ICP-MS ChemStation is required. "Printed" reports, representing the human-read- able form of electronic records, can be stored as PDF files which can be made available for review without the source application being installed on the client machine. An MS viewer is available to view the original electronic record without the originating application.
11.10c	Yes	Are the records protected to enable the accurate and ready retrieval throughout the record retention period?	Raw, meta and result data generated by the Agilent ICP-MS ChemStation are stored and managed in OpenLAB ECM. It resides in a protected storage location and/or archive media. When archived, the media may be on-line, near-line or off- line. Regardless of the physical location of the data, it remains searchable to all users with appropriate privilege. The individual users do not need access to the physical stor- age location of the files.
11.10d	Yes	Is system access limited to authorized individuals?	Each user is identified by a unique user-ID/password combi- nation. Logging on to the ICP-MS ChemStation requires the entry of both identification components. The Windows user management is used for setting up authorized users and assigning user privileges. Users gain access to the ICP-MS ChemStation by being added to one of the three ChemSta- tion user groups that are created upon installation. Access

11.10 Controls for Closed Systems

Section	Result	Question	Answer
			to data maintained in OpenLAB ECM is controlled through a user name, password, and account login. The user management can be integrated with the Windows user management. Once inside the OpenLAB ECM reposi- tory, all file and software functionality access is con- trolled through privileges and roles assigned to individual users or groups of users.
			Setting up and maintaining users as well as specifying security policies and determining the level of access is the responsibility of the system administrator using stan- dard Windows administration functionality (e.g. forcing users to specify a new password during the first logon, defining password expiration periods, number of logon attempts until account is locked etc.)
			The Windows lock function locks the session after a specified period of inactivity. The timeframe is con- figurable by the system administrator. Both identification components are required to start a new session when the inactivity timeout locks the session in order to regain access to the system.
11.10e	Yes	Is there a secure, computer-gen- erated, time-stamped audit trail that independently records the date and time of operator entries and actions that create, modify, or delete electronic records?	Logon and logoff to the ICP-MS ChemStation is docu- mented in the Windows application log. All OpenLAB ECM activities (e.g. data storage, versioning, electronic signatures) are recorded in a secure, computer generated, time stamped audit trail. Entries in the OpenLAB ECM audit trail are non-editable, non-deletable. The user has no influence on the audit trail, e.g. it cannot be switched off or altered nor deleted. Removing records from the database does not affect existing entries in the audit trail. The audit trail lists all modifications with date/time stamp and the user name of the user doing the change.
11.10e	Yes	When records are changed, is previously recorded informa- tion left unchanged?	All entries in the OpenLAB ECM audit trail are non- editable, non-deletable. Removing records from the data- base does not affect existing entries in the audit trail.
11.10e	Yes	Are electronic audit trails kept for a period at least as long as their subject electronic records and available for agency review and copying?	Audit trail entries are stored in the OpenLAB ECM reposi- tory as part of a file's metadata and are kept throughout the electronic records retention period. The audit trail may be reviewed and printed. The Windows application log can also be stored in the same OpenLAB ECM reposi-
11.10(f)	Yes	Are operational system checks used to enforce permitted sequencing of steps and events?	tory. In all ICP-MS ChemStation and OpenLAB ECM functions, when a sequencing of events is required, system checks enforce it.

Section	Result	Question	Answer
11.10(g)	Yes	Are authority checks in place to ensure that only authorized individuals can use the system, electronically sign a record, access the operation or comput- er system input or output device, alter a record, or per-	Each user must logon to the software with a user ID / pass- word combination before they can use the system. This applies at initial program start and after every inactivity timeout on the computer program. Only a successful logon to the system offers access to the chromatographic software functions such as data acquisition, data review, result approval or archiving functionality.
		form the operation at hand?	Users cannot gain access to OpenLAB ECM without a valid user name, password and account. Once logged in, that user's access to files and software functionality (including but not limited to signing a file, inputting values, or altering a record) are determined by the privileges assigned.
11.10(h)	Yes	Are device checks used to deter- mine, as appropriate, the valid- ity of the source of data or oper- ational instruction?	User entry fields provide feedback to the user about the entry types and ranges that are valid for that field.
11.10(i)	Yes	Do the persons who develop, maintain, or use electronic records/signature systems have the education, training, and experience to perform their assigned tasks?	Records of the educational and employment history of employees are verified and kept with personnel records that can be made available during an on-site audit. In addition, all Agilent Technologies employees who deal with regula- tions have attended training workshops on regulatory requirements.
			For system users Agilent provides a basic familiarization during the installation of the product. Training courses for administrators as well as users are available.
11.10(j)	N/A	Have written policies been established, and adhered to, that hold individuals accountable and responsible for actions ini- tiated under their e-signatures in order to deter record and sig- nature falsification?	It is the responsibility of the organization implementing electronic signatures to develop written policies that ensure that individuals responsible for signing documents under- stand that their electronic signature is as equally binding as their handwritten signature.
11.10(k) (1)	N/A	Are there adequate controls over the distribution of, access to, and use of documentation for system operation and mainte- nance?	While documentation is available for ICP-MS ChemStation and OpenLAB ECM users and administrators; controls over the storage and distribution of this material are the respon- sibility of the organization that implements and uses the sys- tem.
11.10(k) (2)	N/A	Are there formal revision and change control procedures to maintain an audit trail that documents time-sequenced development and modification of systems documentation?	The quality process includes formal written revision and change control procedures for system documentation. OpenLAB ECM can be used for development and mainte- nance of system documentation. All revisions to the documents are kept time stamped and audit trailed.

11.30 Controls for Open Systems

A closed system is required. Agilent ICP-MS ChemStation is not designed to operate as an open system.

Section	Result	Question	Answer
11.30	N/A	Are there procedures and con- trols used to protect the authen- ticity, integrity and confiden- tiality of the electronic records from their creation point to the point of their receipt?	When a file is transferred to or within OpenLAB ECM, a byte order dependant checksum is calculated on the file in its source location. A copy of the file is made of the file in the destination location where a second checksum is calculated. The two values are compared and only if they are identical, is the transfer complete. If the values do not match, an error message is generated.
11.30	N/A	Are additional measures used to ensure the confidentiality of the electronic records from the point of their creation to the point of their receipt?	OpenLAB ECM supports the use of Secure Socket Layer (SSL) encryption for security during data transmission. SSL breaks a single file into very small data packets. These data packets are individually encrypted with con- figurable 64-bit or 128-bit encryption before being trans- mitted. On the receiving side the data packets are decrypted and reassembled.

11.50 Signature Manifestation

Section	Result	Question	Answer
11.50 (a)	Yes	Do the signed electronic records contain information associated with the signing that clearly indicates the following: 1. Printed name of signer, 2. Date and time that the signature was executed 3. The meaning associated with the signature?	 OpenLAB ECM 's electronic signature manifestation includes: 1. User name in addition to the full name of the signer 2. Signer's title 3. Date and time that the signature was applied 4. Location where the signing occurred 5. User configurable meaning associated with the signature
11.50 (b)	Yes	Are these items part of any human readable form of the electronic record?	The eSignature Plug-in for Adobe Acrobat places a visible signature manifestation on all human readable forms of the document, electronic display and printed form.

11.70 Signature / Record Linking

Section	Result	Question	Answer
11.70	Yes	Is the electronic signature linked to its respective electron- ic record to ensure that the sig- nature cannot be excised, copied or otherwise transferred to falsi- fy an electronic record by ordi- nary means?	The eSignature Plug-in for Adobe Acrobat encrypts the signature within the document to prevent the signature from being excised or copied to another document. Open-LAB ECM will not recognize a signature that was applied outside its own electronic signature plug-ins.

Section	Result	Question	Answer
11.100 (a)	Yes	Is each electronic signature unique to one individual and not reused by, or reassigned to, anyone else?	OpenLAB ECM uses the user name / password combination (unique to each user) in the electronic signature feature. User names within OpenLAB ECM are required to be unique and cannot be reused or reassigned to another individual.
11.100 (b)	N/A	Are the identities of the individ- uals verified prior to the estab- lishment, assignment, and cer- tification or otherwise sanction- ing an individual's electronic signature or any element of an electronic signature?	This is the responsibility of the organization, which plans, implements and operates the system. Such a verification process is a system requirement that is set before imple- menting electronic signature procedures and / or assigning electronic signature privileges to an individual.
11.100 (c)	N/A	Has the company delivered its corporate electronic signature certification letter to FDA?	It is the company's responsibility, before submitting elec- tronically signed documentation to the FDA, to register their intent to use electronic signatures. In addition, training pro- grams must be in place to ensure that users signing docu- ments electronically understand the legal significance of their electronic signature.
11.100 (c)(1)	N/A	Is it in paper form with a tradi- tional handwritten signature?	This is the responsibility of the organization, which operates the system. See $11.100(c)$.
11.100 (c)(2)	N/A	Can additional certification or testimony be provided that a specific electronic signature is the legally binding equivalent of the signer's handwritten signa- ture?	This is the responsibility of the organization, which operates the system. See $11.100(c)$.

11.100 Electronic Signatures: General Requirements

11.200 Electronic Signature Components and Controls.

Section	Result	Question	Answer
11.200 (a)(1)	Yes	Does the e-signature employ at least two distinct identification components such as user ID and password?	The OpenLAB ECM electronic signature tools consist of two components, user ID (unique) and password.
11.200 (a)(1)(i)	Yes	When an individual executes a series of signings during a sin- gle, continuous period of con- trolled system access, is the first signing executed using all the electronic signature compo- nents?	When an individual signs the first of a series of documents during a single period of controlled access the user is required to enter both signature components; user ID/ password.
11.200 (a)(1)(i)	Yes	When an individual executes a series of signings during a sin- gle, continuous period of con- trolled system access, is each subsequent signing executed	When a OpenLAB ECM user executes a series of continuous electronic signatures (defined as signatures executed within a system administrator determined period of time) they are required to enter user ID, password and reason on the first signature only. Each subsequent signature requires

Section	Result	Question	Answer
		using at least one electronic sig- nature component that is only executable by, and designed to be used by, the individual?	only the user's password, which is known only to the user.
11.200 (a)(1)(ii)	Yes	When an individual executes a series of signings not performed during a single, continuous period of controlled system access; does each signing exe- cuted require all signature com- ponents?	When a Cerity ECM user executes a series of non-contin- uous electronic signatures (defined as signatures execut- ed outside of a system administrator determined period of time) they are required to enter user ID, password and reason on each signature.
11.200 (a)(2)	Yes	Are controls in place to ensure that only their genuine owners can use the electronic signa- ture?	Cerity ECM can be configured such that an administrator can assign an initial password to a user for new account or forgotten password, but the user is required to change that password on their first login. In this manner the user ID/password combination is known only to the individual.
11.200 (a)(3)	Yes	Are the electronic signatures to be administered and executed to ensure that the attempted use of an individual's electronic signa- ture by anyone other than its genuine owner requires the col- laboration of two or more indi- viduals?	Cerity ECM uses the user's user ID and password to initi- ate the electronic signature. An Cerity ECM user's pass- word is stored encrypted within the database and is dis- played as asterisks in all location within the software. The system administrator should only know user IDs as he sets up the users. During setup he can force a pass- word change during the first logon. Then the password is only known to the individual users as it is defined at each user's individual first logon. See also 11.200(a)(2). The enforcement is the responsibility of the organization, which operates the system. Thus it requires active collab- oration with the purpose of sharing passwords to enable irregular use of somebody else' identification.
11.200 (b)	N/A	Are electronic signatures based on biometrics designed to ensure that only their genuine owners can used them?	Cerity ECM does not support signatures based on bio- metrics at this time.

11.300 Controls for Identification Codes / Passwords

Section	Result	Question	Answer
11.300 (a)	Yes	Are controls in place to ensure the uniqueness of each com- bined identification code and password maintained, such that no two individuals have the same combination of identifica- tion code and password?	The company's Windows™ logins are used to validate users, so no two users can have the same user ID/password combination.

Section	Result	Question	Answer
11.300 (b)	Yes	Are controls in place to ensure that the identification code and password issuance is periodi- cally checked, recalled, and revised?	Password renewal interval is configured as part of the Windows password policy setup. The administrator can define that passwords are automatically, periodically revised and users are prevented from reusing passwords.
11.300 (c)	N/A	Are there loss management pro- cedures in place to electronically disable lost, stolen, missing, or otherwise potentially compro- mised tokens, cards, and other devices that bear or generate identification code or password information?	Neither the ICP-MS ChemStation nor Cerity ECM support devices that bear or generate identification codes (such as tokens or cards) at this time.
11.300 (d)	Yes	Are transaction safeguards in use to prevent unauthorized use of passwords and/or identifica- tion codes?	The system can be configured such that only the user knows their user ID/password identification code. Pass- words are always displayed as asterisks and are stored encrypted within the database so that even an administra- tor cannot see them.
11.300 (d)	Yes	Are transaction safeguards in use to detect and report in an immediate and urgent manner, any attempts at their unautho- rized use to the system security unit, and, as appropriate, to organizational management?	The Windows security policy can be configured such that a user defined number of unauthorized access attempts lock out the user account and send email notification to a system administrator. The system audit trail documents general events such as logon attempts to the computer as well as the application or user changes, in the Windows Event logs as a central audit repository for all security information. This includes the system and computer ID along with the operator name and application identifica- tion, allowing for an immediate check of the potential security leak.
11.300 (e)	N/A	Are there controls in place to initially test devices that bear or generate identification code or password information to ensure that they function prop- erly and have not been altered in an unauthorized manner?	Neither the ICP-MS ChemStation nor Cerity ECM support devices that bear or generate identification codes (such as tokens or cards) at this time.

References

Improving 21 CFR Part 11 Compliance with Cerity ECM, Agilent Application Note, publication number 5989-1750EN, 2005.

Copyright © 2006-2007 Agilent Technologies, Inc. All Rights Reserved. Reproduction, adaptation or translation without prior written permission is prohibited, except as allowed under the copyright laws.

Published November 1, 2007 Publication Number 5989-4850EN



Agilent Technologies