

Improved LC Analysis of Human Serum Albumin Tryptic Digest Using 1.8-µm **ZORBAX Rapid Resolution HT SB-C18** Column

Application

Proteomics

William E. Barber

Introduction

The liquid chromatography (LC) analysis of complex biological fluids is significantly improved using fine particle columns. This is demonstrated by a comparison of two columns, differing in length and particle size, on the analysis of a tryptic digest of human serum albumin (HSA). The chromatographic comparison is shown in Figure 1.





Conditions

Instrument	Agilent 1100LC
Column	1A: ZORBAX SB-C18 (2.1 × 150 mm, 5 μm) p/n 883700-922 1B and 1C: ZORBAX Rapid Resolution HT SB-C18 (2.1 × 50 mm, 1.8 μm) p/n 822700-902 SB-C18 was 80 Å for all columns
Temperature	50 °C
Sample concentration Injection	15 pmol/μL 8 μL (120 pmol total on-column)

The Agilent 1100 LC was optimized to minimize peak broadening by using 0.12-mm id tubing from the autosampler to the diode array detector (DAD) and a 1.7-µL flowcell. Peaks were measured using the 214-nm absorption band.

Highlights

- High resolution and rapid analysis using Agilent Rapid Resolution High Throughput (RRHT) columns is demonstrated for a sample of HSA tryptic digest
- ZORBAX SB-C18, 2.1 × 50 mm, 1.8-µm columns provide better resolution than, and comparable sensitivity to 150-mm columns, but in one-sixth the time.

Mobile phase		
A:	0.1% TFA in H ₂ O	
B:	0.1% TFA in ACN	

Gradient timetable				
Flow		Gradient		
Run	mL/min	2% to 50% B		
1A	0.2	70 min		
1B	0.5	10 min		
1C	0.5	30 min		



Run	X/Y Peaks*	Window, min	Peak widths
1A (5 µm)	90/120	60	12–18 s
1B (1.8 µm)	102/125	10	3–6 s
1C (1.8 µm)	134/156	25	6–12 s

* There were X total peaks for the window shown out of Y peaks for the entire chromatogram

Results

Figure 1A shows the results obtained for a column $(2.1 \times 150 \text{ mm}, 5 \mu\text{m})$ and gradient conditions in common use for the gradient HPLC analysis of HSA tryptic digests. Using a shorter 50-mm column of the same phase, but with smaller $(1.8 \mu\text{m})$ particles, higher resolution is obtained in one-sixth the time, or 500% faster, (Figure 1B) due to the combined effects of using a much smaller particle, reducing the column length by a third, increasing the flow rate by a factor of 2.5 and reducing the gradient time from 70 min to 10 min. The resolution can be increased further (Figure 1C) by increasing the gradient time (shallower gradient) from 10 min to 30 min while the analysis time is still less than half that of the original conditions.

Conclusion

LC separation of HSA tryptic digest components is efficiently performed using short RRHT columns with 1.8-µm particles.

Run	Column type	Dimensions	Part number
1A	ZORBAX SB-C18	2.1×150 mm, 5 μm	883700-922
1B	ZORBAX SB-C18	2.1×50 mm, 1.8 μm	822700-902
1C	ZORBAX SB-C18	2.1×50 mm, 1.8 μm	822700-902



For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

The author, William E. Barber, is a LC Applications Specialist based at Agilent Technologies, Wilmington, Delaware.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2004

Printed in the USA January 21, 2004 5989-0542EN

