


Investigation of the Unique Selectivity and Stability of Agilent GS-OxyPLOT Columns



Application

Gas Chromatography

Authors

Yun Zou and Min Cai
Agilent Technologies (Shanghai) Co. Ltd.
412 Ying Lun Road
Waigaoqiao Free Trade Zone
Shanghai 200131
P.R. China

Abstract

The stationary phase of a GS-OxyPLOT column is a proprietary, salt deactivated adsorbent. GS-OxyPLOT columns show unique selectivity to oxygenated hydrocarbons, excellent stability and reproducibility, long column lifetime, and a wide application range.

Introduction

The determination of oxygenated hydrocarbons in different sample matrices is very important for the petrochemical industry, because oxygenates directly influence product quality. Presence of such oxygenates may cause the catalysts to be poisoned and deactivated, resulting in more downtime and higher costs. ASTM has developed several methods for analysis of oxygenates, such as ASTM D7059, D4815, and D5599. The oxygenates include ethers, esters, ketones, alcohols, and aldehydes.

Methanol is one of the oxygenates that often present in light hydrocarbon streams. For example, it is added to natural gas and production of crude oil to prevent hydration of hydrocarbons during transportation via pipelines. Therefore, it is important

to accurately measure the content of methanol from light hydrocarbons at different concentrations, including at trace levels.

To achieve this, a new porous layer open tubular (PLOT) capillary column, the GS-OxyPLOT column, was used. The stationary phase of the GS-OxyPLOT is a proprietary, salt deactivated adsorbent with a high chromatographic selectivity for low molecular weight oxygenated hydrocarbons, while having virtually no interactions with saturated hydrocarbon solutes [1].

Using Capillary Flow Technology, such as back-flush or Deans switch, GS-OxyPLOT columns can provide a turnkey solution for the analysis of trace level oxygenate impurities in complex matrices, such as motor fuels, crude oil, and gaseous hydrocarbon [2]. Meanwhile, a GS-OxyPLOT column can be used as a single analytical column to separate oxygenates for some samples. In this application, methanol was set as an example to investigate the performance of the GS-OxyPLOT column.

Experimental

The experiments were performed on an Agilent 7890A GC system and a 6890N GC system equipped with split/splitless capillary inlet, flame ionization detector (FID), and Agilent 7683 Automatic Liquid Sampler (ALS). The split/splitless inlets were fitted with long-lifetime septa (Agilent p/n 5183-4761) and split/splitless injection liners (Agilent p/n 5183-4711). Injections were done using 10- μ L syringes (Agilent p/n 9301-0714). A glass indicating moisture trap (Agilent p/n LGMT-2-HP), an oxygen trap (Agilent p/n BOT-2), and a



Agilent Technologies

hydrocarbon trap (Agilent p/n 5060-9096) were installed. Agilent ChemStation was used for all instrument control, data acquisition, and data analysis.

Results and Discussion

Analysis of Normal Hydrocarbons and Methanol

A mixture of normal hydrocarbons and methanol was prepared with the following approximate concentrations %(w/w): 34.8% n-pentane, 12.8% n-hexane, 1.8% n-heptane, 1.9% n-octane, 2.1% n-nonane, 3.9% n-decane, 2.1% n-undecane, 9.8% n-dodecane, 11.8% n-tridecane, 4.7% n-tetradecane, 2.4% n-pentadecane, 4.5% n-hexadecane, 2.4% n-heptadecane, 1.0% n-octadecane, 0.9% n-eicosane, 0.9% n-docosane, 1.1% n-tetracosane, and 0.8% methanol.

The analytical conditions are summarized in Table 1. The normal hydrocarbons and methanol analysis was performed on a GS-OxyPLOT column (Agilent p/n 115-4912). The GC chromatogram is shown in Figure 1.

Table 1. Conditions for Normal Hydrocarbons and Methanol Analysis

Column	GS-OxyPLOT, 10 m × 0.53 mm × 10 µm (Agilent p/n 115-4912)
Carrier gas	Helium, constant flow mode, 40 cm/s @ 50 °C
Inlet	Split/splitless at 325 °C
Split ratio	80:1
Oven temperature	50 °C (2 min); 10 °C/min to 290 °C (2 min)
Post-run	300 °C (2 min)
Detector	FID at 325 °C
Injection size	0.2 µL

In Figure 1, the GS-OxyPLOT column shows unique retention characteristics for methanol. The lower boiling point hydrocarbons were not strongly retained on the stationary phase and eluted through the FID very rapidly. The methanol eluted after n-C14, allowing it to be quantified without any interference from the hydrocarbon matrix, and making it feasible for trace-level methanol analysis in a range of hydrocarbon streams.




Figure1. Analysis of methanol and normal hydrocarbons on a GS-OxyPLOT column, 10 m × 0.53 mm × 10 µm.

In addition, the baseline was quite smooth, even when the oven temperature was up to 290 °C. GS-OxyPLOT has an upper temperature limit of 350 °C and exhibits virtually no bleed, making it widely applicable for long-term reliable analysis.

Analysis of Alcohols

A mixture containing a range of primary alcohols from methanol to lauryl alcohol was analyzed on a GS-OxyPLOT column using a temperature-programmed method. Table 2 lists conditions for alcohols separation, and the resulting chromatogram is shown in Figure 2.

Sample

The sample had an approximate concentration (v/v) of 1% methanol, ethanol, propanol, butanol, amyl-alcohol, heptanol, octanol, nonanol, decyl alcohol, and lauryl alcohol in toluene.

As can be seen in Figure 2, all of the alcohols are separated and eluted with good peak shape within

Table 2. Conditions for Alcohols Analysis

Column	GS-OxyPLOT, 10 m × 0.53 mm × 10 µm
Carrier Gas	Helium, constant flow mode, 40 cm/s at 150 °C
Inlet	Split/splitless at 325 °C
Split ratio	50:1
Oven temperature	150 °C (0 min); 10 °C/min to 300 °C (5 min)
Detector	FID at 325 °C
Injection size	0.2 µL

an analysis time of 15 min. In this experiment, oven temperature was set up to 300 °C. Thanks to its advanced dynamic coating process, Agilent's GS-OxyPLOT stationary phase exhibits virtually no detector spiking due to particle generation from the phase coating [3].

Due to the high viscosity of alcohols, especially decyl alcohol and lauryl alcohol, it is necessary to wash the needle after each injection in case of carryover problems.




Figure 2. Separation of alcohols using GS-OxyPLOT, 10 m × 0.53 mm × 10 µm.

Influence of Temperature on the Selectivity of GS-OxyPLOT

To polar stationary phases, the temperature has a direct influence on the selectivity. GS-OxyPLOT offers extremely high polarity. The analysis of normal hydrocarbons and methanol demonstrated that methanol elutes after n-C14. Using a mixture containing methanol, n-tetradecane, and n-pentadecane, isothermal Kovats retention indices were tested at isothermal oven temperatures of 150, 200, 220 and 250 °C, respectively (Table 3). The relationship between Kovats retention indices and oven temperature is shown in Table 4.

Table 3. Conditions for Kovats Retention Indices Test

Column	GS-OxyPLOT, 10 m × 0.53 mm × 10 µm
Carrier gas	Helium, constant flow mode, 30 cm/s at 150 °C
Inlet	Split/splitless at 250 °C 100:1 split ratio
Oven temperature	150, 200, 220, and 250 °C, respectively; isothermal
Detector	FID at 250 °C
Injection size	0.2 µL

Table 4. Kovats Retention Indices and Oven Temperature (n > 3)

Oven temp.	150 °C	200 °C	220 °C	250 °C
LOT1	1419	1418	1418	1413
LOT2	1420	1421	1419	1417

Retention index, I_x , was calculated using the following equation:

$$I_x = 100n + 100[\log(t_x) - \log(t_n)] / [\log(t_{n+1}) - \log(t_n)]$$

Where t_n and t_{n+1} are retention times of the reference n-alkane hydrocarbons eluting immediately before and after chemical compound X; t_x is the retention time of compound X. Here compound X is methanol, the reference n-alkane hydrocarbons are n-tetradecane and n-pentadecane, respectively.

Table 4 shows good repeatability of Kovats retention indices for two different lots of GS-OxyPLOT columns. The retention index for methanol only changed by less than 10 index units over 100 °C temperature difference. Therefore, when the oven temperature changes from 150 to 250 °C, it has little influence on the selectivity of GS-OxyPLOT.

Influence of Moisture on GS-OxyPLOT

Some PLOT columns can adsorb water, which can lead to changes in retention times and selectivity

for analytes. Therefore, column performance will be influenced greatly in the presence of water. Although cumbersome solvent-extraction procedures can be performed before injection, injecting sample that contains water is, in some cases, unavoidable.

From a GC point of view, water is a less-than-ideal solvent. The problems associated with water include large vapor expansion volume, poor wet ability and solubility in many stationary phases, detector problems, and perceived chemical damage to the stationary phase. In order to test the effect of water, a GS-OxyPLOT column that had gone through about 1,500 runs was tested before and after injecting 100% aqueous samples.

Water has a large vapor expansion volume; the vapor volume of water (assuming a 1-µL injection) can easily exceed the physical volume of the injection liner (typically 200 to 900 µL). The volume for the liner used in this experiment (Agilent p/n 5183-4711) is 870 µL, so the injection volume was set as 0.2 µL. Table 5 lists the conditions for the moisture testing, and the resulting chromatograms are shown in Figure 3.

Table 5. Conditions for Moisture Test

Column	GS-OxyPLOT, 10 m × 0.53 mm × 10 µm
Carrier gas	Helium, constant flow mode, 38 cm/s at 150 °C
Inlet	Split/splitless at 300 °C 15:1 split ratio
Oven temperature	150 °C isothermal, post-run: 300 °C (5 min)
Detector	FID at 300 °C, H2:45mL/min, air: 400 mL/min, makeup: 30 mL/min
Injection size	0.2 µL
Sample	0.1% n-Dodecane, Methyl tert-butyl ether, n-Tridecane, Iso-Butyraldehyde, n-Tetradecane, Methanol, Acetone, and n-Pentadecane

As shown in Figure 3, the area of n-pentadecane remained the same before and after 100 injections of water. However, compared with the area before injecting water, the area of methanol (peak 6) decreased by 50%, and the area of acetone (peak 7) decreased by 14.4% after 100 injections of water (see Table 6). It demonstrated that water can affect the activity of GS-OxyPLOT, especially for the analysis of those relatively low molecular weight oxygenated compounds, such as methanol and acetone.




Figure 3. Comparison of test mixture separation before (A) and after (B) 100 injections of water.

As for retention times and column efficiency, they are not strongly influenced. After 100 injections of water, the retention time of C15 changed from 9.689 min to 9.384 min, and the column efficiency of C15 changed from 14,792 to 14,781.

Condition the column at 300 °C for two hours, followed by 12 hours at 250 °C. As shown in Figure 4 and Table 6, it is obvious that GS-OxyPLOT phase can be regenerated by conditioning.




Figure 4. Expanded view shows comparison of test mixture separation on GS-OxyPLOT.
4A. Before injection of water. 4B. After 100 injections of water. 4C. After conditioning the column.


Table 6. Comparison of Test Mixture Separation

	Methanol			Acetone			n-Pentadecane		
	Before injection of water	After 100 injections of water	After conditioning column	Before injection of water	After 100 injections of water	After conditioning column	Before injection of water	After 100 injections of water	After conditioning column
RT (min)	6.022	5.835	5.915	6.429	6.160	6.305	9.689	9.384	9.658
Area	20.23	9.18	20.88	94.53	80.92	98.07	277.79	287.7	287.9
Plates	11887	12920	11616	9532	10357	9573	14792	14781	15100

After conditioning the GS-OxyPLOT column, the peak area and retention time reproducibility were determined. Figure 5 and Table 7 show excellent RT precision, lower than 0.6% over five test mixture runs on this GS-OxyPLOT column. The peak area has a relative standard deviation (RSD%) below 2.5%. It proved that column performance can be restored via conditioning.

Determination of Methanol

The following analysis of methanol followed ASTM D7059 [4]: “Standard Test Method for Determination of Methanol in Crude Oils by Multidimensional Gas Chromatography.” Methanol was determined by gas chromatography with FID using internal standard method with GS-OxyPLOT column.

**Figure 5. Fifth run overlaid using GS-OxyPLOT (after conditioning column).****Table 7. Peak Area Reproducibility and Retention Time Reproducibility on GS-OxyPLOT (after conditioning column)**

Compound (by eluted order)	Dodecane	MTBE	Tridecane	Iso- Butyraldehyde	Tetradecane	MeOH	Acetone	n-C15
Area RSD% (N = 5)	1.18	1.58	1.59	2.49	1.15	2.12	1.98	1.82
RT RSD% (N = 5)	0.18	0.12	0.26	0.55	0.29	0.16	0.19	0.33

Reagents and Materials

Carrier gas, Helium, > 99.95% purity

Methanol, > 99.9% purity

1-propanol, > 99.9% purity, and containing
< 500 ppm methanol

Toluene, > 99.9% purity, and containing < 0.5 ppm
methanol

A set of calibration standards 5, 25, 125, 250, 500,
1,000 and 1,500 ppm (m/m) of methanol, and each
containing 500 ppm (m/m) of 1-propanol internal
standard, were prepared in toluene.

The calibration standard solutions should be
stored in tightly sealed bottles in a dark place
below 5 °C.

Linearity

Under the conditions listed in Table 8, the
methanol calibration standards were analyzed. The
linearity is shown by plotting the response ratio of
methanol and internal standard 1-propanol against

their amount ratio (see Figure 6). For methanol,
good linearity was gained ranging from 5 to
1,500 ppm. The correlation r^2 value for the calibra-
tion curve is higher than 0.999.

Figure 7 and Figure 8 are chromatograms of
methanol at a level of 5 ppm and 1500 ppm,
respectively. At a relatively high concentration of
1500 ppm, methanol still could get a sharp peak.
The limit of quantification (LOQ) was calculated to
be 1 ppm using the chromatogram of 5 ppm
methanol.

Table 8. System Settings for the Calibration Curve

Column	GS-OxyPLOT, 10 m × 0.53 mm × 10 µm
Carrier gas	Helium, constant flow mode, 50 cm/s at 150 °C
Inlet	Split/splitless at 250 °C 10:1 split ratio
Oven temperature	150 °C (3 min); 20/min to 300 °C (5 min)
Detector	FID at 325 °C
Injection size	1 µL




Figure 6. The calibration curve of methanol in toluene.




Figure 7. Test mixture of 5 ppm methanol in toluene.




Figure 8. Test mixture of 1,500 ppm methanol in toluene.

Repeatability

The reproducibility of the GS-OxyPLOT is given in Table 9. Those values were obtained by the replicate analysis of different methanol levels (25, 125, and 1,500 ppm) in different days. The injection was done by ALS with RSD no less than 3% either intraday or interday analysis, which was very low for this type of determination.

Life Span

Under the conditions in Table 5, a mixture was analyzed with a GS-OxyPLOT column which went through 1,500 injections of methanol. It shows that the column has a long lifetime. The GS-OxyPLOT column still has good resolution for each compound and high efficiency of 1,482 plates per meter for n-pentadecane (see Figure 9).

Table 9. Relative Standard Deviations Intraday and Interday at Different Levels (25, 125, and 1,500 ppm) of Methanol

Day	25 ppm (average)	RSD (%)	125 ppm (average)	RSD (%)	1,500 ppm (average)	RSD (%)
D 1	25.2	0.46	123.9	0.45	1507.3	0.55
D 2	25.3	1.53	123.2	0.79	1494.4	0.45
D 3	24.4	0.36	125.4	1.71	1523.5	0.35
D 4	25.9	1.06	123.0	0.90	1537.8	0.51
D 5	23.9	0.44	121.1	0.76	1502.4	1.03
Stand. dev.	0.7		1.70		17.4	
Average	24.97		123.6		1513.1	
RSD (%)	2.8		1.37		1.15	




Figure 9. Chromatogram of performance mixture after 1,500 injections.

Conclusions

GS-OxyPLOT provides good retention and selectivity for oxygenated compounds. Normal alkanes up to C24 and primary alcohols up to lauryl alcohol can elute from GS-OxyPLOT within its program temperature maximum limit of 350 °C. Methanol elutes after n-C14 with retention index higher than 1,400; the retention index is quite stable from 150 to 250 °C, allowing methanol to be measured at low levels in a wide range of hydrocarbon streams.

Methanol has to be measured usually at specs as low as 5 ppm. From 5 to 1,500 ppm, it shows good linearity on GS-OxyPLOT. And the column has proven extremely stable with long lifetime.

GS-OxyPLOT can tolerate a little amount of water in samples, and column performance can be restored via conditioning.

GS-OxyPLOT can be used for a single-column system or in multidimensional GC systems. It offers a unique solution for the analysis of oxygenates in the chemical and petrochemical industries.

References

1. A. K. Vickers, "GS-OxyPLOT: A PLOT Column for the GC Analysis of Oxygenated Hydrocarbons," Agilent Technologies publication 5989-6447EN, March 2007
2. "Agilent J&W GS-OxyPLOT Capillary GC Columns," Agilent Technologies publication 5989-6489EN
3. A. K. Vickers, "A 'Solid' Alternative for Analyzing Oxygenated Hydrocarbons — Agilent's New Capillary GC PLOT Column," Agilent Technologies publication 5989-6323EN, February 2006
4. Standard test method for the determination of methanol in crude oils by multidimensional gas chromatography, ASTM D7059-04, July 2004

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2008

Printed in the USA
June 17, 2008
5989-8771EN



Agilent Technologies