




Agilent Lab Advisor

User Manual



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2011

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

M8550-90001

Edition

11/2011

Printed in Germany

Agilent Technologies
Hewlett-Packard-Strasse 8
76337 Waldbronn

This product may be used as a component of an in vitro diagnostic system if the system is registered with the appropriate authorities and complies with the relevant regulations. Otherwise, it is intended only for general laboratory use.

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will

receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

1	Lab Advisor B.02.02 Overview	5
	Lab Advisor Basic	6
	Lab Advisor Advanced	7
2	Installation	9
	Prerequisites	10
	Deployment Modes	12
	Hard drive installation	13
	Server Installation	18
	Running Lab Advisor from the CD-ROM	19
	Updating Lab Advisor	20
	Migrating from Lab Advisor B.01.xx	21
3	Using Lab Advisor	23
	Navigation	24
	System Overview	28
	Configuration	34
	Apps	39
	Firmware Update	43
	Logs and Results	44
	Service and Diagnostics	47
	Instrument Control	49
	EMFs	54
	Status Report	56



1

Lab Advisor B.02.02 Overview

Lab Advisor Basic 6

Lab Advisor Advanced 7

This chapter describes the two versions of Lab Advisor B.02.02: Lab Advisor Basic and Lab Advisor Advanced.

With its advanced diagnostic and maintenance capabilities, Agilent Lab Advisor helps you to keep your Agilent LC and CE instruments in top condition and thereby achieve high quality chromatographic results. Agilent Lab Advisor is an application-independent tool: it can support Agilent LC and CE instrumentation regardless of whether you are using Agilent or non-Agilent software to control the instruments. The Agilent Lab Advisor is available in two flavors: Lab Advisor Basic and Lab Advisor Advanced.



Lab Advisor Basic

Lab Advisor Basic provides state of the art tests, tools and calibrations to support you in the daily operation, maintenance and basic troubleshooting tasks. The Basic version comes with a full set of diagnostic and maintenance capabilities, allowing you to perform troubleshooting and maintenance tasks efficiently and with little effort.

Lab Advisor Advanced

Lab Advisor Advanced has been designed to support users who need the highest quality data and the utmost reliability from the Agilent LC and CE instrumentation. This is provided with additional functionality that includes tools and features that allow you to carry out sophisticated diagnostic, usage-based maintenance, and generate traceable results, including features such as user log-on with password, traceable result files, advanced EMF functionality and data sharing facilities.

1 Lab Advisor B.02.02 Overview

Lab Advisor Advanced



2 Installation

Prerequisites	10
Deployment Modes	12
Hard drive installation	13
Server Installation	18
Running Lab Advisor from the CD-ROM	19
Updating Lab Advisor	20
Migrating from Lab Advisor B.01.xx	21

This chapter contains instructions for installing Lab Advisor.B.02.01.

Under construction



Prerequisites

Agilent Lab Advisor can run on any Microsoft Windows XP or Windows 7 PC with Microsoft .NET 2.0 or higher installed. The software has been extensively tested with the following software packages:

Table 1 Supported Software List

Element	Revision
Windows XP Pro	SP3
Windows 7	32 bit
Windows 7	64 bit
Windows Server 2008 R2	
Internet Explorer	7.0
Internet Explorer	8.0
.Net Framework	2.0 and higher
Symantec Antivirus	11.0
Adobe Acrobat	9.0

For optimum performance of your Agilent Lab Advisor software, the following minimum requirements should be fulfilled. The minimum supported configuration is based on the installation on a Netbook, but for larger installations, or for higher performance, the recommended configuration should be used.

Table 2 PC Hardware Configurations

	Minimum (Netbook) Configuration	Recommended Configuration
Processor	Intel Atom processor	Pentium D or higher, Intel Dual-Core 3.4GHz or higher
RAM	1.0 GB or more	≥ 1 GB (Windows XP SP3 Pro) ≥ 2 GB (Windows 7)

Table 2 PC Hardware Configurations

	Minimum (Netbook) Configuration	Recommended Configuration
Hard disc	1 GB	160 GB or bigger
Video	1024 × 600 resolution	1280 × 1024 resolution
Removable media	(external) CD-Rom drive	DVD drive
Mouse	MS Windows compatible pointing device	MS Windows compatible pointing device
LAN	10/100baseT	10/100baseT
Operating Systems	Windows 7 Starter	Windows XP Pro Service Pack 3 Windows 7 Windows Server 2008 R2 for server installations
Printer	all printers supported by the operating system in use	all printers supported by the operating system in use

Deployment Modes


Lab Advisor can be deployed in four different modes:

Table 3 Lab Advisor Deployment Modes


Mode	Installation	Start	Data location
Standard	Installed on the local hard drive using setup.exe, see “ Hard drive installation ” on page 13.	Run from desktop icon or Windows Start menu.	Standard: ..\Agilent Lab Advisor\AgilentLabAdvisorData\...
Server	Installed on Windows Server 2008 using setup.exe and published as a shareable application, see “ Server Installation ” on page 18.	Run using RDP link at clients' desktops.	..\Agilent Lab Advisor\AgilentLabAdvisorData-[TerminalHostName]\...
Instant	No installation required.	Run AgilnstrDiag.exe from the installation CD, see “ Running Lab Advisor from the CD-ROM ” on page 19	Standard: ..\Agilent Lab Advisor\AgilentLabAdvisorData\...

Hard drive installation

- 1 Double-click the **Setup.exe** icon to start the installation.




The **Choose Language** dialog box is displayed.



- 2 Click the down-arrow and select the language from the list and click **OK**. Lab Advisor supports English, Chinese and Japanese.

The first page of the installation wizard is displayed.



NOTE


You can change the language after installation, see [“Language”](#) on page 34

2 Installation

Hard drive installation


3 Click **Install**.

The License Agreement is displayed.




4 Select to accept the terms of the license and click **Next**.

The **Custom Setup** page of the installation wizard is displayed.



5 If you want to change the default location, click **Browse** to select a new location; otherwise, click **Next**.


The last page of the installation wizard is displayed.



- 6 Mark the check boxes against the options you want to activate and click **Next**.

Restore old data restores old data from either B.01.xx installations or from previous B.02.xx installations. The data is not removed from the hard drive during de-installation, but if **Restore old Data** is cleared, the old data is permanently removed from the hard drive.

If the **Launch the Agilent Installation Qualification Tool** check box was marked, the Agilent IQT tool launches automatically. To start the Installation Qualification click **Qualify**.



The Installation Qualification Report is opened in a browser window and can be printed. To access the report at a later time, go to \Program Files\Agilent Technologies\Lab Advisor\IQTTool\IQProducts\Agilent Lab Advisor LC & CE\Reports.

2 Installation

Hard drive installation

Installation Qualification Report				
Date:	31, May 2011	Time:	10:05:26 [UTC +02:00]	Host Name: CZC9030ZHW
Windows User Name:	je1208	Base Revision Number:	02.01.013.0	Product Name: Agilent Lab Advisor LC/CE
Install Type:	N/A			

Base Reference File Name : [igtfile32.xml](#)

Summary

Overall Evaluation of Installation Check: : PASS

File Report Summary

- No missing files or invalid files found
- No system file differences found

Files Registration Report Summary


- No Registerable Files found for Qualification

Registry Report Summary

- No registry entries found for Qualification.

To rerun the Installation Qualification at a later date, navigate to Start\All programs\Agilent Lab Advisor\Installation Qualification and start the program.

When you close the IQ tool, the installation is complete.



If the **Enter License** check box was marked in the installation wizard, the software starts and automatically navigates to the Configuration\License screen, where you can enter your the license keys.

The screenshot displays the 'Configuration\License' window. On the left, a vertical sidebar contains four icons with labels: 'Software' (gear icon), 'Licenses' (key icon), 'Traceability' (person icon), and 'About' (globe icon). The main content area is divided into a table and a form section. The table has four columns: 'Key', 'Product', 'Level', and 'Modules'. The table body is currently empty. Below the table, there is a 'License Key:' label, a text input field, and an 'Add' button. A 'Remove' link is located in the bottom right corner of the main area.

Key	Product	Level	Modules
-----	---------	-------	---------

License Key:

[Remove](#)

Server Installation

You can install Lab Advisor on a Windows Server 2008 server and make it available for use on client systems.

- 1** Install Lab Advisor on the Windows Server 2008 server as described in [“Hard drive installation”](#) on page 13.
- 2** Install the role **Remote Desktop Services** on the Windows Server 2008 server.
- 3** Publish the Lab Advisor installation on the Windows Server 2008 server.
- 4** Create an RDP (Remote Desktop Protocol) file for the Lab Advisor installation.
- 5** Distribute the RDP file to all client systems.

Running Lab Advisor from the CD-ROM

When you run Lab Advisor from the CD-ROM, no program files are copied to the local hard drive. However, any data that you generate will be saved in a folder ..\Agilent Lab Advisor\ AgilentLabAdvisorData\ on the local hard drive.

- 1 Insert the Lab Advisor installation CD-ROM into the drive.
- 2 Double-click the **Instant Diagnostic** item in the root directory of the CD-ROM.

Updating Lab Advisor


A fast update process has been developed to update the Agilent Lab Advisor software to a new revision. The fast update process uses .DLA files, and is launched by double-clicking the file itself. An update using the fast update process requires files specific to the update required. The file name carries the name of the origin and the destination revisions, and must be selected accordingly. For example, to update from Lab Advisor B.02.01 [012] to B.02.02 [002], double-click the following .DLA file:

```
LA_B.02.01[028]_B.02.02[042].dla
```

The fast update using .DLA files cannot be used to update portable installations of the software. To upgrade these installations, the upgraded software of a hard drive installation needs to be copied to the portable media. It is highly recommended that you rename the old application folder before copying the updated software folder. After finishing the copy of the application folder, the data folder can be copied from the “old” application folder.

Migrating from Lab Advisor B.01.xx


It is possible to migrate data and licenses from B.01.xx installations to B.02.xx. To do this, the B.01.xx has to be uninstalled, since the uninstallation process creates a data repository that can be used by the B.02.xx installer. To access the old data, mark the **Restore old Data** check box during the installation process.



If B.01.xx is uninstalled after the installation of the B.02.xx, then the restore can be activated by navigating to the application folder, typically `..\Program Files\Agilent Technologies\Lab Advisor` and double-clicking the `Agilent.InstrDiag.DataMigration.exe` program.

2 Installation

Migrating from Lab Advisor B.01.xx



3 Using Lab Advisor


Navigation	24
System Overview	28
Adding a new system	28
Changing system properties	29
Removing a system	30
Adding a new system group	31
Deleting a system group	31
Moving systems between groups	31
Copying Device Details to the Clipboard	32
Fast Connect	32
Configuration	34
Configuring the Software	34
Licenses	35
Traceability	37
Apps	39
Diagnostic Catalog	39
Data Sharing	40
Firmware Update	43
Logs and Results	44
Service and Diagnostics	47
Instrument Control	49
EMFs	54
Status Report	56

This chapter describes the Lab Advisor user interface and provides details about the available features.



Navigation

The Lab Advisor User Interface is divided into six main areas. The content of these areas changes depending on the screen selected in the Agilent Lab Advisor software.




Title Bar



The Title Bar shows which of the configured systems in the System Overview is currently selected. It also hosts the buttons for maximizing, minimizing and closing the Agilent Lab Advisor application.

Navigation Panel




The primary navigation between the different screens of the Lab Advisor software is done in the Navigation Panel. The Navigation Panel consists of three areas:

- The Global Tasks lists system-independent screens that access information or configurations regardless of the configured systems and their current state.
- The System Tasks are System-specific and change with the selected system. The name of the selected system is displayed in the header of the System Tasks and in the Title Bar.

- At the bottom of the Navigation Panel, the Help topics provide information about the software and about the individual screens. Context-sensitive help topics can also be accessed at any time by clicking **F1**.

The Navigation Panel can be minimized by clicking the minimize icon in the top Task bar.

Information Panel




The Information Panel contains information about the currently selected screen and the Agilent Lab Advisor software version. If the Traceability feature is in use, the Information Panel also includes information about the current logged-in user.

Application Panel

This is the primary area where the different screens selected in the Navigation Panel of the Agilent Lab Advisor software are displayed. Refer to the individual applications for more detailed information.

Group Controls



If system grouping has been activated in software configuration, these controls allow you to switch between system groups, and to add a new system group or delete an existing group.

Action Panel



Additional buttons or actions that are applicable to the selected screen may be displayed in the Action Panel. Such controls are applicable to the entire selected screen; buttons or actions that are applicable to individual items within the screen are displayed in the screen itself.


Status Bar



The left side of the Status Bar contains information about the connection used; details about the Agilent Lab Advisor revision, license level and license usage are shown on the right side.


System Overview

The System Overview screen gives a fast overview of the state of all configured and connected systems. The System Overview is also the main selection screen for the System Tasks.



Adding a new system

- 1 In the Action Panel of the **System Overview**, click **Add System**.
The **Add System** dialog box is displayed.



- 2 In the **Add System** dialog box, enter a name for the system in the **System Name** field and the connection details in the **Connection Address** field.

NOTE

The connection details can be an IP address, the host name or, if you are connecting using a serial cable, the COM port.

- 3 Click the **Type** down-arrow and select the type of system you are adding from the list.
- 4 Click **OK** to finish adding the system and close the **Add System** dialog box.
The system becomes visible in the **System Overview**, and Lab Advisor tries to connect to it.

Changing system properties


You can change the name or connection address of an existing system, add additional information or activate the automatic **Reconnect** feature for the system.

- 1 Click on the system in the **System Overview** screen to select it.
- 2 Click **System Properties** in the Action Panel.

OR

Right-click on the system and select **Properties** from the context menu.

The **System Properties** dialog box is displayed.



Note that the **System Group** field is present only when the **Activate Grouping** check box in the **Configuration - Software** screen is marked.

- 3 Add or modify any of the parameters in the **System Properties** dialog box.

NOTE

When you mark the Reconnect check box, Agilent Lab Advisor automatically tries to connect to the system whenever the application is launched. This feature can be activated for all systems configured in the **System Overview** simultaneously.

- 4 Click **Apply** to register the changes and close the **System Properties** dialog box.

Removing a system

- 1 Click on the system in the **System Overview** screen to select it.
- 2 Click **Remove System** in the Action Panel.

The system is removed from the **System Overview**.

NOTE


The data collected for the system will still be available in the **Logs and Results** application, but will be listed as unassigned systems.

Systems that have been removed from the **System Overview** are still counted toward the module limit of the installed license. To permanently delete a module, see [“Permanently deleting a hardware module”](#) on page 36.

Adding a new system group

NOTE

The system group controls are available only when the **Activate Grouping** check box in the **Configuration - Software** screen is marked.

- 1 Click  at the right of the group controls.
A new system group tab is added with a default name.
- 2 Right-click on the new tab, select **Rename** from the context menu and overwrite the default name with a new name.
OR
Double-click on the new tab and overwrite the default name with a new name.
- 3 Click **Add System** to display the **Add System** dialog box, which allows you add a system into the new system group.

Deleting a system group

NOTE

You cannot delete a system group that contains systems. Before deleting a group, move the systems into another group (see [“Moving systems between groups”](#) on page 31).

- 1 Right-click on the tab of the system group that you want to delete.
- 2 Select **Delete** from the context menu.
If the system group is empty, it is deleted; if the group contains systems, a message is displayed.

Moving systems between groups

- 1 In the System Overview, select the system that you want to move and click **System Properties**.
OR

In the System Overview, right-click the system that you want to move and select **Properties** from the context menu.


The **System Properties** dialog box is displayed.

- 2** Click the **System Group** down-arrow and select the target group that you want to move the system to.
- 3** Click **Apply**.

The system is moved from the existing group to the new group and the **System Properties** dialog box is closed.

Copying Device Details to the Clipboard

Sometimes, it can be helpful to have the details presented in the **System Information** section of a device available for copying and pasting into other applications.

- 1** Click on the system in the **System Overview** screen to select it.
- 2** If the system modules are not listed, click  to display them.
- 3** Right-click on the module whose details you want to copy and select **Copy details to Clipboard** from the context menu.

The device information is copied to the clipboard, and can be pasted into another application such as Notepad, Wordpad or a Microsoft Office application.

Fast Connect

If you are using mobile laptop computers for servicing systems, a fast connection can be established using a serial cable (RS232).

- 1** Connect the serial cable between the system and the PC (a USB-to-Serial adapter, p/n 8121-1013, might be necessary). The serial cable should be connected to the module that is providing the most data, usually the detector.
- 2** Click **Connect via Serial Cable** on the Action Panel.

Lab Advisor searches all available COM ports for installed systems and adds them automatically to the **System Overview** screen.

Lab Advisor can also be setup to do this automatically on startup by marking the **On startup** check box in the Action Panel.

NOTE

This feature is especially helpful for connecting systems with no LAN access, because it provides easy access to data such as LAN card configuration, MAC address and IP address without having to reconfigure the Laptop internal IP address or set up a BootP server.

Configuration

Application-wide settings, information and tools are accessible from the **Configuration** screen.

Configuring the Software

The **Software** configuration specifies the **Path** to the data generated by the Lab Advisor software. This is a non-configurable path, and is dependent on the operating system used and the type of installation (USB stick or hard drive).

Import/Export Data

To distribute configured systems and their corresponding data, it is possible to export data from one Lab Advisor installation and import it into other installations. This feature can also be used as a backup solution, where the exported .ZIP file is stored in a safe location.

Tracing

If unexpected behavior is observed from the Lab Advisor software, a trace file can be generated to help the Agilent Technologies technicians locate the problem.

Language

The Lab Advisor software supports English, Chinese and Japanese languages. The language is usually selected during installation of the software. However, there is a possibility to change the language later by selecting the appropriate language in the **Software** configuration screen. After the language has been changed, the software needs to be restarted for the new settings to take effect.

System Groups

Lab Advisor B.02.02 supports the grouping of systems, for example, by laboratory. Up to 25 systems can be assembled into a group; each system can contain up to 50 devices. An unlimited number of groups can be defined.

The grouping of systems is switched on by marking the **Activate Grouping** check box. When the check box is marked, the group controls (see “[Group Controls](#)” on page 26) are appended to the application panel in the System Overview and the Review Client.

Licenses

The licensing scheme of Lab Advisor B.02.xx has changed compared to previous versions. The unique combinations of *Type* and *Serial number* for each configured device are counted and tracked in the license module of the software. For each configured device, a license is deducted from the total number of eligible devices; the license status can be tracked in the Status Bar.

Connection Address: localhost Version B.02.01 [008] **Advanced | Licenses 7/10 ...**

Lab Advisor licenses acquired for previous versions of Lab Advisor and Lab Monitor and Diagnostic software are still eligible for Lab Advisor B.02.xx and are transformed according to [Table 4](#) on page 35.

Table 4 Lab Advisor Licenses

Product Number	Description	HW Modules	Replaces
M8550A	Agilent Lab Advisor Advanced	10	G4800AA, G4809AA
M8551A	5 add-on HW modules *	5	G4801AA
M8552A	25 add-on HW modules *	25	G4802AA
M8553A	50 add-on HW modules *	50	G4803AA
M8554A	100 add-on HW modules *	100	New

Table 4 Lab Advisor Licenses


Product Number	Description	HW Modules	Replaces
M8555A	Agilent Lab Advisor Basic	10	Agilent Instrument Utilities LC/CE

* Requires M8550A installed

Any combination of the listed Product Numbers is possible, and it is also possible to install a Product Number multiple times. Each unique license number adds to the number of modules supported.

Permanently deleting a hardware module

If a hardware module is retired or removed from a system, it can be permanently deleted from the license count in the Agilent Lab Advisor installation by clicking **Permanently delete HW module** in the **Licenses** tab in the **Configuration** screen. **Permanently delete HW module** is active only if all systems are disconnected in the **System Overview** screen. The **Remove HW Module** window that opens lists all modules contributing to the license limit; selecting the appropriate module and pressing **OK** permanently deletes the selected device.



NOTE

When a device is permanently deleted, all data belonging to the device is also permanently deleted.

Traceability

The **Traceability** feature of the software logs the use of Lab Advisor and keeps track of which user did what procedure and when. This information is written into the **Logs & Results** and is included in the printed results. If the **Data Sharing** feature is used, the traceability data is also uploaded.

Traceability Level

Lab Advisor offers three levels of traceability. Change of traceability level requires administrator rights. For the initial setting, the user *Admin* is set up using the password *Admin*. After the first login, this password should be changed to prevent unauthorized access.

The default level is **No Traceability**, which allows any user to access all parts of the Lab Advisor software, depending on Licenses. At this level, no user names are added to **Logs & Results** or printouts.

Medium Traceability requires that the user be selected from a drop-down list. No password is required, and users can register and enable themselves. The selected user name is added to the printouts and **Logs & Results**.


Full Traceability requires that the user logs in with a unique password. Users must be granted access by an Administrator before they can access the Lab Advisor software.

Setting Up a New User

To enable the Lab Advisor to print the user names in the reports and **Log & Results**, the users need to be set up. This is done in the **Traceability** tab of the **Configuration** screen.

- 1 Click **Add** at the bottom right of the **Traceability** screen.


The **Add Contact** dialog box is displayed.

A screenshot of the 'Add Contact' dialog box. It has a title bar with a close button (X). Inside, there are three text input fields labeled 'Name', 'E-Mail', and 'Phone'. At the bottom, there are two buttons: 'Ok' and 'Cancel'.

- 2 Enter the user name as it will be printed on reports, and optionally the user's email address and telephone number.

This information will be included in the **User** section of the **Status Report**.

Alternatively, a new user can be added by clicking **Register as new User** in the **Login** box:



The screenshot shows a 'Login' dialog box with a title bar containing a close button. Inside, there are two input fields: 'Username' and 'Password'. Below these fields are three buttons: 'Ok', 'Cancel', and 'Register as new User'. The 'Register as new User' button is highlighted with a red rectangular border.

If **Full Traceability** is active, an Administrator must activate the user before the user can use Lab Advisor:

Configured Users				
Name	E-Mail	Phone	Active	Admin
Admin			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
John Doe	John_Doe@Company....	0123-456789	<input type="checkbox"/>	<input type="checkbox"/>

Apps

Apps are small applications designed to help you perform specific non-system-related tasks.

Diagnostic Catalog

Diagnostic Catalog

Diagnostic Catalog

Device Class

Agilent LC Pumps

Device Type

G4220A "1290 Bin Pump"

Product Level

Advanced

Print

Results

Tests	Controls	Actuals
- Pump Head Leak Test	- A1	- Pressure
- Pump Leak Rate Test	- A2	- Flow
- System Pressure Test	- B1	- High Pressure Limit
Tools	- B2	- Pump Ripple
- Diagnostic Buffers	- Degasser On	Statemachines
- Module Infos	- Degasser Off	- Outlet Valve
- Purge Pump	- Start Purge Process	- Degasser
- Remove/install pump head	- Stop Purge Process	- Prime
	- Prime On	- Purge
	- Prime Off	- Pump
	- Conditioning On	- Flow Reduction
	- Conditioning Off	- Seal Wash Pump
	- Identify Module	Signals
	- Clear Error	- Pressure [bar]

Details

Name: Pump Head Leak Test

Approx. Time: 5 min

Description: The test determines the leakage of the individual pump heads.

Double-Click the selected results line in order to review online help

Back

The Diagnostic Catalog is a catalog of all tests, calibrations, tools, instrument controls and EMF counters for each module at each Lab Advisor product level. The list is filtered by **Device Class** (that is, instrument or module type), **Device**

Type and **Product Level**. The diagnostic catalog for the selected module at the selected product level is displayed in the **Results** table three columns:

- **Tests, Calibrations** and **Tools** available in Lab Advisor at the selected product level
- **Controls** provided in the **Instrument Control** screen of Lab Advisor at the selected product level
- **EMF Counters** shown the **EMFs** screen of Lab Advisor at the selected product level

For the tests, calibrations and tools, a short description is provided in the **Details** panel below the **Results** table. You can retrieve more information on the test, calibration or tool by double-clicking its name, which displays the online help.

You can print the diagnostic catalog for the current selection by clicking **Print**.

Data Sharing


The **Data Sharing** feature of Lab Advisor allows multiple Lab Advisor installations to upload and synchronize the collected instrument information and data to an upload folder. When activating the synchronization, all data and information for the systems and devices configured in the Lab Advisor installation are uploaded to a specified upload folder. If data and information for these systems and devices has already been uploaded to the upload folder, the uploaded data is merged with the existing data. Additionally, the delta between the two data sets is downloaded to the Lab Advisor installation to complement the data. This ensures that all participating Lab Advisor installations connecting to specific systems and devices have all data available regardless of which installation generated the data. The Lab Advisor installation receives data only from systems and devices configured in that installation. No other data is downloaded.

Review Client

The complete set of data uploaded to the data share folder can be accessed by the built-in Review Client, which is started from the **Data Sharing** app. The Review Client additionally supports the combination of any uploaded devices from any system to allow cross-system comparison of data. This might be

helpful in finding problematic systems or devices, or systems not being used efficiently.

The Review Client requires the number of licenses corresponding to the number of modules hosted on the upload share. Deleting a system from the **System Overview** will not free usable licenses. The number of licenses entered in the normal Lab Advisor installation is reused in the Review Client, so that if you have a 10-module license installed, you will only be able to start the Review Client if there are less than 10 modules uploaded to the upload share.



There are several ways to use the Review Client:

- All portable Lab Advisor installations upload their data to the same folder. Each portable Lab Advisor is then up to date, and a Review Client has the opportunity to observe the entire data pool.
- Each portable Lab Advisor has its own folder, which is used for back-up. By changing the share folder in the Review Client, it is possible to look into each separate user's data.

- Lab Advisor installations on local PCs connected to a single Instrument can use the synchronize function to upload data to a system-specific folder. This can be used as a backup solution, and by changing the share folder the Review Client can be used to look at each system separately.

In each case, the installation starting the Review Client needs enough licenses for all uploaded modules.

If system grouping has been activated in software configuration, you can set up groups of systems; the group controls allow you to switch between the groups. The groups you set up in the Review Client are independent of those set up in the System Overview.

Firmware Update

Lab Advisor can be used to update the device-internal software (called firmware – FW).

The **Firmware Update** screen lists all systems configured in the **System Overview** screen. The devices of these systems can be updated individually, or the entire system can be updated at one time. It is also possible to update multiple systems at one time. To start the FW update procedure, select the **Firmware Update** tab in the global screens section of the Navigation Panel.

Instrument NameProgress

Test System

Update

LockedUnlock

Component	Current Version	Update	Main System Firmware	Update	Resident System Firmware	Progress	Status
G4220A:DEBAA00157	B.06.41[0002], B.06.40	<input checked="" type="checkbox"/>	B.06.42[0001]	<input type="checkbox"/>	B.06.40	<div>0%</div>	
G4226A:DE93000560	A.06.34[008], A.06.32	<input checked="" type="checkbox"/>	A.06.36[003]	<input type="checkbox"/>	A.06.32	<div>0%</div>	
G4208A:PP55055002	B.02.12[0001], B.02.12	<input type="checkbox"/>	B.02.12[0001]	<input type="checkbox"/>	B.02.12	<div>0%</div>	
G4212A:DEBAF00163	B.06.41[0002], B.06.40	<input checked="" type="checkbox"/>	B.06.42[0001]	<input type="checkbox"/>	B.06.40	<div>0%</div>	

Firmware files are currently taken from: \\wadnas1\WAD-firmware\Public\Firmware\Current-Test

Change folder

Lock all

Default Firmware Folder: \\wadnas1\WAD-firmware\Public\Firmware\Current-Test

Set Default Firmware Folder

Logs and Results

The **Logs & Results** screen presents data collected from the configured devices, and helps to review the status of the systems or devices.

The **Logs & Results** data includes:

- Test results
- Error information
- FW revision and updates
- EMF changes
- Maintenance log entries

Each line in the log shows the module identifier (type and serial number), type of information, description and a time stamp. If the **Traceability** feature is in use, user-generated data is logged with user name in the **Message** field.

Filter

Systems: Laboratory 1
Test System 2
test
<Not System assigned>

Devices:
All
G1367C:DE00055018
G1329A:DE60556235
G4212A:DEBAF00163
G1313A:DE72004035
G5667A:LP00055727
G1367B:DE64556177

Source:
All
Result
Error
Maintenance
FW Revision
EMF

Time:
All
Last Hour
Last 24 Hours
Last 7 Days
Last 30 Days

Source Message Time

Source	Message	Time
G1367C	HiP ALS SL	10/20/2011 13:12:44
Serial #	DE00055018	
Result	Module Infos - Passed	10/20/2011 12:50:45
Error	Temperature out of range[EE 4114;:227400]	10/20/2011 10:56:24
Error	Temperature out of range[EE 4114;:227400]	10/20/2011 10:50:14
Error	Leak detected[EE64,0]	8/31/2011 15:22:33
Error	Leak detected[EE64,0]	8/31/2011 15:19:05
Error	Shutdown[EE63,0]	8/31/2011 12:58:33
Error	Shutdown[EE63,0]	5/30/2011 15:11:19
Error	Shutdown[EE63,0]	5/30/2011 15:03:06
Error	Shutdown[EE63,0]	4/29/2011 15:41:48
Error	Leak detected[EE64,0]	4/29/2011 15:39:45
Error	Transport initialization failed - Connection to G2254A fa...	4/27/2011 15:18:56
Error	Transport initialization failed - Connection to G2254A fa...	4/27/2011 15:14:53
Error	Shutdown[EE63,0]	4/15/2011 15:21:17
Error	Shutdown[EE63,0]	4/15/2011 14:13:15
Error	Shutdown[EE63,0]	4/14/2011 19:07:37
Error	Shutdown[EE63,0]	4/14/2011 15:52:09

Load external data Export filtered data Add Log Entry Print Result Review Result

For easy overview, you can filter the data by **Instruments**, **Devices**, **Source** or **Time**. Multiple selections are supported for **Devices** and **Source**, and can be selected by keeping **Ctrl** pressed while clicking the data required in the filter.

Filter

Systems: Laboratory 1
Test System 2
test
<Not System assigned>

Devices:
All
G1367C:DE00055018
G1329A:DE60556235
G4212A:DEBAF00163
G1313A:DE72004035
G5667A:LP00055727
G1367B:DE64556177

Source:
All
Result
Error
Maintenance
FW Revision
EMF

Time:
All
Last Hour
Last 24 Hours
Last 7 Days
Last 30 Days


Time View

The **Logs & Results** screen offers two modes of operation. The default mode is the **Module View**, which lists the devices by system, and presents the information per device. In the alternative **Time View**, the data is presented and sorted by time stamp. This allows for a system-wide overview of the sequence of the data.

3 Using Lab Advisor

Logs and Results

When you click **Add Log Entry**, you can type a log entry to the device, which is stored on the device main board. The information written to the device is limited to 50 characters; this is typically maintenance log data.



The data filtered by the built-in filter can be exported in ZIP format to enable it to be distributed. This is typically helpful if remote engineers need to assess data from the system. Clicking **Load external data** allows you to navigate to the storage location of the ZIP file and load the exported data.

NOTE

You can also use **Load external data** to load LRS files.

Service and Diagnostics

The **Service & Diagnostics** screen hosts the procedures (tests, calibrations and tools) of the Lab Advisor software. To select a procedure, select the device and then select the procedure from the list.

Filter

☒ Tests
 ☒ Calibrations
 ☒ Tools

Test System

	G4220A	1290 Bin Pump
	Serial #	DEBAA00157
	G4226A	1290 ALS
	Serial #	DE93000560
	G4212A	1290 DAD
	Serial #	DEBAF00163
	G4208A	1200 Instant Pilot
	Serial #	PP55055002

Pump Head Leak Test

Pump Leak Rate Test

System Pressure Test

Diagnostic Buffers

Module Infos

Purge Pump

Remove/install pump head

Name:

Pump Head Leak Test

Approx. Time:

5 min





Description:

The test determines the leakage of the individual pump heads.

more...

Run

For a better overview, you can filter the type of procedure you want to use.

Filter		
<input checked="" type="checkbox"/> Tests	<input type="checkbox"/> Calibrations	<input type="checkbox"/> Tools
Test system		
	G4220A 1290 Bin Pump	
	Serial # DEBAA00157	
	G4226A 1290 ALS	
	Serial # DE93000560	
	G4212A 1290 DAD	
	Serial # DEBAF00163	
	G4208A 1200 Instant Pilot	
	Serial # PP55055002	
Pump Head Leak Test		
Pump Leak Rate Test		
System Pressure Test		

Tests

Tests are procedures that result in a Passed/Failed statement, so the results of the test are compared with predefined limits. The Pump Pressure Test is an example of a test.

Calibrations

If internal calculations in the devices need to be corrected, calibration procedures normally take care of this. An example of a calibration is the Detector Wavelength Calibration. If you are operating in a controlled environment, this type of procedure might need to be verified. This could typically be done using a System Suitability Test.


Tools

Tools are procedures that have a supporting function and that do not produce a Passed/Failed statement when finished.

Instrument Control

The **Instrument Control** screen allows you to control a connected system without having to have a CDS running. This might be helpful in complex diagnostic situations, where the built-in diagnostic tests do not give a definitive answer.


Instrument Control in the Basic version provides only limited functionality, whereas the Advanced version provides a full set of controls and a freely configurable Signal Plot.



The Reply Panel of the **Instrument Control** screen displays any replies generated from the device, to verify that the control used was accepted; it shows only the last three replies. In order to get a complete history of replies, click **Save Session Results** in the Action Panel. The reply history is saved and can be viewed in the **Logs & Results** screen.

Actual Status Information

Each device is displayed separately in the Control Panel, and provides information about actual values. If a device has several actual values to display, the **more** link gives access to these values.



G4220A **1290 Bin Pump**
Serial # DEBAA00157
Firmware: B.06.41

Flow	0.000
Pressure	2.83
Pump Valve	0

[more](#)

▶ Controls

Controls

When the **Controls** section is expanded, a complete set of buttons becomes available providing extensive control of the device. This includes setpoints, controls, special commands and module information. When the buttons are clicked, a reply for the action (reply accepted or reply error) is displayed in the Reply panel. Initially, the setpoints display the value already loaded in the device; the displayed value changes on changing the setpoint. When a setpoint is changed, the change must be confirmed by clicking **Send**.

G4220A 1290 Bin Pump	
Serial #	DEBAA00286
Firmware	B.06.32

G4220A 1290 Bin Pump	
Serial #	DEBAA00286
Firmware	B.06.32
Flow [mL/min]	0.000
Pressure [Bar]	0.00
Pump Ripple [X]	0.00

Control

Purge + Prime

Set Purge Flow [mL/min] 10.00 Send
Set Purge Time [minutes] 1.00

Set Purge Composition [KB] 0 Send

Purge Process: ☐ On ☒ Off
Prime: ☐ On ☒ Off

Degasser Control

Set Vacuum Level [mBar] 100 Send

Degasser On/Off: ☒ On ☐ Off

Pump On/Off: ☒ On ☐ Off ☐ Standby
Conditioning: ☐ On ☒ Off

Method Parameters

Set Flow [mL/min] 0.000 Send
Set High Pressure Limit [Bar] 1000

Set Low Pressure Limit [Bar] 0 Send

Solvent Composition

Set Composition [%]

A B


90 10 Send

Signal Plot


The Signal Plot is used for monitoring specific function(s) of a system in real time. Combined with the Controls, it can provide very valuable troubleshooting information for experienced users. It can also be used to monitor the progress of certain tasks, and check when they are complete, saving time.

The Signals that you want to monitor are set up by clicking **Signal Configuration** in the Action Panel. The **Signal Config** dialog box that is displayed contains all available signals for the system. To select a signal, mark its check box and click **OK**.


3 Using Lab Advisor Instrument Control



The selected signals are visible in the Signal Plot starting in “lanes” mode, which divides the available area of the window between the number of configured signals. This gives each signal a limited size in the window, but all are easy to differentiate and each scale is shown on the left side of the plot.




Other scaling parameters can be selected by right-clicking the signal window and selecting the **Auto Scaling** from the context menu. The available scaling options are presented in the submenu.



Using the mouse pointer, it is also possible to scroll the scales directly. Click the scale you want to change and use the scroll wheel to change the scale. You can also change the placement of the scale by pressing the scroll wheel while moving the mouse forward and back.



Alternatively, you can specify a fixed scale window. Double-click the scale to open the **Scale** dialog box and enter the scale range and/or the lower starting value of the scale.



EMFs





Agilent Technologies LC Instrumentation has supported the Early Maintenance Feedback (EMF) feature since the introduction of the 1100 system in 1995, and continues to support this feature. EMF helps to keep the usage of devices in focus, and facilitates usage-based maintenance, which minimizes maintenance costs.

The EMF counters can be read and reset with both the Basic and the Advanced versions of the software, but the Advanced version allows limits to be activated and set in addition. Lab Advisor provides Agilent-recommended EMF limits. These limits have been determined by measurements under standard laboratory conditions, and do not take into account any application-, user- or site-specific conditions; to maximize the lifetimes of system components, they might need to be adjusted based on experience.

▼ Filter

☒ All Counters

☐ Counters with Limit

		Title	Value	Unit	Limit	Progress		
▼ Test system								
 G4220A Serial # DEBAA00157	1290 Bin Pump		0	Hour	3000	0%	★ ↺	
		3.45	Liter	50		6%	★ ↺	
		8.64	Liter	50		17%	★ ↺	
		631	Count	15000		4%	★ ↺	
	Liquimeter (A+B)	12.09	Liter	0		0%	★ ↺	
 G4226A Serial # DE93000560	1290 ALS		0	Count	1000	0%	★ ↺	
	Needle into seat counter	1191	Count	1500		79%	★ ↺	
		1.53	Hour	3000		0%	★ ↺	
	Valve switching counter	2418	Count	60000		4%	★ ↺	
 G4212A Serial # DEBAF00163	1290 DAD							
	Accumulated UV lamp on-time	2519.65	Hour	2000		100%	★ ↺	
	UV lamp ignition counter	28	Count	1500		1%	★ ↺	
	UV lamp on-time	360.65	Hour	0		0%	★ ↺	
 G4208A Serial # PP55055002	1200 Instant Pilot							

Activate EMF

Deactivate EMF

Refresh Counters

The EMF screen can be used to view all possible counters or, for better overview, to filter only those counters that have an activated limit.

▼ Filter

☒ All Counters

☐ Counters with Limit

Status Report

The Lab Advisor **Status Report** screen provides a system-wide overview of the devices in the system.

The information in the **Status Report** includes:

- Lab Advisor software information
- Contact information
- PC information (optional)
- System configuration
- Logbook
- EMF counters
- Test Results
- Instrument actuals (optional)

The information included in the **Status Report** can be used to document the system or to share diagnostic information with a remote engineer when troubleshooting the system.

When you start the **Status Report** screen, you can include initial information for Report Name, Contact Information, Company, Logs and results, PC information and Instrument Actuals.

The screenshot shows a web form for generating a Status Report. The form is titled 'Report Name' with a text box containing 'Agilent System Report'. Below this is a section for 'Contact Information' with fields for 'Name' (John Doe), 'Phone', and 'Email'. A 'Company' field contains 'Big Company Inc.'. The 'Included Information' section has a dropdown menu for 'Logs and results' set to 'Include all data', and two checkboxes: 'Include PC information' (checked) and 'Include instrument actuals' (unchecked). A 'Comments' section with a large text area is at the bottom. A 'Create Report' button is located in the bottom right corner of the form.

Report Name	Agilent System Report
Contact Information	
Name	John Doe
Phone	
Email	
Company	Big Company Inc.
Included Information	
Logs and results	Include all data
	<input checked="" type="checkbox"/> Include PC information
	<input type="checkbox"/> Include instrument actuals
Comments	

Create Report

Contact and Company information can be helpful for easy and precise identification if the report is sent to a Remote engineer during troubleshooting.

The Logs & Result information stored by each individual device might be extensive, so to reduce the amount of data, you can filter the data based on time.

If you mark the **Include PC information** check box, a list of Agilent programs installed on the PC is generated for the report. This includes all programs starting with *Agilent*.

The instrument actuals are the setpoints currently loaded in the system at the time of generation of the report. If a method has been loaded in the CDS (and has not been changed), then Lab Advisor can report these settings. Note, however, that the receiver of the **Status Report** will be able to see method information.

3 Using Lab Advisor Status Report

The screenshot shows a web browser window displaying the 'Agilent System Report' for a 'Test system'. The report is organized into three sections: Software Information, Contact Information, and PC Information. The browser's address bar shows 'Page Width' and search options. A 'Back' button is visible at the bottom right of the report area.

Agilent System Report		Test system
Software Information		
Lab Advisor Version:	Version B.02.01 [010]	
Data updated on:	5/20/2011 15:26:26	
Contact Information		
Name	John Doe	
Phone		
E-Mail		
Company	Big Company Inc.	
Comment		
PC Information		
Name	CZC9030ZHW	
Free Disk Space	C:\ 91.92 GB	
Installed Agilent Programs	Agilent Windows XP Security Template - AD Agilent Lab Advisor Agilent Installation Qualification Tool A.03.03.008 Agilent Internet Explorer8 Settings Agilent OLECM BPM Client	

Back

Index

A

action panel 26
 activate user 38
 add contact 37
 add log entry 46
 add system 28
 advanced version 7
 application panel 26
 applications 39
 apps 39
 automatic reconnect 29

B

backup 34
 basic version 6

C

calibrations 48
 catalog 39
 COM port 32
 COM port 29
 commands 50
 configuration 34
 connect via serial cable 32
 connection address 29
 contact information 56
 controls 50
 copy device details 32

D

data sharing 40
 delete hardware module 36
 device details 32

diagnostic catalog 39
 diagnostics 47

E

early maintenance feedback 54
 EMF changes 44
 EMF counters 56
 EMF 54
 error information 44
 export data 46

F

fast connect 32
 filter 45, 47, 55
 firmware revision 44
 firmware update 43
 full traceability 37

G

global tasks 25
 group controls 26
 groups 35, 42

H

hard drive installation 13
 hardware requirements 10
 help 26
 host name 29

I

import/export data 34
 information panel 26
 installation qualification 15

installation 13
 instrument control 49
 IP address 29, 33

L

LAN configuration 33
 language 13, 34
 license key 17
 license agreement 14
 licenses 35
 load external data 46
 location 14
 log entries 44
 log entry 46
 logbook 56
 logs 44
 LRS files 46

M

MAC address 33
 maintenance log entries 44
 medium traceability 37
 migrate from B.01.xx 21
 minimize 26
 module information 50
 module view 45

N

navigation panel 25
 navigation 24
 new user 37
 new system 28
 no traceability 37

Index

O

old data 15
on startup 33
online help 26
options 15

P

path 34
PC information 56

R

reconnect 29
remove hardware module 36
remove system 30
reply panel 49
report name 56
report 15, 56
requirements 10
restore old data 15
results 44
review client 40
RS232 32

S

scaling 52
serial cable 32
service 47
setpoints 50
signal plot 51
signal configuration 51
signal scaling 52
signals 51
software configuration 34
software requirements 10
software revision 26
status bar 27
status information 50

status report 56
synchronize data 40
system properties 29
system configuration 56
system groups 35
system name 29
system overview 28
system properties 29
system tasks 25
system type 29

T

tasks 25, 28
test results 44, 56
tests 48
time view 45
title bar 24
tools 48
trace file 34
traceability 37
tracing 34
type 29

U

updating 20
user interface 24

V

version 5, 26, 49

In This Book

This manual describes the two versions of Lab Advisor B.02.02: Lab Advisor Basic and Lab Advisor Advanced.

The manual contains the following information:

- Lab Advisor B.02.02 Overview
- Installation
- Using Lab Advisor

© Agilent Technologies 2011

Printed in Germany
11/2011



M8550-90001



Agilent Technologies