

Labware Stacker and Rack-Handling

User Guide

Original Instructions

Notices

© Agilent Technologies, Inc. 2011

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

User Guide Part Number

G5407-90001

Edition

Revision 01, December 2011

Updates: title page, Safety Notices, compliance declarations, general safety warnings, and ActiveX.

Revision 00, November 2003

Contact Information

Agilent Technologies Inc.
Automation Solutions
5301 Stevens Creek Blvd.
Santa Clara, CA 95051
USA

Technical Support: 1.800.979.4811
or +1.408.345.8011
service.automation@agilent.com

Customer Service: 1.866.428.9811
or +1.408.345.8356
orders.automation@agilent.com

European Service: +44 (0)1763208826
euroservice.automation@agilent.com

Documentation feedback:
documentation.automation@agilent.com

Web:
www.agilent.com/lifesciences/automation

Acknowledgements

Microsoft® and Windows® are either registered trademarks or trademarks of the Microsoft Corporation in the United States and other countries.

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or sub-contract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14

(June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

 **A WARNING or !!!INJURY HAZARD!!** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

A **CAUTION** or **!!DAMAGE HAZARD!!** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

Letter to our Customers

Dear Customer,

The Agilent Technologies acquisition of Velocity11 resulted in the following changes:

- Creation of Agilent Technologies Automation Solutions, formerly Velocity11
- Renaming of some Velocity11 products
- New Customer Service and Technical Support contact information
- New website address for product information

Please make a note of the following changes as they impact this user guide.

Velocity11 product name changes

Velocity11 product name	Changes to ...
Access2 Automated Microplate Loader	Automated Centrifuge Loader
Element Automation System	BioCel 900 System
IWorks Device Driver Programming Interface	VWorks DCL Interface
PlatePierce Seal Piercing Station	Microplate Seal Piercer
VCode Barcode Print and Apply Station	Microplate Barcode Labeler
Velocity11 Robot	3-Axis Robot
VHooks Integration Interface	VWorks Hooks Interface
VPrep Pipetting System	Vertical Pipetting Station
VSpin Microplate Centrifuge	Microplate Centrifuge
VStack Labware Stacker	Labware Stacker

New contact information

Documentation feedback: documentation.automation@agilent.com

Technical Support: 1.800.979.4811 or +1.408.345.8011
service.automation@agilent.com

Customer Service: 1.866.428.9811 or +1.408.345.8356
orders.automation@agilent.com

European Service: +44 (0) 1763208826
euroservice.automation@agilent.com

Web: <http://www.agilent.com>

Labware Rack Handling Guide


This guide explains how to handle the labware racks safely. The topics are:

- About the labware racks
- Carrying the racks
- Lifting the racks
- Loading labware into the racks
- Accessing product user documentation
- Contacting Velocity11

About the labware racks

The BenchCel® Microplate Handling Workstation and the VStack® Labware Stacker use labware racks to store the stacks of labware (microplates, tipboxes, and tube racks) that are processed during a protocol run.

Figure BenchCel Workstation with front-load labware racks



For a description of the BenchCel device or VStack device, see the *BenchCel User Guide* or *VStack User Guide*.

The labware racks are available in three models: standard rack, top-load rack, and front-load rack.


Figure Three models of labware racks



All the rack models have the following basic parts:

- *Carrying handle.* The standard rack has a polished top bar that can be used as a carrying handle. The top-load and front-load racks have fold-down carrying handles.
- *Tabs.* A pair of vertical tabs are located at the bottom sides of the rack. The tabs insert into slots on the device when you mount the rack.
- *Stacker grippers.* A gripper is located on the interior bottom of each tab. The pair of grippers hold a microplate during the labware loading, unloading, downstacking, and upstacking processes. A clamp in the device opens and closes the grippers.

Figure Standard rack containing labware



Carrying the racks


WARNING Do not hold a rack by the interior edges. The interior edges can have sharp surfaces that can cause cuts if handled improperly.

CAUTION A rack that is fully loaded with labware can be heavy. Grasp the rack handle firmly to prevent the rack from slipping or tilting.

To carry a rack:

Firmly grasp the rack by the handle.

Figure Carrying a front-load rack



Lifting the racks

Make sure to use the proper lifting technique when mounting a rack on a device or removing a rack from a device.


WARNING Avoid touching the interior edges of a rack when lifting the rack. The interior edges can have sharp surfaces.

IMPORTANT See your *BenchCel User Guide* or *VStack User Guide* for the procedure to mount a rack on the device or to release a rack from the device.

To lift a rack:

Use both hands to grasp the rack securely around the four corners near the base, as the following figure shows.

Figure Lifting a standard rack



Loading labware into the racks

Before loading the labware into a rack:

- Position the rack so that the opening is facing you.
- Determine how the microplates should be oriented in the rack.

For example, if the BenchCel orientation-sensing feature is enabled, make sure the A1 wells are oriented in the rack as specified.

Loading labware in standard and top-load racks

IMPORTANT See your *BenchCel User Guide* or *VStack User Guide* for the details on how to release a rack for removal or to prepare for loading a mounted rack.


To load labware into a standard or top-load rack:

- 1 If possible, remove the rack from the device, and place the rack on a flat, level surface.
- 2 Using both hands, carefully slide a small stack of labware down through the top of the rack.

You can use one hand to support underneath the labware stack, while the other hand holds the top of the labware to keep it level. See the following figure.


WARNING Use care to avoid sliding your hand on the interior edges in the rack. The edges can have sharp surfaces.

Figure Loading a standard rack



- 3** (Standard racks only) When you reach the bottom of the open slot, transfer your hand positions so that you continue supporting the labware through the bottom slot.

Figure Supporting labware through the bottom slot



- 4** Ensure that the bottom labware in the stack rests on the rack stacker grippers.


To unload labware from a standard or top-load rack:

With your hands positioned as shown in the previous figures, carefully slide the labware in small stacks up and out of the top of the rack.

Loading labware in a front-load rack

The doors on the front-load rack provide easy access for loading labware into the front of a rack that is mounted on a device.

Figure Door mechanism on the front-load rack




To load labware in a front-load rack:

IMPORTANT Before you attempt to load the labware in a mounted rack, ensure the device is ready for loading. For example, the clamps in the BenchCel stacker head must be closed (extended). See your device user guide for details.

- 1** On each side of the rack, slide the Door-release (black) buttons forward, while pushing outward on the thumb tabs. The rack doors open.

Figure Opening the front-load rack



- 2** Place the labware directly through the open rack doors so that the bottom labware rests on the rack stacker grippers. Ensure the labware is level and securely in the rack.
- 3** To close the doors, press the thumb tabs inward until the doors snap shut.

To unload labware from a front-load rack:

IMPORTANT Before you attempt to unload the labware from a mounted rack, ensure the device is ready for unloading. For details, see your device user guide.

- 1** To open the rack doors, slide the black Door-release buttons forward on each side of the rack, while pressing outward on the thumb tabs.
- 2** Carefully, lift the labware out through the front of the rack.

Accessing product user documentation

Velocity11 product user documentation is available in the following formats:

- Online help available within the software
- PDF files on the software CD
- Printed books

You can also search the online help or download the latest version of any PDF file from the Velocity11 website at:

http://www.velocity11.com/support/knowledge_base

Contacting Velocity11

- Technical Support: 1.800.979.4811 or +1.650.846.6611
- Customer Service: 1.866.428.9811 or +1.650.846.6601
- Email: info@velocity11.com
- Web: <http://www.velocity11.com>

Contents

Chapter 1. Introduction

About this Guide	2
Computer and Laboratory Requirements	4
General VStack Hazard Warnings	5
VStack Overview	6

Chapter 2. Installation

Unpacking the VStack	10
VStack Parts	12
Mounting the VStack	15
Connecting To Power and Air	17
About VStack Networking Communications	18
Installing the Network Card	20
Installing the DHCP Server	22
Connecting the VStack To the Network	24
Turning On the VStack	25
Installing Software	25
Loading and Releasing a Rack	28
Placing a Plate	30
Verifying Network Communication	31
Identifying VStacks	32
Changing the Device Id	33
Deleting Old VStacks	34

Chapter 3. Configuration and Manual Operation

Opening Diagnostics for a Specific VStack.	36
Selecting a Plate Type	38
Adding a Plate Type To the Labware Database.	39
Testing Plate Handling	40
Adjusting the Grip Pressure	41
Adjusting the Plate Specific Grip Offset	42
Viewing Sensor Readings	43
Configuring Plate Orientation Checking	45
Operating the VStack Manually	52

Displaying a Message.	54
About Teachpoints	54
Changing the Presentation Teachpoint.	56
About the Motor Page.	57

Chapter 4. Maintenance and Troubleshooting

VStack Error Messages.	60
Maintenance Schedule	66

Chapter 5. ActiveX Commands

About VStack ActiveX	69
--------------------------------	----

Introduction

1

This chapter introduces the VStack®, lists the system requirements, gives general safety information and provides an important overview of the product.

About this Guide

Intended Audience

This guide is intended to be read by:

- Everyone who uses the VStack
- Lab personnel who install and configure the VStack
- Developers who integrate the VStack into third-party systems


VStack Software Version

This version of the *VStack User Guide* guide is only for use with VStack ActiveX software version 4.0. You can find the version number after installation.

To find the version number of your VStack ActiveX software:

1. In your Windows® operating system, search for the file named *VStackBioNet.ocx*.
2. Right-click the file and select **Properties**.
3. Click the **Version** tab.

The version number is displayed. You only need to consider the first two parts of the version number.



**Lab Automation
System
Documentation**

The VStack is never used as a stand-alone device; it is always incorporated into a larger lab automation system that requires its own documentation.

If your lab automation platform is the Velocity11[®] BioCel[®], read this guide in combination with the *BioCel User Guide*. If your lab automation platform is a third-party system, read this guide in combination with the documentation that is available for your system.

**Other
Documentation**

In addition to this guide, read *Getting Started* and parts of the *Database and Security User Guide*.

Getting Started gives Velocity11 support information, explains how to use the online help versions of Velocity11 documentation and provides general safety guidelines.

The *Database and Security User Guide* gives information about the labware editor.

Content Summary

This guide covers:

- How to install the VStack hardware
 - How to connect the VStack to the network
 - How to install and use VStack ActiveX software
 - Maintenance and troubleshooting
 - ActiveX commands, which allow control of the VStack by third-party systems
-

Computer and Laboratory Requirements

Computer Requirements

- PC running Windows 2000 or XP
- At least a Pentium 166 processor
- 32 MB RAM
- 10 BaseT Ethernet port

Laboratory Requirements

-
- Bench space: 8.5 inches x 8.0 inches x 7.5 inches (21.60 cm x 20.30 cm x 19.05 cm) (L x W x H)
 - Electrical supply: 110–240 VAC, 50–60 Hz
 - Air supply: 50–80 psi at <1 cfm (minimum 15 psi)

!! DAMAGE HAZARD !! Operating the VStack at the wrong voltage may seriously damage the instrument.

!! IMPORTANT !! All labware (including deepwell and PCR plates) used in the VStack must conform to the SBS microplate standards. Please visit www.sbsonline.org/msdc/crevs.html for current SBS standards, or contact your labware's manufacturer.

General VStack Hazard Warnings

Before installing and using the VStack

Before installing and using the VStack, make sure you are aware of the potential hazards and understand how to avoid being exposed to them. You must be properly trained in the correct and safe installation and operation of the device.

EU installations only. Be aware that the VStack is considered partly completed machinery that is required to be installed with other equipment such that the fully assembled machinery complies with the essential health and safety requirements (EHSRs) of the Machinery Directive 2006/42/EC.

Intended product use

!! INJURY HAZARD !! Do not remove the VStack exterior covers or otherwise disassemble the robot. Doing so can expose you to hazards that could cause serious injury and damage the Labware Stacker.

!! INJURY HAZARD !! Using controls, making adjustments, or performing procedures other than those specified in the user guide can expose you to moving parts and hazardous voltage.

Agilent Technologies products must only be used in the manner described in the Agilent Technologies product user guides. Any other use may result in damage to the product or personal injury. Agilent Technologies is not responsible for any damages caused, in whole or in part, by improper use of the products, unauthorized alterations, adjustments or modifications to the products, failure to comply with procedures in Agilent Technologies product user guides, or use of the products in violation of applicable laws, rules or regulations. Except as otherwise expressly provided in Agilent Technologies product user guides, any alteration, adjustment, or modification to the products will void the product warranty.

The VStack is not intended or approved for diagnosis of disease in humans or animals. You assume full responsibility for obtaining any regulatory approvals required for such use and assume all liability in connection therewith.

Injury hazards

!! INJURY HAZARD !! Never lean against the VStack.

!! INJURY HAZARD !! Keep far enough away from the VStack while it is in motion so you do not catch your hands, clothing, jewelry and so on, in the mechanism.

!! INJURY HAZARD!! Do not start the VStack controlling software while you are near the VStack's plate stage. The stage moves when you turn on the VStack.

Damage hazards

!! DAMAGE HAZARD !! Always make sure that the plate stage is empty before you start to load or release a rack.

!! DAMAGE HAZARD !! Do not use harsh abrasives, corrosive cleaning agents or metal brushes to clean VStack components or accessories. Do not use any concentration of household bleach or sodium hypochlorite on the VStack. Do not allow cleaning agents to contact electrical, or sensitive mechanical components.

Safety labels

Warnings in the user documentation or on the device must be observed during all phases of operation, service, and repair of this device. Failure to comply with these precautions violates safety standards of design and the intended use of the product. Agilent Technologies assumes no liability for the customer's failure to comply with these requirements.

The following table lists the common symbols you might find on the device. The symbol on the label indicates the risk of danger. A description of the warning and information that will help you avoid the safety hazard are provided in this guide.

Symbol	Description
	Indicates that you must read the accompanying instructions (for example, the safety guide) for more information before proceeding.
	Indicates hazardous voltages.
	Indicates pinch, crush, and cut hazard.
	Indicates laser hazard.
	Indicates hot surface hazard.
	Indicates protective conductor terminal, which is bonded to conductive parts of an equipment for safety purposes.
	Indicates that you must not discard this electrical/ electronic product in domestic household waste.

Safety and regulatory compliance

The VStack complies with the applicable EU Directives. See the Declaration of Conformity or Declaration of Incorporation, as applicable, for details. The Labware Stacker is designed to comply with the standards listed in the following table.

Regulatory Compliance	Standard
EMC	
European Union	EMC Directive 2004/108/EC
	IEC 61326-1:2005 / EN 61326-1:2006
Canada	ICES/NMB-001:2004
Australia/New Zealand	AS/NZS CISPR 11:2004
Safety	
European Union	Machinery Directive 2006/42/EC
	Low Voltage Directive 2006/95/EC
	IEC 61010-1:2001 / EN61010-1:2001
Canada	CAN/CSA-C22.2 No. 61010-1-04
USA	ANSI/UL 61010-1:2004

Electromagnetic compatibility

If the VStack causes interference with radio or television reception, which can be determined by turning the device off and on, try one or more of the following measures:

- Relocate the radio or television antenna.
- Move the device away from the radio or television.
- Plug the device into a different electrical outlet, so that the device and the radio or television are on separate electrical circuits.
- Make sure that all peripheral devices are also certified.
- Make sure that appropriate cables are used to connect the device to peripheral equipment.
- Consult your equipment dealer, Agilent Technologies, or an experienced technician for assistance.

Changes or modifications not expressly approved by Agilent Technologies could void the user's authority to operate the equipment.

Sound emission declaration

Sound pressure: $L_p < 70$ dB according to EN 27779:1991.

Schalldruckpegel: $LP < 70$ dB nach EN 27779:1991.

VStack Overview

Product Description The VStack is a microplate stacker used in lab automation systems to dispense and receive plates. It holds deepwell and shallow-well plates, pintools, tipboxes and lids. Typically, several VStacks are used in a single lab automation system.

VStack Network The VStack communicates with the lab automation system using an Ethernet network and DHCP server, which dynamically recognizes VStacks as they are added to, and removed from, the system.

Software Overview The VStack has few hardware controls; it is controlled almost entirely by ActiveX commands issued by software installed on a computer.

You can interact with the VStack by:

- Writing a protocol that is then run automatically on a lab automation platform, according to a programmed schedule controlled by platform software
- Sending individual commands, in real time, directly to the VStack from software supplied with your VStack

The VStack is intended to be routinely used for automated runs, but you also need to control it directly.

Automated Control

When using the VStack during an automated protocol run, it is controlled by platform software. If you are using Velocity11's BioCel, the platform software is VWorks™. If you are running the VStack in a system developed by your own organization, the platform software will have been custom-written using a program such as C++ or Visual Basic.

Manual Control


You need to control the VStack outside a protocol run to:

- Recover from run errors
- Test Ethernet communications
- Troubleshoot problems
- Set plate stage positions

To do this, Velocity11 provides software with a user interface, called VStackBioNet Testcontainer software.

In addition, Velocity11 provides a labware editor and database that you use to define labware dimensions. The labware editor is accessed through the VStackBioNet Testcontainer software or VWorks.

Interaction of Software Components



ActiveX Software

ActiveX is a set of technologies that allows software components to interact with each other regardless of the language in which the components were created. Using the VStack ActiveX, simple commands can be used to initiate complex operations in networked VStacks.

The VStack's ActiveX software includes:

- VStack Diagnostics software
- Methods, which control individual operations
- Properties, which set the values used in methods, such as speed = fast
- Events, which are notifications that methods are complete or have errors

To allow you to develop your own controlling software, the chapter called [“ActiveX Commands” on page 67](#), provides the ActiveX methods, events and properties.

VStack Diagnostics

VStack Diagnostics is a window that is part of the ActiveX software. From the VStack Diagnostics window, you can:

- Perform simple VStack operations, such as Downstack, Upstack and Release Stack
- Enter properties, such as the plate specific grip offset, the type of plate to use and the text to show on the VStack display
- View sensor readings
- Access the labware editor
- Set teachpoint offsets, which provide positional information to the robot

Labware Database

The labware editor is used to enter the dimensions of each type of plate used with the VStack. This information, required for proper plate handling, is stored in the labware database, which is maintained in the Windows registry.

The same labware database is also used with other Velocity11 products. For this reason, the details of the labware database are given in a separate publication, the *Database and Security User Guide*.

Installation

2

This chapter explains how to install the VStack®. The installation process includes:

Step	Procedure
1.	“Unpacking the VStack” on page 10
2.	“Mounting the VStack” on page 15
3.	“Connecting To Power and Air” on page 17
4.	“Installing the Network Card” on page 20
5.	“Installing the DHCP Server” on page 22
6.	“Connecting the VStack To the Network” on page 24
7.	“Turning On the VStack” on page 25
8.	“Installing Software” on page 25
9.	“Loading and Releasing a Rack” on page 28
10.	“Verifying Network Communication” on page 31

Unpacking the VStack

Parts Inspection

As you unpack the VStack, inspect all parts for possible shipping damage. If anything is missing or appears to be damaged, contact the Velocity11 Service Center as soon as possible.

!! IMPORTANT !! Save the carton and packing materials in case you ever need to ship the VStack.

Please note that Velocity11 is not responsible for damage if the VStack is incorrectly packaged and moved by someone other than a Velocity11 employee.

North American VStack Package

The North American VStack package contains:

- A VStack and rack
 - A Power cord
 - An Ethernet cable, 8 feet
 - Polyurethane tubing, 10 feet, 1/4-inch diameter
 - Polyurethane tubing, 10 feet, 1/8-inch diameter
 - The following fittings:
 - ◆ 1/4-inch quick connect–1/4-inch NPT
 - ◆ 1/4-inch quick connect–1/8-inch NPT
 - ◆ 1/4-inch quick connect–1/8-inch quick connect
 - ◆ 1/4-inch quick connect T
 - ◆ 1/8-inch quick connect T
 - An ActiveX CD-ROM
 - The *VStack User Guide*
 - Getting Started*
 - The *Database and Security User Guide*
-


**European VStack
Package**

The European VStack package contains:

- A VStack and rack
 - A country-specific power cord
 - An Ethernet cable, 8 feet
 - Polyurethane tubing, 3 meters, 6-mm diameter
 - Polyurethane tubing, 3 meters, 4-mm diameter
 - Polyurethane tubing, 3 meters, 1/8-inch diameter
 - The following fittings:
 - ◆ 6-mm quick connect–4-mm quick connect
 - ◆ 4-mm quick connect–1/8-inch quick connect
 - ◆ 6-mm quick connect T
 - ◆ 4-mm quick connect T
 - ◆ 1/8-inch quick connect T
 - An ActiveX CD-ROM
 - The *VStack User Guide*
 - Getting Started*
 - The *Database and Security User Guide*
-

VStack Parts

Main Parts: Diagram The main parts of the VStack are shown in the following diagram. It will help you to be familiar with these terms as you read this guide.



**Main Parts:
Description**

The main parts of the VStack are described in the following table:


Part	Description
Air pressure regulator	<p>Controls the pressure of the air that is used to close the grippers.</p> <p>The two similar controls to the left and right are flow controls that affect the speed with which the grippers move in and out. You should not need to adjust the flow controls.</p>
Base	<p>Base of the VStack, which has dowel pin holes and screw holes for attaching the VStack to a custom-prepared lab bench or base plate.</p>
Display	<p>Displays a text string when the void DisplayMessage ActiveX method runs.</p> <p>You can run this method manually to test the display through the VStack Diagnostics software. For more information, see “Displaying a Message” on page 54.</p> <p>This method can also be called during a run by the run-controlling software. For example, on the Velocity11[®] BioCel[®], VWorks[™] uses this method to deliver the name of the type of plate that should be loaded into each VStack.</p>
Gripper	<p>When loaded, grips the bottom microplate, holding up the rest of the stack of plates. During operation it releases the bottom plate to the plate stage when dispensing plates and grips the plate on the plate stage when receiving plates from the robot.</p>
Load/Release button	<p>Pressing the button loads or releases a plate. When the plates are released, you can remove the rack.</p>
Plate stage	<p>Moves up and down to move plates between the rack and Presentation teachpoint.</p> <p>For more information about the Presentation teachpoint, see “About Teachpoints” on page 54.</p>
Rack	<p>Stores microplates. Can be removed from the VStack to add, remove and transport plates.</p>

VStack External Connections

The VStack has the following external connections:

- Quick-connect air input port, 1/8-inch
- Serial port
- Power receptacle
- Ethernet port

These connections are shown in the following diagram.



The serial port is used by Velocity11 to install firmware on the VStack.

Mounting the VStack

About Mounting the VStack

For the robot to move plates to and from the plate stage, the three-dimensional position of the plate stage must be fixed and known to the controlling software of the lab automation system. This requires that the VStack is mounted on a stable base that does not move in relation to the robot.

You can mount the VStack directly on a benchtop or on a mounting plate that also holds other components of your system. You should have the benchtop or mounting plate accurately drilled to create dowel pin holes, and tapped to create screw threads, according to the dimensions given in this section.

When you are planning where to install your VStack, make sure you choose a place that allows access on all sides for cleaning and maintenance.

Procedure



This procedure assumes that you have already prepared your bench area and have had the necessary holes drilled, and screw threads tapped, in your benchtop or mounting plate.

Note: Velocity11 does not provide dowel pins or screws. You will need to obtain pins and screws that are the correct length for your particular lab setup.

To mount the VStack:

1. Measure your prepared mounting holes to make sure that they will accommodate the VStack.

Dimensions in the diagram below are in inches.



2. Place two 1/4-inch dowel pins in the dowel pin holes of your benchtop.
You do not need to use all three dowel pin positions.
 3. Pick up the VStack and align the dowel pin holes in its base with the dowel pins of your benchtop.
 4. Place the VStack on the benchtop, making sure that the dowel pins fit properly into each hole.
 5. Screw two M6 flathead screws into two of the screw holes, making sure they screw completely into the base.
You do not need to use all screw positions.
 6. Check that the VStack is firmly mounted to your benchtop by trying to gently rock it.
-

Connecting To Power and Air

Introduction

This section explains how to connect your VStack to the power and air supplies. Connect the power cord first, followed by the air line.

Connecting the Power Cord

Review electrical requirements in “[Laboratory Requirements](#)” on page 4 before proceeding.

To connect the power cord:

1. Plug one end of the power cord into the power port at the back of the VStack.
2. Plug the other end of the power cord into an appropriately grounded electrical socket.

Connecting the Air Line

!! INJURY HAZARD !! Working with open, charged air lines can result in injury if the end of the line whips around rapidly. Switch the compressed air line off before installing the VStack. If you have questions about setting up your air line, contact your facilities department or the Velocity11 Customer Service Center.

!! DAMAGE HAZARD !! Make sure that air coming into the VStack is properly filtered for moisture and aerosol impurities; there is no air filter on the VStack. Dirt in unfiltered air can build up in air valves eventually causing a malfunction.

To connect the air line:

1. Cut the supplied 1/8-inch tubing to the length you need for the air line.
Cut the tubing with a sharp craft knife, and make sure that the edge of the tubing is square and clean.
2. Plug one end of the 1/8-inch tubing into your lab’s air supply.
You can also use 1/4-inch tubing from your lab’s air supply with an adapter to switch to 1/8-inch tubing. If your house air uses a threaded nozzle, attach the 1/8-inch NPT fitting to your house air.
3. Push the other end of the 1/8-inch tubing into the air input port on the back of the VStack.
4. Tug gently on the air line.
If you feel resistance, the line is in place.
5. Turn on the air and make sure that the air pressure does not exceed 120 psi.

The air pressure at the gripper is internally regulated to 40 psi, although it may be changed for particular plate types. For more information about adjusting the pressure, see “[Adjusting the Grip Pressure](#)” on page 41.

About VStack Networking Communications

Introduction

A VStack receives its operational commands from a computer to which it is networked. The network used for communication is a standard Ethernet network, consisting of a network card in the controlling computer and an Ethernet cable between the computer and VStack. If you want to have more than one VStack in your system, the network must include an Ethernet switch to serve as an adapter between the single Ethernet cable attached to your computer and the Ethernet cables of the multiple VStacks.

!! INJURY HAZARD !! Run your VStacks on a dedicated lab automation network and not a general organizational network. It is unsafe to allow control of VStacks from computers on a general network because remote computer operators may accidentally initiate an operation that causes the plate stage to move unexpectedly when lab personnel are nearby.

DHCP Server Kit

If you do not have an existing network for your VStacks, you will need to set one up, even if you are only using one VStack.

Note: The exception to this is if you are using one VStack and the controlling computer does not need to be connected to your organization's network. In this case, you can use the computer's existing network card and a crossover cable to connect to the VStack.

To set up the network, you can buy the DHCP Server Kit from Velocity11, which includes:

- A networking card
- An Ethernet cable
- An Ethernet switch
- DHCP server software

About IP Addresses

IP addresses allow computers and other devices to be identified on a network so information can be specifically routed to them. When setting up the VStack network, there are two implementations of IP addresses to consider:

- The automatically assigned IP address of each VStack
- The manually assigned IP address of the network card on the controlling computer that interfaces with the VStack network

Automatically Assigned IP Addresses

Each VStack on a network must have a unique IP address to allow the controlling computer to specifically send commands to it. IP addresses are assigned automatically by a DHCP server that you install on the controlling computer.

**Manually Assigned
IP Addresses**

When the DHCP server starts, it reads an initialization file that points the server to the network card through which the computer will serve addresses. When you install the network card, you must assign the IP address and subnet mask so that it matches the information in the initialization file.

**Overview of Setting
Up the Network**

To set up a network for your VStack or VStacks, perform the following procedures:

Step	Procedure
1.	“Installing the Network Card” on page 20
2.	“Installing the DHCP Server” on page 22
3.	“Connecting the VStack To the Network” on page 24
4.	“Turning On the VStack” on page 25
5.	“Installing Software” on page 25
6.	“Verifying Network Communication” on page 31

Installing the Network Card

Introduction

This section assumes that you have purchased the DHCP Server Kit from Velocity11.


The network card supplied with your DHCP Server Kit must be installed for your VStack to be able to communicate with the controlling computer, assuming that you want to use your existing network card to communicate with your organization's network.

Installation Procedure


Note: You may need additional administrator privileges on the controlling computer to install the network card. Please check with your network administrator if you think you may not have sufficient privileges.

To install and configure the network card:


1. Install a new network interface card:
 - a. Locate the PCI Adapter (network card) included in the kit.
 - b. Follow the instructions included with the card to install it.
2. Set the IP address of the card:
 - a. Click **Start > Settings > Network and Dial-up Connections**.
 - b. Right-click the icon for the new local area connection (probably **Local Area Connection 2**), and click **Properties**.



- c. Double-click **Internet Protocol (TCP/IP)**.



- d. Enter the values shown in the following screenshot.



- e. Click **OK**.
f. Close the **Network and Dial-up Connections** window.

For more information about the IP address of the network card, see [“About IP Addresses”](#) on page 18.

Installing the DHCP Server

Introduction

This section assumes that you have purchased the DHCP Server Kit from Velocity11.

The DHCP server automatically assigns IP addresses to VStacks when they are turned on so they can be identified on the network.

The server will run on Windows® XP, NT and 2000.

Installation Procedure

Note: You may need additional administrator privileges on the controlling computer to install the DHCP server. Please check with your network administrator if you think you may not have sufficient privileges.


To install and start the DHCP Server service:

1. Create the following directory on the controlling computer:
C:\Program Files\Velocity11\DHCP server
2. Extract the files from the DHCP Server CD-ROM into this new folder.
3. Double-click the file named dhcprsv.exe.

The **DHCP Server** window opens.



4. Click **Yes**.
This installs the DHCP server as a service on your computer.
The DHCP Server window opens and displays the status as **Running**.



5. Click the close box to close the window.

The DHCP Server service will now automatically run every time you turn on the controlling computer.

Stopping and Starting the Service

You can stop and start the DHCP Server service from within Windows by navigating to **Start > Programs > Administrative Tools > Services**. Consult your Windows documentation for more information.

You can also start and stop the service using the executable.

To stop and start the DHCP Server service:

1. Navigate to C:\Program Files\Velocity11\DHCP server.
2. Double-click the file named dhcpsrv.exe.
3. Click **Stop** to stop the service, and **Start** to start the service.

Connecting the VStack To the Network

Introduction

This section assumes that you have purchased the DHCP Server Kit from Velocity11.

Connecting a VStack To the Network

To connect the VStack to the network, you must connect the:

- Ethernet switch to the PC
- VStack to the Ethernet switch

To connect the Ethernet switch to the PC:

1. Locate the 8-port Ethernet switch included in the DHCP Server Kit.
2. Connect the included power cord to the switch.
3. Locate the Ethernet cable included in the DHCP Server Kit.
4. Plug one end of the cable into any of the eight ports in the Ethernet switch.
5. Plug the other end into the new network card in the PC.

To connect the VStack to the Ethernet switch:

1. Locate another Ethernet cable included with the Accessories.
2. Plug one end of the cable into any of the eight ports on the Ethernet switch.
3. Plug the other end into the Ethernet port of the VStack.

Connecting Multiple VStacks

You can connect up to six more VStacks to the network using the 8-port switch. Up to 25 VStacks can be used on this network if you add more 8-port switches.

Turning On the VStack

Introduction

Turn on the VStack after setting up the network and before installing the ActiveX software.

Turning on the VStack

To turn on the VStack:

1. Toggle the power switch to the on position.
The motor moves to the home position and the VStack processor obtains an IP address. The text on the display reads “Finding Server.”

Installing Software


Introduction

For general information about the software that controls the VStack, see “VStack Overview” on page 6.

Installation Procedure

To install the ActiveX software:

1. Insert the CD-ROM labelled *VStackBioNet ActiveX* into your PC.
2. In Windows, navigate to the CD-ROM drive and open it to view the files on the drive.
3. Double-click the Setup.exe file.
The InstallShield Wizard starts.




4. Click **Next**.
The **Customer Information** dialog box opens.



5. Enter your name, and your company's name.
6. Click **Next**.

The **Choose Destination Location** dialog box opens.



Unless you specify a different location, the installation program will install the ActiveX in a folder it creates, which is C: Program Files\Velocity11\VStack.

7. If you want to choose a different destination folder:
 - a. Click **Browse**.
 - b. Navigate to and select the folder.

- c. Click **Next**.

The **Select Program Folder** dialog box opens.



- d. In the Program Folders text box, type a name for the shortcut you want to be created in the **Start > Programs** menu, or leave the default name.
 8. Click **Next**.
 9. The **Setup Status** screen opens.
 10. When the text **InstallShield Wizard Complete** opens, click **Finish**.
-

Loading and Releasing a Rack

Loading and Releasing Defined

When you stack plates into a rack, the bottom plate sits on the gripper. When you install a rack into a VStack, the bottom plate drops, controlled by the plate stage, to be held on the side by the gripper. The rack is also held so it will not rock. This process is called loading.

When the plate stage moves the bottom plate up so it rests on the gripper the rack becomes free to move. This process is called releasing.

Note: The VStack will not operate until the rack is loaded.


!! INJURY HAZARD !! Keep far enough away from the VStack while it is in motion so that you do not catch your hands, clothing, jewelry, and so on, in the mechanism.



Loading a Rack

To load a rack:

1. Place a rack on the base unit, lowering it so that the rectangular tabs drop into the slots in the gripper unit.

The rack can be placed with the open corner at the front or back.





2. Either:
 - ◆ Press the load/release button on top of the base unit.
 - ◆ Click **Release Stack** in the **VStackBioNet Diagnostics** window (see [“Opening Diagnostics for a Specific VStack”](#) on page 36).The VStack loads the rack.

Releasing a Rack

To release a rack:

1. Either:
 - ◆ Press the load/release button on top of the base unit.
 - ◆ Click **Release Stack** in the **VStackBioNet Diagnostics** window.


!! DAMAGE HAZARD !! Always remove the rack from the base unit before attempting to move the VStack.
 2. Lift the rack up and off the base.
-

Placing a Plate

Plate Orientation

To place a plate on the plate stage:

1. Place the plate so that the A1 well is at the far left as viewed by the robot.




Verifying Network Communication

Procedure


After installing the ActiveX software, follow this procedure to verify that the computer and VStack are communicating correctly.

To verify that the Computer and VStack are communicating:

1. Navigate to **Start > Programs > Velocity11** and select **VStackBioNet**. The **VStackBioNet Testcontainer** software opens.



2. Click **Display Diagnostics**.



If the **Discovered BioNet Devices** window opens and the VStack has a status of **Matched**, you know that your computer is communicating with the VStack.

For more information about identifying VStacks on the network, see [“Identifying VStacks” on page 32](#).

3. Click **OK**.
4. If you want to verify that you are communicating with a particular VStack click **Move Here** at the Orientation Sensor teachpoint of the **VStackBioNet Diagnostics** dialog box. The plate stage will move to the Orientation teachpoint. For more information, see [“Opening Diagnostics for a Specific VStack” on page 36s](#).

Identifying VStacks

Introduction

A VStack's IP address may change as you turn the VStack on and off. To allow you to select a particular VStack to communicate with, there must be another way to identify a specific VStack. The Device Id, which does not change when the VStack is turned on and off, is used for this purpose.

Device Ids are assigned in the order that VStacks are first added to the network.


You should either keep a record of which Device Id corresponds to which VStack, or rename the Device Ids to values that are easily recognizable, see [“Changing the Device Id” on page 33](#).

Procedure

To identify the VStacks on your network:

1. Turn off all VStacks at their on/off switches.
2. Turn on one VStack.
3. Navigate to **Start > Programs > Velocity11** and select **VStackBioNet**.
4. Click **Display Diagnostics**.

One or more VStacks appear in the dialog box, but only one has a **Status** of **Matched**.



5. Note the Device Id of the VStack.
6. Turn on another VStack, click **Cancel** in the **Discovered BioNet Devices** dialog box, and click **Display Diagnostics** again.
Note the **Device Id** of the new VStack.
7. Repeat [step 6](#) on other VStacks in turn, recording their Device Ids.

Changing the Device Id

Introduction


The Device Id identifies a VStack on the network. You can change the Device Id to an easily recognizable name or number, such as the VStack's serial number, or a simple identifying number that you have applied to the outside of the VStack.

Procedure

Note: Device Ids are specific to a computer. If you change a Device Id, the change will not be transferred to other computers used to control the VStacks on the network. You must change the Device Id on each computer.

To change the Device Id:

1. Navigate to **Start > Programs > Velocity11** and select **VStackBioNet**. The VStackBioNet Testcontainer software opens.



2. Click **Display Diagnostics**.
3. Click to select the VStack from the list.
4. Click the number in the **Device Id** column. The cursor appears.
5. Type in the new Device Id.

Device Id	Device Type	IP Address	MAC Address	Status
500000	VStack	10.10.10.173	00-90-C2-C0-33-2C	Matched

Deleting Old VStacks


Introduction

If you remove a VStack from the network it will still show up in the Discovered BioNet Devices window (see [“Identifying VStacks” on page 32](#)). The entry will be reused if you add the VStack back later. If you do not plan to add the VStack back into the network, you may want to delete the entry.

Procedure

To delete an unused device:

1. Navigate to **Start > Programs > Velocity11** and select **VStackBioNet**.
2. Click **Display Diagnostics**.
3. Click in the **Status** column to reveal a drop-down list.
4. From the list box, select **Delete**.



Note: Other commands in this list box are for service use only.

Configuration and Manual Operation

3

This chapter shows you how to set up a VStack® and perform manual, real-time operations.

If your lab automation system is already set up and you only need to check that a particular plate works with a particular stacker, follow the procedures in this table.

Step	Procedure
1.	“Opening Diagnostics for a Specific VStack” on page 36
2.	“Selecting a Plate Type” on page 38
3.	“Testing Plate Handling” on page 40

If you have just installed the VStack you will need to completely configure it, using the procedures in the above table, as well one or more of the following procedures, depending on your situation:

- [“Adding a Plate Type To the Labware Database” on page 39](#)
- [“Adjusting the Grip Pressure” on page 41](#)
- [“Adjusting the Plate Specific Grip Offset” on page 42](#)
- [“Configuring Plate Orientation Checking” on page 45](#)
- [“Displaying a Message” on page 54](#)
- [“Changing the Presentation Teachpoint” on page 56](#)

In addition, you may need to perform some procedures given in the *Database and Security User Guide*. This guide will direct you to those procedures at the appropriate places.

After your VStack is set up and routinely running protocols, you may occasionally need to make minor adjustments using these same procedures.

Opening Diagnostics for a Specific VStack

Introduction


Velocity11 provides the VStackBioNet Diagnostics software for configuration and manual control of the VStack. Many configuration procedures first require you to open the VStackBioNet Diagnostics software for a specific VStack.

Procedure

Note: If you are using the VStack on the BioCel, refer to the *BioCel User Guide* for how to open the VStack Diagnostics software from within VWorks.

To open the VStack Diagnostics software:

1. Navigate to **Start > Programs > Velocity11** and select **VStackBioNet**. This opens the VStackBioNet Testcontainer software.




The Cycle Test button in the VStackBioNet Testcontainer dialog box is used during manufacturing to test the VStack. It has no purpose during normal use.

2. Click **Display Diagnostics**.

The **Discovered BioNet Devices** window appears. For more information, see [“Verifying Network Communication” on page 31](#).

3. Select your VStack from the list.
4. Click **OK**.

The **VStackBioNet Diagnostics** window opens.



Selecting a Plate Type

Introduction


For the VStack to correctly handle a plate during manual, real-time operation, you must select the type of plate you intend to use from the labware database.

If the type of plate you intend to use is not already defined in the labware database, you will need to create an entry for the plate and then configure the entry's properties before returning to this section and selecting the plate type. For more information, see [“Adding a Plate Type To the Labware Database”](#) on page 39.

Procedure

To select a plate for a manual operation:

1. Open the VStackBioNet Diagnostics window.
For directions, see [“Opening Diagnostics for a Specific VStack”](#) on page 36.
2. In the **Plate Selection** group box of the **Diagnostics** page, select a plate from the list box.



3. Select one of the following options.

Option	Description
Normal	The plates you plan to use are not sealed and have no lids.
Lidded	The plates have lids.
Sealed	The plates are sealed.

4. If you want the VStack to sense the orientation of the plate, select the **Orientation Checking** check box.

If this option is not selected, the plate stage still moves to the Orientation teachpoint during upstack and downstack operations, but the orientation is registered as correct regardless of the actual orientation of the plate.

5. Click **Apply** to save the changes.

What To Do Next

If you are setting up the VStack for the first time, read through the remainder of this chapter and check all configuration settings before testing plate handling.

If your VStack is already set up and you are just checking plate handling after making an adjustment or selecting a new plate type, follow the procedure in [“Testing Plate Handling”](#) on page 40 next.

Adding a Plate Type To the Labware Database

Background

For correct operation of the VStack the software controlling the robot must be given the dimensions of the type of plate you intend to use. This information is held in the labware database. If the type of plate you want to use is already defined in the database, you can ignore this section. If the type of plate you want to use is not already defined start by following the procedure in this section.

Opening the Labware Editor


Note: If you are using the VStack on the BioCel, refer to the *BioCel User Guide* for how to open the labware editor from within VWorks. Then refer to the *Database and Security User Guide* for how to use the labware editor.

To open the labware editor:

1. Open VStackBioNet Diagnostics window.

For directions, see [“Opening Diagnostics for a Specific VStack” on page 36](#).

2. In the Plate Selection group box, click **Edit**.



3. Refer to the *Database and Security User Guide* for directions for adding a plate type.

The *Database and Security User Guide* will direct you back to this guide for some of the procedures you need to perform.

Testing Plate Handling

Introduction

After selecting the plate type you want to use, perform some downstacking and upstacking operations with the plate to make sure that the plate is handled correctly by the VStack.

If the VStack performs without error, it is ready for use with an automated protocol.

If you receive errors:

- Check the labware editor parameters for the labware you are using. For more information, see the *Database and Security User Guide*.
- Check the troubleshooting information given in “[Maintenance and Troubleshooting](#)” on page 59.
- Modify the grip pressure, as described in “[Adjusting the Grip Pressure](#)” on page 41.

Downstacking Plates

To downstack a plate:

1. Open VStackBioNet Diagnostics.
For directions, see “[Opening Diagnostics for a Specific VStack](#)” on page 36.
2. Select your plate type from the list box, and click **Apply**.
3. Load a rack containing plates.
Make sure that your plates are properly positioned in the rack. For more information, see “[Placing a Plate](#)” on page 30.
4. Click **Release Stack** to lock the rack into place.
5. Click **Downstack**.
6. Remove the plate from the stage.
7. Click **Downstack**, and repeat until you are satisfied that the downstacking operation is working correctly.

Upstacking Plates

To upstack a plate:

1. Place a plate on the stage.
 2. Click **Upstack**.
 3. Repeat until you are satisfied that the upstacking operation is working correctly.
-

Adjusting the Grip Pressure

Definition

The grip pressure determines the force with which the VStack grips the plates. It is a property of a specific VStack and not the type of plate.

You can see the sensor reading for grip pressure in the VStackBioNet Diagnostics window. For more information about sensor readings, see [“Viewing Sensor Readings” on page 43](#).

When To Adjust the Grip Pressure

If the plates you are using are flimsy, they may distort when held with the default pressure of 40 psi so you should lower the grip pressure.

If the stack of plates you are using is heavy because there are many plates in the stack, and perhaps because they are made of a heavy material and contain a lot of liquid, the bottom plate may fall out of the gripper. In this case, increase the grip pressure. The grip pressure, however, should not exceed the maximum pressure specified in [“Laboratory Requirements” on page 4](#).

Grip Pressure Adjustment

To adjust the grip pressure:

1. Locate the air pressure regulator.
For more information, see [“VStack Parts” on page 12](#).
2. Using your fingers, or a pair of pliers, unscrew the locking ring by turning it counterclockwise.
3. Adjust the grip pressure by turning the regulator clockwise to increase the pressure and counterclockwise to decrease the pressure.

!! DAMAGE HAZARD !! Do not set the grip pressure lower than 15 psi.

4. Toggle the grippers 2–3 times to allow the pressure to adjust.
 5. Lock the locking ring by turning it clockwise.
-

Adjusting the Plate Specific Grip Offset

Definition	The VStackBioNet Diagnostics software's Plate specific grip offset specifies the distance below the standard Grip teachpoint.
When To Adjust	Adjust this value only if the VStack is not gripping the plates at the correct position.
The Plate Specific Grip Offset Is a Proxy	<p>The Plate specific grip offset is a property of the plate and not a property of the VStack. It is actually a proxy of the labware editor's Stacker gripper offset property, meaning that if you change the value of the plate specific grip offset, the same change is made in the labware editor for that particular type of plate. The reverse is also true.</p> <p>For more information about the labware editor, see the <i>Database and Security User Guide</i>.</p>
Procedure	<p>To change the plate specific grip offset:</p> <ol style="list-style-type: none"> Place at least two plates in the rack. Load the rack onto the VStack. For more information, see “Loading and Releasing a Rack” on page 28. Open the VStackBioNet Diagnostics software. For more information, see “Opening Diagnostics for a Specific VStack” on page 36. Select the appropriate plate type. Load the rack by clicking Release Stack. Look at the grippers from the side through the rack. If the grip position is where you want it, do not change the gripper offset. If it is not, continue with step 7. Click Release Stack to move the plate out of the grippers. Change the Plate specific grip offset value, which is in millimeters, and click Apply. <div data-bbox="522 1564 836 1612" style="border: 1px solid #ccc; padding: 2px; margin: 5px 0;"> Plate specific grip offset <input style="width: 60px;" type="text" value="0.4"/> </div> <p>Enter a number higher than the set value to grip the plate above the current grip point.</p> <p>Enter a number lower than the set value to grip the plate below the current grip point.</p> <ol style="list-style-type: none"> Click Release Stack to load the plate into the grippers.

- Repeat [step 7](#) to [step 9](#) to change this value until the VStack grippers hold the labware at the position you want.

Viewing Sensor Readings

Procedure

To view sensor readings:

- Open the VStack Diagnostics software.
For more information, see [“Opening Diagnostics for a Specific VStack”](#) on page 36.
- View the sensor display.



Sensor Types

The following sensors are used on the VStack.

Sensor	Description
Left Notch	Detects the presence of a notch in the left corner of the plate when the plate stage is at the Orientation Sensor teachpoint. A white number indicates the absence of a notch. An orange number indicates the presence of a notch. For more information about the Orientation Sensor teachpoint, see “About Teachpoints” on page 54.
Right Notch	Detects the presence of a notch in the right corner of the plate when the plate stage is at the Orientation Sensor or Grip teachpoint. A white number indicates the absence of a notch. An orange number indicates the presence of a notch.
Grip Pressure	Indicates the grip pressure when the grippers are closed, in pounds per square inch. A positive number indicates that the grippers are closed. When the grippers are open, the value always reads zero. For more information, see “Adjusting the Grip Pressure” on page 41.

Sensor	Description
Rack Present	<p>A white number indicates that a rack is loaded on the VStack.</p> <p>An orange number of less than 20 indicates that there is no rack loaded on the VStack.</p>
Plate Present	<p>A white number indicates that a plate is in the VStack grippers.</p> <p>An orange number indicates that no plate is present in the grippers.</p>
Stage Position	<p>This sensor is the large, unlabeled value. It indicates the plate stage position in relation to the home position at the bottom. The measurements are in millimeters, with positive values meaning above the home position.</p>

Configuring Plate Orientation Checking

About Orientation Sensing

The VStack has two sensors that sense the presence of a notch in the left or right far corners of a plate (from the robot's perspective). Although there are only two notch sensors, the software is aware of the existence of four plate corners, so that if a plate with two notches is on the stage, and no notches are detected, you will receive an error message telling you to turn the plate around 180°.

The controls that affect orientation checking are listed in the following table. The default values will not support accurate orientation checking so they must be correctly set for your VStack and plates.


Control	Description	Where Set
Orientation checking check box	Turns orientation checking on and off for both manual operations and protocol runs using the VStack. <i>Note:</i> There is no relationship between this check box and the Check orientation check box in the labware editor, which is used only for Velocity11's BenchCel. For more information, see the <i>Database and Security User Guide</i> .	VStack BioNet Diagnostics
Plate specific orientation sensor offset	Adjusts the height of the plate stage for the selected plate type, relative to the Orientation sensor teachpoint. <i>Note:</i> This property is a proxy for the Orientation sensor offset in the labware editor. When you change the value here, the value of the Orientation sensor offset also changes. The reverse is also true. For convenience you may decide to only adjust this value in the labware editor. Adjusting this value avoids having to adjust the actual Orientation Sensor teachpoint. !! IMPORTANT !! Do not change the value of this teachpoint.	Either: <input type="checkbox"/> VStack BioNet Diagnostics <input type="checkbox"/> Labware editor equivalent property, Orientation sensor offset
Sensor intensity	Adjusts the relative brightness of <i>all</i> sensor lamps. The default value is 50%. For more information, see the <i>Database and Security User Guide</i> .	Labware editor
Orientation Sensor threshold	Sets the Orientation sensor reading that defines the boundary between the indicated presence and absence of a plate notch. A sensor value below the threshold indicates the presence of a notch, and a sensor value above the threshold indicates the absence of a notch. The default threshold is 0 units, which must be changed for correct orientation checking. For more information, see the <i>Database and Security User Guide</i> .	Labware editor

You may need to adjust one or more of these settings to optimize orientation checking accuracy.

!! IMPORTANT !! The Orientation sensor teachpoint also affects orientation checking, but does not change its value.

Overview of Setting Up Orientation Checking

The overall process of setting up orientation checking is shown in the following diagram.



Use this diagram in combination with the following procedures to set up orientation checking.

- [❑ “Turning On Orientation Checking” on page 48](#)
- [❑ “Changing the Default Orientation Sensor Value” on page 49](#)
- [❑ “Placing a Plate and Examining Sensor Readings” on page 50](#)
- [❑ “Adjusting the Sensor Intensity” on page 50](#)
- [❑ “Adjusting the Orientation Sensor Threshold” on page 50](#)
- [❑ “Adjusting the Orientation Sensor Offset” on page 50](#)


The overall goal is to obtain Left Notch and Right Notch sensor values that are widely separated, and an Orientation sensor threshold that lies roughly half way between those values.

Turning On Orientation Checking

Note: If orientation checking is turned off, the plate stage still moves to the Orientation teachpoint during upstack and downstack operations, but the orientation is registered as correct regardless of the actual orientation of the plate.

To turn on plate orientation checking:

1. Open the VStack Diagnostics software.
For more information, see [“Opening Diagnostics for a Specific VStack” on page 36](#).
2. In the **Plate Selection** group box, select the **Orientation Checking** check box.




3. Click **Apply** to save the changes.
-

Changing the Default Orientation Sensor Value

The default setting for the plate specific orientation sensor offset (and therefore also the labware editor's orientation sensor offset) is zero. You must change this value for orientation sensing to work.

You can calculate an initial value by taking the average of the height of the plate plus the height of the plate's skirt $(x+y)/2$ where x and y are dimensions shown in the following diagram).



This value may not need to be further adjusted.

To enter the initial plate specific orientation sensor offset value:

1. Open the VStack Diagnostics software.
For more information, see [“Opening Diagnostics for a Specific VStack” on page 36](#).
2. Enter the value for the **Plate specific orientation offset**.
3. Click **Apply** to save the changes.

Placing a Plate and Examining Sensor Readings

Before making further adjustments to sensor settings, place the plate you want to use on the VStack and examine the current readings.

To examine the sensor readings:

1. Open the VStackBioNet Diagnostics software.
For directions, see [“Opening Diagnostics for a Specific VStack” on page 36](#).
2. Make sure that the plate you plan to use has only one notch, positioned in either of the far corners, and is selected in the **Plate Selection** list box.
If the plate type does not appear in the **Plate Selection** list box, add definitions for the plate in the labware editor. For more information, see [“Adding a Plate Type To the Labware Database” on page 39](#).
3. Click **Apply**.
4. Place the plate on the plate stage.
For information about plate orientation, see [“Placing a Plate” on page 30](#).
5. In the **Teachpoints** section next to the **Orientation Sensor** value, click **Move Here** to move the stage to the Orientation sensor position.



6. Examine the sensor readings and refer to the flowchart in [“Overview of Setting Up Orientation Checking” on page 47](#) for what to do next.

Adjusting the Sensor Intensity

The sensor intensity is adjusted using the labware editor. See the definition of this property in the *Database and Security User Guide*. Also, see the flowchart in [“Overview of Setting Up Orientation Checking” on page 47](#) for when to adjust this property.

Adjusting the Orientation Sensor Threshold

The Orientation sensor threshold is adjusted using the labware editor. See the definition of this property in the *Database and Security User Guide*. Also, see the flowchart in [“Overview of Setting Up Orientation Checking” on page 47](#) for when to adjust this property.

Adjusting the Orientation Sensor Offset

When you added a labware entry, the plate specific orientation sensor offset was set to zero. Earlier, you entered a default value calculated from plate measurements ([“Changing the Default Orientation Sensor Value” on page 49](#)). If this value is not accurate enough for your specific plate type and VStack, you will have been directed here from the overview flowchart. Follow this procedure to obtain a more accurate offset value.

Note that this property can be changed in either of two places. The plate specific orientation sensor offset is a proxy of the labware editor's

orientation sensor offset. When you make a change in one place, the change is automatically copied to the other place. The procedure in this section explains how to change the plate specific orientation sensor offset. For more information about the orientation sensor offset, see the *Database and Security User Guide*.

To adjust the plate specific orientation sensor offset:

1. Remove the rack from the VStack.
For more information see, [“Loading and Releasing a Rack” on page 28](#).
2. Place a plate that has a single notch on the plate stage.
For more information see, [“Placing a Plate” on page 30](#).
3. In the **VStackBioNet Diagnostics** dialog box, click **Move Here** for the Orientation sensor teachpoint.
The plate stage moves to the Orientation sensor position.
4. In the **Jog Controls** group box, set the slider to a value of one (1 mm), and click **Down** until the left and right sensor readings are both zero.
At this point, the light beam shines above the plate. The Rack Present sensor should also read zero.
5. Click **Up** until the sensor corresponding to the corner with the notch has a white value.
This position indicates that the top of the plate is just aligned with the top of the sensor’s vertical range.
6. Write down the value of the stage position, the large number in the sensor display.
7. Continue to click **Up** until both left and right notch values are large and white.
At this point, the light beams on both sides are hitting the top of the plate skirt and reflecting back.
8. Record the stage position value.
9. Calculate the average of the two values you recorded.
This represents the mid point of the valid sensor range.
10. Click **Move Here** for the **Orientation Sensor** teachpoint to move the plate stage to the teachpoint.
11. Compare the stage position value of the teachpoint with the value that you calculated in [step 9](#).
The difference represents the Plate specific orientation offset you should apply to achieve an optimal reading position.
12. Enter the value for the **Plate specific orientation offset**.
For example, if the Orientation Sensor teachpoint stage position value is 65 mm and the calculated value is 70 mm, enter an offset value of 5 mm.
13. Click **Apply** to save the changes.

Operating the VStack Manually

Introduction

Once your VStack is configured, you can operate it manually using VStackBioNetDiagnostics software. If you configured the VStack you will have already used some of these controls.

Procedure


To perform a manual operation:

1. Open the VStackBioNetDiagnostics software.
For more information, see [“Opening Diagnostics for a Specific VStack” on page 36](#).
2. Click one of the following buttons.

Command	Description
Toggle Gripper	Opens or closes the grippers. Do not use the Gripper Delay slider without proper instruction from a Velocity11 representative. The default value is 250. This value synchronizes the closing of the gripper with the movement of the plate stage and you should not have to change it.
Operations Group Box	
Downstack	Transfers plates from the rack to the plate stage.
Upstack	Transfers plates from the plate stage to the rack.
Release Stack	Allows you to lock or unlock the rack from the base.
Home	Moves the plate stage to the home position and then the presentation position.
Jog Controls Group Box	
Up	Moves the plate stage up by the distance in millimeters set in the adjacent slider control.
Down	Moves the plate stage down by the distance in millimeters set in the adjacent slider control.

The **In Motion** light indicates that the plate stage is moving.

The **Busy** light indicates that the VStack is performing an operation and cannot perform another operation until the current one is complete.



Displaying a Message

Introduction

The VStack display shows text messages. Long messages are scrolled across the display.

You can display your own message on the VStack display.

Procedure

To display a message:

1. Open the VStack Diagnostics software.
For more information, see [“Opening Diagnostics for a Specific VStack” on page 36](#).
2. Type your message in the **VFD Display** text box.
3. Click **Set**.

About Teachpoints

Teachpoints Defined

A teachpoint is a position to which the plate stage moves when commanded. As the plate stage only moves in a vertical plane, teachpoints are defined by a single number, which is the height of the stage above the home position, in millimeters. The home position lies at the bottom of the stage’s range of movement and has a positional value of 0 mm.

Teachpoints are specific to a particular VStack.

Teachpoint Locations


There are four teachpoints, described in the following table.

Teachpoint	Definition	Default Value (mm)
Stack Release	The height at which the plate stage releases the rack.	100.00 !! IMPORTANT !! Do not change this value
Grip	The height at which the grippers grip the plate during an upstack or downstack operation.	84.00 !! IMPORTANT !! Do not change this value for VStacks used on the BioCel.
Orientation Sensor	The height at which the VStack senses orientation notches on a plate.	68.00 !! IMPORTANT !! Do not change this value for VStacks used on the BioCel.

Teachpoint	Definition	Default Value (mm)
Presentation	The height at which the stage presents a plate for pick-up by a robot. See “Changing the Presentation Teachpoint” on page 56 for information about changing this value.	2.00 !! IMPORTANT !! Do not change this value for VStacks used on the BioCel.

Teachpoint Values

Teachpoint values are displayed in the Teach Points group box of the VStackBioNet Diagnostics window.



Changing the Presentation Teachpoint

Introduction

The Presentation teachpoint specifies the height above the home position at which the VStack presents plates. This teachpoint should only be changed if a particular configuration requires a Presentation teachpoint that is different from the 2.00 mm default.

!! IMPORTANT !! Do not change this value for VStacks used on the BioCel.

Note that the labware editor includes a property called Presentation offset. This property is used by the BioCel's VWorks. For more information, see the *Database and Security User Guide*.

Procedure

Note: When controlling the VStack from other computers, the default teachpoint will appear, unless it is changed on each computer.

To change the Presentation teachpoint:

1. Open the VStack Diagnostics software.

For more information, see [“Opening Diagnostics for a Specific VStack” on page 36](#).


2. Click **Move Here** for the **Presentation** teachpoint.

This moves the plate stage to the presentation position.



3. Move the plate stage up or down using the jog controls.


The increment control allows you to move the plate stage up or down in large or small moves. Adjust the increment control to the desired increment and click **Up** or **Down** to move the plate stage.



4. When the plate stage is at the desired height, click **Teach Here**.



5. At the **Modify Teach Point** dialog box, click **OK** to accept the changes or **Cancel** to cancel the operation and return to the previous screen.



The Presentation teachpoint is now saved.

About the Motor Page

Motor Page Screenshot

!! IMPORTANT !! The Motor page is intended for use during manufacture and you should not use it without supervision by a Velocity11® service representative.

If you accidentally change a value, refer to the screenshot below for the default value to change it back to.

	Empty	Full		
Kp:	60	60	High Travel Limit	100
Ki:	2	2	Low Travel Limit	0
Kd:	40	40		
Integrator limit:	10	10		
Output limit:	253	253		
Deadband compensation:	0	0		
Velocity (mm/s):	200	150		
Acceleration (mm/s ²):	3000	1000		
Servo rate:	1	1		
Position error limit (mm):	1000	1000		

Maintenance and Troubleshooting

4

This chapter:

- ❑ Describes the common VStack® errors, and their possible causes and solutions
- ❑ Gives a maintenance schedule that you can use to keep the VStack in excellent working condition

Note: The error messages described in this chapter are generated by the VStack ActiveX. If you are using the VStack with the BioCel's VWorks, the error messages you receive may be slightly different.

VStack Error Messages

Introduction


In this section, the text of the error message is written as a heading, which is followed by a description of the error, possible causes, and solutions.

Labware Definition Errors

If the labware definitions are not set appropriately for the plate you are using, you may get an error of some kind, depending on the specific problem. If you receive a VStack error, make sure you consider the labware settings in the labware editor.

Cannot Perform Operation Without First Loading Rack


The VStack failed to sense a rack when you clicked the Upstack button, even though a rack was present.



Possible Cause	Solution
The rack is not properly positioned in the VStack.	Check that the rack is positioned correctly in the rack grippers and click OK . For more information, see “Loading and Releasing a Rack” on page 28.
The Rack Present sensor has failed.	In the VStack Diagnostics window, make sure that the Rack Present sensor reads a value greater than 20 when a rack is loaded. For more information, see “Opening Diagnostics for a Specific VStack” on page 36. If the sensor reads a value of less than 20, and the rack is properly loaded, the sensor has malfunctioned. Contact the Velocity11 Service Center for help.

Could Not Perform Operation With Plate on Stage

The VStack sensed a plate on the stage when you clicked the Release Stack or Downstack button.



Possible Cause	Solution
There is a plate on the stage.	Remove the plate and click OK .
The sensor threshold is set too low.	In the labware editor, set the sensor threshold to the appropriate value for that plate type and attempt the operation again. For more information, see the <i>Database and Security User Guide</i> .
One or both of the notch sensors has failed.	In the VStack Diagnostics window, both notch sensors should read values close to zero when there is nothing in front of them. If not, the sensor has malfunctioned. Contact the Velocity11 Service Center for help. For more information, see “Opening Diagnostics for a Specific VStack” on page 36.

No Plate on Stage

The VStack failed to sense a plate when you clicked the Upstack button.



Possible Cause	Solution
There is no plate available to upstack.	Place a plate on the stage, click OK , and attempt the command again.
The sensor threshold is set too high.	Adjust the sensor threshold value. For more information, see the <i>Database and Security User Guide</i> .
One or both of the notch sensors has failed.	If the sensor intensity is set correctly, and the notch sensor values do not change when a plate is on the stage in the orientation sensor position, one or both of the notch sensors have failed. Contact the Velocity11 Service Center for help.

No Plate Available

The VStack failed to sense a plate when you clicked the Downstack button.

Possible Cause	Solution
A rack is present (loaded or unloaded) but it contains no more plates.	Add more plates.
A rack is present and it contains plates, but it is not loaded at the time you press the Downstack button.	Load the rack.
The VStack is not correctly configured. For example, you are using a black plate and the sensor intensity is set too low.	Follow the directions in “Configuration and Manual Operation” on page 35 to configure the VStack.
The plate present sensor has failed.	Contact the Velocity11 Service Center for help.

**Plate Is Rotated
 180°**

The VStack sensed a notch in the incorrect position.



Possible Cause	Solution
<p>The plate is rotated by 180° so that the notch is not in the specified position.</p>	<p>Rotate the plate to the specified position and click OK. For more information, see the <i>Database and Security User Guide</i>.</p>
<p>The notches are set incorrectly in the labware editor.</p>	<p>In the labware editor, make sure that the correct notch location is checked, then make sure that the plate is oriented correctly. <i>Note:</i> If the problem persists and the run must continue, turn off orientation checking. For more information, see “Selecting a Plate Type” on page 38.</p>
<p>The orientation sensors are not configured properly.</p>	<p>Make sure that the sensor intensity and sensor threshold configurations are correct. For more information, see the <i>Database and Security User Guide</i>.</p>
<p>One or both of the notch sensors has failed.</p>	<p>If the sensor values do not change with the presence or absence of a plate in the orientation sensor position, then one or more sensors may have failed. Please contact Velocity11’s Service Center for assistance.</p>

Plate on Stage Is the Wrong Type


The VStack sensed a type of plate that is different from the type of plate selected.



Possible Cause	Solution
The incorrect plate was selected.	In the Plate Selection group box, select the plate type you are using.
The notches are set incorrectly in the labware editor.	Choose the correct plate in the Plate Selection menu. Edit the notch location or locations of that plate. For more information, see the <i>Database and Security User Guide</i> .
The plate is rotated 180° so that the notch is not in the specified position.	Rotate the plate to the specified position and click OK . For more information, see the <i>Database and Security User Guide</i> .
The orientation sensors are not configured properly.	Make sure that the sensor intensity and sensor threshold configurations are correct. For more information, see the <i>Database and Security User Guide</i> .
One or both of the notch sensors has failed.	If the sensor values do not change with the presence or absence of a plate in the orientation sensor position, one or more sensors may have failed. Contact the Velocity11 Service Center for help.

Rack Not Present

The VStack failed to sense a rack when you clicked the Release Stack, Upstack or Downstack buttons.



Possible Cause	Solution
The rack is not properly positioned in the VStack, or there is no rack.	Check that the rack is positioned correctly in the rack grippers and click OK . For more information, see “Loading and Releasing a Rack” on page 28.
The Rack Present sensor has failed.	In the VStack Diagnostics window, make sure that the Rack Sensor reads a value greater than 20 when a rack is loaded. If it does not, the sensor has failed.

Stage Position Error

The VStack stage failed to move to the commanded position.



Possible Cause	Solution
Something is obstructing the VStack stage and preventing it from reaching its commanded location.	Check for debris that may be blocking the VStack stage.
The stacking thickness is not set correctly in the labware editor.	Measure the thickness of the plate and make sure that the correct value is entered. For more information, see the <i>Database and Security User Guide</i> .
The plate is not sitting properly on the plate stage.	Reset the plate and click OK .
A mechanical or communications failure has occurred.	Click OK and see if the VStack will recover.

Maintenance Schedule

Schedule

Maintenance tasks for the VStack are listed in the following table.

Maintenance Task	Time Schedule	Symptoms
Wipe grippers with Ethanol	As needed	Grippers appear dirty or are sticking to plates
Replace gripper pads	As needed	Gripper pads are damaged or worn

If you have any questions about maintenance, please call the Velocity11 Service Center.

ActiveX Commands

5

The ActiveX Template Engine (ATL) is used to operate the VStack when integrated with a third-party system. The ATL is located in the following directory:
C:\Program Files\Microsoft Visual Studio\VC98\ATL

This page is intentionally left blank.

Index

A

- ActiveX, 7
 - events, 73
 - finding the version, 2
 - methods, 70
 - properties, 68
- air
 - connecting to, 17
 - default pressure, 17
 - input port, 14
 - maximum pressure, 17
 - operating pressure, 4
- air pressure
 - regulator, 13

B

- base, 13
- boolean Busy ActiveX property, 69
- boolean CheckPlateSensor ActiveX method, 72
- boolean GripperState ActiveX property, 68
- boolean IsRackLoaded ActiveX method, 72
- boolean IsRackPresent ActiveX method, 72
- boolean IsStageEmpty ActiveX method, 71
- boolean Open ActiveX method, 70
- boolean OpenIP ActiveX method, 71
- boolean OrientationChecking ActiveX property, 69
- boolean Vacuum ActiveX property, 69
- BSTR PlateName ActiveX property, 68
- Busy light, 52

C

- communication verification, 31
- computer requirements, 4

D

- damage hazards, 5
- Device Id
 - about, 33
 - changing, 33
- DHCP server
 - installing, 22
 - IP address, 19
 - subnet mask, 19
- Diagnostics software, 52, 54, 56
 - opening, 36
- display, 13
- documentation to read, 3
- dowel pins, 15
- Down command, 52

- Downstack command, 52
- downstacking a plate, 40

E

- error messages, 60
- Ethernet port, 14
- Ethernet switch, 24
- events, ActiveX, 73

G

- grip pressure, 41
- Grip Pressure sensor, 43
- Grip teachpoint, 54
- gripper location, 13

H

- hazards. *see* safety
- Home command, 52

I

- In Motion light, 52
- IP address, 18
- IPictureDisp* ControlPicture ActiveX property, 68

J

- jog controls, 52

L

- laboratory requirements, 4
- labware database, 8
- labware editor, opening, 39
- Left Notch sensor, 43
- load/release button, 13
- loading a rack, 28
- long SensorFrequency ActiveX property, 68
- longGripDelay ActiveXProperty, 68

M

- maintenance schedule, 66
- manual operation, 52
- Motor tabbed page, 57
- mounting the VStack, 15

N

- network
 - connecting VStacks, 24
 - identifying VStacks, 32, 34
 - verifying communication, 31
- network card installation, 20

networking kit, 18

O

Orientation Sensor teachpoint, 54

P

Plate Present sensor, 44

plate stage, 13

plate testing, 40

plate type selection, 38

power

connecting to, 17

operating voltage and current, 4

receptacle, 14

Presentation teachpoint, 55

changing, 56

preventive maintenance, 66

R

rack, 13

loading and releasing, 28

Rack Present sensor, 44

Release Stack command, 52

releasing a rack, 28

Right Notch sensor, 43

S

safety

general hazards, 5

injury hazards, 5

screws, 15

selecting a plate type, 38

sensor

Grip Pressure, 43

Left Notch, 43

Plate Present, 44

Rack Present, 44

Right Notch, 43

Stage Position, 44

types, 43

viewing readings, 43

serial port, 14

short Downstack ActiveX method, 70

short Home ActiveX method, 71

short Jog ActiveX method, 71

short PlateDimensionSelect ActiveX property, 68

short ReleaseStack ActiveX method, 70

short Upstack ActiveX method, 70

software

Diagnostics, 8

installing, 25

overview, 6

Stack Release teachpoint, 54

Stage Position sensor, 44

starting. *see* VStack, 25

switch, 24

system requirements. *see* computer requirements, 4

T

teachpoint

about, 54

Grip, 54

Orientation Sensor, 54

Presentation, 55

changing, 56

Stack Release, 54

testing plates, 40

Toggle Gripper command, 52

U

Up command, 52

Upstack command, 52

upstacking a plate, 40

V

void AboutBox ActiveX method, 72

void Close ActiveX method, 70

void DisplayMessage ActiveX method, 71

void SetButtonMode ActiveX method, 71

void ShowDiagsDialog ActiveX method, 72

VStack

basic description, 6

external connections, 14

identifying on network, 32, 33

manual operation, 52

overview, 6

parts description, 13

parts diagram, 12

parts inspection, 10

turning on, 25

unpacking, 10

VStack error messages, 60

VStackBioNet.ocx. *see* ActiveX, 2



User Guide

G5407-90001

Revision 01, December 2011