

# Agilent Technologies

# Latest Enhancements to the GC/MS Analysis of Gasoline by ASTM Method D-5769

# **Application Note**

#### Author

Roger J. Leibrand Agilent Technologies

#### Abstract

The 6890/5973 GC/MSD system and its latest ChemStation software each contain enhancements that refine the analysis of aromatics in gasoline by ASTM Method D-5769. The improvements in the MSD hardware are the addition of a high energy dynode electron multiplier and increased scanning speed. The improvements in the software are easy entry of weights, easy signal to noise calculation, weight and volume percent calculated in the customized report, use of API index number, and expanded Help files. The improvements in the GC/MS method are fast temperature programming, carrier gas programming, high split ratio, narrow mass range, and electron multiplier programming. The performance of the 6890/5973 GC/MSD system with respect to ASTM Method D-5769 has been characterized. It was found that excellent toluene linearity, fast 20 minute analysis times, and excellent precision on duplicates could be achieved.

#### **Keywords**

Benzene, toluene, aromatics, gasoline, reformulated gasoline, RFG, gas chromatography/mass spectrometry, GC/MS, ASTM D-5769

#### Introduction

ASTM Method D-5769 is for analysts who have to certify the aromatic content in finished motor gasoline. These analysts are gasoline producers, manufacturers, blenders, distributors, independent testing laboratories, and regulatory bodies. A recent review by Vince Giarrocco of Agilent Technologies details the background of the

#### Agilent Technologies 6890/5973 GC/MSD System

reformulated gasoline analytical methodology.<sup>1</sup> An earlier Agilent Technologies Applications Brief describes the use of the 5972 MSD as applied to the first version of this ASTM method.<sup>2</sup>

The latest enhancements afforded by using the 6890/5973 GC/MSD here provide an improved means to validate compliance with government regulations and enhance productivity, and are more specific to the latest version of D-5769. The hardware and software recommendations are listed in Table 1.

Table 1. Hardware and software recommendations

GC:	6890
GC Consumables:	60 m × 0.25 mm ID, 1.0μm HP-1 P/N 19091Z-236
	Split/Splitless liner with glass wool P/N 19251-60540
	Merlin Microseal P/N 5181-8833
ALS:	7673 with nanoliter adapter
MS:	5973A Mass Selective Detector
MSD ChemStation:	Aromatics in Gasoline mode MSD Productivity ChemStation Software

# **Software Enhancements**

The improvements in the software are easy entry of sample weights, internal standard weights, density, easy signal to noise calculation, inclusion of uncalibrated compounds in the calculations, weight and volume percent calculated in the customized report, use of API index number, and expanded Help files that are specific to this analysis.

When recalibration utilizes a new batch of calibration standards, the individual amounts defining the calibration curve will need to be updated. This process is automated if the calibration standards are purchased from a vendor that also supplies an update disk. The individual amounts will be inserted for each compound that has a name and a CAS number that matches the one in the file.

The gasoline analysis method ASTM D5769 calls for a signal to noise check on the extracted ion chromatogram from the GC/MS analysis of a 0.01 weight % solution of 1,4-diethylbenzene. With a few simple click and drag operations the check is made and a special report is printed out.

The calibrated components do not account for all the aromatics in gasoline. These components are identified by the existence of peaks with characteristic ions in specified retention time ranges. The quantitation of the uncalibrated components utilizes the calibration curve of an appropriate calibrated component. The ASTM method also calls for a report of the weight % and volume % for each aromatic compound in the sample. These are provided by the software via a standard internal standard report that is enhanced with sample weight, sample density, ISTD weight, weight %, and volume %. Calculation of these quantities requires sample specific information. The sample weight and sample density need to be input by the operator.

There are also several editing features that are enhancements beyond the standard quantitation software, such as use of sample density or API index number and the y-intercept test for each calibration curve.

# **GC/MS Enhancements**

There are two approaches to GC analysis of gasoline. One requires long analysis times, high resolution columns, and universal detectors. The other uses mass selective detectors and can result in shorter analysis times with less stringent chromatographic resolution requirements. Detailed hydrocarbon analysis (DHA) is an example of the former approach while D-5769 and its enhancements described here is an example of the later. This ASTM method has only a few mandated conditions. The refinement of the final method is at the individual laboratory discretion. The major requirements are that a methylsilicone column is used and the mass spectrometer is operated in full scan mode generating electron impact ions at 70 electron volts. There are also requirements of

minimum column resolution, sensitivity, and data collection rate; these will be discussed later.

# Experimental

In this work several chromatographic method enhancements were employed to enhance throughput. In the past for Method D-5769 run times of 30 to 50 minutes have been common in some laboratories. In order to attain a more desirable 20 minute run time, a starting oven temperature of 120° was used with a relatively fast two-level temperature program.3 Holding the helium carrier gas flow at a constant 35 cm/sec during the run was also employed. This is 40% above the flow used for maximum efficiency, i.e., maximum chromatographic resolution. The details are shown in Table 2.

Even with 0.1 microliter injected and split the suggested 250 to 1, column overloading and toluene nonlinearity were observed. By raising the split ratio to 1100 to 1, the GC peaks shapes improved markedly. But at 15% toluene levels a nonlinear curve was still observed. Method D-5769 calls for the use of m/z 91 for toluene quantitation. This is the base peak. However, use of the carbon 13 isotope peak of m/z 92, which is 54% of m/z 91 produces a linear calibration curve.<sup>4</sup> Use of m/z 65, which is 10% of m/z 91 also, gives a linear calibration curve. Of the three possible mass peaks, the use of m/z 92 gives the best signal to noise along with a linear calibration curve. The high use of helium for the 1100 to 1 split ratio is

#### Table 2. Recommended Instrumental Parameters

GC Injection Port:	250°C Split 1100:1
GC Oven Ramp:	120°C for 0.5 minutes 3°C/minute to 140°C 10°C/minute to 250°C 0 minutes hold at 250°C
Carrier Gas:	Helium
GC Pressure Program:	0.3 mL/minute for 1.0 minutes 1.0 mL/min at 1.0 minutes 35 cm/second constant flow Vacuum compensation ON Gas saver ON after 1 minute
ALS:	5μL syringe Nanoliter adapter, stop = 1 (0.1μL injection size)
MSD:	Transfer Line: 280°C Ion Source: 200°C Quadrupole: 150°C Solvent Delay: 3 minutes Electron Energy: 70 eV Emission current: 35 µamps Mass Range: 60 to 170 Daltons Scan Speed: 2 <sup>3</sup> A/D samples Autotune: Standard Multiplier Voltage: Autotune value (1153 volts) +200 volts after 6.0 minutes

somewhat wasteful and so the "gas saver" mode is used to reduce the split flow at times other than injection.

The 5973 uses a high energy dynode electron multiplier that incorporates enhanced ion optics to reduce noise and increase sensitivity. Special ion focusing elements concentrate the ion beam exiting the quadrupole on the 10 kV high energy dynode where a conversion to electrons takes place. This gives higher sensitivity and longer multiplier life. The high energy dynode is  $90^{\circ}$  off-axis to minimize noise. Further detector noise reduction is accomplished by careful shielding of the electron multiplier and high energy dynode. This type of detector has increasing sensitivity with mass, which improves the detection of the heavier gasoline components.

The electron multiplier voltage obtained with Autotune is used for benzene and toluene. It is then raised programmatically 200 volts higher for the remaining components. This assures that the very highest toluene levels that are in some gasoline blending stocks will not overload the detector electronics and that there is sufficient signal for the low level components containing ten carbons and above.

Method D-5769 requires at least five mass spectral scans to be taken over the top half of the GC peak. According to Mathews and Hayes<sup>5</sup> if the worst case error is to be kept below 1%, at least nine scans must be taken. There is a real quantitation benefit to having as many scans as possible over the GC peak. However, the signal to noise ratio is adversely affected. Therefore, a compromise is needed. The 5973 scans twice as fast as its predecessor, the 5972. In real terms this means that with an A/D setting of 2^3 selected, 6.5 scans per second are taken. For the 1,4-diethylbenzene peak this represents 27 scans, reducing the sampling error to 0.1%, which is quite acceptable. These conditions yield a signal to noise ratio of 71 to 1 for the 0.01 wt % sample of 1,4-diethylbenzene.

Additionally, the mass range to use was explored. While Method D-5769 suggests that a mass range of 50 to 350 could be used, it is important to note that the highest mass required in the quantitation is m/z 148. Some of the higher uncalibrated compounds that may need to be quantitated separately have ions no higher than m/z 168.<sup>6</sup> Therefore, a narrow mass range of m/z 60 to 170 was used to maximize the number of scans.

### Results

Figure 1 is the total ion chromatogram of a gasoline test sample, RFA003 QC. Note all compounds elute between 3 and 20 minutes. Figure 2 is a typical report for an RFA sample using the Aromatics in Gasoline software. Figure 3 shows the final linear calibration curve for toluene utilizing the high split ratio, electron multiplier voltage programming, and change of quant ion. The analysis of gasoline by Method D-5769 is not simple and requires careful attention to operational parameters as well as sample preparation and further method refinements.<sup>4</sup> It was found that excellent toluene linearity, fast 20 minute analysis times, and excellent precision on duplicates could be achieved with the 6890/5973 GC/MSD.



Figure 1. Total ion chromatogram of a gasoline test sample, RFA003 QC.

Quant Time: Sep 2 18:53 1996       Quant Results File: RFG6890B.R         SampWt: 1.2818       SampDens: 0.7501       IstdWE: 0.0569         Quant Method : C:\HPCHEM\2\METHODS\RFG6890B.M (RTE Integrator)       Title         Title : Gasoline Aromatics System       Last Update : Mon Sep 02 18:45:01 1996         Response via : Initial Calibration       DataAcq Meth : RFG6890B         Internal Standards       R.T. QIon Response Conc Units Dev(Min)         1) benzene-d6       3.94       84 143289       22.76 mg       0.02         3) ethylbenzene-d10       6.32 116       191828       2.76 mg       0.02         3) ethylbenzene-d10       6.32 116       191828       2.76 mg       0.02         System Monitoring Compounds       Wt% Vol%       Vol%         Aromatic Compounds       Wt% Vol%       1.33       1.13         4) toluene       4.98       92 492998       9.21       7.92         5) ethylbenzene       6.57       106       978127       6.13       5.29         7) o-xylene       7.72       120       9665       0.08       0.07         9) n-proylbenzene       8.37<120       64838       0.52       0.45 #         11<1-methyl-4-ethylbenzene       8.51<120       1.57941       0.44       0.73	Data File : C:\HPCHEM\2\DATA\92\1801006.D Acq On : 2 Sep 96 17:27 Sample : rfa3 Misc : 5973 MS Integration Params: RTEINT P				Vial: 18 Operator: R Leibrand Inst : 5973 Multiplr: 1.00		
Quant Method : C:\HPCHEM\2\METHODS\RFG6890B.M (RTE Integrator)         Title : Gasoline Aromatics System         Last Update : Mon Sep 02 18:45:01 1996         Response via : Initial Calibration         DataAcq Meth : RFG6890B         Internal Standards R.T. QIon Response Conc Units Dev(Min)         1) benzene-d6 3.94 64 143289 22.76 mg 0.02         3) ethylbenzene-d10 6.32 116 191828 22.76 mg 0.02         3) ethylbenzene-d8 13.52 136 519787 11.38 mg 0.02         System Monitoring Compounds         Aromatic Compounds         Aromatic Compounds         Aromatic Compounds         At toluene 4.98 92 492998 9.21 7.92         5) ethylbenzene 6.43 106 3321130 2.45 2.100         6) máp xylene 6.57 106 978127 6.13 5.29         7) o-xylene 7.72 120 9665 0.088 0.07         9) n-propylbenzene 7.72 120 9665 0.088 0.07         9) n-propylbenzene 8.37 120 64838 0.52 0.45 #         10) 1-methyl-3-ethylbenzene 8.51 120 157941 0.64 0.73         13) 1-methyl-4-ethylbenzene 8.67 120 157941 0.44 0.73         13) 1-methyl-2-2         15) 11, 2, 3-trimethylbenzene 8.67 120 17941 0.44 0.73         13) 1-methyl-2-2         15) 12, 2, 4-trimethylbenzene 10.65 119 197255 0.73 0.63 #         10) 1.2, 4-trimethylbenzene 10.62 107 700 0.30 0.23         11) 1, 2, 4, 5-terramethylbenzene 10.62 107 700 0.30 0.23         11,	Quant Time: Sep 2 18:53 1996 SampWt: 1.2818 SampDens: 0	.7501	Istd	Quant Res Wt: 0.05	ults Fil 69	e: RFG6	890B.RF
Internal Standards         R.T. QIon         Response         Conc Units Dev(Min)           1) benzene-d6         3.94         84         143289         22.76 mg         0.02           3) ethylbenzene-d10         6.32         116         191828         22.76 mg         0.03           17) naphthalene-d8         13.52         136         519787         11.38 mg         0.02           System Monitoring Compounds	Quant Method : C:\HPCHEM\2\MET. Title : Gasoline Aromat Last Update : Mon Sep 02 18:4 Response via : Initial Calibra DataAcq Meth : RFG6890B	HODS\RF ics Sys 5:01 19 tion	G6890 tem 96	B.M (RTE I	ntegrato	r)	
1) benzene-d6       3.94       84       143289       22.76 mg       0.02         3) ethylbenzene-d10       6.32       116       191828       22.76 mg       0.03         17) naphthalene-d8       13.52       136       519787       11.38 mg       0.02         System Monitoring Compounds       3.96       78       106819       1.33       1.13         4) toluene       4.98       92       492998       9.21       7.92         5) ethylbenzene       6.43       106       332130       2.45       2.10         6) m&p xylene       6.57       106       978127       6.13       5.29         7) o-xylene       7.02       106       335231       2.17       1.84         10) 1-methyl-3-ethylbenzene       8.37       120       64838       0.52       0.45 #         11) 1-methyl-4-ethylbenzene       8.67       120       157941       0.84       0.73         13) 1-methyl-2-ethylbenzene       8.67       120       157941       0.49       0.42         14) 1,2,4-trimethylbenzene       10.00       120       79801       0.47       0.40         14) 1.2,3-trimethylbenzene       10.00       120       77801       0.47       0.40 </th <th>Internal Standards</th> <th colspan="2">R.T. QIon Response</th> <th>Response</th> <th colspan="3">Conc Units Dev(Min)</th>	Internal Standards	R.T. QIon Response		Response	Conc Units Dev(Min)		
System Monitoring Compounds       Wt%       Vol%         Aromatic Compounds       3.96       78       106819       1.33       1.13         4) toluene       4.98       92       492998       9.21       7.92         5) ethylbenzene       6.43       106       321130       2.45       2.10         6) m&p xylene       6.57       106       978127       6.13       5.29         7) o-xylene       7.72       120       9665       0.08       0.07         9) n-propylbenzene       8.37       120       64838       0.52       0.45       #         10) 1-methyl-3-ethylbenzene       8.51       120       237310       1.67       1.45         12) 1,3,5-trimethylbenzene       8.67       120       157941       0.84       0.73         13) 1-methyl-2-ethylbenzene       9.0120       434304       2.49       2.12         15) 1,2,3-trimethylbenzene       10.00       120       79801       0.47       0.40         16) indan       10.35       117       110913       0.30       0.23         18) 1,4-diethylbenzene       10.68       91       36861       0.09       0.07         20) 1,2-diethylbenzene       10.68       91 <td>1) benzene-d6 3) ethylbenzene-d10 17) naphthalene-d8</td> <td>3.94 6.32 13.52</td> <td>84 116 136</td> <td>143289 191828 519787</td> <td>22.76 22.76 11.38</td> <td>mg mg mg</td> <td>0.02 0.03 0.02</td>	1) benzene-d6 3) ethylbenzene-d10 17) naphthalene-d8	3.94 6.32 13.52	84 116 136	143289 191828 519787	22.76 22.76 11.38	mg mg mg	0.02 0.03 0.02
Aromatic CompoundsWt%Vol%2) benzene3.96781068191.331.134) toluene4.98924929989.217.925) ethylbenzene6.431063211302.452.106) m&p xylene6.571069781276.135.297) o-xylene7.081063352312.171.848) Isopropylbenzene7.7212096650.080.079) n-propylbenzene8.37120648380.520.45#10) 1-methyl-3-ethylbenzene8.37120648380.520.45#11) 1-methyl-4-ethylbenzene8.511202373101.671.4512) 1,3,5-trimethylbenzene8.671201579410.840.7313) 1-methyl-2-ethylbenzene9.301204343042.492.1215) 1,2,3-trimethylbenzene10.00120798010.470.4016) indan10.351171109130.300.2318) 1,4-diethylbenzene10.65119172550.730.63#20) 1,2-diethylbenzene12.07134436240.210.1821) 1,2,4,5-tetramethylbenzene12.591282857960.470.4022) 1,2,3,5-tetramethylbenzene15.3614837580.020.0223) naphthalene15.97142710050.180.1524) pentamethylbenzene15.3614837580.02<	System Monitoring Compounds						
	Aromatic Compounds 2) benzene 4) toluene 5) ethylbenzene 6) m&p xylene 7) o-xylene 8) Isopropylbenzene 10) 1-methyl-3-ethylbenzene 11) 1-methyl-4-ethylbenzene 12) 1,3,5-trimethylbenzene 13) 1-methyl-2-ethylbenzene 14) 1,2,4-trimethylbenzene 15) 1,2,3-trimethylbenzene 16) indan 18) 1,4-diethylbenzene 20) 1,2-diethylbenzene 21) 1,2,4,5-tetramethylbenzene 22) 1,2,3,5-tetramethylbenzene 23) naphthalene 24) pentamethylbenzene 25) 2-methylnaphthalene 26) 1-methylnaphthalene 27) C10 benzenes 28) C11benzenes 29) C12 benzenes 30) Alkyl Indans	3.96 4.98 6.43 6.57 7.08 7.72 8.37 8.51 8.67 8.97 9.300 10.005 10.65 10.65 10.65 10.65 12.07 12.15 13.59 15.66 15.97 11.000 13.000 14.64 10.39	$\begin{array}{c} 78\\ 92\\ 106\\ 106\\ 120\\ 120\\ 120\\ 120\\ 120\\ 120\\ 120\\ 120$	106819 492998 321130 978127 335231 9665 64838 64838 237310 157941 67186 434304 79801 110913 197255 36861 6707 43624 59607 285796 3758 176980 71005 238008 117323 27702 299343	Wt% 1.33 9.21 2.45 6.13 2.17 0.08 0.52 0.52 1.67 0.84 0.49 2.49 0.47 0.30 0.73 0.09 0.03 0.21 0.30 0.47 0.30 0.47 0.32 0.42 0.42 0.42 0.42 0.42 0.56 0.56 0.13 0.78	Vol% 1.13 7.92 2.10 5.29 1.84 0.07 0.45 0.45 1.45 0.73 0.42 2.12 0.40 0.23 0.63 0.63 0.07 0.03 0.18 0.25 0.40 0.02 0.35 0.15 1.01 0.46 0.41 0.60	# # # #







Figure 3. Final linear calibration curve for toluene utilizing the high split ratio, electron multiplier voltage programming and change of quant ion.

# References

- 1. V. Giarrocco, Fuel Technology and Management, 57, (March/April 1996).
- 2. Analyzing aromatics in reformulated gasoline by GC/MS, MS Application Brief MS 95-1, Agilent Technologies publication number (23) 5964-0116E.
- 3. Personal communication with Fred Feyerherm, Agilent Technologies
- 4. Personal communication with Michael Mathiesen, Marathon Oil.
- 5. D. E. Mathews and J. M. Hayes, Anal. Chem., 48 1375 (1976).
- 6. Personal communication with Axel Lubeck, Marathon Oil.

Agilent Technologies shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this material.

Information, descriptions and specifications in this publication are subject to change without notice.

Copyright © 1999 Agilent Technologies Company All rights reserved. Reproduction and adaptation is prohibited.

Printed in the U.S.A. December 1999 (23) 5966-2798E