



# Measurement of Chromium in Environmental Waters by Zeeman Corrected Graphite Tube Atomization

## Application Note

Atomic Absorption

### Author

Jonathan H. Moffett

### Introduction

The measurement of trace metal levels in environmental waters is ideally performed by atomic absorption spectroscopy, and the literature on the techniques used is very extensive [1–7]. Chromium, along with other elements such as cobalt, copper, lead, manganese and nickel, requires monitoring in the environment.

The levels encountered are often too low to be measured by direct aspiration into the flame of an atomic absorption spectrometer, thus two methods are commonly employed to overcome this problem. The first is to form a metal complex using various ligands (pentan-2, 4-dione, acac; ammonium 1-pyrrolidinedicarbodithioic acid, APDC) and extracting the complex into an organic solvent (4-methylpentan-2-one, MIBK; 2,5-methylheptan-4-one, DISK) [2]. By careful selection of the relative volumes, the analyte can be preconcentrated before aspiration into a flame. This method is relatively inexpensive but tedious and prone to contamination. Solvent extraction can also be used for furnace work [1–6].

The other method is by direct injection of the water sample into a graphite furnace [1,3–7]. Sample preparation and the possibility of contamination are therefore minimized [4]. The use of a background corrector is strongly recommended so that the effect of background interferences is reduced [7].



**Agilent Technologies**

## Practical

Water samples collected under Environmental Protection Agency recommended conditions were used for method development. The method was then checked by measuring the amount of chromium in certified water samples. The first water sample was supplied by the US Environmental Protection Agency [8]. The other sample was a standard reference material obtained from the National Bureau of Standards [9].

## Experimental

An Agilent SpectrAA-30 Zeeman atomic absorption spectrometer with a programmable sample dispenser and a chromium hollow cathode lamp-were used for this study. An Epson MX-80 printer was used for all printouts and signal graphics traces. A partitioned graphite tube was used in the furnace. The instrument parameters used for the final studies are reproduced in Table 1.


The chromium standard used in this study was prepared from BDH Spectrosol stock solution of 1000 mg Cr/L in 1 M nitric acid and diluted using deionized distilled water to give a final working concentration of 20 µg Cr/L.

Table 1. Instrument Parameters


PROGRAM 7 Cr EFFLUENT WATER					
INSTRUMENT MODE			ABSORBANCE		
CALIBRATION MODE			CONCENTRATION		
MEASUREMENT MODE			PEAK HEIGHT		
LAMP POSITION			1		
LAMP CURRENT (mA)			7		
SLIT WIDTH (nm)			0.2		
SLIT HEIGHT			REDUCED		
WAVELENGTH (nm)			357.9		
SAMPLE INTRODUCTION			SAMPLER AUTOMIXING		
TIME CONSTANT			0.05		
MEASUREMENT TIME (sec)			2.0		
REPLICATES			2		
BACKGROUND CORRECTION			ON		
MAXIMUM ABSORBANCE			2.00		
FURNACE PARAMETERS					
STEP NO.	TEMPERATURE (C)	TIME (sec)	GAS FLOW (L/min)	GAS TYPE	READ COMMAND
1	85	5.0	3.0	NORMAL	NO
2	95	60.0	3.0	NORMAL	NO
3	120	10.0	3.0	NORMAL	NO
4	1100	20.0	3.0	NORMAL	NO
5	1100	20.0	3.0	NORMAL	NO
6	1100	2.0	0.0	NORMAL	NO
7	2500	0.7	0.0	NORMAL	YES
8	2500	2.0	0.0	NORMAL	YES
9	2500	2.0	3.0	NORMAL	NO
SAMPLER PARAMETERS					
VOLUMES (µL)		MODIFIER			
SOLUTION	BLANK				
BLANK	--	15			
STANDARD 1	2	13			
STANDARD 2	5	10			
STANDARD 3	10	5			
SAMPLE	10	5			
RECALIBRATION RATE			0		
RESLOPE RATE			0		
MULTIPLE INJECT NO	HOT INJECT	NO	PRE INJECT	NO	

## Results and Discussions

The collection of the environmental water samples and the preparation of all the standards meant that the chromium ions were present in a dilute nitric acid matrix. This particular matrix stabilizes the analyte in solution and allows the use of an ashing temperature which is high enough to remove all organic material. Signal graphics traces of a standard and a sample are reproduced together in Figure 1. Close examination reveals that there is very little difference in the respective background signals.



(i) standard



(ii) sample

Figure 1. Signal graphics traces.

A study was performed to see if there are any matrix effects. Because all organic matter is apparently destroyed, other metal ions present in the matrix could be a potential source of interferences. A standard additions calibration slope was compared to a normal concentration calibration slope. The data points and slopes are shown in Figure 2. As can be seen, the slopes are virtually identical. The concentration calibration mode was used for convenience in subsequent studies.

The validity of the method was checked by analyzing samples of known concentration. A concentration calibration curve was determined followed by analysis of samples of USEPA sample 4 and NBS SRM 1643b respectively. The autorun report is reproduced in Figure 3. The solutions being dilute aqueous solutions, their density would be the same as pure water to a very good approximation. Hence the numerical value for ng/L would be equal to the value in ng/g. The results obtained are compared in Table 2, and are seen to be in very good agreement.





Figure 2. Comparison of slopes.

OPERATOR		JONATHAN MOFFETT			
DATE		10 SEPTEMBER 1986			
BATCH		EPA#4/NBS			
PROGRAM 7		Cr EFFLUENT WATER			
SAMPLE	CONC ug/L	%RSD	MEAN ABS	READINGS	
BLANK	0.00		0.008	0.008	0.007
STANDARD 1	4.00	3.1	0.118	0.115	0.120
STANDARD 2	10.00	0.2	0.282	0.281	0.282
STANDARD 3	20.00	0.3	0.543	0.544	0.542



USEPA #4	9.75	1.6	0.273	0.278	0.272
NBS SRM1643B	18.65	1.3	0.507	0.512	0.503

Figure 3. Autorun report.

Table 2. Results Comparison

Comparison of experimentally determined results with quoted results

Standard	Quoted (ng/g)	Found (ng/g)
US EPA Sample 4	10.2 ± 1.1	9.75 ± 0.16
NBS SRM 1643b	18.6 ± 0.4	18.65 ± 0.24

## Conclusion

The determination of trace levels of chromium in waters and effluents can be done quickly and accurately and with good precision using a graphite furnace with an effective Zeeman background corrector.

## References

1. Trevor McKenzie, Varian Instruments At Work No AA-8, October 1980.
2. Lucinda Voth, Varian Instruments At Work No AA-17, August 1981.
3. Douglas Shrader, Lucinda Voth, Lawrence Covick, Varian Instruments At Work No AA-30, June 1983.
4. Douglas Shrader, Lucinda Voth, Lawrence Covick, Varian Instruments At Work No AA-31, June 1983.
5. Douglas Shrader, Lucinda Voth, Lawrence Covick, Varian Instruments At Work No AA-33, September 1983.
6. Lucinda Voth, Varian Instruments At Work No AA-34, September 1983.
7. Lucinda Voth, Varian Instruments At Work No AA-35, September 1983.
8. Environmental Monitoring and Support Laboratory, US Environmental Protection Agency., Cincinnati, USA.
9. National Bureau of Standards, US Department of Commerce, Washington DC, USA.

## For More Information

For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem)

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 1987  
Printed in the USA  
November 1, 2010  
AA069



**Agilent Technologies**