



Agilent MSD
Productivity ChemStation
G1701 & G1710

**EasyID Resolution, Tailing, and
Degradation**



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G1701-90071

Edition

First Edition

May 14, 2012

Printed in USA

Agilent Technologies, Inc.
5301 Stevens Creek Blvd.
Santa Clara, CA 95051

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data -Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Introduction


This document describes the calculations used by EasyID for finding Peak Resolution and Tailing in the Agilent MSD Productivity ChemStation. EasyID is a tool for updating the retention times in a quantitation database on a *compound-by-compound* basis.

With a quantitation method and data file loaded, EasyID is started by selecting EasyID from the View menu (Tools menu if Environmental Data Analysis mode). This displays the Quick EasyID window. In the Quick EasyID window, double click on a compound entry to display a chromatogram of that compound in window #6.

This resolution calculation has its origin in the EPA method 8280B which uses a resolution suitability requirement to determine the separation of critical isomer pairs. This requirement specifies that the % valley as calculated here must be at or below 25%.

Evaluate Resolution

To display the resolution between adjacent peaks in the chromatogram, select Evaluate Resolution from the ChromEval menu. The software asks for you to select the two peaks used for calculating resolution by manually integrating them. Integrate these peaks manually by drawing a baseline (right-click and drag) under each peak. The software calculates the resolution and displays it as a valley percentage next to the peak.




$$\% \text{ Valley} = \frac{P_{hv}}{P_{ha}} * 100$$

Where:

P_{hv} = baseline height of the valley between the peaks;
 P_{ha} = baseline peak height to the apex.

Evaluate Tailing

In EasyID, Tailing uses the calculation for the Asymmetry factor (see definition below) and not the more common tailing calculation as used in a USP System Suitability method. It is a measure of peak symmetry. For a perfectly symmetric peak, the peak tailing factor = 1. In EasyID, to display the Tailing for a peak, select Evaluate Tailing from the ChromEval menu. The software asks for you to select the peak by manually integrating it. Integrate this peak manually by drawing a baseline under the peak. The software calculates the tailing and displays it next to the peak. The values used for calculating tailing are displayed in the status line at the bottom of the screen.



$$Tailing = \frac{BC}{AB}$$

Where:

B = Retention time at the peak apex.

C = RT of peak end at 10% of peak height;

A = RT of peak start at 10% of peak height;

Asymmetry factor

The asymmetry factor is a measure of peak tailing. It is calculated when a user selects Evaluate Tailing in EasyID. It is defined as the distance from the center line of the peak to the back slope divided by the distance from the center line of the peak to the front slope, with all measurements made at 10% of the maximum peak height. The asymmetry factor of a peak will typically be similar to the tailing factor for the same peak, but the two values cannot be directly converted.

Evaluate Degradation

To display the Degradation between a primary peak and a peak containing a degradation compound of the same primary peak compound, select Evaluate Degradation from the ChromEval menu. The software asks for you to select the two peaks used for calculating degradation by manually integrating them. Integrate these peaks manually by drawing a baseline (right-click and drag) under each peak, starting with the Primary peak. The software calculates Degradation and displays it as a Degradation percentage next to the primary peak.

Degradation is calculated as a percent area ratio of the degradation peak to the primary peak.

$$\text{Degradation} = \frac{D}{P} \times 100$$

Where: D = Degradation peak area
 P = Primary peak area

