New Hybridization Apparatus For In Situ Synthesized Microarrays Which Improves Experimental Performance. Petula N. D'Andrade, Paul K. Wolber, Jacqueline Tso, Allen Thompson, Arthur Schleifer, Jenny Xiao, and Jeffrey M. McMillan Agilent Technologies, 3500 Deer Creek Rd. Palo Alto, 94304 CA

ABSTRACT The standard 1" x 3" microscope slide is the most widely used microarray format. However, sample hybridization to the microarray slide has traditionally been a poorly optimized step using home-made hybridization apparatus that lead to lowered microarray performance. We have developed and characterized a new apparatus and method which improves the following properties of A2 A1 our existing system: Quality of hybridization, EXPERIMENTAL SET-UP Consistency of hybridization, and Ease of use. All experiments were performed with Agilent's Yeast Oligo microarray which contain Hybridization quality has been improved by providing a 60-mer yeast probes. The total number of features on the microarray is 10,807 hybridization chamber made from inert materials that representing 6,256 yeast ORFs. enable easy hybridization buffer mixing. Consistency has been improved via material quality and high-quality Labeled targets were prepared using the Agilent Linear Amplification Kit with S. manufacturing. Finally, ease-of-use has been improved cerevisiae poly A+ RNA. A final concentration of 0.5ug/ml of cRNA for each dve in by protocols that employ standard lab equipment in 200ul of hybridization buffer. A statistically significant number of microarrays for each conjunction with a straightforward hybridization chamber condition was hybridized according to Agilent's standard procedure. assembly. Homologous comparisons (A1) were performed using starting material from Clontech The apparatus has been validated for use with Agilent's (Cat. No. 6999-1). Differential comparisons (A2) were carried out in fluor reversed *in-situ* synthesized oligonucleotide microarrays by pairs using cRNA obtained from yeast grown in synthetic complete (SC) medium comparison to Agilent's current standard hybridization versus yeast grown in sporulation (Spo) medium. method. It exhibited equivalent or superior experimental performance in these studies. Arrays were randomly assigned to the different conditions, loaded in random order, washed in random order and scanned in random order, to reduce any systematic biases. Microarrays were scanned using the Agilent dual laser DNA microarray scanner (Cat. Number G2565AA) and data was extracted using Agilent's Feature Extraction software (Cat. Number G2566AA). The data was analyzed using Rosetta Resolver®, and Microsoft® Access in conjunction with Microsoft Excel. Green Background Subtracted Intensity Red Background Subtracted Intensity New Technique vs. Legacy Inter-Array CV% (6 Arrays) New Technique vs. Legacy ·Green_Legacy - Red_New ·Red_Legacy New = 1.3403 *Legacy - 12.53 New = 1.3763* Legacy - 3.33 R² = 0.9939 50000 25000 30000 Legacy Legacy Figure D1 Figure D2 40% 50% 60% 70% 80% 90% 100% CV% (Std. Dev./Avg) INCREASED SIGNAL Figure E1 In a 12 microarray self vs. self experiment, a background subtracted signal **INTER-MICROARRAY REPRODUCIBILITY** comparison was perform to identify if the new hybridization apparatus offered any improvement in sensitivity (figures D1 and D2). The following table summarizes the findings. Overall, the signal increased an average of 36%. In the same 12 microarray self vs. self experiment in which 6 microarrays were hybridized using the legacy hybridization tool, and the other 6 were Local Background Std. Dev hybridized using the new tool, the CV%s (SD/Avg.) were compared on a Signal Increase Factor Channel feature-by-feature basis. Each feature on the microarray was compared Red_Legacy 1.38 Cy3 to the corresponding feature on the 5 other microarrays in the group. Any Cy5 1.34 Green_New saturated feature, or any feature that was identified by the Feature Extraction software program as a 'Feature Non-Uniformity Outlier' was removed from the calculation. The CV% of the Dye-Normalized Signal was calculated and plotted as a histogram (frequency of occurrence) plot in figure E1, and versus Dye-Normalized Signal in Figure E2.

The plots show a clear difference between the legacy and the new apparatus indicating a 2 fold reduction in the mode and median CV%s.

	Mode CV%	Median CV%
Legacy	15%	14.5%
New	7%	7.4%



UNCHANGED LOCAL BACKGROUND

Figure D3

The standard deviation of all the non-flagged local backgrounds were plotted in a histogram frequency distribution to detect any shift in the background std dev. that might indicate a change in the background noise (Figure D3). The histogram curves nearly overlap indicating no significant change in the local background.

5989-0354EN November 1, 2003

Rosetta Resolver is a U.S. registered trademark of Rosetta Inpharmatics. Microsoft is a U.S. registered trademark of Microsoft Corporation. Information, descriptions and specifications in this publication are subject to change without notice. © Agilent Technologies, Inc. 2003







COMPARISON OF LOG RATIOS

The log ratios for the differentially expressed ORFs, using eight microarrays for each hybridization tool, were compared using the Rosetta Resolver system (Figure F1). At the 99% confidence level, 5,452 ORFs for the new hybridization tool vs. 5,215 ORFs for the legacy tool were detected. Of these, 5,136 ORFs were common between the two tools with zero anticorrelated suggesting no change in the log ratio for these ORFs. This results in 316 unique signature ORFs (green) vs. 79 unique ORFs (red) for the new and legacy hybridization tools respectively.

process more microarrays in the same period of time.

An improvement in hybridization reproducibility is realized through lower microarray-to-microarray variability as demonstrated by the self-self experiment in which the signal CV% (mode method) were reduced by more than half from 14-16% to 6-7%.

An improvement in hybridization sensitivity is achieved by using this new tool. An increase between 34-38% in background subtracted signal is observed in conjunction with no change in the background noise as measured by the local background standard deviations. As a result, the signal-to-noise, the most common indicator of sensitivity reports a 34-38% increase in sensitivity.

Finally, the differential expression study shows that while most of the ORFs in the experiment were equally differentially expressed in both data sets, more unique signature ORFs were detected while using the new hybridization apparatus.



Agilent Technologies