**Agilent P3211 and
P3311
pH Combination
Electrodes
pH 电极**

**Operating Guide
用户手册**

Agilent Technologies

Overview

The P3211 and 3311 Refillable pH Combination Electrodes measure pH in aqueous solutions. The refillable reference electrode helps to prolong the service life of the electrode. They are used with the 3200P pH Meter or similar meters.

P3211 Consists of pH-sensing electrode and reference electrode.

P3311 Consists of pH-sensing electrode and reference electrode, and temperature-sensing probe.

WARNING

Use this electrode according to the operating manual to avoid personal injury.

WARNING

The electrode solution can cause chemical burns or illness if it is taken orally or contacted by human skin. Use protective clothing or gloves to avoid contact. In case of contact, rinse contacted area with tap water or deionized water thoroughly.

CAUTION

The electrode body material is glass. Handle with care to avoid damage to the instrument.

Table 1 P3211 and P3311 pH Combination Electrode specifications

Specification	P3211	P3311
pH range	0 to 14	0 to 14
Zero potential	7.00 ± 0.5	7.00 ± 0.5
Temperature range	0 to 80 °C	0 to 60 °C
Temperature accuracy		± 0.5 °C
Probe impedance	< 300 Mohm	< 300 Mohm
Reference type	Ag/AgCl	Ag/AgCl
Body material	Glass	Glass
Liquid junction material	Ceramics	Ceramics
Reference solution	3 mol/L KCl containing AgCl	3 mol/L KCl containing AgCl
Cable connector	BNC	BNC + miniDIN
Probe diameter	12 mm	12 mm
Probe length	120 mm	120 mm
Cable length	1000 mm	1000 mm


Figure 1 P3211 Combination Electrode assembly


Figure 2 P3311 Combination Electrode assembly

Operation

Preparing the electrode

- 1 Remove the electrode from the storage bottle. Keep the bottle and cap for future use.
- 2 Rinse the measuring tip with distilled or deionized water. With the measuring tip downward, swing the electrode several times to eliminate air bubbles near the sensitive bulb.
- 3 Unplug the filling hole.

Measurement

- 1 Connect the electrode to the meter.
- 2 *P3211 only*. Enter the temperature of the calibrating solution.
- 3 Place the measuring tip in the calibrating solutions, beginning with the most dilute. Calibrate as described in the meter manual.
- 4 After calibration, place the electrode in the sample solutions. When the readings become stable, record the meter readings.
- 5 After measurement, rinse the measuring tip with distilled or deionized water. Place the electrode in the storage bottle and finger tighten the cap. To protect the measuring tip, keep 5 to 10 mm between the tip and the bottom of the bottle.

Operating hints

- The main material of the measuring tip is glass. Ensure that the sample solution will not damage the measuring tip before measurement. Do not immerse the measuring tip in any material containing sand or solid.
- During measurement, keep the level of the reference filling solution higher than the sample solution.
- Keep the calibration solution and sample solution flowing at a uniform speed around the electrode or the accuracy may be severely impacted.

- It is normal for a small amount of reference filling solution to leak through the reference junction during electrode use. The reference filling solution contains silver ions and is not suitable for samples that contain protein. The silver ion will precipitate near the liquid junction and damage the electrode.
- White powder or crystals that collect on the electrode during measurement or transportation do not affect the performance of the electrode. Rinse the residue with tap or deionized water.
- Do not soak the electrode in sample solution for a long time. After measurement, rinse the electrode carefully.
- The reading of a solution is affected by its temperature. The pH reading done with a sample solution at 10 °C will be different than a reading with the sample temperature to 25 °C.
- For acidic sample solutions, use pH 4.00 and pH 7.00 standards for calibration. For alkaline sample solutions, use pH 7.00 and pH 10.01 standards.
- Do not apply force to the electrode cap, the cable, or the cable interface.
- Keep the cable interface dry.

Maintenance

Exchange of the Reference filling solution

- 1 Siphon the electrode reference filling solution and add fresh solution (5190-0545) through the filling hole until the level is 5 mm below the filling hole.
- 2 Repeat this procedure several times.

Cleaning of inorganics

Soak the measuring tip in 0.1 mol/L HCl or EDTA solution for 15 minutes.

Cleaning of organics

Soak the measuring tip in absolute alcohol or organic solvent that can dissolve the material for 15 minutes.

Cleaning of grease

Soak the measuring tip in warm weakly alkaline detergent for 15 minutes.

Cleaning of protein precipitation

Soak the measuring tip in 0.1 mol/L HCl solution that contains 1% pepsin for 15 minutes.

Electrode restoration

WARNING

HF is lethal. Refer to the HF material safety data sheet before using it.

- 1 Soak the sensitive membrane in 4% HF solution for 3 to 5 seconds.
- 2 Rinse in 0.1 mol/L HCl for 10 seconds.
- 3 Rinse with distilled or deionized water.
- 4 After cleaning the electrode, exchange the reference filling solution and soak the measuring tip in reference filling solution for two hours.

Troubleshooting

Meters

Refer to the meter operating manual. Check all related parts, such as the electrode, calibration solution and samples.

Electrode

- 1 Connect the electrode to the meter. Set the meter to display mV.
- 2 Place the electrode in pH 4.01 buffer.
- 3 Record the mV reading in pH 4.01 buffer when the reading becomes stable.
- 4 Rinse the electrode with distilled or deionized water and place the electrode in pH 10.01 buffer.

- 5 Record the mV reading in pH 10.01 buffer when the reading becomes stable.
- 6 If the difference between these two mV readings is in the range of 326 to 366 mV, the electrode is performing properly. If not, clean the electrode. Then repeat the measurements.
- 7 If the readings are not stable, exchange the reference solution. Repeat this action a few times and repeat the measurements.

Accuracy

- 1 Choose any two calibration solutions to calibrate the electrode.
- 2 Use the calibrated electrode to measure the pH of another calibration solution.
- 3 Compare the theoretical and measured values. Check if the measuring accuracy of the electrode is satisfactory.

Sample

- The electrode is only suitable for measuring the pH value of conventional solutions. When measuring the pH of solutions with low conductivity, high viscosity, or nonaqueous reagents, the response time of the electrode may be prolonged or the measurement accuracy may be affected.
- Use effective calibration solution. Make sure the distilled or deionized water used to prepare the solution meets all requirements. Make sure the calibration solution is not contaminated or beyond its shelf life.

For any other problems during electrode use, contact your Agilent Technologies customer service representative.

Storage

Place the electrode and storage bottle in the packing box and store it at ambient temperature in a dry location.

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© Agilent Technologies, Inc. 2012