Agilent P3214
pH Combined
Electrode
pH 复合电极

Operating Guide
用户手册

Agilent Technologies

Overview

The P3214 pH Combined Electrode can measure the pH of semi-solid substances. It is used with the 3200P pH meter or similar meters.

WARNING

Use this electrode according to the operating manual to avoid personal injury.

WARNING

The electrode solution can cause chemical burns or illness if it is taken orally or contacted by human skin. Use protective clothing or gloves to avoid contact. In case of contact, rinse contacted area with tap water or deionized water thoroughly.

CAUTION

The electrode body material is glass. Handle with care to avoid damage to the instrument.

Specifications

Table 1 P3214 pH Combined Electrode specifications

Specification	Value
pH range	0 to 14
Zero potential pH	7.00 ± 0.5
Temperature range	0 to 60 °C
Internal impedance	$< 500 \text{ M}\Omega$
Reference type	Ag/AgCl
Reference filling solution	3 mol/L KCl containing AgCl
Liquid junction material	Fiber
Body material	ABS
Measuring terminal probe diameter	14 mm (including stainless steel sheath)
Measuring terminal probe length	55 mm
Cable interface	BNC
Cable length	1000 mm


Figure 1 P3214 Combined Electrode assembly

Operation

Preparing the electrode

- 1 Remove the electrode from the storage bottle. Store the bottle with its mouth facing upward for future use.
- 2 Rinse the measuring tip with distilled water or deionized water. Hold the measuring tip downwards and swing the electrode several times.
- 3 Unscrew the fixed ring, install the protective stainless puncture on the electrode body, and attach the fixed ring until finger tight.

Calibration

- 1 Connect the electrode to the meter. Input the temperature value of calibration solution according to procedures listed in the meter operating manual.
- 2 Soak the measuring tip in calibration solutions in sequence. Calibrate the electrode according to the procedures listed in the meter operating manual.

Measurement

After calibration, slowly insert the measuring tip into the sample. When the reading becomes stable, record the value from the meter.

Operating hints

- The main materials of the measuring tip are stainless steel and ABS plastic. Ensure the sample solution will not damage the measuring tip before measurement.
- Do not test samples with hard, large, or fixed solid particles. This may damage the sensitive glass spear-tip from the side.
- Do not leave the electrode in the sample for a long time. After measurement, rinse the electrode carefully.

- The pH value of a solution is affected by its temperature. Maintain the temperature of the solution at ambient to ensure accuracy of measurement.
- During electrode calibration, use solutions of pH 4 and pH 7 for acidic samples and solutions at pH 7.0 and pH 10.01 for alkaline samples.
- White powder or crystals that collect on the electrode during measurement or transportation do not affect the performance of the electrode. Rinse the residue with tap or deionized water.
- Do not apply force onto the electrode cap, cable, or cable interface.
- Keep the cable interface dry.

Maintenance

Cleaning of inorganics

Soak the measuring tip in 0.1 mol/L HCL or EDTA solution for 15 minutes.

Cleaning of organics

Soak the measuring tip in absolute ethyl alcohol, or other solvent that can dissolve organics for 15 minutes.

Cleaning of grease

Soak the measuring tip in warm weakly alkaline detergent for 15 minutes.

Cleaning of protein precipitation

Soak the measuring tip in 0.1 mol/L HCL solution that contains 1% pepsin for 15 minutes.

Regeneration of sensitive glass membrane

WARNING

HF is lethal. Refer to the HF material safety data sheet before using it.

- 1 Soak the sensitive glass membrane in 4% HF solution for 3–5 seconds.
- 2 Rinse with 0.1 mol/L HCL several times.
- 3 Rinse with distilled water or deionized water.
- 4 After cleaning the electrode, exchange the reference filling solution and soak the measuring tip in reference filling solution for 2 hours.

Troubleshooting

Meters

Refer to the meter operating manual. Check all related parts, such as the electrode, calibration solution, and samples.

Electrode

- 1 Connect the electrode to the meter. Set the meter to display mV.
- 2 Soak the measuring tip in calibration solutions of pH 4.00 and pH 10.01. When the readings become stable, record mV value respectively.
- 3 If the absolute value of the difference between these two mV readings is 326–370 mV in total, it indicates the electrode has a good performance. If not, perform maintenance on the electrode.

Calibration Solution

- Use the correct calibration solution.
- The distilled water or DI water used to prepare calibration solution should meet with corresponding requirements.
- The calibration solution should not be contaminated, degraded, or expired.

Accuracy

- 1 Choose any two calibration solutions to calibrate the electrode.
- 2 Use the calibrated electrode to measure the pH value in another calibration solution.
- 3 Compare the theoretical value of that calibration solution with the measured value.
- 4 Check if the measuring accuracy of the electrode satisfies the user's need.

For any other problems during electrode use, contact your Agilent Technologies customer service representative.

Storage

- 1 Unscrew the fixed ring.
- 2 Remove the stainless protective puncture.
- 3 Brush the measuring tip gently with a fine brush.
- 4 Rinse the fixed ring and protective stainless puncture.
- 5 Install the storage bottle cap and storage bottle onto electrode body in sequence. To protect the measuring tip, keep 5–10 mm between the bottom of the storage bottle and the electrode measuring tip when installing the electrode. Screw the electrode storage bottle cap onto the storage bottle until finger tight.
- 6 Place the electrode in its storage box and store it at ambient temperature and dry conditions.

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© Agilent Technologies, Inc. 2012