

PLOT Columns

Installation and Conditioning



Column installation and conditioning

Thank you for purchasing this Agilent Technologies PLOT column.

Before proceeding, please read the following information regarding the proper installation and conditioning of your Agilent PLOT column. These steps are critical to proper performance and lifetime of the column.

PLOT columns are highly retentive. Trace amounts of water in your carrier gas stream can adversely affect chromatography. The use of a moisture trap on your carrier gas line is highly recommended.

1. Install the column according to the instructions enclosed in the column box.
2. After installing the column into the inlet, set the carrier gas flow rate. For all columns except the GS-GasPro and GS-CarbonPLOT, slowly increase the carrier gas

pressure at a rate of 2 to 3 psi/min until reaching the appropriate flow rate:

| Inner dia. | Standard flow rate |
|------------|--------------------|
| 0.32 mm | 2 to 4 mL/min |
| 0.53 mm | 6 to 9 mL/min |

Slow pressure ramping is not required for GS-GasPro and GS-CarbonPLOT – these columns can be brought to pressure rapidly due to their unique particle bonding technology.

3. Let the carrier gas purge for 3 to 5 minutes before connecting the column to the detector.
4. Condition the column per the guidelines in Table 1.

Do not exceed the maximum operating temperature of the column when setting injector and detector temperatures.

Table 1. Conditioning Guidelines

| Column | Conditioning temperature | Conditioning time |
|--|--------------------------|-------------------|
| HP-PLOT Molesieve | 300-350 °C | 3-4 hours |
| HP-PLOT Al ₂ O ₃ "KCI" | 200 °C | 8 hours |
| HP-PLOT Al ₂ O ₃ "S" | 200 °C | 8 hours |
| HP-PLOT Al ₂ O ₃ "M" | 200 °C | 8 hours |
| GS-Alumina | 200 °C | 8 hours |
| HP-PLOT Q | 270 °C | 3-6 hours |
| GS-Q | 250 °C | 8-10 hours |
| HP-PLOT U | 190 °C | 3-6 hours |
| GS-GasPro | 260 °C | 3-6 hours |
| GS-CarbonPLOT | 300-350 °C | 3-6 hours |

Column maintenance

PLOT columns can become contaminated by water, polar compounds, and/or hydrocarbons in sample or carrier gas streams. When contamination is suspected, bake the column at its conditioning temperature for several hours, or until the baseline is stable.

Contaminant compounds vary depending on the retention and selectivity of the PLOT stationary phase. Consult Table 2, PLOT Column Selection Guide, for recommended applications and suggestions on possible contaminant compounds.

To maintain the installed column between uses, hold the oven temperature at 100 °C to 150 °C with continuous carrier gas flow. This helps to avoid problems associated with accumulation of water, CO₂, or other carrier gas impurities.

Particle traps for use with PLOT columns

Though highly stabilized, it is impossible to guarantee that no particles will dislodge from the column wall with the following phases: HP-PLOT Molesieve, HP-PLOT Al₂O₃ (any phase), GS-Alumina, HP-PLOT Q, HP-PLOT U, and GS-Q. When used in valve-switching applications, consider using a particle trap to prevent scarring of the valve rotors. Also consider a particle trap when interfacing one of the listed PLOT phases to a mass spectrometer.

Particle traps are short pieces of capillary column with a moderate film thickness of a siloxane stationary phase. Any particles eluting from the

PLOT column are trapped on the siloxane phase layer. Agilent sells particle traps for your convenience. Please see Table 3 for dimensions and part numbers.

Table 3. Particle Traps

| Inner dia. (mm) | Length (m) | Part no. |
|-----------------|------------|-----------|
| 0.32 | 2.5 | 5181-3351 |
| 0.53 | 2.5 | 5181-3352 |

Agilent PLOT columns

Agilent offers a comprehensive line of PLOT columns for analysis of fixed gases, low molecular weight hydrocarbon isomers, volatile polar compounds, and reactive analytes such as sulfur gases, amines, and hydrides. Many of our PLOT phases are offered in dimensions from 0.25 to 0.53 mm ID, allowing for easy column selection for various detector and system requirements. For GC/MS systems, we offer several small diameter columns with truly bonded and immobilized stationary phases, eliminating potential detector fouling due to particle elution.

Refer to Table 2 for information on PLOT column selection. This reference table includes all of the Agilent PLOT column phases, common application areas in which they are used, and potential contaminants. Many PLOT column applications require the use of switching valves for a variety of reasons, including the avoidance of contaminant introduction to a PLOT phase.

Agilent Technical Support

If you require more detailed technical information on column selection, or on column switching valves, please contact Agilent Technical Support. Also, for a complete listing of J&W Scientific columns and PerfectFit supplies, as well as ordering information, please visit our Web site at:

www.agilent.com/chem

Or, you can contact your authorized Agilent distributor or local Agilent representative.

Table 2. PLOT Column Selection Guide

| Column | Stationary phase | Typical applications | Typical contaminants |
|--|--|---|--|
| HP-PLOT Molesieve | 5Å zeolite molecular sieve | Permanent and noble gases. Thick and thin films available. Thick film column will resolve argon and oxygen at 35 °C. | Water; carbon dioxide; hydrocarbons larger than C ₃ . |
| HP-PLOT Al ₂ O ₃ "KCI" | Aluminum oxide deactivated with KCl | Least "polar" Al ₂ O ₃ phase. Lowest retention of olefins relative to comparable paraffin. C ₁ to C ₈ hydrocarbon isomers. Column of choice for accurate quantitation of dienes, especially propadiene and butadiene from ethylene and propylene streams. | Water; carbon dioxide; large hydrocarbons and substituted aromatics; oxygen-containing compounds such as alcohols, ethers, ketones; sulfur-containing compounds. |
| HP-PLOT Al ₂ O ₃ "S" | Aluminum oxide deactivated with sodium sulfate | Excellent midpolarity, general use Al ₂ O ₃ column. C ₁ to C ₈ hydrocarbon isomers. Best for resolving acetylene from butane and propylene from isobutane. | Water; carbon dioxide; large hydrocarbons and substituted aromatics; oxygen-containing compounds such as alcohols, ethers, ketones; sulfur-containing compounds. |
| HP-PLOT Al ₂ O ₃ "M" | Aluminum oxide with proprietary deactivation | Most "polar" of the Al ₂ O ₃ columns. Highest retention of olefins relative to comparable paraffin. Excellent general use Al ₂ O ₃ column. C ₁ to C ₈ hydrocarbon isomers. Good for resolving cyclopropane from propylene. | Water; carbon dioxide; large hydrocarbons and substituted aromatics; oxygen-containing compounds such as alcohols, ethers, ketones; sulfur-containing compounds. |
| GS-Alumina | Aluminum oxide with proprietary deactivation | Most "polar" of the Al ₂ O ₃ columns. Highest retention of olefins relative to comparable paraffin. Excellent general use Al ₂ O ₃ column. C ₁ to C ₈ hydrocarbon isomers. Best for resolving cyclopropane from propylene. Good stability and recovery from water saturation. Selectivity slightly different than "M" column. | Water; carbon dioxide; large hydrocarbons and substituted aromatics; oxygen-containing compounds such as alcohols, ethers, ketones; sulfur-containing compounds. |
| HP-PLOT Q | Polystyrene-divinylbenzene | C ₁ to C ₃ isomers, alkanes to C ₁₂ , CO ₂ , methane, air/CO, water, oxygenated compounds, sulfur compounds, solvents. | Large hydrocarbons and substituted aromatics. |
| GS-Q | Polystyrene-divinylbenzene | C ₁ to C ₃ isomers, alkanes to C ₁₀ , CO ₂ , methane, air/CO. Slightly different selectivity than HP-PLOT Q. Not recommended for quantitation of polar compounds. | Large hydrocarbons and substituted aromatics. |
| HP-PLOT U | Divinylbenzene/ethylene glycol dimethacrylate | More polar than HP-PLOT Q and GS-Q. C ₁ to C ₇ hydrocarbons, CO ₂ , methane, air/CO, water, oxygenates, amines, solvents, alcohols, ketones, aldehydes. | Large hydrocarbons and substituted aromatics. |
| GS-GasPro | Proprietary, bonded silica-based | C ₁ to C ₁₂ hydrocarbons, CO ₂ , trace-level sulfurs, hydride gases, inorganic gases, halocarbons, SF ₆ , oxygen/nitrogen separation at -80 °C. | Large hydrocarbons and substituted aromatics. |
| GS-CarbonPLOT | Bonded, monolithic carbon-layer | C ₁ to C ₅ hydrocarbons, CO ₂ , air/CO, trace acetylene in ethylene, methane. | Large hydrocarbons and substituted aromatics. |



J&W GC Column Performance Summary

| Performance Results |
|---------------------|
| |

| Compound Identification | Retention Time (t _r) | Partition Ratio (k) | Peak Width (w _{1/2}) |
|-------------------------|----------------------------------|---------------------|--------------------------------|
| | | | |

| Test Conditions |
|-----------------|
| |

