


# Process control of polystyrenes

## Application

**Heinz Goetz and  
Angelika Gratzfeld-  
Huesgen**

Polystyrenes are widely used for the production of packaging materials, household goods, cases of electronic equipment, toys, and insulation materials. Polymerization can be performed with pure styrene or by copolymerization with butadiene, acrylonitrile, rubber and methylstyrene. The properties of the product strongly depends on the monomers used, the molecular weights and the molecular weight distribution.

Figure 1 shows an overlay of 3 chromatograms of a technical polystyrene – the original granulate, one after 1st injection moulding and one after second injection moulding. After the first injection moulding there is almost no change in the chromatogram and therefore the molecular



**Figure 1**  
Overlay of three chromatograms of a technical polystyrene

## Conditions

### Sample preparation

Sample was dissolved in THF.  
Polystyrene EasyCal Vial standards (Agilent p/n 79911-60500 and 79911-60501) were used for narrow standard calibration.

### Column

3 x PLgel mixed B, 7.5 x 300 mm, 10 µm (Agilent p/n 79911GP-MXB) in series

### Mobile phase

Tetrahydrofuran

### Flow rate

1.0 mL/min

### Column compartment temperature

20° C

### Injection volume

10 µL

### Detector

Refractive index detector, alternatively VWD, 254 nm



**Agilent Technologies**

Innovating the HP Way



weight distribution. After grinding the chips, and injection moulding a second time there is a significant change which will have an effect on the properties. The visual information is supported by the number average molecular weight,  $M_n$ , as calculated by the ChemStation data analysis software:

$M_n$  (original granulate): 59000

$M_n$  (after second process): 55000

To characterize such small differences in polymers reliably a GPC-SEC instrument with excellent precision, such as the Agilent 1100 Series GPC-SEC system, is required. Further information on the precision of this system is given in application brief "Precision in GPC-SEC analysis", Agilent publication number 5988-0109EN.

## HPLC performance

RSD of  $M_w$  < 1%

RSD of  $M_n$  < 2%

## Equipment

### Agilent 1100 Series GPC-SEC system

consisting of

- vacuum degasser for efficient degassing of the mobile phase
- isocratic pump with large solvent cabinet
- autosampler with single valve design
- thermostatted column compartment for precise column temperatures
- refractive index detector with automatic recycle valve
- ChemStation Plus with GPC-SEC data analysis software

Heinz Goetz and Angelika Gratzfeld-Huesgen are application chemists at Agilent Technologies, Waldbronn, Germany

For more information on our products and services, visit our worldwide website at

<http://www.agilent.com/chem>

© Copyright 2000 Agilent Technologies

Released 09/2000

Publication Number 5988-0112EN



**Agilent Technologies**

Innovating the HP Way