

Agilent ProTune Probe Tuning Accessory

Data Sheet



Introduction

Agilent ProTune allows laboratories to automatically adjust an NMR probe's tune and match to achieve optimal results with every sample. The Agilent ProTune driver modules offer an automated solution to the time consuming and often challenging requirement of manually tuning probes when changing solvents, sample temperature, or nuclei. In conjunction with an automatable NMR probe, the ProTune driver modules yield a completely integrated system that offers significant time savings and ease of use over the traditional manual approach for adjusting the probe's tune and match.



Agilent Technologies

Superior Results on Every Sample

The ProTune accessory automatically adjusts the probe's tune and match to achieve the best possible results in the shortest amount of time, on every sample. Differences in the dielectric properties between samples (especially when different solvents are used) as well as changes in temperature affect the tune and match of a probe. Without optimization, pulse widths and sensitivity can be drastically degraded, compromising the quality of results. The example in Figure 1 shows the effect of a 60 °C temperature change on a sample without (B) and with (C) automatic tuning. Experiments requiring accurate 90 ° and 180 ° pulses would fail, while experiments using a 90 ° pulse would have poor sensitivity.

Automatic

When using ProTune, it is easy to customize automatic tune and match to the needs of any laboratory. VnmrJ offers the option of specifying how and when the system tunes a sample, giving the user complete control. By setting the preference for the software to always use ProTune, automatic tune and match will occur unless the user decides to override that option by clicking a single checkbox (Figure 2).

Integrated

The ProTune driver modules connect to an automatable probe's tune and match sticks with high performance flexible drive shafts, and they connect to the console through an Ethernet port. The driver modules are controlled by software embedded in VnmrJ. The tune reflection can be observed (Figure 3) when the system is tuning in order to confirm the speed and quality of the tune optimization. It is also possible to manually adjust tune and match of the probe using the knobs on the front panel.



Figure 1. ¹H 90 ° pulse width (pw) calibration as a function of sample temperature: tuned and calibrated at 30 °C (A); tuned at 30 °C, but calibrated at 90 °C (B); and tuned and calibrated at 90 °C (C).

| Sample Nam SampleDir: | e: test_sample (test_sample_20120124_01) | Insert Eject |
|--------------------------|---|--|
| Notebook: | 101111 Page: 10 | After EXP 🖌 Autoplot |
| Concentration | n: 0 Eaddr imr1@agilent.com | plotter |
| Solvent | Acetone | |
| | DMSO CDCI3 D2O Other | Before First Evn |
| Comment | | Lock2 Vos (alockson) |
| | | Cooker 165 (alock whit) |
| | | |
| | | Select shimmap: 2H gradient (read probe file at runtime) |

Figure 2. Automatic tuning and matching can be selected or unselected with a single click next to the **Tune** check box.



Figure 3. This ProTune window can be shown during probe optimization to view tuning in the automated mode or when running the system manually.

Automation Friendly

ProTune delivers increased throughput and improved data quality, allowing the NMR user to focus on results rather than the details of data collection. By having optimal tuning for each sample, sensitivity is improved and access to multiple nuclei is possible in a hands-off, high-throughput manner (Figure 4).

Modular Design

The ProTune driver modules are available in various configurations for automatic tune and match of 1–3 probe channels. Each system has a master driver module and may have a combination of dual and single driver modules, depending on the number of tune and match knobs on the probe. The dual and single driver modules are optional, allowing the greatest flexibility in customizing the tuning configuration of a probe. When only the high band is required for automatic tuning, just the master drive module is necessary.

Versatile

ProTune is compatible with a large number of probe styles and can be interchanged between different Agilent probes, providing ultimate flexibility and maximum productivity for any laboratory.

OneNMR probes

The OneNMR probe uses a hybrid approach to deliver the performance advantages of both the classic carbon probe (direct detect probe) and the highly sensitive proton probe (indirect detect probe). It is simultaneously optimized for both high and low band frequencies in a single unit.



Figure 4. A series of data collected automatically with ProTune, using a 500 MHz AutoX Dual Broadband probe. The sample was diethyl(2,2,2-trifluoro-1-hydroxyethyl)-phosphonate in CDCI_a.

AutoX Dual Broadband probes

These probes are optimized for single and double resonance experiments requiring high performance in both high frequency (¹H–¹⁹F) and low frequency (¹⁵N–³¹P) observe and decouple experiments.

AutoX Indirect Detection probes

These probes are optimized for single and double resonance experiments requiring high ¹H sensitivity and overall performance with pulsing and decoupling capability at low frequency ($^{15}N-^{31}P$).

Automatable Triple Resonance probes

These probes are optimized for single, double, or triple resonance experiments requiring high ¹H sensitivity and overall performance with simultaneous pulsing and decoupling capability at two low frequencies.

Cold probes

A ProTune probe adapter module (PAM) is used to provide a universal interface between the ProTune drive module and the Cold probe.

Maintenance

There is no required maintenance for the individual units.

| Minimum requirements | | | | | |
|--|---|---|--|--|--|
| Agilent NMR system | DD2 MR, 400-MR DD2, VNMRS, and 400-MR NMR | | | | |
| Probe | OneNMR, AutoX, automatable room temperature probes, automatable Triple Resonance Cold probe, automatable Dual Cold probe, | | | | |
| Vnmr.l | Vnmr.1.2.2C or newer | | | | |
| Host workstation | Dell 390N to currently released NMR workstation | | | | |
| Host OS | Linux RH Enterprise WS v.3, update 4, 64-bit or newer NMR OS | | | | |
| ProTune software | User software ver. 1.1 and probe-specific tuning files | | | | |
| Specifications | Master driver module | Dual/Single drive module | | | |
| Power | On master driver only | N/A | | | |
| Supply voltage or voltage range, current | 100–120 V, 2 A | N/A | | | |
| | 220–240 V, 1 A | | | | |
| Frequency range | 50/60 Hz | N/A | | | |
| Fuse rating | 115 V fuse type rated 250 V, 2A | N/A | | | |
| | 220 V fuse type rated 250 V, 1A | | | | |
| Connectors | | | | | |
| Input | A/C power cord, Control RJ45 to DIO board, Ethernet | N/A | | | |
| Output | Two flexible drive shafts | Dual: two flexible drive shafts | | | |
| | | Single: one flexible drive shaft | | | |
| Dimensions W x H x D in (cm) | 7.5 (19.05) × 10.94 (27.79) × 8.7 (22.14) | 7.5 (19.05) × 10.94 (27.79) × 8.7 (22.14) | | | |
| Weight lb (kg) | 10 (4.5) | 9 (4) | | | |

Ordering Information

| Product name | Part number | Includes |
|--|---|--|
| ProTune accessory for RT and Cold probes | G5140A | ProTune master driver module Two dual drive modules |
| Compatible with AutoX probes, ProTune style OneN option G5140A-001: ProTune Adapter Module for au | Six flexible shafts for tuning and matching Interface cables | |
| ProTune accessory for Cold probes only | G5141A | |
| Compatible with automatable Cold probes with the | | |
| Triple Resonance Cold probe: only high band channel | | |
| ProTune will tune and match high band and tune lo | w band only. | |
| ProTune accessory 700 MHz magnets | 191740500 | ProTune master driver module only |
| | | • Two flexible shafts for tuning and matching a single channel |
| ProTune accessory for 800–900 MHz magnets | 192250900 | RF front end accessory kit |
| | | Interface cables |

Note:

1. If you are planning to add a room temperature as well as a Cold probe in the future, we advise purchasing p/n: G5140A.

2. The ProTune driver modules require a special module when operating with an ultra high field (UHF), unshielded magnet.

3. HFX RT probe: Requires an additional single drive module. ProTune will tune and match ¹H, tune only ¹⁹F and tune and match low band on applicable probe. In most cases, an HF filter/combiner kit will be required.

www.agilent.com/chem/nmr

This information is subject to change without notice.

© Agilent Technologies, Inc., 2012 Published in the USA, March 22, 2012 5991-0043EN



Agilent Technologies