

# Purification of Four Anabolic Steroids using Agilent's SampliQ C8 and Amino SPE Tubes

## Technical Overview

### Introduction

Urine analysis is the matrix of choice for detecting the illegal use of anabolic steroids in food-producing livestock. Since these steroids act as growth promoters, concern about human health after consumption of illegally exposed animals is the driving force for monitoring their use.

Free steroids get conjugated in biological systems. Conjugated steroids in the urine are enzymatically hydrolyzed to the free steroids. The free steroids are isolated, purified, and concentrated using an Agilent SampliQ C8 solid phase extraction (SPE) tube. Further sample cleanup is achieved by passing the eluent from the C8 SPE through an Agilent SampliQ Amino SPE tube to remove acids, like uric acid. The amino functional group absorbs the acids, while allowing the steroids to pass through. After evaporation of the eluent from the amino SPE tube, the residue is taken up in the appropriate solvent for analysis.

The four representative anabolic steroids used were: prednisolone, dexamethasone, 1,4-androstadiene-3,17-dione and norgestrel, with five replicates each. After initially spiking urine with the steroids, following the above procedure, and using HPLC analysis, recoveries for three of the steroids were 95% to 99%, with norgestrel's recovery at 86%, and RSDs 1.0% or lower. The combined use of Agilent's SampliQ C8 and Amino SPE tubes provided high recoveries with excellent reproducibility.

### Study Purpose and Methodology

The purpose of this study was to identify SPE cleanup conditions for anabolic steroids. Prednisolone, dexamethasone, 1, 4-androstadiene-3,17-dione, and norgestrel were used as representative compounds. See Figure 1 for the chemical structures of the compounds.



**Agilent Technologies**




Figure 1. Steroid structures.

An aliquot (5 mL) of human urine was adjusted to pH 5.2 by the addition of 2 M sodium acetate buffer (2 mL, pH 5.2). The steroids were spiked into the sample, corresponding to an individual 4 µg/mL steroid level in urine. The steroids were subjected to enzymatic hydrolysis, conducted by adding 50 µL of β-glucuronidase type HP-2 from *Helix pomatia* juice (Sigma-Aldrich, Cat No: G7017) and incubating the samples at 55 °C for 3 hrs. Then the samples were centrifuged at 2000 rpm for 10 minutes at 5 °C [1].

An Agilent SampliQ C8 SPE tube (3 mL tube, 500 mg, p/n 5982-1035) was used for cleanup and concentration of the enzymatic hydrolyzed urine sample. The hydrolyzed sample was applied to the C8 SPE tube as illustrated in Figure 2.




Figure 2. Sample cleanup scheme using SampliQ C8 SPE

An Agilent SampliQ Amino tube (3 mL tube, 500 mg, p/n 5982-1835) was used to further purify the sample by removing uric and other acids typically found in urine. The eluted sample from the C8 SPE tube was applied to the amino SPE tube as illustrated in Figure 3. The combined effluent was evaporated to dryness at ambient temperature under a stream of nitrogen. The residue was dissolved in methanol-water (50:50 (v/v), 400  $\mu$ L), and analyzed at 254 nm using external standards.




Figure 3. Sample cleanup scheme using SampliQ Amino SPE

## Results and Discussion

The average of five recoveries for three of the steroids ranged from 95% to 99%, with the recovery of the fourth steroid at 86% (Table 1). Percent RSD for the five replicate recoveries was 1.0% or less. A representative chromatogram for the four steroids after SPE is shown in Figure 4.

Table 1. Recovery Data

	Prednisolone % Recovery	Dexamethasone % Recovery	Andro* % Recovery	Norgestrel % Recovery
Recovery #1	98	98	94	86
Recovery #2	98	97	94	86
Recovery #3	100	98	95	86
Recovery #4	100	98	95	87
Recovery #5	99	99	96	88
Average	99	98	95	86
%RSD	0.9	0.8	1.0	1.0

\* 1,4-Androstadiene-3,17-Dione




Figure 4. Representative sample chromatogram, 5 µg/mL 1. Prednisolone; 2. Dexamethasone; 3. 1,4-Androstadiene-3,17-Dione; 4. Norgestrel

### HPLC Analysis

Column:	Agilent ZORBAX Rapid Resolution Eclipse Plus C18, 4.6 mm × 150 mm, 3.5 µm, (p/n 959963-902)	
Mobile phase:	A: H <sub>2</sub> O (0.1% Formic Acid) B: ACN (0.1% Formic Acid)	
Gradient profile:	Time (min.)	%B
	0.0	40
	5.0	40
	16.0	75
	17.0	75
	17.1	40
Flow rate:	1 mL/min	
Column temperature:	25 °C	
Detection DAD:	254 nm	
Injection Volume:	30 µL	

## Conclusion

The combined use of Agilent's SampliQ C8 and Amino SPE tubes provided high recoveries with excellent reproducibility for the four tested anabolic steroids. This convenient procedure using Agilent's SPE tubes provided a significantly more concentrated and cleaner sample for analysis. One would expect this procedure to be applicable to additional steroids.

## Reference

1. S.A. Hewitt, M. Kearney, J.W. Currie, P.B. Young, D.G. Kennedy, *Analytica Chimica Acta* 473 (2002) 99-109

## Agilent SPE part numbers

Description	Part number
SampliQ C8 SPE tube, 3 mL tube, 500 mg	5982-1035
SampliQ Amino tube, 3 mL tube, 500 mg	5982-1835

## For More Information

For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

## Author/Contact

Joan Stevens is a sample preparation application chemist based at Agilent Technologies, Inc., Wilmington DE, USA.

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2010

Printed in the USA

April 6, 2010

5990-5543EN



**Agilent Technologies**