

# Quantitation of Nitrofurán Metabolites in Shrimp and Poultry by LC/MS/MS Using the Agilent LC/MSD Trap XCT

## Application

Food

### Authors

Bernhard Wüst, Christian Sauber  
Agilent Technologies GmbH,  
Hewlett-Packardstr. 8, 76337  
Waldbronn  
Germany

Hans (J.) A. van Rhijn  
RIKILT -Institute of Food Safety  
Bornsesteeg 45,6700 AE Wageningen  
the Netherlands

shrimp, and poultry. In 1995 the four drugs Furalolidone, Furaltadone, Nitrofurantoin, and Nitrofurazone were banned by the European Union (EU) for their usage in food-producing animals. Due to their rapid metabolism nitrofurán parent substances are not suitable for monitoring and typically their metabolites are analyzed. A liquid chromatography /mass spectrometry/mass spectrometry (LC/MS/MS) method was developed for the sensitive, qualitative, and quantitative analysis of four derivatized nitrofurán metabolites found in shrimp and poultry. See Figure 1.

### Abstract

**An LC/MS/MS method was developed for the qualitative and quantitative measurement of nitrofurán metabolites in chicken and shrimp using the Agilent LC/MSD XCT Ion Trap. The limit of quantitation for all four nitrofuráns investigated easily met the specified European Union Minimum Required Performance Level of 1 µg/kg and ranged from 0.125 µg/kg to 0.5 µg/kg.**

### Introduction

Nitrofurán antibiotics are widely used for the treatment of infectious diseases in cattle, pigs,

### Experimental

All liquid chromatography/mass spectrometry (LC/MS) experiments were performed using an Agilent 1100 Series LC system coupled to a mass selective detector (MSD) Ion Trap XCT mass spectrometer. The Ion Trap was operated with an orthogonal electrospray source in positive ion mode using multiple reaction monitoring (MRM). A gradient method was used for chromatography. Deuterated NBA-AMÖZ was used as the internal standard (ISTD) for NBA-AMÖZ and deuterated NBA-AÖZ was used as the ISTD for the other metabolites.



Agilent Technologies



---

**Chemicals**

Methanol	(Biosolve, 13683502)
Hydrochloric acid, HCL 32%	(Merck, 100319)
2-Nitrobenzaldehyde (2-NBA), C <sub>7</sub> H <sub>5</sub> NO <sub>3</sub>	(Sigma, N6001)
Tri-sodiumphosphate-dodecahydrate, Na <sub>3</sub> PO <sub>4</sub> 12(H <sub>2</sub> O)	(Merck, 106578)
Sodium hydroxide, NaOH	(Merck, 106498)
Ethyl acetate, CH <sub>3</sub> COOC <sub>2</sub> H <sub>5</sub>	(Biosolve 05402602)
Acetic acid, CH <sub>3</sub> COOH	(Merck 10063)
Acetonitrile, CH <sub>3</sub> CN	(Biosolve, 01203502)
Methanol-d, 99.5% D	(Aldrich, 15.193-9)

---

**LC/MS/MS Method Details**

HPLC:	Agilent 1100
Flow rate:	0.4 mL/min
Column:	Zorbax XDB-C18, 2.1 mm × 150 mm, 3.5 µm
Mobile phases:	A: Water + 0.1% acetic acid B: Acetonitrile + 0.1% acetic acid Gradient: 0–14 min: 10% A - 45% A; 14–16 min: 45% A - 90% A
Injection:	50 µL out of 500 µL

---

**MS 1100 LCMSD XCT Ion Trap**

Ionization mode:	Positive ESI
Nebulizer pressure:	45 psi
Drying gas flow:	12 L/min
Drying gas temperature:	350 °C
Skimmer:	20 V
Capillary exit:	55 V
Trap drive:	55
ICC:	On

---

**MRM mode**

4 Segments:	0–2.2 min, 2.2–7.4 min, 7.4–10.6 min, 10.6–13.9 min			
Segment 1:	Divert valve: to waste, no spectra			
Segment 2:	MS/MS of <i>m/z</i> 335,	Isolation width: 2.0,	Cut-off: 140;	Amplitude: 1.16
	MS/MS of <i>m/z</i> 340,	Isolation width: 2.0,	Cut-off: 140;	Amplitude: 1.16
Segment 3:	MS/MS of <i>m/z</i> 209,	Isolation width: 2.0,	Cut-off: 120;	Amplitude: 1.28
	MS/MS of <i>m/z</i> 249,	Isolation width: 2.0,	Cut-off: 100;	Amplitude: 1.25
Segment 4:	MS/MS of <i>m/z</i> 236,	Isolation width: 2.0,	Cut-off: 100;	Amplitude: 1.25
	MS/MS of <i>m/z</i> 240,	Isolation width: 2.0,	Cut-off: 100;	Amplitude: 1.25

---

**Quantitation**


NBA-AMZ:	EIC of 261 + 291 (MS/MS of 335),	Ret. Time: 4.5 min
NBA-dAMZ:	EIC of 266 + 296 (MS/MS of 340),	Ret. Time: 4.5 min
NBA-SEM:	EIC of 166 + 192 (MS/MS of 209),	Ret. Time: 9.9 min
NBA-AHD:	EIC of 134 (MS/MS of 249),	Ret. Time: 10.0 min
NBA-AOZ:	EIC of 134 (MS/MS of 236),	Ret. Time: 10.8 min
NBA-dAOZ:	EIC of 134 (MS/MS of 240),	Ret. Time: 10.8 min

Maximum accumulation time:	150 ms
Smart target:	100.000
Scan:	100–350

## Results and Discussion


Very low limits of detection (LOD) are required for nitrofuran metabolites and the derivatization method increased the ionization efficiency, as well as improving the chromatographic separation. A liquid-liquid extraction procedure was used which resulted in a relatively high concentration factor to further improve LOD.

The ion trap mass spectrometer was operated in MRM mode. In this mode, only precursor ions are chosen and full-scan MS/MS-spectra of the corresponding analytes are acquired. These full scan-MS/MS spectra are then used for identification by comparing them with MS/MS-spectra stored in a library. No further qualifier ion has to be monitored. See Figure 2.




**Figure 2.** Representative chromatograms (2A) and MS/MS spectra (2B) for all analytes plus ISTDs (1 µg/kg).

Quantitation is performed by selecting one or more product ions to create extracted ion chromatograms for each analyte and ISTD. The product ions used for quantitation were selected for best signal-to-noise (S/N) ratio post-acquisition. See Figures 3 and 4.




**Figure 3. Limit of quantitation (LOQ) for NBA-AM0Z, 0.125 µg/kg in shrimp matrix.**



**Figure 4. Calibration curve for NBA-AM0Z, 0.125 µg/kg – 2 µg/kg poultry matrix, three replicates.**

NBA-AMTZ and NBA-SEM were quantified using the sum of two product ions, while NBA-AHD and NBA-AOT were quantified using one product ion.

The European Union has set a Minimum Required Performance Level (MRPL) of 1 µg/kg for nitrofurantol metabolites. These detection limits are easily reached using this method with LOQs ranging from 0.125 µg/kg for NBA-AMTZ to 0.25 µg/kg for NBA-AOT and NBA-SEM, and 0.5 µg/kg for NBA-AHD. See Figure 5.



**Figure 5. Representative chromatogram of a positive shrimp sample at a level of 0.25 µg/kg.**

Linearity of the method was evaluated up to twice the MRPL (2 µg/kg) and showed a linear weighted regression (1/×) with coefficients of correlation of 0.99 or better. Intraday relative standard deviations (RSDs) were below 10% for all analytes at all concentrations. See Table 1.

**Table 1. Method Reproducibility and Accuracy for the Four Target Derivatized Metabolites**

Standard (µg/kg)	NBA-SEM		NBA-AOZ		NBA-AMOZ		NBA-AHD	
	SD %	Accuracy % (average)	SD %	Accuracy % (average)	SD %	Accuracy % (average)	SD %	Accuracy % (average)
0.125	3.77	98.81	2.29	98.94	4.34	101.10		
0.25	2.26	102.72	2.52	101.76	5.87	95.55	6.05	98.31
0.5	3.40	100.72	3.35	100.58	6.21	105.19	4.59	103.87
1	3.56	96.62	3.01	101.90	5.46	99.11	6.53	100.22
2	2.94	101.12	2.13	96.82	5.77	99.05	6.97	97.60

All calibration curves linear weighted 1/x      n = 6

## Conclusions

An LC/MS/MS method was developed for the qualitative and quantitative measurement of nitrofurans metabolites in chicken and shrimp using the Agilent XCT Ion Trap. The LOQ for all four nitrofurans investigated easily met the specified EU MRPL of 1 µg/kg and ranged from 0.125 µg/kg to 0.5 µg/kg.

## For More Information

For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2004

Printed in the USA  
March 25, 2004  
5989-0738EN



**Agilent Technologies**