Determination of Trace levels of PAHs in Food Grade Oils using Heart-cutting 2D-GC and Triple Quadrupole MS

INTRODUCTION

Polycyclic Aromatic Hydrocarbons (PAH) can potentially be present in food grade oils (white mineral oil), that are used in food production. Quality control of these oils at ppm level can be done by spectroscopic techniques (UV), but for speciation and for accurate determination of the most toxic high molecular weight PAHs, a combination of chromatographic and spectroscopic techniques is needed.

For the determination of 4-5 ring PAHs at ppb levels, labor intensive sample preparation, often including column chromatography, is needed. Hereby the aromatics are isolated from the bulk of the apolar alkanes. Column chromatography has the advantage of high sample capacity, but solvent consumption and process time are high.

Can ppb levels of PAHs in mineral oil be detected if direct injection is used?

In this work, a comparison was made between GC-MS using a single quadrupole system, GC-MS/MS using a triple quadrupole system and GC-GC-MS/MS using heart-cutting twodimensional GC.

METHOD 1: GC-MSD

Direct injection of a 10% dilution, using 1/10 split ratio with GC-MSD in SIM mode

Analytical parameters

Instrument:	Agilent 7890A GC & 5975 MSD
Inlet:	SSL at 280 °C, 1 μL, split ratio 10:1
Column:	30 m x 0.25 mm ID x 0.25 μm HP-5MS
Oven:	70 °C – 1 min – 10 °C/min – 325 °C – 8.5 min
MSD:	SIM mode, ions 228, 252, 276 (50 ms dwell time)

Results

Below the extracted ion chromatograms (from SIM) are given for both oils. In each case, overlays are shown for respectievely the blank oil (blue trace), the oil spiked at 100 ppb $(= 10 \text{ pg}/\mu\text{L} \text{ in solution, red trace})$ and a 10 pg/ μL reference standard (black trace). The PAHs cannot be detected in the oils due to the very high background.



F. David, C. Devos, K. Jacq, P. Sandra

Research Institute for Chromatography, Kennedypark 26, B-8500 Kortrijk, Belgium

METHOD 2: GC-MS/MS

Direct injection of a 10% dilution, using 1/10 split ratio with GC-triple quadrupole MS in MRM mode

Analytical parameters

Instrument: Inlet: Column: Oven: MSD:

Agilent 7890A GC & 7000A QQQ SSL at 280 °C, 1 µL, split ratio 10:1 30 m x 0.25 mm ID x 0.25 µm HP-5MS 70 °C – 1 min – 10 °C/min – 325 °C – 8.5 min MRM mode, 2 transitions per target compound

Results

Below the extracted ion chromatograms (from MRM) are given for respectievely oil 1 and oil 2 spiked at 100 ppb. All PAHs are detected. A higher background is observed for BaA/Chrys in oil 1 and Ipyr/BghiP in oil 2. LOD (S/N = 3) \approx 15 ppb. After several analyses of oil 2, source contamination was observed as well.



SAMPLES

Two food grade oils with different boiling point/MW range were selected. The oils were diluted at 10% in iso-octane (1 g in 10 mL). For each oil, two samples spiked with PAHs at respectively 100 and 10 ppb (μ g/kg) were prepared too: OIL 1: C18-C32 range (elution window in method 1: 14-26 min) OIL 2: C24-C40 range (elution window in method 1: 20-34 min)

For easy comparison, only the results for benzo(a)anthracene (BaA), chrysene (Chrys), benzo(b)fluoranthene (BbF), benzo(k)fluoranthene (BkF), benzo(a)pyrene (BaP), indeno(123cd)pyrene (Ipyr) and benzoghi)perylene (BghiP) are discussed.

METHOD 3: GC-GC-MS/MS



Direct injection of a 10% dilution, using 1/10 split ratio with 2D-GC and MS/MS in MRM mode

Analytical parameters

Instrument: Inlet:

¹D column: ¹D gas flow: GC oven: ²D column: ^{2}D gas flow: LTM oven: Heart-cuts: Detection:

Agilent 7890A GC & 7000A QQQ with capillary flow Deans switch 1 μL, MMI in split mode 100 °C – 720 °C/min – 300 °C (5min) HP-5MS, 30 m x 0.25 mm x 0.25 μ m – **boiling point separation** 1.4 mL/min, helium, constant flow 80 °C (1 min) – 10 °C /min – 325 °C (30 min) DB-17MS, 30 m x 0.25 mm x 0.25 µm - **∆RI ≈ 500** 146 kPa, helium, constant pressure 40° °C (27.8 min) – 20 °C/min – 250 °C – 7 °C /min – 320 °C – 5 °C /min – 330 °C windows set at retention times measured at FID, +/- 0.1 min FID (monitor detector) and MS (in MRM mode) MRM: 2 transitions per target compound

Results

Below the extracted ion chromatograms (from MRM) are given for oil 1 and oil 2 spiked at **10 ppb**. Using 2D-GC-MS/MS, both the selectivity of the separation and detection are enhanced. LOD (S/N = 3) < 5 ppb. No source contamination was observed.







Research Institute for Chromatography

