



SEC Analysis of Sodium Polystyrene Sulfonate

Application Note

Authors

Greg Saunders, Ben MacCreath
Agilent Technologies, Inc.

Introduction

Sodium polystyrene sulfonate is an ion-exchange resin used to treat high levels of blood potassium by substituting potassium in the diet and body for sodium in the resin. The potassium is then excreted and levels in the body are reduced. A sample of sodium polystyrene sulfonate was analyzed by aqueous SEC using Agilent PL aquagel-OH columns. These columns combine high pore volume and high column efficiency (>35,000 plates/meter) for maximum resolution. As the polymers are both ionic and relatively hydrophobic, eluent conditions are chosen to minimize sample to column interaction, which would otherwise result in late elution times.



Agilent Technologies

Conditions

Samples: Sodium polystyrene sulfonates
Columns: 2 x PL aquagel-OH 40 8 μm ,
300 x 7.5 mm (p/n PL1149-6840)
Eluent: 80 % 0.3 M NaNO_3 + 0.01 M
 NaH_2PO_4 at pH 9 + 20 % Methanol
Flow Rate: 1.0 mL/min
Detection: RI

Conclusion

SEC using PL aquagel-OH columns successfully analyzed samples of sodium polystyrene sulfonate. Aqueous SEC not only provides molecular weight data but also provides information on the polydispersity and the shape of the molecular weight distribution. The excellent chemical and mechanical stability of these columns offer high performance with good repeatability and column lifetime.

This product is not cleared or approved by the United States FDA or other international regulatory agency for use in clinical diagnostic procedures. No use of this product as a diagnostic medical device is implied or should be inferred.




Figure 1. Raw data chromatograms of three samples of sodium polystyrene sulfonate

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 2, 2010

SI-01574



Agilent Technologies