

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



Agilent Technologies

SPECIFICATION FOR A 7.0 TESLA/400MM **ROOM TEMPERATURE BORE MAGNET SYSTEM**

*Prepared by:-
Magnex Scientific Limited
The Magnet Technology Centre
6 Mead Road
Oxford Industrial Park
Yarnton, Oxford OX5 1QU, UK*

*Tel : +44 (0)1865 853800
Fax : +44 (0)1865 842466
E-mail : sales@magnex.com
www : magnex.com*

Document Ref : TS1047E

Date : January 2005

CONTENTS

1. Description of System
2. The Superconducting Magnet
 - i General Description
 - ii Specifications
 - iii Superconducting Shim Coils
3. The Cryostat
 - i General Description
 - ii Specifications
4. Scope of Supply

1. DESCRIPTION OF SYSTEM

The MRBR 7.0/400 System is a complete superconducting magnet system intended primarily for Research Studies on the Clinical/Biological applications of NMR Imaging (MRI) and NMR Spectroscopy (MRS).

The system essentially consists of a highly homogeneous superconducting magnet (300MHz ^1P , 7.05 Tesla) housed in a horizontal room temperature bore (400mm), low-loss helium cryostat. Field shimming is normally accomplished using both superconducting and room temperature shim coils. For Imaging applications X, Y and Z gradient coils can also be provided with the system.

Normally the room temperature shims and gradients are mounted on independent non-conducting formers, and are positioned in the room temperature bore of the system.

The system can be complemented with the supply of electronic consoles housing the superconducting magnet power supply with integral switching unit for energisation of the superconducting shim coils, multi channel room temperature shim coil power supply and cryomonitors for helium and nitrogen. An emergency quench heater control unit is also provided.

2. THE SUPERCONDUCTING MAGNET

i. General Description

The magnet is wound from multi-filamentary NbTi conductor with a high percentage of copper to superconductor. The windings are placed on a precision machined aluminium alloy former and then fully vacuum impregnated for robustness and long-term reliability.

The field homogeneity is defined over a 20cm diameter spherical volume and all orders of impurity up to and including 12th order are theoretically cancelled within this volume. Inevitably winding tolerances and small amounts of environmental influence will distort the central field. Corrections for these distortions are made in the first instance by superconducting shim coils loaded on a former surrounding the main coil. Final corrections are made by room-temperature correction coils placed in the bore of the system.

The magnet coils are fully protected from accidental damage due to a quench by a cold diode network located within the helium reservoir.

In the event of the need to activate an emergency discharge of the magnet a quench heater circuit is incorporated within the windings.

The magnet is designed to conservative levels of stress and mechanical stability to ensure reliable and stable operation. In addition the use of high quality superconducting wire ensures that a highly stable magnet system is achieved.

ii. Specifications

Magnet type : Multi-coil superconducting

Central field : 7.05 Tesla (300MHz ¹P)

Field stability measured a minimum of 72 hours after energisation : Less than 0.05 ppm/hour

Operating current : 178 Amps (nominal)

Field homogeneity values

Superconducting only shimmed : Less than 20ppm over 20cm dsv*

Fully shimmed using RT shims : Less than 5ppm over 20cm*
0.1ppm hhlw over 13cm dsv**

Estimate of helium consumption during ramping to full field : 150 litres

Fringe field (position of 5 gauss contour in unshielded state) : See Figure 1

Axially from magnet centre line : 10.9 metres

Radially from magnet centre line : 8.6 metres

* Defined as the peak to peak variations of points plotted over a seven plane plot on the surface of the stated spherical volume.

** hhlw measurement

iii. Superconducting Shim Coils

These coils are positioned on a non-conducting former surrounding the main coil in the helium reservoir. Each coil set is fitted with a superconducting switch for persistent mode operation.

Coil details:-

Shims provided : Z1, Z2, Z3, Z4, X, Y, ZX, ZY, XY,
X2-Y2, Z2X, Z2Y, ZXY & Z(X2-Y2)


Maximum recommended current : 25 amps

Coupling : All shims are decoupled from main coil
except Z1, Z2, Z3 and Z4

Typical shim strengths over 20cm diameter:

Shim	Strength at 10 cm radius (ppm/amp)
Z1	36.6
Z2	16.3
Z3	2.39
Z4	1.20
X	5.90
Y	5.90
ZX	1.21
ZY	1.21
XY	0.47
X2-Y2	0.47
Z2X	0.18
Z2Y	0.18
ZXY	0.032
Z(X2-Y2)	0.032

Figure 1 Stray Field (Unshielded)



3. THE CRYOSTAT

i. General Description

The cryostat is of conventional design, consisting of a central all-welded stainless steel helium vessel which is surrounded by an aluminium gas-cooled radiation shield and liquid nitrogen reservoir. The complete assembly is contained in a stainless steel outer vacuum vessel with a vertical service turret located centrally on top of the cryostat. The turret provides access to the helium reservoir for the demountable magnet leads, helium level probe, and helium transfer siphon. The outer vessel has end-flange closures constructed from aluminium which are sealed to the main body and bore-tube by compressed rubber 'O' ring seals. The room-temperature bore-tube is constructed from stainless steel.

The cryostat is supplied with a support stand that consists of load-spreading plates which have provision for fixing to the floor of the installation room. The helium reservoir contains in total approximately 1030 litres of liquid helium of which approximately 190 litres volume is above the minimum operating level. Details of refill intervals are given below.

Cryogen level monitors are incorporated into both the liquid helium and liquid nitrogen vessels and the associated electronics provide liquid level display and low level alarms. A back-up liquid helium level probe is included for use in the event of failure of the primary probe. The probes will monitor helium levels continuously from empty to full conditions.

ii. Specifications

The cryostat specifications (see figure 2) are as follows:-

Dimensions:-

Length of cryostat	: 1840mm
Height with support frame	: 2210mm
Room temperature clear bore (without shims and gradients)	: 400mm
Room temperature bore-tube material	: Stainless steel
Centre of field to base of stand	: 875mm

Cryostat end flange to centre of field : 920mm
Outside diameter : 1720mm
Minimum ceiling height for service tools : 3033mm
Weight of cryostat (excluding cryogens) : 5200kg (approx)


Liquid helium cryogen details :-


Volume for initial installation (includes cooling the magnet from 77K to 4.2K, volume required to completely fill helium reservoir and top-up after magnet energisation) : 2500 litres
Recommended refill volume during normal operation : 190 litres
Maximum volume of reservoir : 1030 litres
Hold-time during normal operation (static magnetic field, leads withdrawn) : More than 45 days


Liquid nitrogen cryogen details :-


Volume for initial installation (includes pre-cool of magnet to 77K and volume required to completely fill LN2 reservoir) : 4000 litres
Volume of reservoir : 170 litres nominal
Refill volume : 170 litres nominal excluding transfer losses
Hold-time in static condition : More than 12 days

Figure 2 Drawing of cryostat









				CUSTOMER INTERFACE MRBR 7.0T/400 AHZ323396E				SHEET 4 OF 5	
drawn BY:- N.COLEMAN checked BY:- approved BY:- DATE:- 20.03.01				TITLE:- OUT TRAY				DRAWING NUMBER	
ALL WEIGHTS AND DIMENSIONS ARE APPROXIMATE AND CAN CHANGE WITHOUT NOTICE.				MATERIAL				FINISH	
REVISION				PAINT RAL 9010				SCALE	
BY				NO				1:40	
CHKD.				ID				23396	
REV. LEVEL				STOP				NO PART IS TO BE USED WITHOUT PRIOR PERMISSION.	
REV. DATE				REVISION				DETAILS	
MP				E				07.01.05	
MP				D				20.10.04	
PJ				C				25.08.04	
PJ				B				17.03.04	
--				A				20.03.01	
1ST ISSUE				N/A				1ST ISSUE	
STAND SIZE CORRECTED, FOOT SHIMS ADDED				N/A				N/A	
JACKING POINTS AND STAND ADDED, GROUND FILL KIT ADDED				N/A				N/A	
SHEET 4 ADDED WITH UPDATED HEATSINK ASSEMBLY				N/A				N/A	
SERVICE TURRET UPDATED & FOOTPRINT DETAIL ADDED				N/A				N/A	


4. SYSTEM COMPONENTS

1 off	7T/400mm magnet system with integral s/c shims, housed in a low loss cryostat.
1 off	Stand
1 off	De-mountable main current lead
1 off	De-mountable s/c shim current lead
1 off	E5011 helium level monitor
1 off	E5031 nitrogen level monitor
1 off	Head oscillator
1 off	E7002 emergency discharge unit
1 off	Set of service cables (8.5m)
1 off	Helium monitor cable
1 off	Nitrogen monitor cable
1 off	Flexible siphon (2.0m)
1 off	Nitrogen blow-out tube
1 off	Nitrogen ground level fill kit
1 off	Spares kit
1 off	De-mountable helium level probe
1 off	De-mountable nitrogen level probe
1 off	System manual

APPENDIX

7T/400 Passive Room Shield Guidelines

Typical shield dimensions 3 m wide x 3 m tall x 5 m long



Magnex Scientific Ltd
Doc ref :7T 400 shield.doc