NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to **www.agilent.com/chem.** 

Agilent Technologies

## **TECHNICAL SPECIFICATIONS**

# FOR AN MRBR 4.7 TESLA/160mm ACTIVELY SHIELDED ROOM TEMPERATURE BORE MAGNET SYSTEM

Prepared by:-Magnex Scientific Limited The Magnet Technology Centre 6 Mead Road Oxford Industrial Park Yarnton, Oxford OX5 1QU, UK

| Tel    | : | +44 (0)1865 853800      |
|--------|---|-------------------------|
| Fax    | : | +44 (0)1865 842466      |
| E-mail | : | <u>sales@magnex.com</u> |
| www    | : | <u>www.magnex.com</u>   |

Document Ref : TS1348A Date : March 2004

## **CONTENTS**

- 1. Description of System
- 2. The Superconducting Magnet
  - i. General Description
  - ii. Specifications
  - iii. Superconducting Shim Coils
- 3. The Cryostat
  - i. General Description
  - ii. Specifications
- 4. System Components
  - i. Superconducting Magnet System Components
  - ii. Standard Ancillary Parts

## **<u>1. DESCRIPTION OF THE SYSTEM</u>**

The system essentially consists of a highly homogeneous superconducting magnet (4.7 Tesla) housed in a horizontal room temperature bore (160mm), low-loss helium cryostat. Field shimming is accomplished using a combination of superconducting shim coils and passive shim pieces.

Passive shims can be mounted on a tube in the bore of the magnet leaving an overall system bore size of 155mm diameter. However, if a gradient coil is fitted then the passive shims are mounted on the gradient.

Cryomonitors for helium and nitrogen and an emergency quench heater control unit are also provided.

## 2. THE SUPERCONDUCTING MAGNET

#### i. <u>General Description</u>

The magnet is wound from multi-filamentary NbTi conductor with a high percentage of copper to superconductor. The windings are placed on precision machined aluminium alloy and GRP formers and then fully vacuum impregnated for robustness and long-term reliability.

The field homogeneity is defined over two cylinders and a sphere. All orders of impurity up to 10th order are theoretically cancelled within these volumes. Inevitably winding tolerances and small amounts of environmental influence will distort the central field. Corrections for these distortions are made in the first instance by superconducting shim coils. Final corrections are made by passive shim pieces placed in the bore of the system.

The magnet coils are fully protected from accidental damage due to a quench by a diode resistor network located within the helium reservoir.

In the event of the need to activate an emergency discharge of the magnet a quench heater circuit is incorporated within the windings.

The magnet is designed to conservative levels of stress and mechanical stability to ensure reliable and stable operation. In addition the use of high quality superconducting wire ensures that a highly stable magnet system is achieved. Long term field stability is enhanced by the use of an internal superconducting field lock coil.

#### ii. Specifications

| Magnet type   | : | Multi-coil superconducting |
|---|---|----------------------------|
| Central field   | : | 4.7 Tesla                  |
| Field stability measured a minimum of 72 hours after energisation | : | Less than 0.05 ppm/hour    |

The above values assume that the room temperature is maintained at  $+/-1^{\circ}C$ 

| Operating current | : 250 Amps (nominal) |
|-------------------|----------------------|
|-------------------|----------------------|

Field homogeneity values

S/C and passively shimmed

: 10 ppm p-p over 8.5cm  $dsv^1$ 

Residual Gradients over 8.5 cm dsv<sup>1</sup>

| $Z^3$       | < 1ppm |
|-------------|--------|
| ZX          | < 1ppm |
| ZY          | < 1ppm |
| $X^2 - Y^2$ | < 1ppm |
| XY          | < 1ppm |

|   | : | 5ppm p-p over 8.0 cm dsv <sup>2</sup><br>10ppm p-p over cylinder 60mm diameter<br>by 60mm long <sup>3</sup> , and<br>4 ppm p-p over cylinder 30mm diameter<br>by 60mm long <sup>3</sup> |
|---|---|---|
| Fully shimmed using RT shims  | : | 0.1ppm FWHH over 4cm dsv <sup>4</sup>   |
| Time to energise magnet to full field                               | : | Less than 40 minutes  |
| Estimate of helium consumption during ramping to full field         | : | 10-15 litres  |
| Fringe field  | : | See Figure 1  |
| Position of 5 gauss contour <sup>5</sup>                            |   |   |
| Axially from magnet centre line<br>Radially from magnet centre line | : | 1.9 metres<br>1.1 metres  |

<sup>&</sup>lt;sup>1</sup> Defined as the peak to peak variations of points plotted over a seven plane plot on the surface of the stated spherical volume.

<sup>&</sup>lt;sup>2</sup> Projected over a seven plane plot.

<sup>&</sup>lt;sup>3</sup> Defined as the peak to peak variation of points plotted over a five plane plot on the surface of the stated cylinder.

<sup>&</sup>lt;sup>4</sup> If using Magnex standard RT shim set.

<sup>&</sup>lt;sup>5</sup> **Safety Note:** In the event of a quench it is possible for the magnetic field to momentarily bloom beyond this limit. For further details please consult the Magnex site planning guide for this magnet.



Figure 1 Fringe field plot of 4.7 T 160 mm actively shielded magnet

#### iii. <u>Superconducting Shim Coils</u>

These coils are positioned on a former in the helium reservoir. Each coil set is fitted with a superconducting switch for persistent mode operation.

Coil details:-

| Shims provided              | : | Z1, Z2, X, Y, ZX, ZY, XY, X2-Y2                        |
|-----------------------------|---|--|
| Maximum recommended current | : | 25 Amps  |
| Coupling                    | : | All shims are designed to be de-coupled from main coil |

Typical shim strength over 60mm diameter

| <u>Strength</u>         |
|-------------------------|
| (ppm/amp of main field) |
| 15.0                    |
| 6.5                     |
| 5.0                     |
| 0.42                    |
| 0.40                    |
|                         |

## 3. THE CRYOSTAT

#### i. General Description

The cryostat is of conventional layout, consisting of a central all-welded stainless steel helium vessel which is surrounded by an aluminium gas-cooled radiation shield and liquid nitrogen reservoir. The complete assembly is contained in a stainless steel outer vacuum vessel. The outer vacuum vessel is fitted with a vertical service turret located centrally on top of the cryostat. The turret provides access to the helium reservoir for the demountable magnet leads, helium level probe, and helium transfer siphon. The outer vessel is an all welded stainless steel construction with a room-temperature boretube constructed from stainless steel.

The cryostat is supplied with an adjustable height support stand. The helium reservoir contains in total approximately 400 litres of liquid helium of which approximately 200 litres volume is above the minimum operating level. Details of refill intervals are given below.

Cryogen level monitors are incorporated into both the liquid helium and liquid nitrogen vessels and the associated electronics provide liquid level display and low level alarms. A back-up liquid helium level probe is included for use in the event of failure of the primary probe. The probes will monitor helium levels continuously from empty to full conditions.

#### ii. Specifications

| The cryostat is shown in drawing no.             | AHZ328473 |                  |
|--|-----------|------------------|
| Dimensions:-                                     |           |                  |
| Length of cryostat                               | :         | 1012mm           |
| Room temperature clear bore (with passive shims) | :         | 155mm            |
| Room temperature bore-tube material              | :         | Stainless steel  |
| Centre of field to base of stand                 | :         | 1075-1145mm      |
| Cryostat end-flange to centre of field           | :         | 506mm            |
| Minimum ceiling height for siphon                | :         | 3130mm           |
| Weight of cryostat (excluding cryogens)          | :         | 1300kg (approx.) |

## Liquid helium cryogen details:-

| Volume for initial installation (includes<br>cooling the magnet from 77K to 4.2K,<br>volume required to completely fill helium<br>reservoir and to top-up helium reservoir<br>after magnet energisation) | : | 800 litres         |
|--|---|--------------------|
| Recommended refill volume during normal operation  | : | 200 litres         |
| Maximum volume of reservoir  | : | 400 litres         |
| Hold-time during normal operation<br>(static magnetic field, leads withdrawn)  | : | More than 150 days |
| Liquid nitrogen cryogen details:-  |   |                    |
| Volume for initial installation (includes<br>pre-cool of magnet to 77K and volume<br>required to completely fill LN2 reservoir)  | : | 800 litres         |
| Volume of reservoir  | : | 200 litres         |
| Refill volume  | : | 190 litres         |
| Hold-time in static condition  | : | More than 14 days  |







Tech Spec: TS1348A March 2004



## **<u>4 SYSTEM COMPONENTS</u>**

#### i. <u>Superconducting Magnet System Components</u>

- 1 off 4.7T 160 mm actively shielded MRBR magnet system with integral superconducting shims and lock coils housed in a low loss cryostat.
- 1 off Stand

### ii. <u>Standard Ancillary Parts</u>

| 1 off         | Helium level monitor              |                | E5011                               |
|---------------|-----------------------------------|----------------|-------------------------------------|
| 1 off         | Nitrogen level monitor            |                | E5031                               |
| 1 off         | Head oscillator                   |                | E5030                               |
| 1 off         | Emergency discharge unit          |                | E7007                               |
| 1 off         | Service cable                     |                | C0091085                            |
| 1 off         | Helium monitor cable              |                | C0090003                            |
| 1 off         | Nitrogen monitor cable            |                | C0113085                            |
| 1 off         | Flexible siphon (2.0m)            |                | P222000005                          |
| 1 off         | Nitrogen blow out tube            |                | AUC400198                           |
| 1 off         | Nitrogen fill tube                |                | AHU327799                           |
| 1 off         | Spares kit                        |                | AKZ509324                           |
| 1 off         | De-mountable Helium Level Probe   |                | AUE113867                           |
| 1 off         | De-mountable Nitrogen Level Probe |                | DUE100194                           |
| 1 off         | System manual                     |                |                                     |
| 1 off         | 3m nitrogen transfer line         |                | ATU327865                           |
| iii. <u>(</u> | Optional Ancillary Parts          |                |                                     |
|               | Quench duct assembly.             | 4"<br>6"<br>8" | AHC327457<br>AHC327456<br>AHC327455 |