

NOTICE: Varian, Inc. was acquired by Agilent Technologies in May 2010. This document is provided as a courtesy but is no longer kept current and thus will contain historical references to Varian. For more information, go to [www.agilent.com/chem](http://www.agilent.com/chem).



## Application Note SI-01998

# Synthesis of ampR1 Primer using StratoSpheres™ DNA

Andrew Coffey  
Varian, Inc.

### Introduction

The synthesis of primer sequences (for PCR, cloning and sequencing applications) is perhaps the most commonly used application for oligonucleotide synthesis. The ampR1 (or SR1) reverse primer sequence 5' CTT TCT GCT ATG GAG GTC AGG TAT G 3' is recommended for use in sequencing and amplification with the pSMART cloning vector. The StratoSpheres DNA dG support matches the performance of the leading commercial macroporous PS/DVB in the synthesis of this primer 20-mer sequence, and exceeds the capability of a controlled pore glass support (CPG).


### Oligonucleotide Synthesis

The oligonucleotides were prepared on 0.2  $\mu$ M scale using an Applied Biosystems 392 DNA/RNA Synthesizer and standard

chemistry. Following synthesis, the oligonucleotides were cleaved using ammonium hydroxide. Deprotection of the side chain protecting groups (A, C and G) was accomplished by heating to 55 °C overnight as required. The "DMT on" oligonucleotides were diluted with water in order to give an on-scale response during HPLC analysis.

### HPLC Analysis

Column: Pellicular SAX, 4 x 250 mm  
Eluent: A = 25 mM Tris-HCl, 0.5% ACN, pH 8.0; B = 25 mM Tris-HCl, 0.8 M Ammonium Chloride, 0.5% ACN, pH 8.0  
Gradient: 0 -100% B in 26 min  
Flow Rate: 1.5 mL/min  
Temp: 60 °C  
Detector: UV 260 nm



E. J. Summer *et al.* (2004) *Burkholderia cenocepacia* phage BcepMu and a family of Mu-like phages encoding potential pathogenesis factors. *J. Mol. Biol.*, 340, 49-65.

These data represent typical results.  
For further information, contact your local Varian Sales Office.

StratoSpheres, Varian and the Varian Logo are trademarks or registered trademarks of Varian, Inc. in the U.S. and other countries.

© Varian, Inc. 2009

Varian, Inc.  
[www.varianinc.com](http://www.varianinc.com)  
North America: 800.926.3000 – 925.939.2400  
Europe: *The Netherlands*: 31.118.67.1000  
Asia Pacific: *Australia*: 613.9560.7133  
Latin America: *Brazil*: 55.11.3238.0400

