

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



Agilent Technologies

TECHNICAL SPECIFICATIONS
FOR A WATER COOLED ACTIVELY SHIELDED
GRADIENT SYSTEM WITH RT SHIM SET

SGRAD 205/120/HD

Prepared by:-
Magnex Scientific Limited
The Magnet Technology Centre
6 Mead Road
Oxford Industrial Park
Yarnton, Oxford OX5 1QU, UK

Tel : +44 (0)1865 853800
Fax : +44 (0)1865 842466
E-mail : sales@magnex.com
www : magnex.com

Document Ref : TS1442D

Date : March 2005

CONTENTS

Description

1. Mechanical

- 1.1 Dimensions and weight
- 1.2 Finish
- 1.3 Mounting
- 1.4 Electrical connectors
- 1.5 Temperature sensors
- 1.6 Water cooling system

2. Gradient Coils

- 2.1 Strength
- 2.2 Linearity
- 2.3 Residual eddy currents
- 2.4 Inductance
- 2.5 Resistance
- 2.6 Safe operating conditions and duty cycle
- 2.7 Peak strength and rise-time (estimated)
- 2.8 Orthogonality
- 2.9 Insulation

3. Shim System

- 3.1 Shim characteristics

4. System Components

- 4.1 Gradient
- 4.2 Standard ancillaries
- 4.3 Optional extras

GENERAL DESCRIPTION

The SGRAD 205/120/HD is a fully self shielded gradient system designed to suit $\geq 210\text{mm}$ room temperature bore superconducting magnets and gradient systems.

The design incorporates fully optimised X, Y and Z coil configurations. The X and Y coils are made from the highest quality copper plates machined with CNC technology. The Z coil is wound from heavy duty copper strip. The gradient has been optimized to allow for high duty cycle experiments.

The room temperature shim set has been specially designed to minimise coupling between gradients and shims during pulsing. Finally the gradient set is fully vacuum impregnated to minimise mechanical vibration and noise.

1. MECHANICAL

1.1	Dimensions and weight		
	Total length	:	To suit magnet
	External diameter	:	205mm $\pm 3\text{mm}$
	Internal diameter	:	120mm $\pm 0.5\text{mm}$
	Approximate weight	:	TBD
1.2	Finish		
	Bore tube	:	Natural GRP
	End plates	:	Plated aluminium or natural composite
1.3	Mounting		
	Method at access end	:	TBD
	Method at service end	:	TBD
	Adjustment	:	$\pm 5\text{mm}$ axially
	Loading	:	Service or Access End
1.4	Electrical connectors ¹		
	Gradients	:	5 way bayonet
	Temperature sensors	:	15 way bayonet
	RT shims	:	41 way bayonet
1.5	Temperature sensors		
	Type	:	Type T thermocouples
	Number on inner section	:	4 typical
	Number of outer section	:	2 typical
1.6	Water cooling system		
	Volumetric flow rate	:	3.4 litres/min
	Supply pressure (typical)	:	4 bar
	Internal pressure drop	:	2-3 bar
	Heat extraction	:	2.55 kW for $\Delta T=10\text{K}$ and $T_{\text{inlet}}=10^0\text{C}$
	Gradient connectors	:	Double shut off connector
	Fitting for supply hose	:	3/8" ID barbed hose fitting
	Recommended water supply	:	Recirculating water chiller

¹ Also available with Lemo connectors.

2. GRADIENT COILS

- | | | | |
|-----|---|---|---|
| 2.1 | Strength | | |
| | X/Y/Z axis | : | 2.0mT/m/A $\pm 5\%$ |
| 2.2 | Linearity per axis (see figure 1) over 80mm d.s.v. | | |
| | X/Y/Z | : | $\leq 5\%$ |
| | Definition | : | % Linearity = Max spatial deviation as a percentage of the plotting radius. |
| 2.3 | Residual eddy currents (before pre-emphasis) ² | | |
| | X/Y/Z axis | : | $< 1\%$ |
| 2.4 | Inductance | | |
| | X axis | : | 155 μ H $\pm 15\%$ |
| | Y axis | : | 120 μ H $\pm 15\%$ |
| | Z axis | : | 140 μ H $\pm 15\%$ |
| 2.5 | DC Resistance | | |
| | X axis | : | 140m Ω $\pm 25\text{m}\Omega$ |
| | Y axis | : | 100m Ω $\pm 25\text{m}\Omega$ |
| | Z axis | : | 150m Ω $\pm 25\text{m}\Omega$ |
| 2.6 | Safe operating conditions | | |
| | Peak voltage | : | $\leq 300\text{V}$ |
| | Peak current | : | $\leq 200\text{A}$ |
| | RMS current | : | $\leq 70\text{A}$ rms indefinitely all axes together |
| | Peak internal temperature | : | $\leq 60^\circ\text{C}$ |
| 2.7 | Typical peak strength and rise-times (estimated) | | |
| | Peak strength @200A X/Y/Z | : | 400mT/m |
| | Rise-time ³ (0-98%) X | : | 110 μ s |
| | @ 200A, 300V Y | : | 85 μ s |
| | Z | : | 102 μ s |
| 2.8 | Orthogonality | | |
| | X to Y | : | 90 $^\circ$ $\pm 1^\circ$ |
| 2.9 | Insulation | | |
| | Between X, Y, Z, shims, sensors, cooling & enclosure. | : | $> 200\text{ M}\Omega$ at 1000 Volts DC |

² Measured 5ms after a 20ms trapezoidal pulse, and extrapolated to 1ms after the pulse.

³ Excluding the effects of lead resistance and filter impedance.




Figure 1
Theoretical image distortion of a 80mm cube phantom. DSV is denoted by the dashed circle

3. SHIM SYSTEM

3.1 Room temperature shims

First order shimming is achieved by DC offsets to the gradient coils. The nominal shim performances are shown below.

Shim	Strength (mG/cm ⁿ /A)	Inductance (mH)	Resistance (Ω)	Peak current (A)
shielded Z^0 (B_0)	530	0.05	0.30	10
Z^2	87.0	2.22	1.25	10
ZX,ZY	41.0	1.55	1.65	10
XY, X^2 - Y^2	12.7	0.80	1.65	10

4. SCOPE OF SUPPLY

4.1 Gradient:-

1 off Actively shielded gradient, type SGRAD MkIII 205/120/HD

4.2 Standard ancillaries:- (Not included on inserts)

1 off Set of X/Y/Z cables, standard length 15m C0082125

1 off RT shim cable, standard length 15m C0202150

1 off Thermometry cable, standard length 15m C0259150

4.3 Optional extras:-

2 off RF doors ARZ325130

1 off Stand alone temperature monitor unit E3515g

1 off Gradient management unit, consisting of:- E3500

(i) Computer controlled X, Y and Z pre-emphasis

(ii) Computer controlled Zo (Bo) pre-emphasis