

NOTICE: This document contains references to Varian.  
Please note that Varian, Inc. is now part of Agilent  
Technologies. For more information, go to  
**[www.agilent.com/chem](http://www.agilent.com/chem)**.



**Agilent Technologies**

**PRELIMINARY**

**TECHNICAL SPECIFICATIONS**

**FOR A WATER COOLED ACTIVELY SHIELDED**

**GRADIENT SYSTEM WITH RT SHIM SET**

**VECTOR 115/60/S**

*Prepared by:-*  
***Magnex Scientific Limited***  
***The Magnet Technology Centre***  
***6 Mead Road***  
***Oxford Industrial Park***  
***Yarnton, Oxford OX5 1QU, UK***

***Tel*** : +44 (0)1865 853800  
***Fax*** : +44 (0)1865 842466  
***E-mail*** : [sales@magnex.com](mailto:sales@magnex.com)  
***www*** : [magnex.com](http://magnex.com)

***Document Ref : TS1475A***

***Date : March 2005***

## **CONTENTS**

### Description

#### 1. Mechanical

- 1.1 Dimensions and weight
- 1.2 Finish
- 1.3 Mounting
- 1.4 Electrical connectors
- 1.5 Temperature sensors
- 1.6 Water cooling system

#### 2. Gradient Coils

- 2.1 Strength
- 2.2 Linearity
- 2.3 Residual eddy currents
- 2.4 Inductance
- 2.5 Resistance
- 2.6 Safe operating conditions and duty cycle
- 2.7 Peak strength and rise-time (estimated)
- 2.8 Orthogonality
- 2.9 Insulation

#### 3. Shim System

- 3.1 Shim characteristics

#### 4. System Components

- 4.1 Gradient
- 4.2 Standard ancillaries
- 4.3 Optional extras

## **GENERAL DESCRIPTION**

The VECTOR 115/60/S is a fully self shielded gradient system designed to suit  $\geq 120\text{mm}$  room temperature bore superconducting magnets and gradient systems.

The design incorporates fully optimised X, Y and Z coil configurations. The X and Y coils are made from the highest quality copper plates, etched with photo-lithographic techniques. The Z coil is wound from heavy duty copper strip.

The VECTOR range of gradients have been engineered to allow for high duty cycle experiments.

The room temperature shim set has been specially designed to minimise coupling between gradients and shims during pulsing. Finally the gradient set is fully vacuum impregnated to minimise mechanical vibration and noise.

## 1. MECHANICAL

1.1	Dimensions and weight		
	Total length	:	To suit magnet
	External diameter	:	115mm $\pm 4\text{mm}$
	Internal diameter	:	60mm $\pm 0.5\text{mm}$
	Approximate weight	:	20kg
1.2	Finish		
	Bore tube	:	Natural GRP
	End plates	:	Plated aluminium
1.3	Mounting		
	Method at access end	:	Flange bolted to larger gradient
	Method at service end	:	'O'-ring clamp
	Adjustment	:	$\pm 5\text{mm}$ axially
	Loading	:	Service end only
1.4	Electrical connectors		
	Gradients	:	Lemo (EGJ-5B-304-CLA)
	Temperature sensors	:	Lemo (EGJ-4B-320-CLA)
	RT shims	:	Lemo (EGJ-4B-324-CLA)
1.5	Temperature sensors		
	Type	:	Type T thermocouples
	Number on inner section	:	4 typical
	Number of outer section	:	2 typical
1.6	Water cooling system		
	Volumetric flow rate	:	1.2 litres/min
	Supply pressure (typical)	:	4 bar
	Internal pressure drop	:	2-3 bar
	Heat extraction	:	900 W for $\Delta T=10\text{K}$ and $T_{\text{inlet}}=10^0\text{C}$
	Gradient connectors	:	Double shut-off connector
	Fitting for supply hose	:	3/8" ID barbed hose fitting
	Recommended water supply	:	Recirculating water chiller

## 2. GRADIENT COILS


2.1	Strength		
	X/Y/Z axis	:	5.0mT/m/A $\pm 5\%$
2.2	Linearity per axis (see figure 1) over 40mm DSV		
	X/Y/Z	:	$\leq 5\%$
	Definition	:	% Linearity = Max spatial deviation as a percentage of the plotting radius.
2.3	Residual eddy currents (before pre-emphasis) <sup>1</sup>		
	X/Y/Z axis	:	< 1%
2.4	Inductance <sup>2</sup>		
	X axis	:	28 $\mu$ H $\pm 15\%$
	Y axis	:	32 $\mu$ H $\pm 15\%$
	Z axis	:	40 $\mu$ H $\pm 15\%$
2.5	DC Resistance <sup>2</sup>		
	X axis	:	55m $\Omega$ $\pm 25$ m $\Omega$
	Y axis	:	65m $\Omega$ $\pm 25$ m $\Omega$
	Z axis	:	80m $\Omega$ $\pm 25$ m $\Omega$
2.6	Safe operating conditions		
	Peak voltage	:	$\leq 300$ V
	Peak current	:	$\leq 200$ A
	RMS current	:	$\leq 50$ A RMS indefinitely all axes together
		:	$\leq 75$ A RMS indefinitely Y axis only
	Peak internal temperature	:	$\leq 60^{\circ}$ C
2.7	Typical peak strength and rise-times (estimated)		
	Peak strength @200A X/Y/Z	:	1000mT/m
	Rise-time <sup>2</sup> (0-98%) X	:	20 $\mu$ s
	@ 200A, 300V Y	:	25 $\mu$ s
	Z	:	30 $\mu$ s

---

<sup>1</sup> Measured 5ms after a 20ms trapezoidal pulse, and extrapolated to 1ms after the pulse.

<sup>2</sup> Excluding the effects of lead resistance and filter impedance.

- |     |   |   |                           |
|-----|---|---|---------------------------|
| 2.8 | Orthogonality<br>X to Y   | : | 90°+/-1°                  |
| 2.9 | Insulation<br>Between X, Y, Z, shims,<br>sensors, cooling & enclosure | : | > 200 MΩ at 1000 Volts DC |



**Figure 1**  
**Theoretical image distortion of a 40mm cube phantom. DSV is denoted by the dashed circle.**

### 3. SHIM SYSTEM

#### 3.1 Room temperature shims

First order shimming is achieved by DC offsets to the gradient coils. The nominal shim performances are shown below.

Shim	Strength (mG/cm <sup>n</sup> /A)	Inductance (mH)	Resistance ( $\Omega$ )	Peak current (A)
shielded $Z^0$ ( $B_0$ )	510	0.006	0.09	5
$Z^2$	157	0.093	0.51	5
ZX,ZY	124	0.247	0.81	5
XY, $X^2$ - $Y^2$	40	0.159	1.00	5



#### **4. SCOPE OF SUPPLY**

##### **4.1 Gradient:-**

1 off Actively shielded gradient, type SGRAD mk. IV 115/60/HD/S

##### **4.2 Standard ancillaries:- (Not included on inserts)**

1 off Set of X/Y/Z cables, standard length 15m C0397150

1 off RT shim cable, standard length 15m C0398150

1 off Thermometry cable, standard length 15m C0399150

##### **4.3 Optional extras:-**

2 off RF doors (Lemo type) ARZ330642

1 off Stand alone temperature monitor unit E3515g

1 off Gradient management unit, consisting of:- E3500

(i) Computer controlled X, Y and Z pre-emphasis

(ii) Computer controlled Zo (Bo) pre-emphasis