

# Amino Acid Neurotransmitter Analysis with Star 9080 EC Detector

# LC


## Varian Application Note

Number 17

Fred Klink  
Varian Chromatography Systems

**Key Words:** Star 9080, Amino Acids, Clinical, Biomedical

### Standard Separation



**Analysis of amino acid neurotransmitters after OPA-sulphite derivatisation.**


The amount injected is 20 pmol (1  $\mu$ mol/L) for each compound.

Peaks are: serine (1), glycine (2), taurine (3), glutamate (4), arginine (5), alanine (6), OPA reagent (7) and GABA (8).

<b>Detector</b>	Varian Star 9080 Amperometric Electrochemical Detector
<b>Column</b>	C18, 250 x 4.6 mm, 5 $\mu$
<b>Flow rate</b>	0.70 mL/min
<b>Mobile phase</b>	100 mM Phosphate buffer, pH 4.5, 0.5 mM EDTA, 25% methanol
<b>Sample</b>	20 $\mu$ L inj., amino acids derivatised with OPA-sulphite
<b>Temperature</b>	ambient
<b>Flowcell</b>	2.74 mm Glassy Carbon working electrode
<b>E-cell</b>	850 mV (vs. Ag/AgCl)

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to [www.agilent.com/chem](http://www.agilent.com/chem).

Analysis of Rat Hippocampal Tissue



Analysis of rat ventral hippocampal dialysate after OPA-sulphite derivatisation.  
Concentrations (amounts) are: 1. glutamate 3.6 nmol/L (72 fmol) and 2. GABA 1.7 nmol/L (34 fmol).