Agilent D6111
Dissolved Oxygen
Probe
溶解氧电极

Operating Guide
用户手册

Agilent Technologies

Overview

The D6111 Dissolved Oxygen (DO) Probe can measure the dissolved oxygen concentrations in aqueous solutions. It is used with the 3200D Dissolved Oxygen Meter or 3200M Multi-Parameter Analyzer.

WARNING

Use this probe according to the operating manual to avoid personal injury.

WARNING

The probe solution can cause chemical burns or illness if it is taken orally or contacted by human skin. Use protective clothing or gloves to avoid contact. In case of contact, rinse contacted area with tap water or deionized water thoroughly.

Specifications

Table 1 D6111 Dissolved Oxygen Probe specifications

Specification	Value
Temperature range	0 to 45 °C
Concentration range	0 to 20 mg/L
Body material	ABS
Probe diameter	12 mm
Probe length	120 mm
Cable interface	miniDin4
Cable length	1000 mm


Figure 1 D6111 Dissolved Oxygen Probe assembly

Operation

Preparing the probe

- 1 Remove the probe membrane cap (5190-0548).
- 2 Rinse the inside and surface of membrane cap with distilled or deionized water and shake gently to dry.
- 3 Rinse the probe cathode and anode with distilled or deionized water and wipe them dry.
- 4 Add DO filling solution (5190-0547) to the membrane cap to 3/4 capacity.
- 5 Keep the probe membrane cap upright and screw it slowly onto the probe body until finger tight.
- 6 Check that the membrane touches the cathode tip and there is no gas bubble or gap between the membrane and the cathode tip. If there is a gap or gas bubble between the membrane and the cathode tip, unscrew the membrane cap, add DO solution to the top, and reinstall the cap. Wipe the outside dry.

Polarizing the probe

- 1 Connect the probe to the meter and keep the meter at measuring state.
- 2 Expose the measuring tip to air or immerse the measuring tip into water saturated with air to polarize the probe.
- 3 For a new DO probe, a DO probe that has been refilled with new DO filling solution, or a DO probe that has been disconnected from the meter for more than 1 hour, polarize the probe for more than 1 hour. For a DO probe that has been disconnected from the meter for less than 1 hour, or if the meter has been idle for less than an hour, then polarize the probe for 25 minutes.

Calibration

Before measuring the sample solution, perform both a full scale calibration and zero oxygen calibration.

Full scale calibration

There are two methods of performing a full scale calibration. One method is with air for regular measurement. The other method is in water saturated with air at constant temperature for accurate measurement.

- To perform a full scale calibration with air:
 - 1 Place the probe in air at constant flow and temperature.
 - 2 Keep the membrane dry.
 - 3 When the temperature reading and DO readings become stable, perform a full scale calibration according to the meter operating manual.
- To perform a full scale calibration in water saturated with air at constant temperature:
 - 1 Soak the measuring tip in water saturated with air and wait until the temperature reading becomes stable.
 - 2 Maintain a constant velocity by agitating the probe with a constant speed of 20–80 cm/s.
 - 3 Ensure there are no air bubbles trapped on the membrane.
 - 4 When the reading becomes stable, perform a full scale calibration according to the meter operating manual.

Zero oxygen calibration

- 1 Prepare zero-oxygen water by dissolving 12.5 g AR grade anhydrous sodium sulfite in 250 mL distilled or deionized water.
- 2 Soak the measuring tip in zero-oxygen water and gently shake the probe to remove bubbles on the gas permeation membrane surface.
- 3 When the DO reading becomes stable, perform a zero-oxygen calibration according to the meter operating manual.

Measurement

- 1 Rinse the measuring tip with distilled water or deionized water.
- 2 Soak the measuring tip completely in sample solution until the temperature reading becomes stable,
- 3 Maintain a constant velocity by agitating the probe with a constant speed of 20–80 cm/s.
- 4 Ensure there are no air bubbles trapped on the membrane.
- 5 Measure the temperature and DO value of sample solution according to the meter operation manual.
- 6 Rinse the probe measuring tip and wipe dry.

Operating hints

- The measuring tip is made of ABS plastic and porous membrane. Ensure the sample solution will not damage the measuring tip before measurement.
- The probe consumes oxygen near the gas permeation membrane during measurement. Gently shake the probe or stir the sample solution to maintain a constant oxygen concentration during measurement.
- Do not touch the gas permeation membrane.
- During measurement or transportation, white powder or crystals may form on the DO probe. This will not affect the performance of the probe. Rinse away the crystals with tap or deionized water.
- Press the gas permeation membrane tightly against the tip of the cathode without bubbles or gaps to maintain measuring accuracy. If the membrane is damaged or dirty, replace it with a new probe membrane cap.
- Replace the DO filling solution frequently, depending upon the service condition of the probe. If the probe is used frequently, replace DO filling solution weekly.
- Do not apply force onto the probe cap, cable, or cable interface.
- Keep the cable interface dry.

Maintenance

- 1** Remove the probe membrane cap and wipe the cathode and anode dry.
- 2** Polish the cathode and anode gently with abrasive cloth.
- 3** Rinse the cathode and anode with distilled or deionized water.

Troubleshooting

Meters

Refer to the related section in the meter operating manual.

Probe

- DO reading is not stable.
 - a** Check if the membrane cap is screwed tight. The membrane should touch the platinum tip and there should be no gas bubble between the membrane and platinum.
 - b** Check whether the flow speed of solution is in the proper range.
 - c** During calibration or measurement, wait until the temperature and DO meter reading reach equilibrium.
 - d** Check whether there are air bubbles in the solution or there is an air bubble trapped on the membrane.
 - e** Check the time since last refilling or maintenance. If necessary, perform maintenance or refill DO filling solution.
- Reading is drifting.
 - a** If this happens for a new probe, just after maintenance, or if it has not been in use for a long time, perform more preconditioning.
 - b** Check whether the temperature of solution is the same as indicated on the meter. If not, wait more time to reach temperature equilibrium.

- Error in reading.
 - a** Not enough preconditioning.
 - b** Check that the calibration of the probe was performed properly (see “[Calibration](#)” on page 7).
 - c** Check the flow speed of the solution.
 - d** Check that the sample solution does not contain too many bubbles. Bubbles affect probe measurement.
- Reading is always zero.
 - a** Check that the probe is connected to the meter and the meter is on.
 - b** Check that the volume of the DO filling solution is not too low.

For any other problems during probe use, contact your Agilent Technologies customer service representative.

Storage

Short-term storage

Soak the measuring tip in distilled water or deionized water.

Long-term storage

- 1** Unscrew the probe membrane cap and empty the DO filling solution.
- 2** Rinse the probe membrane cap, cathode, and anode of the probe with distilled or deionized water.
- 3** Gently shake the membrane cap to dry.
- 4** Wipe the probe cathode and anode to dry.
- 5** Screw the dry membrane cap onto the probe body.
- 6** Place the probe in the storage box and store at ambient temperature and dry conditions.

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© Agilent Technologies, Inc. 2012