

## Eclipse XDB Family of Bonded Phases Designed for Optimum Lifetime at Intermediate pH

Application  
Technical  
Robert Ricker


During use at intermediate and basic pH (<sup>37</sup>), column packings can change. The change is generally believed to be the result of additional hydroxyl groups created as the silica is dissolved, rather than bonded phase hydrolysis. Compounds which react with silanols will exhibit retention time and peak-shape changes.

Below, the stability of two C8 columns is compared. One has a retention time increase for amitriptyline of 55% over a 125-day exposure to a phosphate buffer mobile phase with a pH of 7; whereas retention time on the other column changes by only 10%. There were two main differences between these two columns. The more quickly degraded material was a xerogel type of silica, and the other was a solgel type. The second difference was that the more quickly changing column was singly endcapped; whereas, the other column was doubly endcapped. Clearly, both columns "aged" during use, as exemplified by the change in the relative capacity factor ( $k'$ ) of the analyte.

To use an analogy, steel is painted to keep "corrosive agents" away from the surface. In HPLC, endcapping may be thought of as a coating (like paint), and double endcapping is a more effective protectant from the "corrosive" mobile phase. Thus, Agilent ZORBAX Eclipse XDB is designed for optimal lifetime in the intermediate pH range.

### Highlights

- Agilent ZORBAX Eclipse XDB-C8 outlives C8 bonded Xerogels.
- Eclipse is designed for long life at pH <sup>3</sup> 6-8.



Agilent Technologies

*Robert Ricker is an application chemist  
based at Agilent Technologies, Wilmington,  
Delaware.*

For more information on our products and  
services, visit our website at:  
[www.agilent.com/chem](http://www.agilent.com/chem)

Copyright© 2002 Agilent Technologies, Inc.  
All Rights Reserved. Reproduction,  
adaptation or translation without prior  
written permission is prohibited, except as  
allowed under the copyright laws.

Agilent shall not be liable for errors  
contained herein or for incidental or  
consequential damages in connection with  
the furnishing, performance, or use of this  
material.

Information, descriptions, and specifications  
in this publication are subject to change  
without notice.

Printed in the USA  
April 25, 2002  
5988-6458EN



**Agilent Technologies**