


Effect of Bonded-Phase Functionality on the Recovery of a Synthetic Lipopeptide

Application
Technical
Robert Ricker

This highly hydrophobic lipopeptide requires high temperature for separation. The recovery of the peptide (determined by protein assay) is higher on the more hydrophobic, long-chain alkyl bonded-phases (SB-C8 and SB-C18). The higher recovery is probably due to higher organic solvent concentration in the mobile phase, at elution. The lower polarity of the solvent is better for peptide solubility. Increasing recovery of hydrophobic peptides with a more hydrophobic bonded phase is counter to conventional wisdom.

Highlights

- The low pH, high-temperature compatible StableBond bonded phases allow choice of separation materials, ranging from highly hydrophobic (SB-C18) to polar (SB-CN).
- The broad range of StableBond functionalities permits the choice of separation media that **best suits the sample**.



Conditions: LC: Hewlett Packard 1090 Series II
Columns: 4.6 mm ID x 150 mm
Mobile Phase: 15 μ L (15 μ g) of peptide in 0.1% TFA / DMSO, 1 mL/min, 60°C
10-90% B in 40 min: A – 0.1% TFA / water, B – 0.1% TFA / ACN



Agilent Technologies

*Robert Ricker is an application chemist
based at Agilent Technologies, Wilmington,
Delaware.*

For more information on our products and
services, visit our website at:
www.agilent.com/chem

Copyright© 2002 Agilent Technologies, Inc.
All Rights Reserved. Reproduction,
adaptation or translation without prior
written permission is prohibited, except as
allowed under the copyright laws.

Agilent shall not be liable for errors
contained herein or for incidental or
consequential damages in connection with
the furnishing, performance, or use of this
material.

Information, descriptions, and specifications
in this publication are subject to change
without notice.

Printed in the USA
April 25, 2002
5988-6446EN



Agilent Technologies