Agilent I9121
Ammonia
Combination ISE
氨气敏复合电极

Operating Guide
用户手册

Agilent Technologies

Overview

The I9121 Ammonia Combination ISE can measure ammonia and ammonium concentration in aqueous solutions. It is used with the 3200I Ion Meter or similar meters.

WARNING Use this probe according to the operating manual to avoid personal injury.

WARNING The probe solution can cause chemical burns or illness if it is taken orally or contacted by human skin. Use protective clothing or gloves to avoid contact. In case of contact, rinse contacted area with tap water or deionized water thoroughly.

CAUTION The probe body material is glass. Handle with care to avoid damage to the instrument.

Specifications

Table 1 I9121 Ammonia Combination ISE specifications

Specification	Value
Concentration range	10^{-1} to 1×10^{-6} mol/L
Temperature range	0 to 50 °C
Probe impedance	$\leq 1500 \text{ M}\Omega$
Reference type	Ag/AgCl
Probe diameter	16 mm
Probe length	180 mm
Cable interface	BNC
Cable length	1000 mm


Figure 1 I9121 Ammonia Combination ISE outer body assembly


Figure 2 I9121 Ammonia Combination ISE assembly

Operation

Preparing the probe

- 1** Hold the probe upright and unscrew the probe cap (see [Figure 2](#)).
- 2** Remove the internal probe from the outer body. Rinse the internal probe and outer body with distilled or deionized water.
- 3** Remove the membrane cap off from the outer body (see [Figure 1](#)).
- 4** Remove the original composite membrane and O-ring from the membrane cap. Rinse the membrane cap with distilled water or deionized water.
- 5** Using tweezers, carefully place the membrane into the membrane cap. The membrane should evenly touch against the bottom of the membrane cap.
- 6** Set a new O-ring on the outer body and tightly screw the membrane cap onto the outer body. The composite membrane should be flat, and pressed against the bottom of the membrane cap without wrinkles.
- 7** Add 3.5 mL filling solution (5190-0544) of ammonia combination ISE into the outer body.
- 8** Tilt the probe slightly to insert it into the outer body. Be certain that the inner body is fully inserted to the bottom of the outer body and screw on the probe cap.
- 9** Pull the cable behind the probe cap backward slightly, and tighten the probe cap to the outer body. Allow the sensitive glass membrane at the tip of the internal probe to touch the composite membrane solely by the force of the spring.

Calibration

If the meter is an ion meter:

- 1 Choose more than two NH_4Cl calibration solutions. Calibrate the probe in these solutions in the order from low to high concentration.
- 2 Add appropriate amount of alkalizer (such as AR grade pure NaOH), dropwise into the calibration solution until its pH value is controlled in the range of 11-12.
- 3 Exchange the calibration solution during calibrating, rinse the measuring tip thoroughly and absorb away the water on the measuring tip.

If the meter is an mV meter (such as pH meter):

- 1 Calibrate the probe in more than two NH_4Cl calibration solutions in order from low to the high concentration.
- 2 Add appropriate amount of alkalizer (such as AR grade pure NaOH) dropwise into the calibration solution until its pH value is controlled in the range of 11-12.
- 3 Record mV readings in each calibration solution.
- 4 Plot mV- pNH_4 linear diagram or calculate the linear equation for the probe (pNH_4 is the negative logarithm of ammonium concentration).
- 5 Exchange the calibration solution during calibrating, and rinse the measuring tip thoroughly.
- 6 Absorb away the water on the measuring tip.
- 7 After calibration repeat probe preparation.

Measurement

- 1 Add appropriate amount of alkalizer, such as AR grade NaOH, to the sample solution. Soak the measuring tip in sample solution. Generally the pH value of the sample solution is controlled in the range of 11-12.
- 2 Record the readings of the meter when the readings become stable.
- 3 If the meter is an ion meter, read the concentration value of sample solution directly. If the meter is an mV meter (such as a pH meter), input the mV reading into the mV-p NH_4 linear diagram or the linear equation for the probe to calculate the pNH_4 value of the sample solution.
- 4 Rinse the measuring tip with deionized water. Absorb the water on the measuring tip.

Operating hints

- Before using the probe, refer to the appropriate measurement standards and select the correct measuring method.
- For measurement of samples with low concentration, use distilled water or deionized water that does not contain ammonium for making standard solution.
- Soak the internal probe in ammonia combination ISE filling solution for 24 hours before measurement.
- Once assembled, make sure the composite membrane is flat and sealed, with no leaks when the filling solution is added into the outer body.
- Keep the calibration solution and sample solution at a uniform temperature and flow rate by swinging the probe gently or turning on the magnetic stirrer.
- Keep the probe upright to prevent loss of filling solution. Replenish filling solution immediately in the event of a spill.
- Maintain constant temperature of solution, environment, and probe to obtain maximum accuracy.
- During calibration or measurement, make sure there are no bubbles in the surface of the composite membrane.
- Use the auxiliary device for the Ammonia Gas Sensing Probe, for accurate measurement of ammonia combination ISE.

- Exchange composite membrane and ammonia combination ISE filling solution if the probe response is too slow.
- If the temperature of the sample solution and the ammonium ion concentration of solution are very low, the response time will be prolonged.
- Dilute samples with high concentration (>500 ppm) before measurement.
- Rinse the composite membrane free of any substance, such as precipitation, that could clog the membrane. Do not wipe the composite membrane directly or you will damage the composite membrane.
- Take probe reading as soon as it becomes stable. Use a sealed container to decrease the evaporation of ammonia.
- After measurement, soak the measuring tip in 0.01 mol/L sulfuric acid solution or rinse it in distilled water or deionized water to remove excess ammonia diffused into the filling solution.

Maintenance

Exchange of reference filling solution

- 1 Siphon the probe filling solution and add fresh solution (5190-0545) until the level is 5 mm lower than the filling hole.
- 2 Repeat several times.

Cleaning of inorganics

Soak the measuring tip in 0.1 mol/L HCl or EDTA solution for 15 minutes.

Cleaning of organics

Soak the measuring tip in absolute ethyl alcohol, or other solvent that can dissolve organics for 15 minutes.

Cleaning of grease

Soak the measuring tip in warm, weakly alkaline detergent for 15 minutes.

Cleaning of protein precipitation

Soak the measuring tip in 0.1 mol/L HCl solution that contains 1% pepsin for 15 minutes.

Troubleshooting

Meters

Refer to the operating manual of the meter. Check all relevant parts, such as electrode, calibration solution, and samples.

Probe

- 1 Check the mV value after the calibration solution concentration has been changed ten times. The change of mV value should be greater than 52 mV. If not, clean the probe as described in the maintenance procedures.
- 2 Use effective calibration solution. Make sure that the distilled or deionized water used to prepare the calibration solution meets all requirements. Make sure the calibration solution is not contaminated or beyond its shelf life.

Method

Refer to relative standards to ensure the measuring methods are correct.

For any other problems during probe use, contact your Agilent Technologies customer service representative.

Storage

Short-term

Soak the internal probe in pH 4.01 calibration solution.

Long-term

- 1 Remove the probe cap and the membrane cap in sequence. Drain the filling solution in the ammonia combination ISE. Rinse with distilled water or deionized water.
- 2 Place the internal probe in the outer body and store the probe in dry condition.

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© Agilent Technologies, Inc. 2012