

## Large-Volume Injection for Gas Chromatography **Using COC-SVE**

Gas Chromatography

## Authors

Bill Wilson, Don Nixon, and Matthew S. Klee Agilent Technologies, Inc. 2850 Centerville Road Wilmington, DE 19808-1610 USA

## Abstract

Large-volume injection (LVI) is a valuable tool for trace analysis by gas chromatography (GC). New inlets and injection techniques have been developed in recent years that enable volumes up to 1 mL to be injected into a standard GC.

This application note describes LVI with a cool on-column (COC) inlet and solvent vapor exit (SVE) kit and is illustrated with examples of straight-chain hydrocarbons, polycyclic aromatic hydrocarbons (PAHs), and semivolatiles.

## Key Words

Large-volume injection, LVI, solvent vapor exit, SVE, gas chromatography, GC, PAH analysis, hydrocarbon analysis, semivolatile analysis.

## Introduction

The increasingly lower detection limits required by regulatory agencies have fueled a demand for more sensitive gas chromatography (GC) analysis. Large-volume injection (LVI) techniques are one response to this demand. Other approaches include sample concentration techniques, the use of more sensitive and selective detectors, and decreasing system noise. See "Large-Volume Injection for Gas Chromatography Using PTV" for a discussion of approaches to increase sensitivity.<sup>1</sup>

The typical injection volume for capillary column GC analysis is 0.5 to 2 µL. The Agilent 6890 and 5890 Series GCs allow approximately two times the normal injection volume (up to  $5 \,\mu$ L, depending on the solvent) using "pulsed" splitless injection. Injecting still larger volumes by standard techniques can lead to contamination of the system, irreproducible results, and loss of sample. In LVI, good chromatography is obtained with injection volumes of 5 to 500 uL or more. While putting a large volume of sample into the column is an intuitive and conceptually simple approach, special techniques are required to remove the

excess solvent. To be successful, LVI demands hardware modifications to remove the solvent and careful method development.

To inject large sample volumes, the bulk of the solvent in which the sample is dissolved must be selectively evaporated before transfer of the sample to the analytical column and the start of the analytical separation. There are two primary techniques to eliminate solvent: programmable temperature vaporization (PTV) and cool on-column (COC) injection with solvent vapor exit (SVE). LVI with PTV is recommended for samples with no early-eluting compounds of interest or for dirty samples. This is discussed elsewhere.

## COC-SVE

LVI by COC injection and SVE uses a COC inlet to introduce the sample. A precolumn assembly, consisting of a retention gap/precolumn combination, is placed ahead of the analytical column. A vent line is attached at the juncture of the precolumn and analytical column as shown in figure 1.

Liquid sample is injected into the precolumn. The solvent evaporates and





Figure 1. COC Injection with SVE-Hardware Configuration

exits through the vent. The oven temperature, column head pressure, and tower are set so there is always some liquid solvent in the precolumn. Sample components, including volatiles, are concentrated in this small volume of solvent. The time required to evaporate this solvent can be predicted with the software included in the COC-SVE kit. After evaporating the bulk of the solvent, the vent is closed and the normal oven temperature ramp is started. The remaining sample is then evaporated from the precolumn and moves to the analytical column for separation. This process allows solutes that elute very close to the tail of the solvent to be retained and analyzed.

An SVE kit is available with all the necessary columns, valves, fittings, and software (G2399A). The software in the COC-SVE kit simplifies method development. Enter the oven temperature, inlet pressure, outlet pressure, vent flow, and amount injected. When you select the solvent from the drop down list, the software automatically calculates the delay time. Figure 2 shows the COC-SVE top level screen.

The 6890 Series GC uses a standard automatic sampler fitted with up to 50-µL syringes to introduce LVI samples. Up to 25 µL can be injected with a single injection from the 50-µL syringe. The injection parameters are set in the Agilent ChemStation software.

With COC-SVE, there is an increased probability of contaminating the precolumns because such large volumes of sample are injected directly. Nonvolatile components remain in the precolumn and may influence subsequent injections. The precolumn may have to be cut or replaced frequently. COC-SVE is primarily used for clean samples such as drinking water extracts.

Table 1 lists the advantages and disadvantages of LVI by COC-SVE. See "Large-Volume Cool On-Column Injection using the G2399A Solvent Vapor Exit Kit" for further discussion and other examples of COC-SVE.<sup>2</sup>

### **Experimental**

A 6890 GC with electronic pneumatics control (EPC) was used for these experiments. Samples were injected with the G1916A automatic liquid sampler (ALS) with a G1513A injector. A G1876AA ChemStation (version 4.02) was used for instrument control and data acquisition.

Table 2 shows the hardware and software revisions that support LVI with COC-SVE.

The SVE configuration described in the introduction is shown in figure 1. The experimental conditions for the GC methods are given with each chromatogram.

-	SVE Val∨e time				
SVE Valve time					
Elir	nination Rate	28.1489	ul/min		
Valve Time		0.89	min		
Suggested Valve Time		0.8	min		
Oven Temperature		40	] C		
Inlet Pressure		2	] bar		
Outlet Pressure		1	bar		
Vent Flow		30	] ml/min		
Amount Injected		25	] ul		
Sol	Solvent Hexane 🛓				
Calculate SVE Valve Time					
Close					

Figure 2. COC-SVE Software

#### Table 1. Advantages and Disadvantages of LVI with COC-SVE

Advantages	Disadvantages	
Can be used with early-eluting compounds	Requires clean samples	
Minimizes degradation of unstable components	Less flexible technique	
Minimum sample discrimination	Complex to optimize for new solvents	

#### Table 2. Software, Hardware, and Firmware Versions that Support LVI

ltem	Software/Firmware
ChemStation software	A.04.02 or higher
G1513A injector	A.09.10
G1512A ALS controller	A.01.08
6890 Series GC	A.02.01 or higher

### **Results and Discussion**

Fairly large sample volumes can be injected using COC injection without SVE, as shown in figure 3, as long as a retention gap is used.

Figure 3A shows the chromatogram obtained from 10  $\mu$ L of a PAH sample without SVE, and figure 3B shows 20  $\mu$ L of the same sample. Twice the volume of sample was injected for the chromatogram in 3B as in 3A and yielded a twofold response with no deterioration of peak shape (see table 3 for the method parameters for these experiments). The ability to inject such large volumes by COC requires use of a very non-polar solvent such as isooctane.

Figure 3C shows what happens if additional sample is injected and the retention gap overfill ed. Injecting smaller amounts of a polar solvent, such as methanol, would have a similar effect. Using COC-SVE allows the injection of such large samples without overfilling the retention gap.

Figure 3C shows what happens if additional sample is injected and the retention gap overfilled. Injecting smaller amounts of a polar solvent, such as methanol, would have a similar effect. Using COC-SVE allows the injection of such large samples without overfilling the retention gap.



Figure 3. PAH Analysis by LV-COC injection without SVE peaks

Figure 4 shows a COC injection of 10  $\mu$ L of a C<sub>10</sub> to C<sub>14</sub> mixture without SVE. The C<sub>10</sub> is well separated from the solvent. To compare the recovery of LVI, 25- $\mu$ L of a 1:25 dilution of the same sample using SVE is also shown. The area counts of the C<sub>10</sub> and C<sub>12</sub> peaks of the two chromatograms are similar. This demonstrates the value of large-volume injections of dilute samples, which yields proportionally lower detection limits.

Notice the pronounced solvent impurity peaks with the larger injections. This is a universal weakness of all LVI techniques and highlights the need for meticulous lab technique and ultraclean solvents.

A valuable use of COC-SVE is in the analysis of semivolatiles in drinking water extracts. Figure 5A shows the chromatogram from an injection of  $1-\mu L$  of an undiluted test solution using a standard method (SVE valve closed during whole run). The test solution was diluted 1:25 and reanalyzed using a 25-µL injection and a 0.45-minute solvent vent delay time (see table 4 for other analytical parameters). The chromatogram is shown in figure 5B. By comparing figures 5A and 5B, it is clear that no chromatographic resolution is lost and that the peak heights for each solute are very similar.

# Table 3. PAH Analysis by LVI by COC Injection Without SVE Conditions for Figure 3A, 3B, and 3C

tem	Without SVE
nlet type	Cool on-column
njection volume	Dichloromethane, 10, 20, or 25 $\mu L$
njection speed	Slow, 3 sec postdwell time
Column	30 m x 0.25 mm id x 0.24 $\mu m$ HP-5 MS
Retention gap	5 m x 0.53 mm id
Oven temperature program	50 °C, 5 min; 10 °C/min to 320 °C
Carrier	1.5 mL/min hydrogen; 1.6 CF 66 kPa initial P



Figure 4.  $C_{10}$  to  $C_{14}$  in hexane by COC injection—comparison of 1-µL standard and 25-µL LVI injections



Figure 5. COC-SVE of an extract of semivolatiles in drinking water-for conditions, see table 3

Table 4.	Conditions	for Figure 5
----------	------------	--------------

Sample	Semivolatiles	
Solvent	Acetone	
Solute concentration	4 ppm	
Parameter		
Injection	Cool on-column	
Inlet temperature	Oven track	
Injection volume	25 μL	
Carrier gas	Hydrogen	
Column flow	2 mL/min	
Oven temperature 320 °C, (2 min)	45 °C (1 min), 20 °C/min to	
SVE on time	0 minute	
SVE off time	0.45 minute	
Detector	FID	
Detector temperature	320 °C	

#### **Conclusions**

LVI is a useful technique for lowering detection limits for GC analysis. LVI with COC-SVE is ideal for clean samples with early-eluting compounds. With careful method development, it is a reliable technique with broad applicability for sample applications that require more sensitivity than can be obtained with standard injection volumes. The COC-SVE kit provides the hardware and software necessary to be successful with performing the SVE technique.

#### References

- M. S. Klee, et al, "Large-Volume Injection for Gas Chromatography Using PTV," Agilent Technologies, Inc., Application Note 228-374, Publication (23) 5965-7770E, March 1997.
- F. David, et al, "Large-Volume Cool On-Column Injection Using the G2399A Solvent Vapor Exit Kit," Agilent Technologies, Inc., Application Note 228-378, Publication (23) 5965-7925E, March 1997.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

Copyright©2000 Agilent Technologies, Inc.

Printed in the USA 3/2000 5965-7923E

