

Agilent MassCode PCR Research Solution

cDNA Synthesis, PCR, and PCR Purification

Protocol

Version C, April 2012

For Research Use Only. Not for use in diagnostic procedures.



Notices

© Agilent Technologies, Inc. 2011

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G9200-90004

Edition

Version C, April 2012

Printed in USA

Agilent Technologies, Inc. 5301 Stevens Creek Rd Santa Clara, CA 95051 USA

Technical Support

For technical product support, contact Agilent at (800) 227-9770 or masscode.support@agilent.com.

Limited Use License for End Users

For Research Use Only. Not for use in diagnostic procedures.

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data - Commercial Items) and DFARS 227.720-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

In this Guide...

This document describes how to synthesize cDNA from a sample of nucleic acids, perform MassCode PCR with the cDNA and purify the MassCode tag labeled PCR products using the the MassCode cDNA Synthesis Kit and the MassCode PCR Reagent Kit.

If you have comments about this protocol, send an email to masscode.support@agilent.com.

1 Before You Begin

Consult the Agilent MassCode PCR Research Solution Quick Start Guide for an overview of the MassCode PCR application and a description of how the protocol provided in this document fits into the overall MassCode PCR workflow. Before starting an experiment, make sure you read and understand the information in this chapter and have the necessary equipment and reagents needed.

2 **Procedures**

This chapter provides instructions on synthesizing cDNA from a mixture of nucleic acids, setting up and running the PCR amplification of the MassCode targets, and purifying the MassCode PCR reactions to isolate the tagged products.

3 Quick Reference Protocols

This chapter contains quick reference protocols with streamlined instructions covering the cDNA synthesis, MassCode PCR amplification, and PCR purification procedures. For more detailed instructions, see Chapter 2, "Procedures".

Contents

1	Before You Begin 7
	Overview 8 Purpose of the protocol 8 Protocol procedure 8
	Important Information for Using the ProtocolPrecautions to avoid RNase contaminationContents of the kits10Storage conditions12Catalog information12Required equipment and supplies13
2	Procedures15Synthesizing First-Strand cDNA16Prepare for the calibrator reactions16Dilute the MS2 phage (if needed)17Anneal primers to template17Perform reverse transcription19Amplifying MassCode Targets by PCR21
	Design a PCR plate setup21Prepare the PCR reactions22Run the PCR program23
	Purifying the PCR Products 24 Prepare the reagents 24 Purify the PCR products 25
	Place the plate in the autosampler 28

10 10

Contents

3 Quick Reference Protocols 29

Synthesize first-strand cDNA30Amplify MassCode targets by PCR32Purify the PCR products33



MassCode cDNA Synthesis Kit and MassCode PCR Reagent Kit Protocol

Before You Begin

1

Overview 8 Purpose of the protocol 8 Protocol procedure 8 Important Information for Using the Protocol 10 Precautions to avoid RNase contamination 10 Contents of the kits 10 Storage conditions 12 Catalog information 12 Required equipment and supplies 13

Consult the *Agilent MassCode PCR Research Solution Quick Start Guide* for an overview of the MassCode PCR application and a description of how the protocol provided in this document fits into the overall MassCode PCR workflow.

Before starting an experiment, make sure you read and understand the information in this chapter and have the necessary equipment and reagents needed.



1 Before You Begin Overview

Overview

Agilent's MassCode PCR Research Solution is a high-throughput platform for the detection of target nucleic acids based on multiplex PCR and LC/MS detection. An overview of the MassCode PCR Research Solution is shown in Figure 1.

Purpose of the protocol

This document provides a protocol for performing the cDNA synthesis, PCR amplification and PCR purification steps of the MassCode PCR workflow (see Figure 1) using Agilent's MassCode cDNA Synthesis Kit, the PCR Reagent Kit, and a target panel of MassCode tagged primers.

Protocol procedure

The starting material for this protocol is a nucleic acid mixture extracted from a sample of cells. The nucleic acids (a mixture of RNA and DNA) are first reverse transcribed to convert the RNA into first-strand cDNA using random primers and MassCode reverse transcriptase.

The DNA/cDNA samples are then used as the template in a multiplex PCR reaction with a mix of primers labeled with the MassCode tags.

The labeled PCR products are then purified using a 96-well purification plate that contains a silica-based fiber matrix in each well. The PCR products are loaded onto the purification plate, washed, and then eluted into an autosampler plate that can go directly into the LC/MS. The purification procedure has been optimized to recover MassCode PCR products.





Important Information for Using the Protocol

Precautions to avoid RNase contamination

Ribonucleases, such as RNase A, are very stable enzymes responsible for RNA hydrolysis. RNase A can survive autoclaving and other standard methods of protein inactivation. Since RNases are present in the oils of the skin, wear gloves at all times when working with nucleic acid solutions that contain RNA. Use sterile tubes and micropipet tips handled with gloves only.

Contents of the kits

MassCode cDNA Synthesis Kit, 192 Reactions, Box #5190-3553				
MassCode Reverse Transcriptase 192 reactions				
RNase Block	105 μl			
10× MassCode RT Buffer	420 μl			
dNTP	168 μl			
Random Primers	660 μl			
RNase-Free Water 2 × 1.5 ml				

 Table 1
 Contents of the MassCode cDNA Synthesis Kit, Part Number 5190-3553

MassCode 2× PCR Master Mix Kit, 192 Reactions, Box #5190-3564			
MassCode 2× PCR Master Mix	192 reactions		
RNase-Free Water	2 × 1.5 ml		
MassCode PCR Purification Kit*, 192 Purific	ations, Box #5190-3605		
DNA Binding Solution	75 ml		
PCR 5× Wash Buffer	50 ml		
MassCode 96-well Binding Plates	2 plates		
Adhesive Plate Sealer	15 sheets		
96-well Waste Plates	2 plates		
Autosampler 96-well Plate Covers	2 covers		
Autosampler 96-well Plates	2 plates		
MassCode System Suitability Calibrator, Bo	x #5190-3743 (see Note below)		
20× System Suitability Calibrator 2 × 50 μl			
MassCode Internal Amplification Control, 1	92 Reactions, Box #5190-3744		
MS2 Phage 200 µl			
RNase-Free Water	2 × 1.5 ml		
KNase-Free Water	2 × 1.5 ml		

Table 2 Contents of the MassCode PCR Reagent Kit, Part Number 5190-3745

* You can purchase the PCR Purification Kit separately as a refill. See "Catalog information" on page 12.

NOTE

Use the System Suitability Calibrator for the system suitability test that is performed on the LC/MS. Refer to the *LC/MS Prep and Maintenance Guide* for instructions on preparing the system suitability calibrator vial.

Table 3 MassCode Target Panels and Primers Available from Agilent

MassCode Respiratory Pathogen Panel Kit	Contact Agilent for ordering information
Custom MassCode Oligo, 10 nmol	Part Number 750010

Storage conditions

MassCode cDNA Synthesis Kit

Store all components at -20°C upon receipt.

MassCode PCR Reagent Kit

Store the following boxes at -20°C upon receipt:

• MassCode 2× PCR Master Mix Kit, Box #5190-3564

After the initial thaw, you may store the $2 \times PCR$ Master Mix at $4^{\circ}C$ for up to 6 months or continue storing at $-20^{\circ}C$ for long term storage.

- MassCode System Suitability Calibrator, Box #5190-3743
- MassCode Internal Amplification Control, Box #5190-3744

MassCode PCR Purification Kit

Store the PCR Purification Kit, Box #5190-3605, at room temperature.

Catalog information

The MassCode cDNA Synthesis Kit and MassCode PCR Reagent Kit are part of Agilent's MassCode PCR Research Solution. Catalog information for the reagent kits that are part of this solution is described in Table 4. To find part numbers for the available MassCode target panels, refer to Table 3.

Part Number	Description
5500-0111	MassCode Respiratory Sample Prep Kit
5190-3553	MassCode cDNA Synthesis Kit
5190-3745	MassCode PCR Reagent Kit
5190-3605	MassCode PCR Purification Kit (refill)

 Table 4
 Reorder Part Numbers for MassCode Reagent Kits

Required equipment and supplies

Table 5	Required Equipment and Supplies
Description	
1.5-ml microo	centrifuge tubes
Microcentrifu	ıge
Agilent Sure	Cycler 8800 with 96-well block (p/n G8800A and G8810A)
MassCode ta	rget panel or primers (see Table 3 for list)
PCR reaction e.g. Agilent p	tube strips and cap strips for cDNA synthesis reactions (optional), art number 410082 and 410086
Centrifuge fo	r PCR reaction tube strips, if using
96-well PCR	plates, non-skirted, and cap strips, e.g. Agilent part numbers 401333 and 410096
Centrifuge fo	r 96-well plates
Multichanne	pipettors (10 μl, 20 μl, 200 μl, and 1 ml volumes)
Disposable p	ipetting reservoirs, sterile, e.g. VWR part number 89094-676
Benchtop wo	rking rack for 96-well plates, e.g. Agilent part number 410094
96- well vacu	um manifold with a deep well plate reservoir, e.g. Agilent part number 5185-5776
MS-grade, de	ionized water
100% ethano	

1 Before You Begin

Required equipment and supplies



MassCode cDNA Synthesis Kit and MassCode PCR Reagent Kit Protocol

Procedures

2

Synthesizing First-Strand cDNA 16 Prepare for the calibrator reactions 16 Dilute the MS2 phage (if needed) 17 Anneal primers to template 17 Perform reverse transcription 19 Amplifying MassCode Targets by PCR 21 Design a PCR plate setup 21 Prepare the PCR reactions 22 Run the PCR program 23 Purifying the PCR Products 24 Prepare the reagents 24 Purify the PCR products 25 Place the plate in the autosampler 28

This chapter provides instructions on synthesizing cDNA from a mixture of nucleic acids, setting up and running the PCR amplification of the MassCode targets, and purifying the MassCode PCR reactions to isolate the tagged products.



Synthesizing First-Strand cDNA

The cDNA synthesis step converts any RNA that is present in a nucleic acid mixture into first-strand cDNA.

Prepare for the calibrator reactions

In addition to setting up cDNA synthesis reactions with each of your nucleic acid test samples, set up three no-template calibrator reactions (CAL-NTC) and one internal amplification calibrator reaction (CAL-IACRNA). These reactions are described below.

The No-Template Calibrator Reactions (CAL-NTC)

In the CAL-NTC cDNA synthesis reactions, the nucleic acid sample is replaced with buffer (use the same buffer that your nucleic acid test samples are stored in). Since the reactions contain no nucleic acids, no PCR products are amplified during the subsequent PCR step.

The CAL-NTC reactions are used during LC/MS analysis to help set the cutoff values for target detection in the MassCode PCR software.

The Internal Amplification Calibrator Reaction (CAL-IACRNA)

The CAL-IACRNA cDNA synthesis reaction uses a nucleic acid sample prepared from MS2 phage alone (with no experimental sample added to the extraction reaction). The MS2 phage cDNA serves as the template for a phage-specific target during PCR. The abundance of this target is used during LC/MS analysis to help set the cutoff values for target detection in the MassCode PCR software. If you did not prepare an extraction reaction with MS2 phage alone, as instructed in the Respiratory Sample Preparation Kit Protocol, set up a cDNA synthesis reaction that includes whole MS2 phage particles (provided in the MassCode PCR Reagent Kit). The phage particles will lyse during the cDNA synthesis reaction and the genomic RNA will be accessible for reverse transcription.

Note that amplification of the phage-specific target also occurs in the test samples as an internal amplification control (IAC). The IAC verifies that the extraction, reverse transcription and PCR steps were successful for that sample.

Dilute the MS2 phage (if needed)

If your nucleic acid samples already contain MS2 phage RNA, proceed to "Anneal primers to template" on page 17. If you did not add MS2 phage to your samples prior to extraction, dilute the MS2 phage for addition into the cDNA synthesis reactions.

1 In a 1.5-ml microcentrifuge tube, combine 99 μ l of RNase-free water with 1 μ l of MS2 phage.

MS2 phage and RNase-free water are provided in the PCR Reagent Kit (Box #5190-3744, MassCode Internal Amplification Control).

2 Mix the diluted phage well and keep on ice until needed.

Do not store the diluted phage; prepare it fresh with each use.

When setting up the primer annealing reactions for your test samples, add 2 μl of the diluted MS2 phage to each reaction.

NOTE Adding the MS2 phage to your experimental samples prior to nucleic acid extraction is preferable to adding it to the cDNA synthesis reactions. Adding the phage prior to extraction controls for possible nucleic acid loss that may occur during the extraction procedure.

Anneal primers to template

The primer annealing reactions are incubated in the SureCycler 8800. Set up the reactions in PCR reaction tube strips or on a PCR reaction plate.

1 Set up the primer annealing reaction tubes for the nucleic acid test samples.

For each nucleic acid test sample, combine the components listed in Table 6. The table lists component volumes for reactions containing 4 μ l or 10 μ l of nucleic acid sample. A volume of 4 μ l of nucleic acid sample is the minimum quantity we recommend. A volume of 10 μ l of nucleic acid sample is optimal for maximizing the sensitivity of the assay.

2 Procedures

Anneal primers to template

	Volume per reaction			
Component	If adding 4 μI nucleic acid	If adding 10 μI nucleic acid		
Nucleic acid sample	4 μl	10 µl		
Random primers	3 μl	3 μΙ		
RNase-free water	8.7 μl*	2.7 μl*		
Total volume	15.7 μl	15.7 μl		

Table 6	Primer	Annealing	Reagent	Mixture [•]	for Test	t Sample	Reactions
---------	--------	-----------	---------	----------------------	----------	----------	-----------

 * If needed, add 2 μ I of the diluted MS2 phage to each reaction and reduce the volume of water by 2 μ I.

2 Set up the primer annealing reaction tubes for the three CAL-NTC reactions.

For each CAL-NTC reaction, combine the components listed in Table 7. The table lists component volumes for reactions containing 4 μ l or 10 μ l of nucleic acid storage buffer. The volume of storage buffer in the CAL-NTC reactions needs to be the same as the volume of nucleic acid used for the test sample reactions in step 1.

	Volume per reaction			
Component	If adding 4 μI storage buffer	If adding 10 μI storage buffer		
Nucleic acid storage buffer	4 µl	10 µl		
Random primers	3 µl	3 μl		
RNase-free water	8.7 μl	2.7 μl		
Total volume	15.7 μl	15.7 μl		

 Table 7
 Primer Annealing Reagent Mixture for the CAL-NTC Reactions

3 Set up the primer annealing reaction tubes for the CAL-IACRNA reaction.

Combine the components listed in Table 8. The table lists component volumes for reactions containing MS2 phage nucleic acids or diluted MS2 phage (see "Dilute the MS2 phage (if needed)" on page 17).

	Volume per reaction			
Component	If using MS2 nucleic acids	If using whole MS2 phage		
MS2 phage nucleic acid sample	4 μl	_		
Diluted MS2 phage	_	2 μl		
Random primers	3 μl	3 μl		
RNase-free water	8.7 μl	10.7 μl		
Total volume	15.7 μl	15.7 μl		

 Table 8
 Primer Annealing Reagent Mixture for the CAL-IACRNA Reaction

4 Incubate the primer annealing reactions at 70°C for 5 minutes in the SureCycler 8800.

We recommend bringing the SureCycler 8800 up to temperature and then pausing the program before adding the reaction tubes.

5 Remove the reactions from the SureCycler 8800, spin them briefly in a microcentrifuge, and place them directly on ice.

Perform reverse transcription

1 Prepare the reverse transcription reagent mixture by combining the components listed in Table 9.

The volumes listed are sufficient for one reaction. Prepare enough reagent mixture for all your reactions plus *at least* one reaction volume excess. Keep all reactions and reagents on ice.

Component	Volume per reaction
10× MassCode RT Buffer	2 μl
dNTP Mix	0.8 µl
RNase Block	0.5 μl
MassCode Reverse Transcriptase	1 μl

 Table 9
 Reverse Transcription Reagent Mixture

CAUTION

Perform reverse transcription

- **2** Add 4.3 μ l of the reverse transcription reagent mixture to each primer annealing reaction. Mix the reactions well (do not vortex).
- **3** Incubate the reactions in the SureCycler 8800 according to the program listed in Table 10.

Table 10	cDNA	Synthesis	Incubation	Program
----------	------	-----------	------------	---------

Incubation Step	Temperature	Duration
Primer Annealing	25°C	10 minutes
1st Strand Synthesis	42°C	60 minutes
Heat Inactivation of Reverse Transcriptase	70°C	15 minutes
Storage	4°C	Hold indefinitely

Do not let the cDNA synthesis reactions incubate at 4°C overnight.

4 When the program is complete, remove the reactions from the SureCycler 8800.

You can now perform the MassCode PCR amplification protocol (see "Amplifying MassCode Targets by PCR" on page 21) or store the cDNA at -80° C for later use.

2

Amplifying MassCode Targets by PCR

The PCR amplification step amplifies target sequences using a multiplex primer mix. The primers in this mix are each labeled with a unique chemical moiety known as a MassCode tag that will later be used for target identification by mass spectrometry. The primer mix can be of a custom design or purchased as part of an Agilent MassCode target panel (see Table 3).

Design a PCR plate setup

The PCR is carried out in a 96-well PCR plate (not provided). Perform PCR on your test cDNA samples, the CAL-IACRNA cDNA sample and the three CAL-NTC cDNA samples. Additionally, set up PCR reactions for any positive template calibrator (PTC) samples that were included with your MassCode target panel. The respiratory panel includes 4 PTC samples; run duplicate PCR reactions for each PTC. Your plate needs a total of 12 calibrator reactions in order for the MassCode PCR software to properly analyze the data.

	1	2	3	4	5	6	7	8	9	10	11	12
Α	CAL-NTC	CAL-NTC	CAL- PTC1	CAL- PTC2	CAL- PTC3	CAL- PTC4	CAL- IACRNA	CAL-NTC	CAL- PTC1	CAL- PTC2	CAL- PTC3	CAL- PTC4
в	Unknown reaction 1	Unknown reaction 2	Unknown reaction 3	Unknown reaction 4	Unknown reaction 5	Unknown reaction 6	Unknown reaction 7	Unknown reaction 8	Unknown reaction 19	Unknown reaction 10	Unknown reaction 11	Unknown reaction 12
с	Unknown reaction 13	Unknown reaction 14	Unknown reaction 15	Unknown reaction 16	Unknown reaction 17	Unknown reaction 18	Unknown reaction 19	Unknown reaction 20	Unknown reaction 21	Unknown reaction 22	Unknown reaction 23	Unknown reaction 24
D	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown
	reaction 25	reaction 26	reaction 27	reaction 28	reaction 29	reaction 30	reaction 31	reaction 32	reaction 33	reaction 34	reaction 35	reaction 36
Е	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown
	reaction 37	reaction 38	reaction 39	reaction 40	reaction 41	reaction 42	reaction 43	reaction 44	reaction 45	reaction 46	reaction 47	reaction 48
F	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown
	reaction 49	reaction 50	reaction 51	reaction 52	reaction 53	reaction 54	reaction 55	reaction 56	reaction 57	reaction 58	reaction 59	reaction 60
G	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown
	reaction 61	reaction 62	reaction 63	reaction 64	reaction 65	reaction 66	reaction 67	reaction 68	reaction 69	reaction 70	reaction 71	reaction 72
н	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown	Unknown
	reaction 73	reaction 74	reaction 75	reaction 76	reaction 77	reaction 78	reaction 79	reaction 80	reaction 81	reaction 82	reaction 83	reaction 84

Figure 2 An example of a plate setup for PCR

Prepare the PCR reactions

An example of a recommended plate setup is shown in Figure 2. The calibrator reactions are set up in row A of the plate. In this example, the target panel includes 4 CAL-PTC samples and each is set up in duplicate. The three CAL-NTC reactions are in wells A1, A2 and A8. The PCR reaction with the CAL-IACRNA sample is set up in a single well (A7). The test reactions are set up in rows B through H. With this setup, you can test up to 84 test samples in a single assay.

Prepare the PCR reactions

1 Prepare the PCR reagent mixture by combining the components in Table 11. Keep all reactions and reagents on ice.

The volumes listed are sufficient for one reaction. Prepare enough reagent mixture for *at least* one reaction volume excess.

Component	Volume per reaction
RNase-free water	3.2 μl*
MassCode 2× PCR Master Mix	10 µl
MassCode primer mix	$2.8\mu l$ (250 nM each primer; see Note below)

Table 11PCR Reagent Mixture

* 3.2 μ l of water is appropriate for PCR reactions containing 4 μ l of cDNA. If you need to maximize the sensitivity of the assay, you can add as much as 7.2 μ l of cDNA and reduce the volume of water accordingly.

NOTE

2.8 μ l is the recommended volume for a 28-plex primer mix (e.g. the primer mix provided in the Respiratory Pathogen Panel Kit). If you are designing your own MassCode primers, maintain a concentration of 250 nM for each primer in the final PCR reactions. For primers with a stock concentration of 100 mM, you need to add 0.05 μ l of each primer (or 0.1 μ l of each primer set) per reaction.

Note that you may need to adjust the volume of water in the PCR reagent mixture to obtain a final volume of 20 μ l in each PCR reaction.

2 Pipet the reagent mixture up and down to mix, then distribute 16 μ l to each well on the PCR plate.

3 Add 4 μl of cDNA or PTC to each well.See "Design a PCR plate setup" on page 21 for guidelines.

Run the PCR program

1 Place the PCR plate in the SureCycler 8800. If you are using the respiratory pathogen panel, run the program in Table 12. If you using a custom target panel, run a PCR program suitable for your primer mix.

Segment	Number of Cycles	Temperature	Duration
1	1	95°C	3 minutes
2	10	95°C	20 seconds
		65–56°C (temperature increment of –1°C/cycle)	50 seconds
3	15	95°C	20 seconds
		55°C	50 seconds
4	10	95°C	20 seconds
		55°C	80 seconds
5	10	95°C	20 seconds
		55°C	2 minutes
6	1	72°C	5 minutes

Table 12 PCR Cycling Program

2 Remove the plate from the SureCycler 8800.

You can now perform the PCR purification protocol (see "Purifying the PCR Products" on page 24) or store the plate at -20 °C for later use.

NOTE

The MassCode tags on the PCR products are light sensitive. Minimize their exposure to light.

Purifying the PCR Products

PCR purification removes excess primers, nucleotides, buffer components, and enzymes. It also removes any primer dimers that may have formed during cycling. Eliminating these excess reagents and primer dimers is important for reducing background signals in the mass spectrometry data.

The MassCode PCR Purification Kit is included as part of the PCR Reagent Kit, but you can purchase it separately as a refill (p/n 5190-3745).

Prepare the reagents

DNA Binding Solution (also used as the wash #1 solution)

- 1 Prepare 100 ml of 30% ethanol by combining 30 ml of 100% ethanol with 70 ml of MS-grade dH_20 .
- **2** Mix 75 ml of the provided DNA Binding Solution with 75 ml of 30% ethanol.

Store the prepared DNA Binding Solution tightly sealed at room temperature.

PCR Wash Buffer (the wash #2 solution)

- Add 200 ml of 100% ethanol to the stock bottle of 5× PCR Wash Buffer. The addition of ethanol brings the PCR wash buffer to a 1× concentration.
- **2** Using a permanent marker, mark the "ADD ETHANOL BEFORE USE" check box printed on the bottle.

Store the 1× PCR Wash Buffer tightly sealed at room temperature.

80% Ethanol (the wash #3 solution)

• Mix 160 ml of 100% ethanol with 40 ml of MS-grade dH_20 .

Store the 80% ethanol at room temperature.

Purify the PCR products

NOTE

If you are not using all 96 wells of the binding plate, you can cover the unused portion of the plate with adhesive plate sealer and use those wells in a later experiment.

Immobilize the PCR products on the binding plate matrix

- **1** Add 180 μl of the prepared DNA Binding Solution (refer to "Prepare the reagents" on page 24) to each 20-μl PCR reaction on the PCR plate. Pipet up and down to mix.
- **2** Place a MassCode 96-well binding plate on a vacuum manifold that contains a deep well plate reservoir. Using a multichannel pipettor, transfer the PCR product/DNA Binding Solution mixtures into the wells of the binding plate. Perform the transfer in a way that maintains the original plate setup configuration.
- **3** Perform the vacuum/centrifugation procedure on the binding plate according to the table below:

St	tep	Details/Comments			
a	Apply 400 mbar of vacuum to the binding plate until each well is dry (approximately 1 minute) and then continue for 1 minute more.	 The vacuum force may decrease as the wells become dry. 			
b	Release the vacuum and remove the binding plate from the manifold. Blot the bottom of the binding plate to remove excess liquid.	 To blot, tip the plate at an approximately 45° angle and touch just one corner of the plate to a paper towel. 			
C	Place the binding plate on top of a 96-well waste plate (provided). Cover the top of the binding plate with adhesive plate sealer and then adhere the two plates together.	 To adhere the plates, wrap them together with a rubber band or tape them together along the edges. 			
d	Spin the stacked plates in a room temperature centrifuge for 3 minutes at $1500 \times g$.	 If any liquid remains in the wells at the end of the spin, perform the spin again. 			
e	Disassemble the stacked plates.	 The PCR products are retained in the matrix of the binding plate wells. Set aside the waste plate. You will use it again during the wash procedure. 			

Vacuum/Centrifugation Procedure for Immobilizing the PCR Products

Purify the PCR products

Wash the PCR products on the matrix

1 Remove the plate sealer from the binding plate and add 500 μ l of the prepared DNA Binding Solution to each well.

This step is wash #1.

2 Perform the vacuum/centrifugation procedure on the binding plate according to the table below:

Vacuum/Centrifugation Procedure for Washing the PCR Products

St	tep	Details/Comments			
а	Apply 400 mbar of vacuum to the binding plate until each well is dry (approximately 1 minute) and then continue for 1 minute more.	 The vacuum force may decrease as the wells become dry. 			
b	Release the vacuum and remove the binding plate from the manifold. Blot the bottom of the binding plate to remove excess liquid.	 To blot, tip the plate at an approximately 45° angle and touch just one corner of the plate to a paper towel. 			
C	Place the binding plate on top of the waste plate. Cover the top of the binding plate with adhesive plate sealer and then adhere the two plates together.	 To adhere the plates, wrap them together with a rubber band or tape them together along the edges. 			
d	Spin the stacked plates in a room temperature centrifuge at $1500 \times g$.	 For wash #1 and wash #2, spin for 3 minutes. For wash #3, spin for 5 minutes. 			
e	Open the vacuum manifold and discard the wash solution from the waste plate. Replace the waste plate back on the vacuum manifold	• This step can be performed while the plates are spinning in the centrifuge.			
f	Remove the stacked plates from the centrifuge and disassemble.	• The PCR products are retained in the matrix of the binding plate wells.			

3 Remove the plate sealer from the binding plate and add 700 μ l of 1× **PCR wash buffer** (see "Prepare the reagents" on page 24) to each well.

This step is wash #2.

4 Repeat the vacuum/centrifugation procedure for washing the PCR products outlined in step 2 above.

- 5 Remove the plate sealer from the binding plate and add 700 μl of 80% ethanol (see "Prepare the reagents" on page 24) to each well.
 This step is wash #3.
- **6** Repeat the vacuum/centrifugation procedure for washing the PCR products outlined in step 2 above.

Elute the PCR products into the LC/MS autosampler plate

- 1 Place the binding plate on top of a 96-well plate that is suitable for use in the LC/MS autosampler (provided in the MassCode PCR Purification Kit; p/n 5190-3728). Check that the two plates are oriented in the same direction (i.e. well A1 of the binding plate is sitting on top of well A1 of the autosampler plate).
- **2** Remove the plate sealer from the binding plate and add 70 μl of MS-grade dH₂0 directly on top of the fiber matrix of each well.

If you want to run the plate in the LC/MS twice, increase the volume of MS-grade dH_20 to 140 µl. Note that eluting in a higher volume decreases sensitivity of the assay.

- **3** Replace the plate sealer on the binding plate. Tape the two plates together along the edges or wrap them together with a rubber band. Allow the plates to incubate at room temperature for 2 minutes.
- **4** Spin the stacked plates in a room temperature centrifuge at $1500 \times g$ for 5 minutes.
- **5** Remove the plates from the centrifuge and cover the autosampler plate with an autosampler plate cover (provided).
- **6** Place the plate in the autosampler (instructions below) to begin the LC/MS analysis or store it at 4°C (for up to 1 week) or -20°C (for long term storage).

2 **Procedures**

Place the plate in the autosampler

Place the plate in the autosampler

- Place the autosampler plate in the front position of the LC/MS autosampler (2 positions for plates: one in front and one in back).
- Place the plate so that row A is in the back and row H faces the front.
- If you want to turn off the light inside the autosampler, see the *Online Help* of the MassCode PCR program. (The default setting is On.)

You are now ready to set up the plate and run parameters with the MassCode PCR program.



MassCode cDNA Synthesis Kit and MassCode PCR Reagent Kit Protocol

Quick Reference Protocols

Synthesize first-strand cDNA 30 Amplify MassCode targets by PCR 32 Purify the PCR products 33

3

This chapter contains quick reference protocols with streamlined instructions covering the cDNA synthesis, MassCode PCR amplification, and PCR purification procedures. For more detailed instructions, see Chapter 2, "Procedures".



Synthesize first-strand cDNA

- **1** If you will be adding MS2 phage to the cDNA synthesis reactions, prepare a 1:100 dilution of the provided stock.
 - **a** In 1.5-ml microcentrifuge tube, combine 99 μ l of RNase-free water with 1 μ l of MS2 phage.
 - **b** Mix the diluted phage well and keep on ice until needed.
- **2** Set up the primer annealing reactions in PCR reaction tubes for the nucleic acid test samples, the three CAL-NTC reactions, and the CAL-IACRNA reaction.
 - For each test sample reaction, combine the components in Table 13.

 Table 13
 Primer Annealing Reagent Mixture for Test Sample Reactions

	Volume per reaction				
Component	If adding 4 μI nucleic acid	If adding 10 μI nucleic acid			
Nucleic acid sample	4 μl	10 µl			
Random primers	3 μl	3 μl			
RNase-free water	8.7 μl*	2.7 μl*			
Total volume	15.7 μl	15.7 μl			

* If needed, add 2 µl of the diluted MS2 phage to each reaction and reduce the volume of water by 2 µl.

• For each CAL-NTC reaction, combine the components in Table 14.

Table 14 Primer Annealing Reagent Mixture for the CAL-NTC Reactions

	Volume per reaction			
Component	If adding 4 μI storage buffer	If adding 10 μI storage buffer		
Nucleic acid storage buffer	4 μl	10 µl		
Random primers	3 μl	3 μl		
RNase-free water	8.7 μl	2.7 μl		
Total volume	15.7 μl	15.7 μl		

٠	For	the	CAL-	IACRNA	reaction,	combine	the	components	in	Table	15.
---	-----	-----	------	--------	-----------	---------	-----	------------	----	-------	-----

	Volume per reaction				
Component	If using MS2 nucleic acids	If using whole MS2 phage			
MS2 phage nucleic acid sample	4 μl	—			
Diluted MS2 phage	_	2 μl			
Random primers	3 μl	3 μl			
RNase-free water	8.7 μl	10.7 μl			
Total volume	15.7 μl	15.7 μl			

Table 15 Primer Annealing Reagent Mixture for the CAL-IACRNA Reaction

- **3** Incubate the reactions at 70°C for 5 minutes in the SureCycler 8800.
- **4** Remove the reactions from the cycler, spin them briefly in a microcentrifuge and place directly on ice.
- **5** Prepare the reverse transcription reagent mixture by combining the components listed in Table 16. Keep all reactions and reagents on ice.

 Table 16
 Reverse Transcription Reagent Mixture

Component	Volume per reaction
10× MassCode RT Buffer	2 μl
dNTP Mix	0.8 µl
RNase Block	0.5 μl
MassCode Reverse Transcriptase	1 µl

6 Add 4.3 μ l of the reverse transcription reagent mixture to each primer annealing reaction. Mix the reactions well (do not vortex).

3 Quick Reference Protocols

Amplify MassCode targets by PCR

7 Incubate the reactions in the SureCycler 8800 using the program listed in Table 17.

Incubation Step	Temperature	Duration
Primer Annealing	25°C	10 minutes
1st Strand Synthesis	42°C	60 minutes
Heat Inactivation of Reverse Transcriptase	70°C	15 minutes
Storage	4°C	Hold indefinitely

 Table 17
 cDNA Synthesis Incubation Program

8 When the program is complete, promptly remove the reactions from the cycler. You can now amplify MassCode targets by PCR or store the cDNA at -80°C for later use.

Amplify MassCode targets by PCR

1 Prepare the PCR reagent mixture by combining the components in Table 18. Prepare enough reagent mixture for all reactions plus *at least* one reaction volume excess.

 Table 18
 PCR Reagent Mixture

Component	Volume per reaction
RNase-free water	3.2 μl*
MassCode 2× PCR Master Mix	10 µl
Target Panel Primer Mix	2.8 μl**

 * 3.2 μl of water is appropriate for PCR reactions containing 4 μl of cDNA.

**2.8 μl of primer mix per reaction is appropriate for the respiratory pathogen panel.

- **2** Vortex the mixture, then distribute 16 μ l to each well on the PCR plate.
- **3** Add 4 μ l of cDNA or PTC to each well. See "Design a PCR plate setup" on page 21 for guidelines.

4 Place the plate in the SureCycler 8800. For the respiratory target panel, run the program in Table 19. For a custom target panel, run a program suitable for your primer mix.

Segment	Number of Cycles	Temperature	Duration
1	1	95°C	3 minutes
2	10	95°C	20 seconds
		65–56°C (temperature increment of –1°C/cycle)	50 seconds
3	15	95°C	20 seconds
		55°C	50 seconds
4	10	95°C	20 seconds
		55°C	80 seconds
5	10	95°C	20 seconds
		55°C	2 minutes
6	1	72°C	5 minutes

Table 19 PCR Cycling Program

5 Remove the plate from the SureCycler 8800. You can now purify the PCR products or store the plate at -20°C for later use.

Purify the PCR products

- **1** Add 180 μl of the prepared DNA Binding Solution (refer to "Prepare the reagents" on page 24) to each 20-μl PCR reaction on the PCR plate. Pipet up and down to mix.
- **2** Place a MassCode 96-well binding plate on a vacuum manifold that contains a deep well plate reservoir. Using a multichannel pipettor, transfer the PCR product/DNA Binding Solution mixtures into the wells of the binding plate. Perform the transfer in a way that maintains the original plate setup configuration.

Purify the PCR products

3 Perform the vacuum/centrifugation procedure on the binding plate according to the table below:

Vacuum/Centrifugation Procedure for Immobilizing the PCR Products

Step		Details/Comments	
a	Apply 400 mbar of vacuum to the binding plate until each well is dry (approximately 1 minute) and then continue for 1 minute more.	 The vacuum force may decrease as the wells become dry. 	
b	Release the vacuum and remove the binding plate from the manifold. Blot the bottom of the binding plate to remove excess liquid.	 To blot, tip the plate at an approximately 45° angle and touch just one corner of the plate to a paper towel. 	
C	Place the binding plate on top of a 96-well waste plate (provided). Cover the top of the binding plate with adhesive plate sealer and then adhere the two plates together.	 To adhere the plates, wrap them together with a rubber band or tape them together along the edges. 	
d	Spin the stacked plates in a room temperature centrifuge for 3 minutes at $1500 \times g$.	 If any liquid remains in the wells at the end of the spin, perform the spin again. 	
e	Disassemble the stacked plates.	 The PCR products are retained in the matrix of the binding plate wells. Set aside the waste plate. You will use it again during the wash procedure. 	

4 Remove the plate sealer from the binding plate and add 500 µl of the prepared DNA Binding Solution (refer to "Prepare the reagents" on page 24) to each well.

This step is wash #1.

5 Perform the vacuum/centrifugation procedure on the binding plate according to the table below:

Vacuum/Centrifugation Procedure for Washing the PCR Products

Step		Details/Comments	
a	Apply 400 mbar of vacuum to the binding plate until each well is dry (approximately 1 minute) and then continue for 1 minute more.	 The vacuum force may decrease as the wells become dry. 	
b	Release the vacuum and remove the binding plate from the manifold. Blot the bottom of the binding plate to remove excess liquid.	 To blot, tip the plate at an approximately 45° angle and touch just one corner of the plate to a paper towel. 	
C	Place the binding plate on top of the waste plate. Cover the top of the binding plate with adhesive plate sealer and then adhere the two plates together.	 To adhere the plates, wrap them together with a rubber band or tape them together along the edges. 	
d	Spin the stacked plates in a room temperature centrifuge at $1500 \times g$.	 For wash #1 and wash #2, spin for 3 minutes. For wash #3, spin for 5 minutes. 	
e	Open the vacuum manifold and discard the wash solution from the waste plate. Replace the waste plate back on the vacuum manifold	• This step can be performed while the plates are spinning in the centrifuge.	
f	Remove the stacked plates from the centrifuge and disassemble.	• The PCR products are retained in the matrix of the binding plate wells.	

- 6 Remove the plate sealer from the binding plate and add 700 μl of 1× PCR wash buffer (see "Prepare the reagents" on page 24) to each well. This step is wash #2.
- 7 Repeat the vacuum/centrifugation procedure for washing the PCR products outlined in step 5 above.
- 8 Remove the plate sealer from the binding plate and add 700 μl of 80% ethanol (see "Prepare the reagents" on page 24) to each well.
 This step is wash #3.
- **9** Repeat the vacuum/centrifugation procedure for washing the PCR products outlined in step 5 above.

- **10** Place the binding plate on top of a 96-well autosampler plate (provided). Check that the two plates are oriented in the same direction.
- 11 Remove the plate sealer from the binding plate and add 70 μ l of MS-grade dH₂0 directly on top of the fiber matrix of each well.

If you want to run the plate in the LC/MS twice, increase the volume of MS-grade dH_20 to 140 μl (note that this decreases sensitivity of the assay).

- **12** Replace the plate sealer on the binding plate. Tape the two plates together along the edges or wrap them together with a rubber band. Allow the plates to incubate at room temperature for 2 minutes.
- 13 Spin the stacked plates in a room temperature centrifuge at $1500 \times g$ for 5 minutes.
- **14** Remove the plates from the centrifuge and cover the autosampler plate with an autosampler plate cover (provided).

The purified PCR products are in the bottom of each well of the autosampler plate. You can now place the plate in the LC/MS autosampler for mass spectrometry analysis or store the covered plate at 4° C (for up to 1 week) or -20° C (for long term storage).

- **15** Place the covered autosampler plate in the LC/MS autosampler or store it at 4°C (for up to 1 week) or -20°C (for long term storage).
 - To place the plate in the autosampler, place it in the front position of the autosampler so that row A is in the back and row H faces the front.
 - To store the plate, place it at 4°C (for up to 1 week) or -20°C (for long term storage).

www.agilent.com

In This Book

This document describes how to synthesize cDNA from a sample of nucleic acids, perform MassCode PCR with the cDNA and purify the MassCode tag labeled PCR products using the the MassCode cDNA Synthesis Kit and the MassCode PCR Reagent Kit.

 ${\ensuremath{\mathbb C}}$ Agilent Technologies, Inc. 2011

Version C, April 2012



G9200-90004

