

Agilent MSD Productivity ChemStation G1701 & G1710

Performance Report Parameters



Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G1701-90101

Edition

First Edition

March 30, 2012

Printed in USA

Agilent Technologies, Inc. 5301 Stevens Creek Blvd. Santa Clara, CA 95051

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user cus-tomers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Soft-ware) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data -Com-mercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Contents

Introduction	4
Report Layout	5
Heading	5
RTE Integration Parameters	5
Graphics Area	5
Noise characteristics	5
Performance Report Results	6

1. Introduction

This document describes the calculations used by the Performance Report in the Agilent MSD Productivity ChemStation. These calculations are dependent on the integrator selected for the Data Analysis. The MSD Productivity ChemStation allows use of the Chemstation integrator and the RTE integrator.

The report is created by selecting Performance Report from the Chromatogram menu. The user is then presented with a dialog and asked to enter a retention time range where the noise will be calculated. A Total Ion Chromatogram (TIC) of this specified noise range is displayed in a new window and the report shown below is printed:



2. Report Layout

The report has several sections as described below. Calculations for noise are located in the Noise characteristics section. Calculations relating to individual peaks are shown in the table at the bottom of the report.

Heading

The heading displays the data file and path used by the report, the date the data was acquired, the operator's name, sample name, misc, and vial number. This is followed by the integration parameters file and the name of the integrator used.

RTE Integration Parameters

When the integrator is the RTE integrator, the RTE integration parameter settings are displayed next.

Graphics Area

The top window contains a Total Ion Chromatogram (TIC) of the data file. The lower graphic contains the noise area specified by the user.

Noise characteristics

Noise Start and Noise End

These are the retention time values entered by the user to define the range used for calculating noise.

Number of points

The number of data points found in the specified noise range.

Maximum noise

The abundance value of the data point with the maximum abundance in the noise range.

Minimum noise

The abundance value of the data point with the minimum abundance in the noise range.

Peak-to-Peak noise

PPNoise = Maximum noise - Minimum noise

Root-mean-square noise

RMS Noise =
$$\sqrt{\left(\frac{(N\sum(x^2) - (\sum x)^2)}{N(N-1)}\right)}$$

Where:

N = Number of points (see definition under Noise characteristics above)

x = abundance of an individual data point

Performance Report Results

#:

Peak number.

RetTime:

Retention time of the peak.

Area:

Area of the peak as determined by the selected integrator.

BsIPH:

The peak height as measured from the corrected baseline.

PWidth:

Peak width for the ChemStation Integrator is the width of the peak at $\frac{1}{2}$ the peak height. For the RTE Integrator, peak width is the width of the peak start minus the peak end.

Resol:

ChemStation Integrator Resolution, as determined using the "ChemStation" Integrator.



$$Resol = 1.18 * \frac{RT_2 - RT_1}{PW_{hh1} + PW_{hh2}}$$

Where:

 RT_2 = retention time of the second peak; RT_1 = retention time of the first peak; PW_{hh1} = peak width at ½ height of the first peak; PW_{hh2} = peak width at ½ height of the second peak.

RTE Integrator Resolution, as determined using the "RTE" Integrator.



$$Resol = 1.18 * \frac{RT_2 - RT_1}{PW_{BL1} + PW_{BL2}}$$

Where:

 RT_2 = retention time of the second peak; RT_1 = retention time of the first peak; PW_{BL1} = peak width at baseline of the first peak; PW_{BL2} = peak width at baseline of the second peak.

Please observe that the calculation for *Resolution* when using the RTE integrator does not use the peak width at ½ height. Instead the resolution is calculated using the peak width at the baseline. This is not one of the classical methods for calculating resolution but can be used as a measure to compare resolutions on the Agilent MSD Productivity ChemStation.

Tailing:

Peak Tailing Factor

This is a measure of peak symmetry. For a perfectly symmetric peak, the peak tailing factor = 1.



$$Tailing = \frac{RT_b - RT_f}{2.0 * (RT - RT_f)}$$

Where: RT = retention time of the peak. $RT_b = RT$ of peak back at 10% of peak height;

 $RT_f = RT$ of peak front at 10% of peak height;

PtP S/N:

Peak-to-Peak signal to noise ratio.

PtP S/N = (BslPH)/(PPNoise)

Where: *BslPH* = baseline corrected peak height; PPNoise (see definition under Noise characteristics above)

RMS S/N:

Root Mean Square signal to noise ratio

RMS S/N = (BslPH)/(RMSNoise)

Where: *BslPH* = baseline corrected peak height; RMSNoise (see definition under Noise characteristics above)

Plates:

Number of Theoretical Plates

This is a measure of the column efficiency (the ability to produce narrow peaks). The result is the number of theoretical plates per meter of column length.





© Agilent Technologies, Inc. 3/30/2012 G1701-90101