


Optimized Column Configuration and Operating Parameters for High-Throughput Analyses

Application
Technical
Robert Ricker

Recent advances in chemistry and improvements in LC-MS have led to a need for ultra-rapid HPLC separations. The following chromatograms show optimization of a method in relation to optimal column, time for separation, mobile-phase re-equilibration, and drift in baseline.



Conditions:
Column: ZORBAX Eclipse XDB-C18, 3.5 μ m, Agilent PN 931975-932
Mobile Phase: A: 10/90, MeOH / H₂O + 0.2% H₃PO₄
B: 90/10, MeOH / H₂O + 0.2% H₃PO₄
Inj. Vol.: 20 mM each compound, 5 μ L; 23°C

Highlights

- By switching to an Agilent ZORBAX Eclipse XDB-C18, the same conditions reduced run time from 3 min to just over 1 min. The baseline change was also improved by the Agilent 1100 detector.
- Reducing gradient time from 2 min to 1 min reduced run time to less than 1 min. Use of 268 nm further improved the baseline, for optimal peak quantitation.
- For optimal combinatorial and LC-MS studies, full gradients were run on these short columns, in 1 or 2 minutes.
- Time required for a re-equilibration with 5 column volumes on these short cartridge columns at 4 mL/min is only 22.5 seconds.



Agilent Technologies

*Robert Ricker is an application chemist
based at Agilent Technologies, Wilmington,
Delaware.*

For more information on our products and
services, visit our website at:
www.agilent.com/chem

Copyright© 2002 Agilent Technologies, Inc.
All Rights Reserved. Reproduction,
adaptation or translation without prior
written permission is prohibited, except as
allowed under the copyright laws.

Agilent shall not be liable for errors
contained herein or for incidental or
consequential damages in connection with
the furnishing, performance, or use of this
material.

Information, descriptions, and specifications
in this publication are subject to change
without notice.

Printed in the USA
April 25, 2002
5988-6463EN



Agilent Technologies