Agilent P3111
pH Electrode
pH 电极

Operating Guide
用户手册

Agilent Technologies

Overview

The P3111 pH Electrode can measure pH values. It is used with the 3200P meter or similar meters, with the R8111 or similar reference electrode.

WARNING Use this electrode according to the operating manual to avoid personal injury.

WARNING The electrode solution can cause chemical burns or illness if it is taken orally or contacted by human skin. Use protective clothing or gloves to avoid contact. In case of contact, rinse contacted area with tap water or deionized water thoroughly.

CAUTION The electrode body material is glass. Handle with care to avoid damage to the instrument.

Specifications

Table 1 P3111 pH Electrode specifications

Specification	Value
pH range	0 to 14
Zero potential pH with R8111 reference electrode	7.00 ± 0.5
Temperature range	0 to 60 °C
Electrode impedance	$< 250 \text{ M}\Omega$
Body material	Glass
Electrode diameter	12 mm
Electrode length	120 mm
Cable interface	BNC
Cable length	1000 mm


Figure 1 P3111 pH Electrode assembly

Operation

Preparing the electrode

- 1 Remove the storage bottle and cap from the electrode and store it upright for future use.
- 2 Rinse the measuring tip with distilled or deionized water.
- 3 Hold the electrode measuring tip downwards and swing it several times to remove air bubbles near the sensitive glass bulb.

Calibration

- 1 Connect the electrode and a reference electrode to the meter. Input the temperature of calibration solution according to the meter operating manual.
- 2 Soak the measuring tip in calibration solutions in sequence. Calibrate the electrode according to the procedures listed in the operating manual.

Measurement

- 1 Soak the measuring tip in sample solution. When the reading becomes stable, read the value from the meter.
- 2 Rinse the measuring tip.

Operating hints

- The main material of the measuring tip is glass. Ensure the sample solution will not damage the measuring tip before measurement. Solutions containing solids or crystals will damage the measuring tip.
- Keep calibration solution and sample solution flowing at a uniform speed around the electrode or the accuracy may be impacted.
- Do not soak the electrode in a sample solution for an extended time. After measurement, rinse the electrode carefully.
- White powder or crystals that collect on the electrode during measurement or transportation do not affect the performance of the electrode. Rinse the residue with tap or deionized water.
- The pH of a solution is affected by its temperature.
- For acidic sample solutions, use pH 4.00 and 7.00 standards for calibration. For alkaline sample solutions, use pH 7.00 and 10.01 standards.
- Do not apply force onto the electrode cap, cable or cable interface.
- Keep the cable interface dry.

Maintenance

Cleaning of inorganics

Soak the measuring tip in 0.1 mol/L HCl or EDTA for 15 minutes.

Cleaning of organics

Soak the measuring tip in ethyl alcohol or another solvent that can dissolve organics for 15 minutes.

Cleaning of grease

Soak the measuring tip in weakly alkaline detergent for 15 minutes.

Cleaning of protein precipitation

Soak the measuring tip in 0.1 mol/L HCl solution that contains 1% pepsin for 15 minutes.

Regeneration of the glass sensitive membrane

WARNING

HF is lethal. See the HF material safety data sheet before using it.

- 1 Soak the sensitive membrane in 4% HF solution for 3 to 5 seconds.
- 2 Rinse the membrane in 0.1 mol/L HCl several times.
- 3 Rinse with distilled or deionized water.
- 4 After one or more cleaning procedures, thoroughly rinse the electrode with distilled or deionized water.
- 5 If necessary, empty the reference electrode and add fresh reference solution. Repeat this action 2 or 3 times.
- 6 Soak the electrode in pH 4.01 solution for at least one hour.

Troubleshooting

Meters

Refer to the meter operating manual. Check all related parts, such as the electrode, calibration solution, and samples.

Electrode

- 1 Connect the electrode and a reference electrode to the meter. Set the meter to display mV.
- 2 Soak the measuring tip in pH 4.00 and pH 10.01 calibration solutions. When the readings become stable, record mV values respectively.
- 3 When the electrode is used with a R8111 reference electrode, the absolute value of the difference between these two mV values should be 326–370 mV in total. This indicates that the electrode has a good performance. If not, perform maintenance on the electrode.

Accuracy

- 1 Choose any two calibration solutions to calibrate a electrode. Use effective calibration solution. Make sure the distilled or deionized water used to prepare the solutions meets requirements, and the solution is not contaminated or beyond its shelf life.
- 2 Use this electrode to measure the pH value in another calibration solution.
- 3 Compare the theoretical value of that calibration solution with the measured value.
- 4 Ensure the measuring accuracy of the electrode is sufficient for the application.

Samples

The electrode is only suitable for measuring the pH value of conventional solutions. Nonaqueous reagents or solutions with low conductivity or high viscosity may have long measurement times, which can affect measurement accuracy.

For any other problems during electrode use, contact your Agilent Technologies customer service representative.

Storage

- 1** Install the electrode storage bottle and cap onto the electrode body in sequence.
- 2** To protect the measuring tip, keep 5–10 mm between the bottom of the storage bottle and the measuring tip of the electrode.
- 3** Tighten the storage bottle cap.
- 4** Place the electrode in the storage box and store at ambient temperature and dry conditions.

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© Agilent Technologies, Inc. 2012