**Agilent P3212**  
**pH Combined**  
**Electrode**  
**pH 电极**

**Operating Guide**  
**用户手册**

Agilent Technologies

## Overview

The P3212 pH Combined Electrode can measure the pH in aqueous solutions. It is used with a 3200P Meter or similar meter.

**WARNING** Use this electrode according to the operating manual to avoid personal injury.

---

**WARNING** The electrode solution can cause chemical burns or illness if it is taken orally or contacted by human skin. Use protective clothing or gloves to avoid contact. In case of contact, rinse contacted area with tap water or deionized water thoroughly.

---


**CAUTION** The electrode measuring tip is glass. Handle with care to avoid damage to the instrument.

---

# Specifications

**Table 1** P3212 pH Combined Electrode specifications

Specification	Value
pH range	0 to 14
Zero potential pH	$7.00 \pm 0.5$
Temperature range	0 to 60 °C
Electrode impedance	<300 M $\Omega$
Reference type	Ag/AgCl
Body material	PC (polycarbonate)
Liquid junction material	Porous polymer
Reference filling solution	3 mol/L KCl containing AgCl
Electrode diameter	12 mm
Electrode length	120 mm
Cable interface	BNC
Cable length	1000 mm

**Figure 1** P3212 pH Combined Electrode assembly

## Operation

### Preparing the electrode

- 1 Remove the electrode storage bottle and the storage bottle cap from the electrode. Store the bottle upright for future use.
- 2 Rinse the measuring tip with distilled or deionized water. Hold the measuring tip downwards and shake the electrode gently to dry.
- 3 Unplug the filling hole.

### Measurement

- 1 Connect the electrode to the meter. Enter the temperature of calibration solution according to the meter operating manual.
- 2 Soak the measuring tip in calibration solutions in sequence. Calibrate the electrode according to the meter operating manual.
- 3 Soak the measuring tip in sample solution. When the reading becomes stable, read the value from the meter.

## Operating hints

- The main material the of measuring tip is polycarbonate. Ensure that the sample solution will not damage the measuring tip before measurement.
- Maintain the calibration solution and sample solution at a uniform flow rate, or the accuracy of measurement may be impacted.
- The reference filling solution contains silver ions. Do not use to measure sample solutions that contain protein. The silver ion will precipitate near the liquid junction and damage the electrode.
- Do not soak the electrode in sample solution for an extended period of time. Rinse the electrode carefully after use.
- The pH value of a solution is affected by its temperature.
- White powder or crystals that collect on the electrode during measurement or transportation do not affect the performance of the electrode. Rinse the residue with tap or deionized water.
- For acidic sample solutions, use pH 4.00 and pH 7.00 standards for calibration. For alkaline sample solutions, use pH 7.00 and pH 10.01 standards.
- Do not apply force onto the electrode cap, cable, or cable interface.
- Keep the cable interface dry.

## Maintenance

### Cleaning of inorganics

Soak the measuring tip in 0.1 mol/L HCl or EDTA solution for 15 minutes.

### Cleaning of organics

Soak the measuring tip in absolute ethyl alcohol or other organic solvent that can dissolve organics for 15 minutes.

### Cleaning of grease

Soak the measuring tip in warm, weakly alkaline detergent for 15 minutes.

### Cleaning of protein precipitation

Soak the measuring tip in 0.1 mol/L HCl solution that contains 1% pepsin for 15 minutes.

### Regeneration of the glass sensitive membrane

**WARNING**

**HF is lethal. Refer to the HF material safety data sheet before use.**

---

- 1 Soak the sensitive membrane in 4% HF solution for 3 to 5 seconds.
- 2 Rinse in 0.1 mol/L HCl for 10 seconds.
- 3 Rinse with distilled or deionized water.
- 4 After one or more cleaning procedures, thoroughly rinse the electrode with distilled or deionized water. Soak the measuring tip in reference filling solution for 2 hours.

## Troubleshooting

### Meters

Refer to the meter operating manual. Check all related parts, such as the electrode, calibration solution, and samples.

### Electrode

- 1 Connect the electrode to the meter with mV displayed.
- 2 Soak the measuring tip in pH 4.00 calibration solution. Record the mV when the reading becomes stable.
- 3 Soak the measuring tip in pH 10.01 calibration solution. Record the mV when the reading becomes stable.
- 4 The difference between these two mV readings should be in the range of 326 to 366 mV. If not, clean the electrode, then repeat the measurements.
- 5 Hold the electrode with the measuring tip downwards and swing it several times.
- 6 Check the accuracy of the electrode by repeating steps 2 and 3 using any two other calibration solutions. Compare the measurements against the expected mV responses.

### Calibration Solution

- Use the correct calibration solution.
- The distilled water or DI water used to prepare calibration solution should meet with corresponding requirements.
- The calibration solution should not be contaminated, degraded, or expired.

## Sample

The electrode is only suitable for measuring the pH value of a conventional solution. When measuring the pH value of a solution with low conductivity, high viscosity, or nonaqueous reagents, the response time of electrode may be prolonged or the measurement accuracy may be affected.

For any other problems during electrode use, contact your Agilent Technologies customer service representative.

## Storage

- 1 Rinse the measuring tip.
- 2 Install the storage bottle cap, and the storage bottle onto the electrode body.
- 3 Screw the storage bottle cap onto the storage bottle until finger tight.
- 4 Place the electrode in the storage box and store at ambient temperature.

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© Agilent Technologies, Inc. 2012