

Quality control of dextran

Application

Heinz Goetz

Dextrans are widely used for different applications in pharmaceutical, cosmetics, paint, adhesive and paper industries. Their molecular weight can range from 3×10^4 to 3×10^7 . The properties strongly depend on the molecular weight, the molecular weight distribution and impurities. The chromatogram shows the analysis of a dextran used in cosmetic formulations.

The initial analysis with just one PL aquagel-OH mixed column already indicated the presence of

impurities by showing a shoulder in the main product. By adding the PL aquagel-OH 30 column the resolution could be increased further, thus clearly proving the impurity. Using the ChemStation GPC data analysis software not only the molecular weight averages and the molecular





Conditions

Sample preparation

Sample was dissolved in the mobile phase (concentration 0.1 %) and filtered Column PL aguagel-OH Mixed, 7.5 × 300 mm, 8 µm (Agilent p/n 79911GF-MXA) in series with PL aguageI-OH 30, $7.5 \times$ 300 mm, 8 µm (Agilent p/n 79911GF-083) Mobile phase water Flow rate 1 mL/min **Column compartment temperature** 25 ° C **Injection volume** 100 ul Detector Refractive index detector **Polymer standards** Polyethylene oxide EasyCal standards in vials for calibration (Agilent p/n 5064-8280)



weight distributions can be calculated, but also additional rapid numeric data is provided in the form of the user-defined report subsets. Figure 2 shows a typical GPC report containing information on the sample and the method used, the molecular weight distribution and the molecular weight results. Three subsets in form of weight fractions at the user-defined molecular weights 10000, 280000 and 780000 are reported:

- Subset 1: for molecular weights up to 10000 dalton: 0.16 % contribution to the polymer
- Subset 2: for molecular weights up to 280000 dalton: 80.76% contribution to the polymer
- Subset 3: for molecular weights up to 780000 dalton: 100% contribution to the polymer

The contribution of the impurity peak with molecular weights from 280000 to 780000 dalton is 19.1 %. Report subsets can be also specified for determination of molar masses at given values of the cumulative distribution.



HPLC performance

RSD of M_w	< 1.5%
RSD of M _n	< 5%

Equipment

Agilent 1100 Series GPC-SEC system

consisting of

- vacuum degasser for efficient degassing of the mobile phase
- isocratic pump with large solvent cabinet
- autosampler with single valve design
- thermostatted column compartment for precise column temperatures
- refractive index detector with automatic recycle valve
- ChemStation Plus with GPC-SEC data analysis software

Heinz Goetz is an application chemist at Agilent Technologies, Waldbronn, Germany

For more information on our products and services, visit our worldwide website at <u>http://www.agilent.com/chem</u>

© Copyright 2000 Agilent Technologies Released 09/2000 Publication Number 5988-0118EN



