Agilent 3200C **Conductivity Meter** 电导率仪

Operating Guide 用户手册

Agilent Technologies

Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

5973-1780

Edition

First Edition, March 2012

Printed in China

Agilent Technologies, Inc.
412 Ying Lun Road
Shanghai 200131 P.R. China

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING


A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

1	Installation	7
	Tools and Components needed for Installation	8
	Installation of the 3200C Conductivity Meter	9
	Installation of electrochemical probes	9
	Installing the power adapter	11
	Installing the ground line	12
	Installing optional software	13
2	Operation	15
	Introduction	16
	Where to Find Information	16
	Agilent customer portal	16
	Definition of Terms	16
	Cell constant	16
	TDS factor	16
	Temperature coefficient	16
	Proper Use	17
	Features	18
	Specifications	19
	Physical overview	20
	The Display	21
	The Keypad	22
	Instrument Control	23
	Turn on the Meter	25
	Setup	26
	Measuring Mode	26

Descriptions of the modes	27
Set Manual Temp	29
Set Condition	30
System Setup	31
Descriptions of the System Setup settings	32
Language Select	34
Auto Power Off	34
Set Default	35
Prepare Probes	36
Calibrate	37
To set or calibrate the cell constant (Cell Const)	37
One-Point calibration	37
Two-Point Calibration	40
Manually enter a cell constant	40
To Set or Calibrate the TDS Factor	41
Manually enter the TDS factor	41
Calibrate the TDS factor	41
Temperature Coefficient (Set Temp Coef)	43
Perform a Measurement	44
Choose the correct probe	44
To measure conductivity	45
To measure resistivity	46
To measure temperature	46
To measure TDS	46
Measure TDS from the initial state	46
Measure TDS using the TDS Measuring Mode	46
To measure salinity	47
Measure salinity from the initial screen	47
Measure salinity using the Salinity Measuring Mode	47
To view other parameters while measuring	48
Manage Data	49

To save data	49
To delete data	49
To output data	49
To view data	50
To View Saved Cond	50
View conductivity parameter	51
View TDS parameter	52
To Turn off the 3200C Conductivity Meter	53
Software Operation	54
To use the optional G4390A Electrochemical Data Collecting software	54
To use the Data Printing software	54
3 Troubleshooting and Maintenance	55
Troubleshooting the 3200C Conductivity Meter	56
General Troubleshooting Procedure	58
Check the temperature measurement	58
Check the conductivity measurement	59
Restore factory default settings	59
To restore factory defaults (Set Default)	59
Meter Self-Diagnostics and Messages	60
Severe errors	60
Error warning	61
Maintenance	63
Probe cleaning	63
Probe storage	63
Consumables and Ordering Information	64
接地线的安装	74
数据打印软件界面操作	103

1 Installation

Tools and Components needed for Installation	8
Installation of the 3200C Conductivity Meter	9

Tools and Components needed for Installation

Agilent provides the tools needed for installation. The following accessories are included in the shipping case.

- Electrode Holder (G4389A)
- Conductivity Probe, for example, C5111 Conductivity Probe (5190-3994)
- Power Adaptor (5185-8389)

Installation of the 3200C Conductivity Meter

Open the 3200C Conductivity Meter shipping case. Remove the meter, electrode holder, and other accessories.

Installation of electrochemical probes

- 1 Place the electrode holder near the meter and move the arm into position. Use the thumbscrews shown in [Figure 1](#) to secure the arm in place.


Figure 1 The electrode holder

- 2 Clip the conductivity probe and ATC probe into the electrode fixture shown in [Figure 1](#).

3 Route the electrode cable as shown [Figure 2](#).


Figure 2 Installed meter and probes

- 4 Locate the conductivity and ATC probe sockets on the back of the meter. Plug the C5111 conductivity probe (5190-3994) connector and ATC probe (5190-3998) connector into the sockets on the rear of the meter. See [Figure 3](#).


Figure 3 The back view of the 3200C Conductivity Meter

Installing the power adapter

A universal power adapter is included with the 3200C. The power adapter operates at 100–240 VAC, 1 A, 50/60 Hz. The output from the power adapter is 9 VDC 1 A.

CAUTION

Do not use any other power adapter with the 3200C.

- 1 The power adapter includes several outlet adapters. Choose the appropriate outlet adapter for the power outlets in your country.

- 2 Slide the outlet adapter on to the two metal prongs on the power adapter. A click will be heard when the plug is properly engaged. See [Figure 4](#).

Figure 4 Assembling the power adapter

- 3 Connect the power adapter cord to the power socket on the back of the meter. See [Figure 3](#) on page 11.


Installing the ground line

The meter provides a ground line, but it is not necessary to install it during measurement. Sometimes, equipment such as a constant temperature bath can be a source of interference, leading to unstable readings. In order to shield from an electromagnetic field, the ground line must be installed to avoid electrical leaking. To shield the wire, connect one terminal with the meter, and the other terminal with the source of interference.

Installing optional software

If purchased, install the optional G4390A Electrochemical Data Collecting Software now. The software provides communications between the meter and a computer. Connect the computer to the meter with a USB cable. See the G4390A software documentation for more information.

If not using the G4390A software, Agilent provides downloadable data printing software on the Agilent Customer Portal (see [“Agilent customer portal”](#) on page 16) To use this software, download it from the portal and install it. Then connect the meter to the PC using the USB cable. Refer to the data printing software documentation for more information.

2 Operation

Introduction	16
Instrument Control	23
Turn on the Meter	25
Setup	26
Prepare Probes	36
Calibrate	37
Temperature Coefficient (Set Temp Coef)	43
Perform a Measurement	44
Manage Data	49
To Turn off the 3200C Conductivity Meter	53
Software Operation	54

This section describes the operation of the 3200C Conductivity Meter. There are five major steps to measure conductivity with the meter:

- Meter setup (“[Setup](#)” on page 26)
- Prepare the conductivity probe (“[Prepare Probes](#)” on page 36)
- Set or calibrate the cell constant (“[Calibrate](#)” on page 37)
- Set the temperature coefficient (“[Temperature Coefficient \(Set Temp Coef\)](#)” on page 43)
- Measure conductivity (“[Perform a Measurement](#)” on page 44)

In each phase of operation, the LCD displays the current state of the meter.

Introduction

Where to Find Information

This document describes how to install, operate, maintain, and troubleshoot the 3200C Conductivity Meter.

Before operating your meter, be sure to read the meter installation and operation information.

Agilent customer portal

Agilent also provides customized information for the products you own through a customer portal. This web service provides many customizable services as well as information related directly to your Agilent products and orders. Log onto the portal at <http://www.agilent.com/chem>.

Definition of Terms

Cell constant

The cell constant of a pair of platinum sheets is the ratio of distance between the two sheets to the area of each sheet. The ratio is measured in cm^{-1} .

TDS factor

The Total Dissolved Solids (TDS) factor is the factor to convert conductivity to TDS.

Temperature coefficient

The temperature coefficient is the % change of conductivity with temperature in $\%/^{\circ}\text{C}$.

Proper Use

To avoid common safety issues:

- Ground the meter using the ground line connected to the back of the meter.
- The meter can be used continuously for a long time. After each measurement, soak the conductivity probe in distilled water. If an electrode is out of use for more than 6 hours, rinse and store it in a protective container.
- Improper positioning of the conductivity probe can result in abnormal measurements. When immersing the conductivity probe in a sample solution, place the conductivity probe in a location where the solution can flow freely around the probe.

CAUTION

Do not expose to corrosive gas. Keep the sockets on the back of the meter clean and dry. Do not allow contact with acid, alkaline, or salt solutions.

CAUTION

Only use the power adapter included with the meter.

Features

The Agilent 3200C Conductivity Meter is a state-of-the-art and customer-friendly benchtop analytical instrument, which is designed for highly accurate measurements of conductivity, resistivity, TDS, and salinity in aqueous solutions. It can also measure purity of water and the salinity (NaCl) of sea water in a desalination process. The meter can measure temperature and conductivity simultaneously.

3200C Conductivity Meter features include:

- Measurement of conductivity, resistivity, TDS, salinity, and temperature
- Auto temperature compensation, auto-calibration, auto-range, and auto-frequency switch over the entire measuring range
- User calibrated cell constants or TDS factors
- A clear dot-matrix LCD display
- GLP norm for conductivity, TDS, and salinity including recorded operator numbers, operation history, calibration data, and 200 sets of measurement data
- The ability to view, output, and delete stored data
- Three measuring modes to meet the needs of various users:
 - Continuous Mode
 - Timed Reading (Timed Reading Mode)
 - Auto-Lock Reading (Auto-Lock Reading Mode)
- USB PC connectivity with available communication software
- Power-off protection. When the meter is manually or automatically shut off, the stored measuring data, calibration data, and setting parameters will not be lost.
- A back-lit design that can be used to take measurements in a dark environment
- A durable key pad

Specifications

Measuring range	
Conductivity	0.000 $\mu\text{S}/\text{cm}$ –2000 mS/cm
Resistivity	5.00 Ωcm –100.0 $\text{M}\Omega\text{cm}$
TDS	0.000 mg/L –1000 g/L
Salinity	0.00–8.00%
Temperature	-5.0–110.0 $^{\circ}\text{C}$
Accuracy	
Temperature resolution	0.1 $^{\circ}\text{C}$
Conductivity	$\pm 0.5\%$ (FS)
Resistivity	$\pm 0.5\%$ (FS)
TDS	$\pm 0.5\%$ (FS)
Salinity	$\pm 0.1\%$
Normal working conditions	
Environmental temperature	0–40 $^{\circ}\text{C}$
Relative humidity	$\leq 85\%$
Power supply	Power adaptor UE15WCP1-090100SPA
Power input	100–240 VAC, 1 A
Power output	9 VDC, 1 A
Size	(length \times width \times height, mm): 190 \times 190 \times 105
Weight	about 1 kg

CAUTION

Do not use where nearby vibrations will affect the performance. Do not use if corrosive gas is in the air. Do not use near strong electromagnetic fields.

Physical overview

The Agilent 3200C Conductivity Meter is composed of a meter and a probe system. The probe system contains a conductivity probe and an ATC probe held by the electrode holder. See [Figure 5](#), [Figure 6](#), and [Figure 7](#).


Figure 5 The meter, electrode holder, and probes


Figure 6 The front view of the 3200C Conductivity Meter


Figure 7 The back view of the Agilent 3200C Conductivity Meter

The Display

The left side of the screen shows current system time and the right side shows the current measuring mode, parameter, cell constant, and temperature coefficient.


Figure 8 The display

During operation, the display shows the working condition and current meter settings.

The Keypad

The keypad has 15 keys. All keys have double functions except for **[Enter]**, **[Cancel]**, and **[On/Off]**. Functionality automatically changes based on context. Key functions are listed in the [Table 1](#) on page 23:


Figure 9 Keypad

Uses for the 3200C Conductivity Meter

Conductivity is a measure of a material’s ability to conduct an electric current. The meter measures the conductivity of solutions with high accuracy. When equipped with a probe with the proper cell constant, the conductivity meter can measure the conductivity of ultra pure water.

Users can select different probes each with a different cell constant to match the expected conductivity of samples to be measured. Refer to the final packing list for the probes that are shipped.

Instrument Control

The 3200C Conductivity Meter is usually controlled by the keypad. Most keys have two functions. Functionality automatically changes based on context. Key functions are listed in the [Table 1](#):

Table 1 Keypad functions




Button	Function 1	Function 2
	Type the number 1 .	Output data when viewing stored or calibrated data.
	Type the number 2 .	Move the selection cursor upward when selecting.
	Type the number 3 .	Store measurement data.
	Type the number 4 .	Move the selection cursor left when selecting.
	Type the number 5 .	Opens the Setup menu and also is used as a general "select" button.
	Type the number 6 .	Move the selection cursor right when selecting.
	Type the number 7 .	View stored or calibrated data.
	Type the number 8 .	Move the selection cursor down when selecting.
	Type the number 9 .	Switch the display window or parameter when measuring.

Table 1 Keypad functions

Button	Function 1	Function 2
	Type the number 0 .	Begin measurement.
	Type a decimal point.	Calibrate probe.
	Type a negative number.	Delete the data being viewed.
	Meter power switch.	

Alternately, the 3200C Conductivity Meter can be controlled through a computer. For how to connect to a computer and how to use data collecting software (the software is optional), refer to installation help in the DVD shipped with the meter.

Turn on the Meter

- 1 Press [**On/Off**] to turn on the meter. When the meter is powered on, the display shows the Agilent logo, the meter model, and other information.
- 2 After the self-check is completed, the display shows the initial state. See [Figure 10](#). The left side shows the system time. The right side shows the setting: measuring mode, parameter, cell constant, and temperature coefficient.

Mode&Para	
10:33	Continuous Mode
2010/06/20	Meas COND
	Cell Const:1.000
	Temp Coef:2.00%/t

Figure 10 The initial state

Setup


Figure 11 The setup screen

Measuring Mode

This meter supports three measuring modes including Continuous Mode, Timed Reading Mode, and Auto-Lock Reading Mode.

To set the measuring mode:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#).
- 2 Use the directional arrow keys to navigate to **Measuring Mode**.
- 3 Press [**Enter**] to select **Measuring Mode**. If you need to go back after selecting a menu item, press [**Cancel**] to exit.
- 4 The **Measuring Mode** screen displays **Parameters** on the left and **Modes** on the right. See [Figure 12](#). Check marks indicate the current selections in use.


Figure 12 Select a measuring mode

- 5 Choose a parameter to display on the initial screen during measurement. Use the arrow keys to highlight the desired parameter and press **[Setup]** to select it. A check appears next to the selected parameter. During measurement, you can still view, save and print the other two parameters. See [“To view other parameters while measuring”](#) on page 48. The parameter choices are conductivity, total dissolved solids (TDS), and salinity.
- 6 Use the arrow keys to highlight the desired mode and press **[Setup]** to select it. See [“Descriptions of the modes”](#) on page 27 for more information. A check appears next to the selected mode.
- 7 Press **[Enter]** to save the new setup, exit Setup, and return to the initial state. Press **[Cancel]** to exit Setup without saving changes and return to the initial state.

The new parameters and measuring mode are displayed as shown in [Figure 13](#).


Figure 13 Measuring mode displayed

Descriptions of the modes

Continuous Mode The meter will simply continue measuring, calculating, and displaying the data. It is up to the user to manually save data or end the measurement.

This is the most frequently used measuring mode. When measurement begins, the meter continuously measures, calculates, and displays the results. You can view the calibrated parameter, calibrate the probe, save, or print results at any time during measurement. To end a measurement, press **[Cancel]** followed by **[Enter]**.

Timed Reading The meter will automatically save measurement data periodically while measuring.

When you select Timed Reading, specify a time interval from 1 to 99 minutes. The default time interval is 10 minutes.

When measurement begins, the meter calculates and displays results as it normally does. When the set interval time elapses, the meter automatically stores the data and continues measuring. Data is constantly recorded at the specified time interval.

If the USB interface is connected to a PC with the Data Printing Software installed, the meter can print the data.

Read about Auto Delete before using this mode. See [“Descriptions of the System Setup settings”](#) on page 32.

To end a measurement run, press **[Cancel]** followed by **[Enter]**.

Auto-Lock Mode The meter will stop measuring when recorded values stabilize enough to fall within a set precision and length of time.

Before performing a measurement in Auto-Lock reading mode, set the Auto-Lock conditions (see [“Set Condition”](#) on page 30). When the measurement begins, the meter automatically measures, calculates, and displays the results. Once the measurements meet the preset Auto-Lock condition, the measurement finishes and the screen locks with the final readings.

You can view the calibrated parameter, calibrate the probe, save, or print results at any time during measurement. After measurement, you can save and print the results. Press **[Cancel]** to exit measuring mode or press **[Measure]** to begin the next measurement.

Set Manual Temp

If the meter is connected to an ATC probe, the probe will supply a temperature value for the meter. If an ATC probe is not connected to the meter, the user must enter a temperature value manually.

NOTE

If the measuring parameter is set to **SALI**, the meter uses a temperature of 18 °C.

To set the manual temperature:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Set Manual Temp**.
- 3 Press [**Enter**] to select **Set Manual Temp** or press [**Cancel**] to exit.
- 4 The Set Manual Temp screen appears on the display. See [Figure 14](#).


Figure 14 Set Manual Temp

- 5 With the Manual Temp value highlighted, press [**Setup**] to edit the value.
- 6 Use the numeric keys to enter a temperature.
- 7 Press [**Enter**] to save the new setup, exit, and return to the initial state. Press [**Cancel**] to exit without saving changes.

Set Condition

Use Set Condition to specify Auto-Lock parameters when measuring in the Auto-Lock measuring mode. An Auto-Lock measurement ends when all measured parameters meet the conditions set. For example, if the condition of **Conductivity** is set to 1.0%, when the measured conductivity holds stable within a 1.0% tolerance for the set Auto-Lock time interval, the measurement ends and the screen locks displaying the final reading.

To set the conditions:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Set Condition**.
- 3 Press [**Enter**] to select **Set Condition** or press [**Cancel**] to exit.
- 4 The **Set Condition** screen appears on the display. See [Figure 15](#).


Figure 15 The Set Condition screen

- 5 Use the arrow keys to select a condition to edit. Press [**Setup**] to edit the value.
- 6 Enter the desired value. The Auto-Lock time has a range of 1–200 seconds.
- 7 Press [**Enter**] to save the new setup, exit, and return to the initial state. Press [**Cancel**] to exit without saving changes.

System Setup

The System Setup screen is where the user sets the time on the system clock, sets the calibration interval, and enters the information for the Good Laboratory Practices (GLP) standard.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **System Setup**.
- 3 Press [**Enter**] to select **System Setup** or press [**Cancel**] to exit.
- 4 The **System Setup** screen appears on the display. See [Figure 16](#).

System Setup		
Calib Interval		Time
Cell Const:	0h	10/06/20
TDS Factor:	0h	10:43:27
Operator No		Auto Delete
Operator No	000	Auto Del <input checked="" type="checkbox"/>

Figure 16 The System Setup screen

- 5 Use the arrow keys to select a value to edit. Press [**Setup**] to edit the value.
- 6 Enter the desired value. See “[Descriptions of the System Setup settings](#)” on page 32 for more information.

- 7 If needed, edit the system date and time. To modify the date and time:
 - a Use the arrow keys to highlight the Time field shown in [Figure 16](#) on page 31 and press **[Setup]**.
 - b The **Set Date & Time** screen opens. This screen displays the current year, month, day, hour, minute, and second (see [Figure 17](#)). Press the arrow keys to highlight the proper cell, then press **[Setup]** to select it.


Figure 17 The Set Date & Time screen

- c Enter the correct value using the numeric keys and press **[Enter]**.
 - d Once all changes are made, press **[Enter]** to save the settings and return to the **System Setup** screen.
- 8 Once all of the system settings are made, press **[Enter]** to save the new setup, exit, and return to the initial state. Press **[Cancel]** to exit without saving changes.

Descriptions of the System Setup settings

Calib Interval is the number of hours between the calibration messages that the Meter periodically displays as a reminder to recalibrate the probe. The meter begins counting the interval time from the end of the previous calibration. When the calibration interval has elapsed, the meter displays a popup window to remind the user to recalibrate the probe. A **Calib Interval** value of **0** disables this function.

Operator No is a three-digit number from 000–200 used to identify the person operating the meter. The **Operator No** is recorded every time the meter saves data.

When stored data reaches the meter's memory limit **Auto Delete** allows you to automatically overwrite the oldest data with new data. If **Auto Delete** is not enabled, the meter does not save new data when its memory is full. For example, the meter supports storing up to 200 sets of data. When you want to store the 201st set of data, if **Auto Delete** is on, meter will delete the first set of data and store the 201st data. If **Auto Delete** is disabled, the new data will not be stored.

CAUTION

If **Auto Delete** is disabled, data could be lost.

Language Select

This meter supports Chinese and English languages.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Language Select**.
- 3 Press [**Enter**] to select **Language Select** or press [**Cancel**] to exit.
- 4 The **Language Select** screen appears. Use the arrow keys to choose either **Chinese** or **English** and press [**Enter**].

Auto Power Off

The meter can power off automatically after a set time period. To set the **Auto Power Off** time:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Auto Power Off**.
- 3 Press [**Enter**] to select **Auto Power Off** or press [**Cancel**] to exit.
- 4 The **Auto Power Off** screen appears. Enter a power off time from 10–480 minutes. Enter a time of 0 to disable **Auto Power Off**.
- 5 Press [**Enter**].

When the **Auto Power Off** time elapses after the start of a measurement, the meter shuts down. If the meter is connected to the Electrochemical Data Collecting Software, Auto Power Off is disabled.

Set Default

Use **Set Default** to reset the meter to its original factory settings.

- 1 From the initial state, press **[Setup]**.
- 2 Use the directional arrow keys to navigate to **Set Default**.
- 3 Press **[Enter]** to select **Set Default** or press **[Cancel]** to exit.
- 4 The **Set Default** screen appears on the display. Press **[Enter]** to restore the meter to the original factory settings.

Table 2 Default values for meter parameters

Parameter	Default value
Cell constant calibration data	One standard solution. Standard Value: 1408 $\mu\text{S}/\text{cm}$ Measured Value: 1408 $\mu\text{S}/\text{cm}$ Calibration temperature: 25.0 $^{\circ}\text{C}$ Default cell constant: 1.000
TDS calibration data	One standard solution. Standard Value: 704 mg/L Measured Value: 704 mg/L Calibrated temperature: 25.0 $^{\circ}\text{C}$ Default TDS coefficient: 0.500
Temperature coefficient	2.00%/ $^{\circ}\text{C}$
manual temperature	25.0 $^{\circ}\text{C}$
measuring mode	Continuous
Timed Reading Mode interval	10 minutes
Calib Interval	message disabled (0 hours)
Operator No	000
Auto Delete	Enabled
Auto-Lock Set Condition	conductivity 1.0% TDS 1.0% salinity 1.0%
Auto-Lock time	5 s
Auto Power Off	Disabled (0 minutes)

Prepare Probes

Perform this procedure when using a new probe, or when using a probe that has been in storage for a long time.

- 1 Submerge the measuring tip in absolute ethanol for 1 minute.
- 2 Repeat [step 1](#) twice with fresh absolute ethanol each time.
- 3 Rinse the measuring tip thoroughly with DI water.
- 4 Soak the measuring tip in DI water for 1 hour.
- 5 Connect the probe to the meter. See [“Installation of electrochemical probes”](#) on page 9.
- 6 Begin measuring. After finishing the measurement, clean the measuring tip with the proper method for cleaning the test solution used.

CAUTION

Avoid damaging platinum black coating on the C5111 Conductivity probe during cleaning.

-
- 7 Calibrate the cell constant. See [“To set or calibrate the cell constant \(Cell Const\)”](#) on page 37. The cell constant of each probe is printed on the label. However, transportation, prolonged storage, and extended use may change the cell constant to some degree.
 - 8 Before use, ensure there is no precipitate on the platinum sensing sheets. Ensure there is no corrosion on the platinum sensing sheets or chemical reaction in solution with platinum.
 - 9 If there is unusual performance during measurement or you suspect an error in measurement, see [Chapter 3](#), “Troubleshooting and Maintenance” to diagnose the problem.

Calibrate

To set or calibrate the cell constant (Cell Const)

The cell constant for each probe was calibrated and printed on the probe label before shipping. However, transportation, storage, and continued use may change the cell constant to some degree. Recalibrate the conductivity probe before measurement.

One-Point calibration

- 1 In order to calibrate, you need to make a standard solution. Read the printed cell constant from the label. Use [Table 3](#) to find the suitable KCl concentration of the standard solution.

Table 3 Concentration of Standard Solution based on cell constant

Cell constant (cm^{-1})	0.01	0.1	1	10
KCl concentration (mol/L)	0.001	0.01	0.01 or 0.1	0.1 or 1

- 2 Use the concentration found in [Table 3](#) to find the weight (g) of KCl per liter of DI water ([Table 4](#)) and prepare the standard solution. For example, to prepare a concentration of 0.1 mol/L, add 7.4365 g of KCl to 1 L of DI water.

Table 4 Composition of Standard Solution

Concentration (mol/L)	Concentration of KCl in g/L (20 °C)
1	74.2650
0.1	7.4365
0.01	0.7440
0.001	Dilute 0.01 mol/L solution from 100 mL to 1 liter

- 3 If the meter is connected to a conductivity probe without an ATC probe, the meter will use the manual temperature (see “Set Manual Temp” on page 29) value as current temperature. If the manual temperature is set at 25.0 °C, 20 °C, 18 °C, or 15 °C, conductivity is displayed without temperature compensation.
- 4 Rinse the conductivity probe with distilled or DI water.
- 5 Soak the conductivity probe in the standard solution. The temperature of the standard solution must be either 30 °C, 25 °C, 20 °C, 18 °C, or 15 °C within ± 0.1 °C.
- 6 Press [**Calibrate**].
- 7 Select **Cell Const**.
- 8 The **Calibration** screen appears. (see Figure 18).


Figure 18 The Calibration screen

- 9 Press **[Setup]** to enter a Standard Solution Value (STD Value).
Type the standard value listed in [Table 5](#) for the corresponding solution concentration and temperature.


Figure 19 The Standard Solution Value

Table 5 STD Value (μS/cm)

	15.0 °C	18.0 °C	20.0 °C	25.0 °C	30.0 °C
1 mol/L	92120	97800	101700	11310	131100
0.1 mol/L	10455	11163	11644	12852	15353
0.01 mol/L	1141.1	1220.0	1273.7	1408.3	1687.6
0.001 mol/L	118.5	126.7	132.2	146.5	176.5

- 10 When the reading stabilizes, press **[Enter]**. The meter will calculate and set the new cell constant.
- 11 The calibration is done. Press **[Cancel]** to return to the initial screen.

Two-Point Calibration

When measuring a solution with high conductivity (>10 mS/cm), perform a two-point calibration.

- 1 Prepare two standard solutions. Standard solutions with high conductivity should have similar chemical composition to sample solution.
- 2 Calibrate the cell constant as described in “[One-Point calibration](#)” on page 37.
- 3 Press **[Enter]** to end the first calibration.

NOTE

Do not press **[Cancel]** to avoid exiting calibration.

- 4 Put the conductivity probe into the second standard solution.
- 5 When the reading stabilizes, press **[Enter]**. The meter will automatically calculate and set the two-point calibration.
- 6 The calibration is done. Press **[Cancel]** to return to the initial screen.

Manually enter a cell constant

- 1 While performing a measurement, press **[Setup]**.
- 2 select **Cell Const** and press **[Enter]**.
- 3 The meter displays a popup window. Enter the cell constant given on the probe label.

NOTE

Manually entering a cell constant will overwrite a previous cell constant calibration.

To Set or Calibrate the TDS Factor

Manually enter the TDS factor

- 1 While measuring, press [**Setup**] and select **Set TDS Factor**. Press [**Enter**]
- 2 The meter displays a popup window. Enter the TDS factor according to your actual need.

Calibrate the TDS factor

- 1 Calibrate the meter with a conductivity standard solution. Select a proper standard solution according to the nature and measuring range of the sample. (The relationship between conductivity and the TDS of the standard solution is shown in [Table 6](#)).

Table 6 Relationship Between Conductivity and TDS of the Standard Solution

Conductivity ($\mu\text{S}/\text{cm}$)	TDS standard value		
	KCl (mg/L)	NaCl (mg/L)	442 (mg/L)
23	11.6	10.7	14.74
84	40.38	38.04	50.0
447	225.6	215.5	300
1413	744.7	702.1	1000
1500	757.1	737.1	1050
2070	1045	1041	1500
2764	1382	1414.8	2062.7
8974	5101	4487	7608
12880	7447	7230	11367
15000	8759	8532	13455
80000	52168	48384	79688

NOTE

The value listed in the form is the one at 25°C.

442 indicates 40% Na₂SO₄, 40% NaHCO₃, 20% NaCl.

- 2 While measuring, press [**Setup**] and select **cell constant**.
- 3 Set the cell constant to the value written on the conductivity probe label, or recalibrate the conductivity probe to obtain a new cell constant.
- 4 Press [**Mode**] and select **TDS Measuring mode**.
- 5 Press [**Calibrate**] and select **TDS Factor**. Press [**Enter**] to enter Calibration mode as shown in [Figure 20](#).


Figure 20 Calibration mode

- 6 Rinse the conductivity probe with distilled water or DI water.
- 7 The temperature of the calibration solution must be 25.0 ± 0.1 °C. Soak the conductivity probe in the standard solution
- 8 Press [**Setup**].
- 9 Select **STD Value** and then press [**Enter**].
- 10 Enter an STD value given in [Table 6](#) on page 41.
- 11 When the readings become stable, press [**Enter**].
- 12 The meter automatically calculates and sets a new TDS factor.
- 13 Calibration is finished. Press [**Cancel**] to exit.

Temperature Coefficient (Set Temp Coef)

Temperature affects the accuracy of conductivity measurement. When you set the temperature coefficient, the meter automatically compensates for the temperature while measuring.

- 1 While measuring, press **[Setup]** and select **Temp Coef**. Press **[Enter]**.
- 2 The meter displays a popup window. Input the new temperature coefficient.

To view current temperature coefficient:

- 1 From the initial screen, press **[View]**.
- 2 Select **View COND Para** and press **[Enter]**. The meter displays the **View COND Para** screen. See [Figure 21](#).

The user can view current calibration data including cell constant and temperature coefficient. Usually the user doesn't need to setup the temperature coefficient. The default temperature coefficient is 2.00%/°C.


Figure 21 The View COND Para screen

Perform a Measurement

During measurement, make sure the tip of the probe is completely submersed in the solution. Place the probe in a location where the solution can flow freely around the measuring tip. To measure ultra pure water, install the conductivity probe in a sealed container with ultra pure water flowing through. The flow rate should be not too high.

Choose the correct probe

Before measuring, choose the proper conductivity probe. Refer to [Table 7](#) on page 44 when choosing a conductivity probe. There are two kinds of conductivity probes with cell constant $K=1.0$:

- Platinum black - the platinum sheet has a platinum black coating.
- Bright - the platinum sheet does not have a platinum black coating. The Bright probe is suitable for wider measurement range (2 to 3000 $\mu\text{S}/\text{cm}$).

If conductivity is $\geq 3000 \mu\text{S}/\text{cm}$, the measurement may not be as accurate. To measure high conductivity, use a conductivity probe with a high cell constant.

When conductivity is $\geq 200.00 \text{ mS}/\text{cm}$, use a conductivity probe with cell constant = 5 or 10.

When conductivity is $\geq 500.00 \text{ mS}/\text{cm}$, use a conductivity probe with cell constant = 10.

Table 7 Conductivity Range and Corresponding Cell Constant

Conductivity range	Resistivity range ($\Omega\text{-cm}$)	Recommended cell constant (cm^{-1})
0.000–19.99 $\mu\text{S}/\text{cm}$	20.00 M–50.0 k	0.01
0.20–200 $\mu\text{S}/\text{cm}$	5.00 M–5.00 k	0.1
2.000–20 mS/cm	500–5	1.0
20.00–200 mS/cm	50–5	10

To measure conductivity

- 1 Calibrate the meter if necessary. See “Calibrate” on page 37.
- 2 From the initial state display, press **[Setup]** and select a measuring mode and the COND parameter. See “Measuring Mode” on page 26.
- 3 Press **[Measure]** to start a measurement according to the selected mode. See Figure 22. The upper area of the screen shows the current measuring mode and system time. The central area of the screen shows conductivity, resistivity and temperature.


Figure 22 Measuring conductivity

Measurement methods and displays change with different measuring modes to some extent.

- 4 During a measurement, you can recalibrate the cell constant and setup parameters.
- 5 After a measurement, press **[Save]** to save data. Press **[Output]** to print data. Press **[Cancel]** to end measurement.

NOTE

When measuring conductivity and TDS with an ATC probe connected to the meter, the measurement will be temperature-compensated into the corresponding value at 25.0 °C according to the set temperature coefficient.

When measuring salinity with an ATC probe connected to the meter, the measurement will be temperature-compensated into the corresponding value at 18.0 °C according to the set temperature coefficient.

If no ATC probe is connected, the measurement is given without temperature compensation.

To measure resistivity

In the conductivity measuring mode, the meter displays both conductivity and resistivity. You can measure resistivity directly in the conductivity measuring mode. See “[To measure conductivity](#)” on page 45’.

To measure temperature

In any measuring mode, the meter displays the current temperature of the solution. See “[To measure conductivity](#)” on page 45.

To measure TDS

Measure TDS from the initial state

- 1 Select a proper conductivity probe before measuring TDS. The selection principle is the same as the selection of a conductivity probe.
- 2 Set corresponding cell constant. See “[To set or calibrate the cell constant \(Cell Const\)](#)” on page 37.
- 3 Set the temperature coefficient. See “[Temperature coefficient](#)” on page 16.
- 4 Set the TDS factor. See “[To Set or Calibrate the TDS Factor](#)” on page 41.
- 5 To skip steps 2 through 4, recalibrate conductivity probe.
- 6 Once the TDS parameters have been set, you can measure TDS directly from the initial state.

Measure TDS using the TDS Measuring Mode

- 1 Set the TDS parameter with **Setup Measuring Mode**. See “[Measuring Mode](#)” on page 26.
- 2 For a description of display, measuring procedure, and operation in TDS measurement, see “[To measure conductivity](#)” on page 45.

- 3 After measurement, press [**Save**] to save data,. Press [**Output**] to output measuring data. press [**Cancel**] to end measurement.

To measure salinity

Measure salinity from the initial screen

- 1 Select a proper conductivity probe before measuring Salinity. See “[Choose the correct probe](#)” on page 44.
- 2 Set the corresponding cell constant. See “[To set or calibrate the cell constant \(Cell Const\)](#)” on page 37.
- 3 Set the salinity parameters. See “[Measuring Mode](#)” on page 26. Now you can measure salinity directly from the initial screen.

Measure salinity using the Salinity Measuring Mode

- 1 Set the salinity parameter with **Setup Measuring Mode**. See “[Measuring Mode](#)” on page 26.
- 2 For a description of display, measuring procedure, and operation in salinity measurement, see “[To measure conductivity](#)” on page 45.
- 3 After measurement, press [**Save**] to save data,. Press [**Output**] to output measuring data. press [**Cancel**] to end measurement.

To view other parameters while measuring

During a measurement, the user still can view other parameters not currently in use without changing the measuring mode. For instance, while measuring with the **Cond** parameter, a user can view, save, and print other parameters such as **TDS** and **SALI**.

- 1 While performing a measurement, press **[Mode]** to conveniently view other parameter values. The screen display highlights the measuring window as shown in [Figure 23](#).


Figure 23 Viewing parameters while measuring

- 2 Press **[4/◀]** or **[6/▶]** to switch the view to TDS, salinity, resistivity, conductivity, and so on.
- 3 If you do not press a key for a few seconds, the meter will automatically exit from viewing mode and return to the regular initial screen.

Manage Data

To save data

The Agilent 3200C Conductivity Meter can save 200 sets of conductivity, TDS, and salinity data. The procedure to save data varies with different measuring modes.

- In continuous mode and Auto-Lock reading mode, press [**Save**] to save data when readings become stable.
- In timed reading mode, you can still save data manually, but the meter automatically saves measuring data at a fixed periodic time interval. See “[Measuring Mode](#)” on page 26.

To delete data

You can delete individual data entries or delete all data at once.

- 1 In either the initial screen or while measuring, press [**View**].
- 2 Select the data and press [**Delete**].

To output data

You can output current measuring data, the most recent calibration data, or saved data. There are two ways to output.

There are two ways to output current measurements, calibration data or saved data. One method is to use the optional data collecting software (G4390A) with the meter. The other method is to use the free Data Printing Software provided.

- 1 Connect the meter to the PC with a USB cable, and run the software.
- 2 Press the [**Output**] key during measurement, while viewing calibration or stored data. This allows the user to print the corresponding data. Refer to the software operating instructions for details.

To view data

The Agilent 3200C Conductivity Meter allows you to view parameters such as the last calibration data and current parameter. It also allows you to set the parameter directly and view saved data. The meter stores data according to the parameter measured. All saved data meets GLP standards. The meter can store 200 sets of conductivity, TDS, or salinity data.

- 1 In the initial state, press **[View]**.
- 2 Navigate to the type of data to view and press **[Enter]**. See [Figure 24](#).


Figure 24 View Data

To View Saved Cond

To view conductivity data, select **View Saved COND**.

- 1 In the initial state, press **[View]**.
- 2 Use the arrow keys to highlight **View Saved COND**.
- 3 Press **[Enter]**.
- 4 The display is shown in [Figure 25](#). The upper area displays the current viewing mode and actual storage amount. Each page can display up to 10 saved data. The display format varies with different viewing modes. The data displayed includes time and operator No.

View Saved Cond			Num=009
001 2010/06/20 11:52			OPERATOR:000
2.49kΩ.cm			401μS/cm 27.9℃
Temp Coef:2.00%/℃			Cell Const:1.000
002	06/20	2.48kΩ.cm	404μS/cm 25.7℃
003	06/20	2.48kΩ.cm	404μS/cm 25.7℃
004	06/20	2.48kΩ.cm	404μS/cm 25.6℃
005	06/20	2.48kΩ.cm	404μS/cm 25.6℃
006	06/20	2.48kΩ.cm	404μS/cm 25.6℃
007	06/20	2.48kΩ.cm	404μS/cm 25.6℃
008	06/20	2.48kΩ.cm	404μS/cm 25.6℃
009	06/20	2.48kΩ.cm	404μS/cm 25.6℃

Figure 25 View saved data

- 5 Press the arrow keys to view the saved data.
- 6 Press **[Output]** or **[Delete]** to print and delete data. (First connect meter with PC.)

View conductivity parameter

- 1 In the initial state, press **[View]**.
- 2 Use the arrow keys to highlight **View COND Para**.
- 3 Press **[Enter]**.
The display is shown in Figure 26. The upper area displays last calibration data. The bottom area displays current conductivity parameter including cell constant and temperature coefficient

View COND Para		STD Num:1	<input checked="" type="checkbox"/>
Last Calib		2008/01/09 08:12	
 No1	STD Value		1273pS/cm
	Measured:		1273pS/cm
	Temp:		25.0℃
Cell Const:		Temp Coef:	
1.000		2.00%/℃	

Figure 26 View conductivity parameters

View TDS parameter

- 1 In the initial state, press **[View]**
- 2 Use the arrow keys to highlight **View TDS Para.**
- 3 Press **[Enter]**. The display is shown in [Figure 27](#).


Figure 27 View TDS parameter

To Turn off the 3200C Conductivity Meter

After using the meter, press [**Save**] to save data. Press [**On/Off**] to turn off the meter.

When not in use, all probes should be soaked in distilled water. If the meter is out of use for a long time:

- Disconnect the power adaptor from power cable to avoid damaging power adaptor and the meter.
- Keep the socket of the meter clean and dry. Avoid contact of the socket with acid, alkaline and salt solution.
- Keep the input terminal of the meter (connector of probe) clean and dry. If the probe has been exposed to high humidity or water, wipe the cable interface of the probe with dry gauze or Kimwipe.
- When the conductivity probe will be out of use for a short time, soak it in DI water. If the probe will be out of use for 6 hours or longer, rinse it with DI water and store it in a protective bottle.

Software Operation

See also “[Installing optional software](#)” on page 13.

To use the optional G4390A Electrochemical Data Collecting software

- 1 Connect the meter to the power adapter.
- 2 Press [**On/Off**] to turn on the meter.
- 3 Connect the meter to the computer using the USB cable.
- 4 Run the Electrochemical Data Collecting Software in the computer and the software will recognize the meter model and type. The software provides the functions corresponding to the meter type. Refer to the software documentation for details.

To use the Data Printing software

- 1 Connect the meter to the power adapter.
- 2 Press [**On/Off**] to turn on the meter.
- 3 Connect the meter to the computer using the USB cable.
- 4 Launch the printing software.
- 5 With the meter taking measurements or when viewing calibration data or stored data, press [**Output**]. The software receives the data. Once received, print the data as desired. For details, refer to the Data Printing software documentation.

3 Troubleshooting and Maintenance

Troubleshooting the 3200C Conductivity Meter 56

General Troubleshooting Procedure 58

Meter Self-Diagnostics and Messages 60

Maintenance 63

Consumables and Ordering Information 64

This chapter describes how to verify whether the meter is working properly and how to maintain the meter.

Troubleshooting the 3200C Conductivity Meter

Table 8 lists common problems, possible causes, and how to resolve them.

Note that improper operation of the conductivity probe will lead to abnormal readings. During measurement, soak the probe tip completely in solution.

Table 8 3200C troubleshooting

Number	Failure mode	Failure cause	Solutions
1	No display after meter starts up. The display is not lit up.	1 Power adaptor is not installed correctly. 2 The power supply does not meet with the requirements. 3 The power adaptor is damaged. 4 The power socket has poor contact. 5 The LCD is damaged (After start up for a few minutes, the user can hear a buzzing sound after pressing [On/Off] key, but there is no display on LCD.)	1 Reconnect. 2 Use the required power supply. 3 Replace the power adaptor. 4 Ensure there is a good contact. 5 Contact Customer Service.
2	Display is dim.		1 Check the power adapter.
3	No buzzing when pressing key	1 You may have pressed invalid keys under the current setup. 2 The buzzer has been damaged.	1 The buzzer only sounds when you press a valid key.
4	No response when pressing key	1 You may have pressed invalid keys under the current setup. 2 The key has been damaged.	1 Press the valid key for operation.
5	Meter does not display 25 °C when ATC probe is not connected with meter.	When the meter is not connected with an ATC probe, the meter will automatically switch to manual temperature setting.	Manually input the correct solution temperature.

Table 8 3200C troubleshooting (continued)

Number	Failure mode	Failure cause	Solutions
6	The reading does not stabilize.	1 The electrode has been damaged or aged. 2 There is strong electrical signal interference source nearby (electrical leaking, strong electro-magnetic field, etc.).	1 Replace the electrode. 2 Remove the electrical signal interference source. Lift the beaker to a higher position to avoid electrical leaking. Shield from the electromagnetic field. Use the ground line shipped with the meter. Connect one terminal to the meter and the other terminal to the interference source.
7	Severe drift during conductivity measurement.	Probe tip contaminated.	Clean the probe.

General Troubleshooting Procedure

Many factors may affect measurement, including the electrode, the standard solutions used to calibrate electrode, the sample solution, the temperature during measurement, incorrect operation, and so forth. When the measured results are significantly different from what was expected, first determine whether the meter itself or factors other than the meter caused the error. Follow the suggestions described below to diagnose the problem. Since there are multiple measuring methods, you will need to diagnose the meter based on your application.

Diagnose whether the errors are caused by the meter hardware, electrode calibration, parameters, the sample solution, solution preparation, buffer solutions used for calibration, or other causes.

For factors other than the meter, compare measurement results between known solutions. Put the electrode in different standard solutions. Based on the comparison, judge whether the deviation resulted from the electrode, the solution, or something else.

Check the temperature measurement

Because temperature measurements are required for all measurement, always ensure temperature measuring is functioning well.

- 1** Connect the meter with the ATC temperature diagnostic tool (5185-8390) shipped with the meter.
- 2** Turn on the meter and begin measuring.
- 3** With the ATC temperature diagnostic tool connected, the meter should display a temperature reading between 49.0 to 51.0 °C. If yes, the meter is correctly measuring temperature. If the meter displays a temperature reading significantly different from 50 °C, there is a problem with the meter. Contact Agilent service.

Check the conductivity measurement

- 1 Disconnect the conductivity probe.
- 2 Connect the meter with the conductive diagnostic tool (5185-8391) shipped with the meter.
- 3 Turn on the meter and disconnect the ATC probe.
- 4 At the initial state, set the manual temperature to 25 °C.
- 5 Check the cell constant. If it is about 1.000, then go to next step. If not, set the cell constant to 1.000. (**View Last COND Calib**).
- 6 Begin measuring conductivity. The meter should display a resistance value of 9.00–11.00 k Ω and a conductivity value of 90.00–110.00 μ S/cm. If so, the meter is functioning properly. If the value deviates significantly from these values, there is a hardware problem with the meter. Contact the customer service department.

Restore factory default settings

If other troubleshooting does not solve the problem, you can restore the meter to the factory default settings. This clears all user calibrations and settings so that you can verify whether an erroneous calibration or setting is the cause of the problem. After restoring the factory default settings, retest the meter.

To restore factory defaults (**Set Default**)

Select **Set Default** from the initial setup display to restore all parameters to the factory default settings.

Meter Self-Diagnostics and Messages

The meter supports self-diagnosis, which can find some common errors caused by meter, probe, solution, or operation. It deals with these errors differently according to their severity.

There are two levels of severity: **Severe system errors** hinder further use of the meter or make the meter unable to finish necessary tasks. In this case, replacement or maintenance of the meter is required.

Minor errors, which are prompts or warnings, are caused by various reasons and can be corrected in several ways. However, you must pay attention to these errors to ensure measurement integrity and reliability. If you ignore these error messages the operation can continue.

Severe errors

When severe errors occur in the meter, the meter will shut down and display an error message similar to [Figure 28](#). In this case, the only operation that can be performed is to turn off the meter.


Figure 28 Severe errors

Error warning

Error warnings indicate improper electrode installation, solution problems, and electrode setup during use. Read these error warnings carefully. Proper maintenance and operation will reduce errors.

Figure 29 shows an example error warning that temperature is out of range.


Figure 29 Error warning

is the error warning icon. The icon includes the error code.

NOTE

This table describes all error messages which may occur in use. These error messages have different definition for different meters.

Table 9 Error warning codes

Code	Description	Solution	Memo
00	Conductivity measuring module error	Contact Customer Service.	Serious error
01	DO measuring module error	Contact Customer Service.	Serious error
02	Temperature measuring module error	Contact Customer Service.	Serious error
03	Data storage error	Contact Customer Service.	Serious error
20	Potential is out of range	Replace electrode.	–1999.9 to 1999.9 mV
21	pH/pX is out of range	Replace electrode.	–3.000 to 21.000 pH

Table 9 Error warning codes (continued)

Code	Description	Solution	Memo
22	Conductivity is out of range	Replace probe and solution.	0 to 2000 mS/cm
23	Resistivity is out of range	Replace probe and solution.	0 to 100 MΩ•cm
24	Temperature is out of range	Replace electrode and decrease solution temperature.	–6.0 to 120.0 °C
25	DO electric current is out of range	Replace probe.	0 to 4000 nA
26	pH electrode slope is out of range	Replace electrode and recalibrate it.	80 to 120 %
27	Failed to recognize pH buffers	Replace electrode, setup proper buffer group and replace buffers.	
28	Temperature of pH buffer is out of range	Cool or heat buffer.	
29	Calibrate the same buffer repeatedly	Replace the buffer.	Due to incorrect operation
30	Number of pH buffers exceeds maximum.	Remove one or more buffers from the group.	5 buffers at most
31	pH buffers conflict with each other	Remove one or more buffers with overlapping pH values.	
32	The data storage is full.	Overwrite the previous data and store new data.	Do not delete all data
33	The internal clock has a low battery.	Set time manually.	
34	Number of customer-defined ions exceeds maximum	Delete one or more unnecessary customer-defined ions.	
35	The maximum ion mode number stored	Delete some storage data of certain ion mode.	

Maintenance

This section describes how to maintain and store the probe.

Probe cleaning

General electrode cleaning includes inorganic cleaning, organic cleaning, grease cleaning, and protein precipitation cleaning. The type of cleaning needed depends on the contaminants and the probe. After one or more cleaning procedures, rinse the outside of the probe with distilled water.


Probe storage

Soak the measuring tip with distilled water or deionized water for temporary storage. When the electrode is not used for a long time, put the electrode back in the box and store it a dry place at ambient temperature.

Consumables and Ordering Information

Order No.	Model and Name	Description
G4388A	3200SA Stirrer	Electrode holder and magnetic stirrer : combined to stir solution with stable and precise speed with a large adjustable range.
G4389A	3200EA Electrode Holder	It is used to fix different probes.
G4390A	Software Package	3200C Conductivity Meter can communicate with computer after this software package is installed (including software).
5185-8389	Power Adaptor	100–240VAC, 1A, 50/60Hz
G4388-27000	Stirring Bar	It is used with stirrer.
5190-3994	C5111 Conductivity Probe	Glass, K=1, platinum black, Measuring range: 2–20000 $\mu\text{S}/\text{cm}$
5190-3995	C5112 Conductivity Probe	Glass, K=1, bright, Measuring range: 2–3000 $\mu\text{S}/\text{cm}$
5190-3996	C5113 Conductivity Probe	Glass, K=0.1, bright, Measuring range: 0.05–200 $\mu\text{S}/\text{cm}$
5190-3998	T7111 ATC Probe	Stainless Steel body, Measuring range: 0–100 $^{\circ}\text{C}$
5185-8390	ATC temperature diagnostic tool	Temperature diagnostics
5185-8391	Conductivity diagnostic tool	Conductivity diagnostics

Log onto the portal at <http://www.agilent.com/chem> for detail information. For new purchase, contact distributors of Agilent Technologies or log on the Agilent official website.

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