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
**Agilent 3200D  
Dissolved Oxygen Meter  
溶解氧测定仪**

**Operating Guide  
用户手册**

**Agilent Technologies**

# Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

## Manual Part Number

5973-1779

## Edition

First edition, March 2012

Printed in China

Agilent Technologies, Inc.  
412 Ying Lun Road  
Waigaoqiao Free Trade Zone  
Shanghai 200131 P.R.China

## Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

## Safety Notices

### CAUTION

A CAUTION notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a CAUTION notice until the indicated conditions are fully understood and met.

### WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

# Contents

<b>1 Installation</b>	<b>7</b>
Tools and Components for Installation	8
Installing the 3200D Dissolved Oxygen Meter	9
Installing the electrode holder	9
Installing the probe	12
Installing the power adaptor	12
Installing the optional ground line	13
Installing Optional Software	14
<b>2 Operation</b>	<b>15</b>
Introduction	16
Where to find information	16
Term definitions	16
Proper use	17
Features	18
Specifications	19
Physical overview	20
Instrument Control	23
Turn on the Meter	25
Setup	26
Measuring Mode	26
Set Condition	30
System Setup	31
Language select	34
To set the air pressure (barometric pressure)	34
To set the salinity value (SALI Value)	35

## Contents

Auto power off	35
Set Default	36
Prepare the DO Probe	37
To use the DO probe	37
Calibrate the DO Probe	39
Zero oxygen calibration (Calib Zero)	39
Full scale calibration (Calib Full)	40
Perform another calibration	40
Perform a Measurement	41
To measure DO	41
To measure saturation	42
To measure temperature	42
To view other parameters while measuring	42
Managing Data	44
To save data	44
To delete data	44
To output data	44
To view data	45
Turn off the 3200D Dissolved Oxygen Meter	50
Software Operation	51
To use the optional G4390A Electrochemical Data Collecting software	51
To use the Data Printing software	51
<b>3 Troubleshooting and Maintenance</b>	<b>53</b>
Troubleshooting the 3200D Dissolved Oxygen Meter	54
General Troubleshooting Procedure	56
Check the current response	56
Check the DO probe	56
Check the temperature measurement	57

Check other possible causes	57
Meter Self-Diagnostics and Messages	58
Severe errors	58
Error warning	59
Maintenance	61
Probe maintenance	61
Short term storage	61
Long term storage	61
Consumables and Ordering Information	62

## **Contents**

## 1 **Installation**

- Tools and Components for Installation 8
- Installing the 3200D Dissolved Oxygen Meter 9
- Installing Optional Software 14

The tools and accessories necessary for installation are shipped with the 3200D Dissolved Oxygen Meter.

## **Tools and Components for Installation**

Agilent provides the tools needed for installation. The following accessories are included in the shipping case.


- Multifunction electrode holder (G4389A)
- D6111 DO Probe (5190-3997)
- Power adaptor (5185-8389)

# Installing the 3200D Dissolved Oxygen Meter

Open the 3200D Dissolved Oxygen Meter shipping case. Remove the meter, electrode holder, and other accessories.

## Installing the electrode holder

- 1 Place the electrode holder near the meter and move the arm into position. Use the thumbscrews shown in [Figure 1](#) to secure the arm in place.

**Figure 1** The electrode holder

- 2 Clip the DO (dissolved oxygen) probe into the electrode fixture shown in [Figure 1](#).

**NOTE**

The DO probe includes an ATC probe.


## 1 Installation

3 Route the probe cable as shown Figure 2.

**Figure 2** Installed meter and probes

- 4 Locate the **DO/Temp** probe socket on the back of the meter. Plug the DO probe connector (5190-3997) into the **DO/Temp** socket on the rear of the meter. See [Figure 3](#).

**Figure 3** The back view of the 3200D Dissolved Oxygen Meter

## Installing the probe

- 1 Locate the socket for the DO probe on the back panel of the meter. Insert the probe connector into the socket.
- 2 Feed the probe cable through the hole of the multifunction electrode holder.

**Figure 4** Installing the probe

## Installing the power adaptor

The universal power adapter included with your benchtop meter is the only power adapter recommended for use with this unit.

The external electrical power adapter is rated at 100 to 240 VAC, 1 A, 50/60Hz.

This power adaptor includes a plug converter that supports several kinds of plugs.

- 1 Choose the appropriate power plug.
- 2 Connect the output plug of the power adapter to the power input on the meter. A click will be heard when the plug is properly engaged.

**Figure 5** Installing the power adaptor

## Installing the optional ground line

The meter ships with an optional ground line. During use, if another device (such as a constant temperature bath) causes electrical interference and an unstable reading, install the ground line. Connect one terminal of the ground line to the meter ground socket ( $\oplus$ ) and the other terminal to the interference source.

## Installing Optional Software

If purchased, install the optional G4390A Electrochemical Data Collecting Software now. The software provides communications between the meter and a computer. Connect the computer to the meter with a USB cable. See the G4390A software documentation for more information.

If not using the G4390A software, Agilent provides downloadable data printing software on the Agilent Customer Portal (see “[Agilent customer portal](#)” on page 16). To use this software, download it from the portal and install it. Then connect the meter to the PC using the USB cable. Refer to the data printing software documentation for more information.

## 2 Operation

- Introduction 16
- Instrument Control 23
- Turn on the Meter 25
- Setup 26
- Prepare the DO Probe 37
- Calibrate the DO Probe 39
- Perform a Measurement 41
- Managing Data 44
- Turn off the 3200D Dissolved Oxygen Meter 50
- Software Operation 51

This chapter provides an overview of the individual components that make up the 3200D Dissolved Oxygen Meter.

## Introduction

### Where to find information

This document describes how to install, operate, maintain, and troubleshoot the 3200D Dissolved Oxygen Meter.

Before operating your meter, be sure to read the meter installation and operation information.

#### Agilent customer portal

Agilent also provides customized information for the products you own through a customer portal. This web service provides many customizable services as well as information related directly to your Agilent products and orders. Log onto the portal at <http://www.agilent.com/chem>.

## Term definitions

<b>Dissolved oxygen concentration</b>	The extent of oxygen molecules dissolved in the water solution under certain conditions in mg/L
<b>Dissolved oxygen saturation</b>	The ratio of actual DO concentration versus saturated DO concentration in percent (%)
<b>Barometric pressure</b>	Barometric pressure in kPa
<b>Salinity</b>	The salinity in water in g/L
<b>Zero point calibration</b>	Probe calibration in an oxygen-free solution (freshly prepared 5% sodium sulfite solution)
<b>Full scale calibration</b>	Probe calibration in air-saturated water

## Proper use

The most common safety issues when working on the meter include:

- Ground the meter using the connection on the back.
- The meter can be used continuously for a long time. After each measurement, soak the DO probe in distilled water, which is first boiled and then cooled. Do not soak the DO probe in sodium sulfite solution. If a probe is out of use for more than 6 hours, rinse and store it in a protective container.

### CAUTION

Do not expose to corrosive gas. Keep the sockets on the back of the meter clean and dry. Do not allow contact with acid, alkaline, or salt solutions.

### CAUTION

Only use the power adapter included with the meter.

### CAUTION

The probe must be polarized for at least 60 minutes after installing a new membrane or after replacing the electrolyte.

### WARNING

**Be sure that the meter has a good ground connection.**

## **Features**

The 3200D Dissolved Oxygen Meter is a state-of-the-art and customer-friendly benchtop analytical instrument for the measurement of DO concentration and saturation in aqueous solutions with high accuracy. It can measure temperature simultaneously.

3200D Dissolved Oxygen Meter features include:

- Measurement of DO concentration, saturation and temperature
- Calibration, both Zero Point and Full Scale
- A clear dot-matrix LCD display
- Supports GLP standards:
  - Requires an operator number
  - Records, views, and prints calibration data
  - Stores 200 sets of measuring data
- Can view, print, and delete stored measuring data
- Three measuring modes to meet the needs of various users:
  - Continuous Mode
  - Timed Reading Mode
  - Auto-Lock Reading Mode
- USB PC connectivity with available communication software
- Power-off protection: When the meter is manually or automatically shut off, the stored measuring data, calibration data, and setting parameters will not be lost.
- A back-lit design that can be used in a dark environment
- A durable key pad

## Specifications

<b>Measuring range</b>	
DO concentration	0.00 to 45.00 mg/L
DO saturation	0.0 to 300.0%
Temperature	-5.0 to 110.0 °C
<b>Resolution</b>	
Temperature resolution	0.1 °C
DO concentration	0.01 mg/L
DO saturation	0.1%
<b>Accuracy of electric unit</b>	
DO concentration	±0.10 mg/L
DO saturation	±2.0%
Temperature	±0.1 °C
<b>Normal working conditions</b>	
Environmental temperature	0–40 °C
Relative humidity	<=85%
Power supply	Power adaptor (5185-8389)
Power input	100–240 VAC, 1 A
Power output	9 VDC, 800 mA
<b>Size</b>	(width × depth × height, mm): 190 × 190 × 105
<b>Weight</b>	about 1 kg

### CAUTION


Do not use where nearby vibrations will affect the performance. Do not use if corrosive gas is in the air. Do not use near strong electromagnetic fields.

## Physical overview


The 3200D is composed of a meter and a probe system. The probe system is held in place by the electrode holder. See [Figure 6](#), [Figure 7](#), and [Figure 8](#).

**Figure 6** The meter, electrode holder, and probe


**Figure 7** The front view of the 3200D Dissolved Oxygen Meter

**Figure 8** The back view of the Agilent 3200D Dissolved Oxygen Meter

### The Display


The display shows the working condition and current setting of the 3200D Dissolved Oxygen Meter.

**Figure 9** The initial state display

### The Keypad

The meter has 15 operating keys. The keys are shown in [Figure 10](#) and described in [Table 1](#) on page 23. Most keys have two values, a numeric entry and a function. Which one applies depends on the display; if it is waiting for a number, the numeric entry applies, otherwise, the function applies.





**Figure 10** Operating keys

# Instrument Control




The 3200D Dissolved Oxygen Meter is usually controlled by the keypad.

**Table 1** Keypad functions

Button	Function 1	Function 2
 1 Output	Type the number 1.	Output data when viewing or calibrating stored data.
 2 ▲	Type the number 2.	Move the selection cursor upward when selecting.
 3 Save	Type the number 3.	Store measurement data.
 4 ◀	Type the number 4.	Move the selection cursor left when selecting.
 5 Setup	Type the number 5.	Opens the Setup menu from the initial screen and also is used as a general "select" button.
 6 ▶	Type the number 6.	Move the selection cursor right when selecting.
 7 View	Type the number 7.	View stored or calibrated data.
 8 ▼	Type the number 8.	Move the selection cursor down when selecting.
 9 Mode	Type the number 9.	Switch the display window or parameter when measuring.


## 2 Operation

**Table 1** Keypad functions

Button	Function 1	Function 2
	Type the number <b>0</b> .	Begin measurement.
	Type a decimal point.	Calibrate the DO probe.
	Type a negative number.	Delete the data being viewed.
	Meter power switch	

## Turn on the Meter

- 1 Press [**On/Off**] to turn on the meter. When the meter is powered on, the display shows the Agilent name, the meter model, and other information.
- 2 After the system self-check, the meter goes to the initial state shown in [Figure 11](#). The left of the screen shows the current system time and the right of the screen shows the current measuring mode, parameter and last calibrated result.

**Figure 11** The initial state

## Setup


<b>Measuring Mode</b>	Sets up current measuring mode (Continuous Mode, Timed Reading Mode, Auto-Lock Reading Mode) and measuring parameter.
<b>Set Condition</b>	Use this mode to set condition for Auto-Lock reading mode.
<b>System Setup</b>	Use this mode to set system time and GLP Standard.
<b>Language Select</b>	The meter supports both Chinese and English languages.
<b>Air Pressure</b>	Set the barometric pressure value.
<b>Set SALI Value</b>	Set the salinity value.

### Measuring Mode

This meter supports three measuring modes including Continuous Mode, Timed Reading Mode, and Auto-Lock Reading Mode.


To set the measuring mode:

- 1 From the initial state, press [Setup]. The setup screen appears. See [Figure 12](#).

**Figure 12** The setup menu


- 2 Use the directional arrow keys to navigate to **Measuring Mode**.
- 3 Press [**Enter**] to select **Measuring Mode**. If you need to go back after selecting a menu item, press [**Cancel**] to exit.
- 4 The **Measuring Mode** screen displays **Parameters** on the left and **Modes** on the right. See [Figure 13](#). Check marks indicate the current selections in use.

**Figure 13** Measuring modes

- 5 Choose a parameter to display on the initial screen during measurement. Use the arrow keys to highlight the desired parameter and press [**Setup**] to select it. A check appears next to the selected parameter. During measurement, you can still view, save and print the other parameter. See "[To view other parameters while measuring](#)" on page 42. The parameter choices are dissolved oxygen and saturation.
- 6 Use the arrow keys to highlight the desired mode and press [**Setup**] to select it. See "[Descriptions of the modes](#)" on page 28 for more information. A check appears next to the selected mode.
- 7 Press [**Enter**] to save the new setup, exit Setup, and return to the initial state. Press [**Cancel**] to exit Setup without saving changes and return to the initial state.

The new parameters and measuring mode are displayed as shown in [Figure 14](#).

**Figure 14** The Continuous Mode display

### Descriptions of the modes

**Continuous Mode** The meter will simply continue measuring, calculating, and displaying the data. It is up to the user to manually save data or end the measurement.

This is the most frequently used measuring mode. When measurement begins, the meter continuously measures, calculates, and displays the results. You can view the calibrated parameter, calibrate the probe, save, or print results at any time during measurement. To end a measurement run, press [**Cancel**] followed by [**Enter**].

**Timed Reading** The meter will automatically save measurement data periodically while measuring.

When you select Timed Reading, specify a time interval from 1 to 99 minutes. The default time interval is 10 minutes.

When measurement begins, the meter calculates and displays results as it normally does. When the set interval time elapses, the meter automatically stores the data and continues with the next measurement.

Read about Auto Delete before using this mode. See “[Descriptions of the System Setup settings](#)” on page 32.

To end a measurement run, press [**Cancel**] followed by [**Enter**].

**Auto-Lock Mode** The meter will stop measuring when values stabilize enough to fall within the preset time.

Before performing a measurement in Auto-Lock reading mode, set the Auto-Lock conditions (see “[Set Condition](#)” on page 30). When the measurement begins, the meter automatically measures, calculates, and displays the results. Once the measurements meet the preset Auto-Lock condition, the measurement finishes and the screen displays the final readings.


You can view the calibrated parameter, calibrate the probe, save, or print results at any time during measurement. After measurement, you can save and print the results. Press [**Cancel**] to exit measuring mode or press [**Measure**] to begin the next measurement.

## Set Condition

Use Set Condition to specify Auto-Lock parameters when measuring in the Auto-Lock measuring mode. An Auto-Lock measurement ends when all measured parameters meet the conditions set on this screen. For example, if the condition of **DO** is set to 1.0%, when the measured DO holds stable within a 1.0% tolerance for the set Auto-Lock time interval, the measurement ends and the screen locks displaying the final reading.

To set the conditions:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 12](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Set Condition**.
- 3 Press [**Enter**] to select **Set Condition** or press [**Cancel**] to exit.
- 4 The **Set Condition** screen appears on the display. See [Figure 15](#).

**Figure 15** Setting a condition


- 5 Use the arrow keys to select a condition to edit. Press [**Setup**] to edit the value.
- 6 Enter the desired value. The Auto-Lock time has a range of 1–200 seconds.
- 7 Press [**Enter**] to save the new setup, exit, and return to the initial state. Press [**Cancel**] to exit without saving changes.

## System Setup

The System Setup screen is where the user sets the time on the system clock, sets the calibration interval, and enters the information for the Good Laboratory Practices (GLP) standards.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 12](#) on page 26.
- 2 Use the directional arrow keys to navigate to **System Setup**.
- 3 Press [**Enter**] to select **System Setup** or press [**Cancel**] to exit.


The **System Setup** screen appears on the display. See [Figure 16](#).

**Figure 16** System Setup

- 4 Use the arrow keys to select a value to edit. Press [**Setup**] to edit the value.
- 5 Enter the desired value. See “[Descriptions of the System Setup settings](#)” on page 32 for more information.

- 6** If needed, edit the system date and time. To modify the date and time:
- Use the arrow keys to highlight the Time field shown in [Figure 16](#) on page 31 and press [**Setup**].
  - The **Set Date & Time** screen opens. This screen displays the current year, month, day, hour, minute, and second (see [Figure 17](#)). Press the arrow keys to highlight the proper cell, then press [**Setup**] to select it.

**Figure 17** Set Date & Time

- Enter the correct value using the numeric key pad and press [**Enter**].
  - Once all changes are made, press [**Enter**] to save the settings and return to the **System Setup** screen.
- 7** Once all of the system settings are made, press [**Enter**] to save the new setup, exit, and return to the initial state. Press [**Cancel**] to exit without saving changes.

### Descriptions of the System Setup settings

**Calib Interval** is the number of hours between the calibration messages that the Meter periodically displays as a reminder to recalibrate the probe. The meter begins counting the interval time from the end of the previous calibration. When the calibration interval has elapsed, the meter displays a popup window to remind the user to recalibrate the probe. A **Calib Interval** value of **0** disables the reminder.

**Operator No** is a three-digit number from 000–200 used to identify the person operating the meter. The **Operator No** is recorded every time the meter saves data.

When stored data reaches the meter's memory limit, **Auto Delete** allows you to automatically overwrite old data with new data. If **Auto Delete** is not enabled, the meter does not save new data when its memory is full. For example, the meter permits up to 200 sets of data. When you want to store the 201<sup>st</sup> set of data, if **Auto Delete** is on, meter will delete the first set of data and store the 201<sup>st</sup> data. If **Auto Delete** is disabled, the new data will not be stored.

**CAUTION**

If **Auto Delete** is disabled, data can be lost.

---

## Language select


This meter supports Chinese and English languages.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 12](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Language Select**.
- 3 Press [**Enter**] to select **Language Select** or press [**Cancel**] to exit.
- 4 The **Language Select** screen appears on the display. Use the arrow keys to choose either **Chinese** or **English** and press [**Enter**].

## To set the air pressure (barometric pressure)

You do not usually need to set barometric pressure. The default value is 101.3kPa.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 12](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Air Pressure**.
- 3 Press [**Enter**] to select **Air Pressure** or press [**Cancel**] to exit.
- 4 The **Air Pressure** screen appears on the display. Enter the atmospheric pressure.


**Figure 18** Set Air Pressure (barometric pressure)

## To set the salinity value (**SALI Value**)

You do not usually need to set salinity. The default salinity value is 0.0 g/L.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 12](#) on page 26.
- 2 Use the directional arrow keys to navigate to **SALI Value**.
- 3 Press [**Enter**] to select **SALI Value** or press [**Cancel**] to exit.

The **SALI Value** screen appears on the display. Enter the salinity pressure.

**Figure 19** Set salinity value

## Auto power off

The meter can power off automatically after a set time period. To set the **Auto Power Off** time:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 12](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Auto Power Off**.
- 3 Press [**Enter**] to select **Auto Power Off** or press [**Cancel**] to exit.
- 4 The **Auto Power Off** screen appears on the display. Enter a power down time from 10 to 480 minutes. Enter a time of 0 to disable **Auto Power Off**.
- 5 Press [**Enter**].

When the **Auto Power Off** time elapses after the start of a measurement, the meter shuts down. If the meter is connected to the Electrochemical Data Collecting Software, **Auto Power Off** is disabled.

## Set Default

Use **Set Default** to reset the meter to its original factory settings.

- 1 From the initial state, press [**Setup**].
- 2 Use the directional arrow keys to navigate to **Set Default**.
- 3 Press [**Enter**] to select **Set Default** or press [**Cancel**] to exit.
- 4 The **Set Default** screen appears on the display. Press [**Enter**] to restore parameters to their original settings.

**Table 2** Default values for meter parameters

Parameter	Default value
DO calibration data	Zero Oxygen: 0 nA Full Scale: 1000 nA Calibration temperature: 25.0 °C Calibrated Barometric Pressure: 101.3 kPa Calibrated Salinity: 0.0 g/L
Measuring mode	Continuous
Timed Reading Mode interval	10 minutes
Calib Interval	message disabled (0 hours)
Operator No	000
Auto Delete	Enabled
Auto-Lock Set Condition	DO 1.0% saturation 1.0%
Auto-Lock time	5 s
Auto Power Off	Disabled (0 minutes)

## Prepare the DO Probe

- 1 Unscrew and remove the membrane cap.
- 2 Thoroughly rinse with distilled or deionized water.
- 3 Swing the membrane cap a few times to eliminate water drops.
- 4 Thoroughly rinse the silver and platinum parts of the stem with distilled or deionized water and dry it with a soft tissue.
- 5 Fill the membrane cap with DO filling solution to about 3/4 height.
- 6 Slowly screw the membrane cap onto the probe until finger tight. Ensure the membrane touches the platinum tip and there is no gas bubble between the membrane and the platinum tip.
- 7 Polarize the DO probe before use. Connect the DO probe with a powered-on 3200D DO meter. Polarization takes 60 minutes. The DO probe is always being polarized when connected with the meter, so there is no need to repolarize unless the DO probe has undergone maintenance or it has been disconnected from the meter for over 1 hour. If the DO probe is disconnected from the meter for less than 1 hour, polarize it for only 25 minutes before use.

## To use the DO probe

### Stirring

Some oxygen is consumed during measurement, so oxygen concentration becomes low near the membrane. The sample solution near the membrane must be refreshed, either by swinging the probe horizontally at the rate of 20–80 cm/s, or by using a stirrer to maintain consistent flow of the sample solution to the sensitive membrane.

### DO filling solution

The DO probe should be refilled with new DO filling solution after a certain period of operation, depending upon application conditions. Refill the DO filling solution every 2 weeks to 2 months.

#### CAUTION

There should not be any damage to the membrane. Avoid touching the membrane to any solid or hard materials. Ensure the membrane touches the platinum tip and there is no gas bubble between the membrane and the platinum tip. If there is any damage to the membrane, change to a new membrane cap.

---

## Calibrate the DO Probe

To obtain accurate measurements, calibrate the DO probe before use. The meter can perform a Zero Point Calibration and a Full Scale Calibration.


Figure 20 Calibrating the probe

### Zero oxygen calibration (Calib Zero)

- 1 Rinse the DO probe with distilled water and put it in a zero oxygen solution (5% sodium sulfite solution or as described in the instruction manual for the DO probe).
- 2 From the initial state or while measuring, press [**Calibrate**].
- 3 Select **Calib Zero**. Press [**Enter**] to select Zero Point Calibration as shown in [Figure 21](#). During calibration, press [**Setup**] to set calibration parameters such as barometric pressure and salinity. Press [**Cancel**] to end calibration.


Figure 21 Zero Oxygen calibration

- 4 When reading becomes stable, press [**Enter**] to automatically record the Zero Point value. This ends calibration.
- 5 Press [**Cancel**] to finish calibration.

## Full scale calibration (Calib Full)

- 1 Take the DO probe out of solution and thoroughly rinse it with distilled water.
- 2 Dry it with a soft tissue on the membrane surface.
- 3 To calibrate with air, put the probe in a well-ventilated area. Otherwise, calibrate with water saturated with air.
- 4 Swirl the probe at a speed of 20–80 cm/s or stir the solution to produce a similar linear speed near the membrane. Keep the flow speed constant during calibration. Ensure there is no air bubble trapped under the membrane.
- 5 From the initial state or while measuring, press [**Calibrate**].
- 6 Select **Full Scale**. Press [**Enter**] to open Full Scale Calibration as shown in [Figure 22](#). During calibration, press [**Setup**] to set calibration parameters such as barometric pressure and salinity.

**Figure 22** Calib Full

- 7 When reading becomes stable, press [**Enter**] to automatically record the Full Scale response. This ends Full Scale Calibration.
- 8 Press [**Cancel**] to finish calibration.

## Perform another calibration


After performing a calibration, the meter displays a prompt asking if another calibration is wanted. For instance, after Zero Point Calibration, the meter will ask if a Full Scale Calibration is wanted, and vice versa. During calibration, press [**Mode**] to switch between Zero Point Calibration and Full Scale Calibration.

## Perform a Measurement

During measurement, make sure the tip of the probe is completely submerged in the solution. Place the probe in a location where the solution can flow freely around the measuring tip.

### To measure DO

- 1 When the DO probe is used for the first time, rinse it with distilled water and put it in sample solution.
- 2 Calibrate the meter. See “[Calibrate the DO Probe](#)” on page 39.
- 3 From the initial state display, press [**Setup**] and select a measuring mode and the DO parameter. See “[Measuring Mode](#)” on page 26.
- 4 Put the probe in the sample solution and press [**Measure**] to start taking a measurement according to the selected mode. See [Figure 23](#). The upper area of the screen shows the current measuring mode and system time. The central area of the screen shows DO data, the corresponding DO current value, and temperature. The bottom area shows the most recent calibration data.

**Figure 23** The Measurement Mode display

The measuring and display methods vary slightly with different measuring modes. During measurement, you can calibrate the probe and reset the parameter. After measurement, you can save and output measuring data.

- 5 After a measurement, press [**Save**] to save data. Press [**Output**] to print data. Press [**Cancel**] to end measurement.

## To measure saturation

- 1 When the DO probe is used for the first time, rinse it with distilled water and put it in sample solution.
- 2 Calibrate the meter. See “[Calibrate the DO Probe](#)” on page 39.
- 3 From the initial state display, press [**Setup**] and select a measuring mode and the **SATU** parameter. See “[Measuring Mode](#)” on page 26.
- 4 Put the probe in the sample solution and press [**Measure**] to start taking a measurement according to the selected mode.
- 5 After a measurement, press [**Save**] to save data. Press [**Output**] to print data. Press [**Cancel**] to end measurement.


## To measure temperature

In any measuring state, the meter directly displays current solution temperature.

## To view other parameters while measuring

During a measurement, the user still can view other parameter not currently in use without changing the operational mode. For instance, while measuring with the **DO** parameter, a user can view, save, and print the **SATU** parameter.

- 1 While performing a measurement, press [**Mode**]. The screen display highlights the measuring window as shown in [Figure 24](#).

**Figure 24** Viewing parameters while measuring

- 2 Press [4/◀] or [6/▶] to switch the view to the other parameter.
- 3 If you do not press a key for a few seconds, the meter will automatically exit from viewing mode and return to the regular initial screen.

## Managing Data

### To save data

The Agilent 3200D Dissolved Oxygen Meter can save 200 sets of data. The procedure to save data varies with different measuring modes.

- In continuous mode and Auto-Lock reading mode, press [**Save**] to save data when readings become stable.
- In timed reading mode, you can still save data manually, but the meter automatically saves measuring data at a fixed periodic time interval. See “[Measuring Mode](#)” on page 26.

### To delete data

You can delete individual data entries or delete all data at once.

- 1 In either the initial screen or while measuring, press [**View**].
- 2 Select the data and press [**Delete**].


### To output data

You can output current measuring data, the most recent calibration data, or saved data. Control the data output from the connected software.

## To view data

The Agilent 3200D Dissolved Oxygen Meter allows you to view parameters such as the last calibration data and current parameter. It also allows you to set the parameter directly and view saved data. The meter stores data according to the parameter measured. All saved data meets GLP standards. The meter can store 200 sets of measurement data.

- 1 In the initial state, press [View].
- 2 Navigate to the type of data to view and press [Enter]. See Figure 25.

**Figure 25** View selection

## To View DO Data

The meter stores data according to the parameter measured. All saved data meets GLP standards. The meter can save 200 sets of concentration or saturation data.

- 1 In the initial state, press [View] and select **View Saved DO**. See Figure 26.

**Figure 26** Viewing saved data

The upper area shows the present viewing mode and actual

storage amount. Each page can display up to 10 saved data. The display format depends on the viewing mode.

- 2 The data displayed include **save time** and **operator no**. Press the arrow keys to scroll through the saved data.
- 3 Press [**Delete**] to delete data.
- 4 To output stored data, connect the PC with the meter (see “To output data” on page 44.) and press [**Output**]. See Figure 27 on page 47 for the output format.


```
=====
MODEL          3200D
              DISSOLVED OXYGEN
              METER
SOFTWARE        VER 1.00
PRINT TIME      09:03:28
                  2010/06/20
OPERATOR NO     000
*****
STORED NUM:    002
*****
NO:001
OPERATOR NO:   000
STORED TIME:   10:43:00
                2010/06/20
PRESSURE:       101.3kPa
SALINITY:        0.0g/L
ZERO:           0nA
FULL :          1000nA

CURRENT:         0nA
DO:               0.00mg/L
SATURATION:     0.0%
TEMP:            25.0c
*****
NO:002
OPERATOR NO:   000
STORED TIME:   10:43:00
                2010/06/20
PRESSURE:       101.3kPa
SALINITY:        0.0g/L
ZERO:           0nA
FULL :          1000nA

CURRENT:         0nA
DO:               0.00mg/L
SATURATION:     0.0%
TEMP:            25.0c
=====
```

**Figure 27** Output format

### To View the Last Calib

- 1 In the initial state, press [View], select **Last Calib Data**, and then press [Enter]. The display is shown in [Figure 28](#). The upper area displays last calibration data. The bottom area displays current calibrated parameters.

**Figure 28** Viewing Last Calib

- 2 To print current data, connect the meter with the PC. Press [**Output**] to print data. See "[To output data](#)" on page 44 for more information. See [Figure 29](#) on page 49 for the output format.

```
=====
MODEL          3200D
              DISSOLVED OXYGEN
              METER
SOFTWARE        VER 1.00
PRINT TIME      09:14:03
                  2010/06/20
OPERATOR NO     000
*****
DO CALIB DATA
CALIB TIME:    08:12:00
                2008/01/09
OPERATOR NO:   000

PRESSURE:       101.3kPa
SALINITY:        0.0g/L
ZERO:           0nA
FUL:            1000nA
TEMP:           25.0c
=====
```

**Figure 29** Output format

## Turn off the 3200D Dissolved Oxygen Meter

- 1 Before turning off the meter, save any data that you want to keep. See “[To save data](#)” on page 44.
- 2 Press [**On/Off**] to turn off the meter.

When not in use, all probes should be soaked in distilled water. If the meter is out of use for a long time, check the following:

- Disconnect the power adaptor from the power line to avoid damaging the adaptor and the meter.
- The socket of the meter must be kept clean and dry. Avoid contact of the socket with acid, alkaline and salt solution.
- The input terminal of the meter (connector of the DO/ATC probe) must be kept clean and dry. If the DO probe is used in a strongly humid place, wipe the DO probe plug with clean gauze.
- If the DO probe is not used, store it in distilled water. Do not soak it in sodium sulfite solution. If the sodium sulfite solution leaks into the DO probe chamber, the probe performance will deteriorate.
- Unscrew the membrane cap and empty the filling solution.
- Rinse the cathode and anode components with distilled water and then dry with a soft tissue.
- Store the DO probe in a dry place.
- For a new DO probe, a DO probe refilled with new DO filling solution, or a DO probe with a new membrane cap, polarize the DO probe for more than one hour before use.

## Software Operation

See also “[Installing Optional Software](#)” on page 14.

### To use the optional G4390A Electrochemical Data Collecting software

- 1 Connect the meter to the power adapter.
- 2 Press [**On/Off**] to turn on the meter.
- 3 Connect the meter to the computer using the USB cable.
- 4 Run the Electrochemical Data Collecting Software in the computer and the software will recognize the meter model and type. The software provides the functions corresponding to the meter type. Refer to the software documentation for details.

### To use the Data Printing software

- 1 Connect the meter to the power adapter.
- 2 Press [**On/Off**] to turn on the meter.
- 3 Connect the meter to the computer using the USB cable.
- 4 Launch the printing software.
- 5 With the meter taking measurements or when viewing calibration data or stored data, press [**Output**]. The software receives the data. Once received, print the data as desired. For details, refer to the Data Printing software documentation.

## **2      Operation**

## 3

# Troubleshooting and Maintenance

Troubleshooting the 3200D Dissolved Oxygen Meter 54

General Troubleshooting Procedure 56

Meter Self-Diagnostics and Messages 58

Maintenance 61

Consumables and Ordering Information 62

This chapter describes how to verify whether the meter is working properly and how to maintain the meter.

## Troubleshooting the 3200D Dissolved Oxygen Meter

**Table 3** lists common problems, possible causes, and how to resolve them.

Note that improper operation of the DO probe causes abnormal readings. During measurement, soak the probe tip completely in solution.

**Table 3** 3200D Troubleshooting

Number	Failure mode	Failure cause	Solutions
1	No display after meter starts up. The display is not lit up.	<ul style="list-style-type: none"> <li>1 Power adaptor not installed correctly.</li> <li>2 The power supply does not meet with the requirements.</li> <li>3 The power adaptor is damaged.</li> <li>4 The power socket has poor contact.</li> <li>5 The LCD is damaged (After start up for a few minutes, the user can hear a buzzing sound after pressing [On/Off] key, but there is no display on LCD.)</li> </ul>	<ul style="list-style-type: none"> <li>1 Reconnect.</li> <li>2 Use the required power supply.</li> <li>3 Replace the power adaptor.</li> <li>4 Ensure there is a good contact.</li> <li>5 Contact Customer Service.</li> </ul>
2	Display is dim.		Check the power adapter.
3	No buzzing when pressing key	<ul style="list-style-type: none"> <li>1 You may have pressed invalid keys under the current setup.</li> <li>2 The buzzer has been damaged.</li> </ul>	The buzzer only sounds when you press a valid key.
4	No response when pressing key	<ul style="list-style-type: none"> <li>1 You may have pressed invalid keys under the current setup.</li> <li>2 The key has been damaged.</li> </ul>	Press the valid key for operation.

**Table 3** 3200D Troubleshooting (continued)

<b>Number</b>	<b>Failure mode</b>	<b>Failure cause</b>	<b>Solutions</b>
5	The reading is not stable for a long time.	<p>1 The probe has been damaged or aged.</p> <p>2 There is strong electrical signal interference source nearby (electrical leaking, strong electromagnetic field, and so forth).</p>	<p>1 Replace the probe.</p> <p>2 Remove the electrical signal interference source. Lift or move the beaker away from the source of interference. Shield the meter and beaker from the electromagnetic field. Use a ground wire with one terminal connected to the meter and the other terminal connected to the interference source.</p>

## General Troubleshooting Procedure

During use, there are many factors that may impact measurement results, including the probe, the standard solutions used to calibrate the probe, the sample solution, ambient temperature during measurement, the stirring rate, incorrect operation, and parameter setup. When the measured results are significantly different from what is expected, first determine whether the meter itself or factors other than the meter caused the error. Follow the suggestions described below to diagnose the problem. You will need to diagnose the meter based on your application conditions. For best results, follow the troubleshooting order shown below.

### Check the current response

- 1 Connect the meter to the ATC temperature diagnostic tool (5185-8390) shipped with the meter.
- 2 Turn on the meter.
- 3 Enter the measurement state. At this moment, the meter should display a current between 680 to 720 nA. When the ATC temperature diagnostic tool (5185-8390) is disconnected, the meter should display a current between 0 to 5 nA. If yes, the DO measurement of this meter is functioning well. If the display deviates significantly from these values, there is a problem with the meter. Contact Agilent service.

### Check the DO probe

Generally, after the meter is connected to a DO probe for a period of time, the current for the DO probe in air should be at a stable value ranging from 300 to 1000 nA (refer to DO probe manual for details). If not, there is some problem with the DO probe.

## Check the temperature measurement

- 1 Connect the meter to the ATC temperature diagnostic tool (5185-8390) shipped with the meter.
- 2 Turn on the meter to enter into the measurement state. With the ATC temperature diagnostic tool connected, the meter should display a temperature reading between 49.0 to 51.0 °C. If yes, the meter is correctly measuring temperature. If the meter displays a temperature reading significantly different from 50 °C, there is a problem with the meter. Contact Agilent service.

## Check other possible causes

If the meter is functioning properly, that is, the meter and probe passed the above tests, check for other possible problems by comparing readings from different standard solutions. Put the DO probe in different, fresh standard solutions. Check the current values measured. Based on the comparison, judge whether the problem resulted from the probe, the sample solution, or something else. Also consider the stirring rate used during DO measurement. The stirring rate can impact measurement results.

## Meter Self-Diagnostics and Messages

The meter supports self-diagnosis, which can find some common errors caused by meter, probe, solution, or operation. It deals with these errors differently according to their severity.

There are two levels of severity:

**Severe system errors** hinder further use of the meter or make the meter unable to finish necessary tasks. In this case, replacement or maintenance of the meter is required.

**Minor errors**, which are prompts or warnings, are caused by various reasons and can be corrected in several ways. However, you must pay attention to these errors to ensure measurement integrity and reliability. If you ignore these error messages the operation can continue.

### Severe errors

When severe errors occur in the meter, the meter will shut down and display an error message similar to Figure 30. In this case, the only operation that can be performed is to turn off the meter.



Figure 30   Severe errors

## Error warning

Error warnings indicate improper probe installation, solution problems, and probe setup during use. Read these error warnings carefully. Proper maintenance and operation will reduce errors.

**Figure 31** shows an example error warning that temperature is out of range.

**Figure 31** Error warning

is the error warning icon. The icon includes the error code.

### NOTE

This table describes all error messages which may occur in use. These error messages have different definition for different meters

**Table 4** Error warning codes

Code	Description	Solution	Memo
00	Conductivity measuring module error	Contact Customer Service.	Serious error
01	D0 measuring module error	Contact Customer Service.	Serious error
02	Temperature measuring module error	Contact Customer Service.	Serious error
03	Data storage error	Contact Customer Service.	Serious error
20	Potential is out of range	Replace electrode.	-1999.9 to 1999.9 mV
21	pH/pX is out of range	Replace electrode.	-3.000 to 21.000 pH

### 3 Troubleshooting and Maintenance

**Table 4** Error warning codes (continued)

Code	Description	Solution	Memo
22	Conductivity is out of range	Replace probe and solution.	0 to 2000 mS/cm
23	Resistivity is out of range	Replace probe and solution.	0 to 100 MΩ•cm
24	Temperature is out of range	Replace electrode and decrease solution temperature.	-6.0 to 120.0 °C
25	DO electric current is out of range	Replace probe.	0 to 4000 nA
26	pH electrode slope is out of range	Replace electrode and recalibrate it.	80 to 120 %
27	Failed to recognize pH buffers	Replace electrode, setup proper buffer group and replace buffers.	
28	Temperature of pH buffer is out of range	Cool or heat buffer.	
29	Calibrate the same buffer repeatedly	Replace the buffer.	Due to incorrect operation
30	Number of pH buffers exceeds maximum.	Remove one or more buffers from the group.	5 buffers at most
31	pH buffers conflict with each other	Remove one or more buffers with overlapping pH values.	
32	The data storage is full.	Overwrite the previous data and store new data.	Do not delete all data.
33	The internal clock has a low battery.	Set the time manually.	
34	Number of customer-defined ions exceeds maximum	Delete one or more unnecessary customer-defined ions.	
35	The maximum ion mode number stored	Delete all storage data of certain ion mode.	

# Maintenance

This section describes how to maintain and store the probe.

## Probe maintenance

- 1 Disconnect the probe from the meter.
- 2 Unscrew the membrane cap.
- 3 Empty the membrane cap.
- 4 Rinse thoroughly with distilled or deionized water, then swing the membrane cap for a few times to eliminate water drops.
- 5 Rinse the silver and platinum parts of the stem with distilled or deionized water and dry with a soft tissue.

## Short term storage

Between tests, immerse the measuring tip in distilled water or calibration solution while the probe is connected to the meter.

## Long term storage


- 1 Disconnect the probe from the meter.
- 2 Unscrew the membrane cap.
- 3 Empty the membrane cap.
- 4 Rinse the membrane cap with distilled or deionized water.
- 5 Swing the probe a few times to throw away water drops.
- 6 Dry the silver and platinum parts of the stem. Screw the membrane cap back (without adding any DO filling solution). Store the probe in a dry place at ambient temperature.

## Consumables and Ordering Information

**Table 5** Consumables and replacement parts

Order No.	Model and name	Description
G4388A	3200SA Stirrer	Electrode holder and magnetic stirrer: combined to stir solution with stable and precise speed with a large adjustable range.
G4389A	3200EA Electrode Holder	It is used to support several electrodes.
5185-8389	AC Adaptor	100–240VAC, 1A, 50/60Hz
G4388-27000	Stirring Bar	It is used with stirrer.
5190-3997	D6111 DO Probe	Plastic body, polarographic probe with built-in ATC probe, (with DO filling solution and membrane kit *3), Measuring range: 0–20 mg/L, 0–45 °C
5190-0547	DO Filling Solution	Bottles, 3 x 30 mL
5190-0548	DO Membrane Kit	3 pieces/kit
5185-8390	ATC temperature diagnostic tool	For troubleshooting temperature problems

Log onto the portal at <http://www.agilent.com/chem> for detailed information.

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