

Ultra-Inert chemistry for Trace Level Analysis

Cikui Liang, Ph.D.



Challenges and Needs of Today's Laboratories

- **Challenges**

- Qualification/quantification of trace samples
- Keep instrument up and running

- **Needs**


- Lower detection limits
- Improved stability in GC or GC/MS system

| Lower Detection Limit | |
|---|--|
| Reduce noise | Increase signal |
| Injection system (septa, liners, connections) | Sample concentration |
| Carrier gas and detector gases | Sample size |
| Leaks | Inert injection and detection port sleeves/liner |
| Temperature setting | Gas velocity or temp program rate |
| Stationary phase and column bleed | Column inertness |



What Does Column Activity look like?

column activity → peak tailing and loss of response → reduce sensitivity



1. 1-Propionic acid
2. 1-Octene
3. n-Octane
4. 4-Picoline
5. n-Nonane
6. Trimethyl phosphate
7. 1,2-Pentanediol
8. n-Propylbenzene
9. 1-Heptanol
10. 3-Octanone
11. n-Decane

Column activity is the loss of detectable analytes to the GC column through non-ideal chemical reactions and interactions.

How Important Column Inertness to Overall Flow Path Inertness?

GC Flow Path Surface Areas:

| | l (cm) | d (cm) | pi | Surface Area (cm ²) |
|--------|--------|--------|-------|------------------------------------|
| Liner | 7 | 0.2 | 3.142 | 4.4 |
| Seal | 0.4 | 0.8 | 3.142 | 1.0 |
| Column | 3000 | 0.025 | 3.142 | 235.6 |

What is an Ultra Inert GC Column?

- Not just another column claims better performance
- Not hand picked from standard production based on “exceptional performance”
- New approach to column manufacturing with significantly improved sensitivity and accuracy for trace samples
- Treated with new, proprietary processes for surface treatment and deactivation, resulting in much better peak shape for acidic, basic, and other active compounds
- The deactivated surface coupled with the low-bleed bonded phases lead to much improved analyses and lower detection limits




Test Probes and Column Activity QC Testing


- Test probes are vital to ensure the quality and reproducibility of GC columns
 - Properly deactivated
 - Contain the correct amount of stationary phase
 - consistent column-to-column relative retention time
- Test probes can either highlight or mask the deficiencies of a column
 - An organic acid
 - A base
 - An alcohol
 - Non-active probes (e.g. alkanes)

* Peak tailing or loss response of an acid (or base) indicates the column is basic (or acidic).

Weak Probes vs. Strong Probes




2,6-Dimethylphenol




2,6-Dimethylaniline

Weak Probes

Acidic and basic portion of the molecules are shielded by the methyl groups of the 2,6-dimethyl substituted phenyl ring



1-Propionic acid



4-Picoline


Strong Probes

Active end of each compound is available to interact with any active sites on the columns




Grob-Type Test Mixture - Not Probative

Agilent J&W DB-5ms
Ultra Inert
30m x 0.25mm x 0.25um
(P/N 122-5532UI)



- Elevated oven temperature at 120°C
- Probes sweep past active sites and can mask solute/column interactions.
- Least probative probes for column activity




Competitive 5ms column
30m x 0.25mm x 0.25um

Sampler: Agilent 7683B, 5 µL syringe (Agilent part # 5181-1273), 1.5 µL split injection, 4 ng each component
Carrier: Hydrogen constant pressure 37 cm/s
Inlet: Split/splitless; 250 °C, 1.4 ml/min. column flow, split flow 100 ml/min.
Liner: Deactivated single taper w glass wool (Agilent part # 5183-4647)
Oven: 120 °C isothermal
Detection: FID at 325 °C, 450 ml/min. air, 40 ml/min. hydrogen, 45 ml/min. nitrogen makeup


DB-5ms Test Mix – More Probative

**Agilent J&W DB-5ms
Ultra Inert
30m x 0.25mm x 0.25um
(P/N 122-5532UI)**




1. 2-Ethylhexanoic acid
2. 1,6-Hexanediol
3. 4-Chlorophenol
4. Tridecane
5. 1-Methylnaphthalene
6. 1-Undecanol
7. Tetradecane
8. Dichlorohexylamine


**Competitive 5ms column
30m x 0.25mm x 0.25um**



2-Ethylhexanoic acid



Dicyclohexylamine



Sampler: Agilent 7683B, 5 μ L syringe (Agilent part # 5181-1273), 1.5 μ L split injection, 4 ng each component
Carrier: Hydrogen constant pressure 38 cm/s
Inlet: Split/splitless; 250 $^{\circ}$ C, 1.4 ml/min. column flow, split flow 75 ml/min.
Liner: Deactivated single taper w glass wool (Agilent part # 5183-4647)
Oven: 125 $^{\circ}$ C isothermal
Detection: FID at 320 $^{\circ}$ C, 450 ml/min. air, 40 ml/min. hydrogen, 45 ml/min. nitrogen makeup


Ultra Inert Test Probe Mixture – QC Testing for Today's Demanding Applications

| Probe | (ng on column) | Column functional test |
|------------------------|----------------|------------------------|
| 1. 1-Propionic acid | 1.0 | Basicity |
| 2. 1-Octene | 0.5 | Polarity |
| 3. n-Octane | 0.5 | Hydrocarbon marker |
| 4. 4-Picoline | 1.0 | Acidity |
| 5. n-Nonane | 1.0 | Hydrocarbon marker |
| 6. Trimethyl phosphate | 1.0 | Acidity |
| 7. 1,2-Pentanediol | 1.0 | Silanol |
| 8. n-Propylbenzene | 1.0 | Hydrocarbon marker |
| 9. 1-Heptanol | 1.0 | Silanol |
| 10. 3-Octanone | 1.0 | Polarity |
| 11. n-Decane | 1.0 | Hydrocarbon marker |



Ultra Inert Test Probe Mixture on a Competitive Column

Fully Probative



Sampler: Agilent 7683B, 0.5 μ L syringe (Agilent part # 5188-5246), 0.02 μ L split injection

Carrier: Hydrogen constant pressure, 38 cm/s

Inlet: Split/splitless; 250 $^{\circ}$ C, 1.4 ml/min. column flow, split flow 900 ml/min., gas saver flow 75 ml/min. on at 2.0 min.

Liner: Deactivated single taper w glass wool (Agilent part # 5183-4647)


Oven: 65 $^{\circ}$ C isothermal

Detection: FID at 325 $^{\circ}$ C, 450 ml/min. air, 40 ml/min. hydrogen, 45 ml/min., nitrogen makeup



Agilent Technologies

Ultra Inert Test Probe Mixture on an Agilent J&W DB-5ms Ultra Inert column



Sampler: Agilent 7683B, 0.5 μ L syringe (Agilent part # 5188-5246), 0.02 μ L split injection
Carrier: Hydrogen constant pressure, 38 cm/s
Inlet: Split/splitless; 250 $^{\circ}$ C, 1.4 ml/min. column flow, split flow 900 ml/min., gas saver flow 75 ml/min. on at 2.0 min.
Liner: Deactivated single taper w glass wool (Agilent part # 5183-4647)
Oven: 65 $^{\circ}$ C isothermal
Detection: FID at 325 $^{\circ}$ C, 450 ml/min. air, 40 ml/min. hydrogen, 45 ml/min., nitrogen makeup


Test Probes – What Do We Learn?

- Test probes can highlight or mask the deficiencies of column activity
- Grob-type mix is not probative for inertness
- DB-5ms test mix is a good test for the 90s
- Ultra Inert test mix probes inertness and differentiates an excellent column from a mediocre one
- Well designed test mix uncovers potential adsorption of acid and base analytes and raises the bar in inertness QC

**Column Inertness - Proof in
EVERY GC column box – Performance Summary Sheet.**



Excellent Inertness and Same Selectivity



10 overlaid TIC for semivolatiles, 2 ng on-column (black-DB-5ms Ultra Inert, red DB-5ms)


Table 2. Fifty Semivolatile Compound Separation using DB-5 ms and DB-5 ms Ultra Inert Columns

| Peak number | Compound | Ava. RT of DB-5 ms (n=5) | Ava. RT of DB-5 ms Ultra inert (n=5) | *RSD% of RT | Rs of DB-5 ms | Rs of DB-5 ms Ultra inert | m/z |
|-------------|---------------------------------------|--------------------------|--------------------------------------|-------------|---------------|---------------------------|---------|
| 1 | Isophorone | 5.647 | 5.642 | 0.068 | — | — | 82 |
| 2 | 1,3 Dimethyl-2-nitrobenzene(SS) | 6.047 | 6.044 | 0.061 | 8.4 | 8.41 | 134 |
| 3 | Dichlovos | 6.351 | 6.349 | 0.059 | 9.05 | 9.03 | 109 |
| 4 | Hexachlorocyclo-pentadiene | 7.129 | 7.126 | 0.065 | 21.58 | 21.8 | 237 |
| 5 | Dimethyl phthalate | 8.084 | 8.087 | 0.071 | 23.14 | 24.1 | 163 |
| 6 | 2,6-Dinitrotoluene | 8.195 | 8.196 | 0.074 | 2.51 | 2.52 | 165 |
| 7 | Acenaphthylene | 8.342 | 8.342 | 0.080 | 3.18 | 3.17 | 152 |
| 8 | Acenaphthylene-d10(IS #1) | 8.611 | 8.612 | 0.076 | 5.57 | 5.74 | 164 |
| 9 | 2-Chlorobiphenyl | 8.737 | 8.739 | 0.088 | 2.6 | 2.74 | 188 |
| 10 | 2,4-Dinitrotoluene | 8.993 | 8.996 | 0.085 | 5.44 | 5.6 | 165 |
| 11 | Diethyl phthalate | 9.572 | 9.579 | 0.098 | 12.21 | 12.35 | 149 |
| 12 | Fluorene | 9.804 | 9.808 | 0.088 | 4.65 | 4.6 | 166 |
| 13 | 2,3-Dichlorobiphenyl | 11.153 | 11.159 | 0.095 | 24.87 | 25.63 | 222/152 |
| 14 | Hexachlorobenzene | 11.218 | 11.219 | 0.090 | 1.14 | 1.12 | 284 |
| 15 | Pentachlorophenol | 11.795 | 11.798 | 0.092 | 9.78 | 10.07 | 266 |
| 16 | Phenanthrene-d10(IS #2) | 12.357 | 12.363 | 0.088 | 9.28 | 9.34 | 188 |
| 17 | Phenanthrene | 12.426 | 12.432 | 0.091 | 1.16 | 1.13 | 178 |
| 18 | Anthracene | 12.585 | 12.591 | 0.091 | 2.66 | 2.71 | 178 |
| 19 | 2,4,5-Trichlorobiphenyl | 13.133 | 13.140 | 0.089 | 9.27 | 9.61 | 256 |
| 20 | Heptachlor | 14.001 | 14.008 | 0.090 | 14.36 | 14.80 | 100 |
| 21 | Di-n-butyl phthalate | 14.587 | 14.600 | 0.095 | 9.99 | 10.14 | 149 |
| 22 | 2,2',4,4'-Tetrachlorobiphenyl | 14.733 | 14.741 | 0.085 | 2.51 | 2.42 | 292 |
| 23 | chlorpyrifos | 14.882 | 14.892 | 0.088 | 2.50 | 2.53 | 197/97 |
| 24 | 2,2',3',4',6-Pentachlorobiphenyl | 16.247 | 16.255 | 0.083 | 22.89 | 23.28 | 326 |
| 25 | Butachlor | 17.058 | 17.070 | 0.081 | 13.17 | 13.74 | 176/160 |
| 26 | Pyrene-d10(SS) | 17.153 | 17.163 | 0.078 | 1.49 | 1.41 | 212 |
| 27 | Pyrene | 17.214 | 17.223 | 0.076 | 0.87 | 0.91 | 202 |
| 28 | p,p'-DDE | 17.901 | 17.913 | 0.075 | 11.03 | 11.34 | 246 |
| 29 | 2,2',4,4',5,6'-Hexachlorobiphenyl | 18.077 | 18.088 | 0.074 | 2.79 | 2.74 | 360 |
| 30 | p,p'-DDD | 19.162 | 19.176 | 0.075 | 7.54 | 7.04 | 235/165 |
| 31 | o,p'-DDT | 19.235 | 19.250 | 0.078 | 1.17 | 1.18 | 235/165 |
| 32 | Benzyl butyl phthalate | 20.157 | 20.177 | 0.078 | 14.76 | 15.08 | 149 |
| 33 | p,p'-DDT | 20.301 | 20.316 | 0.073 | 2.39 | 2.35 | 235/165 |
| 34 | bis(2-Ethylhexyl)adipate | 20.731 | 20.755 | 0.084 | 7.17 | 7.25 | 129 |
| 35 | Triphenylphosphate (SS) | 20.851 | 20.873 | 0.082 | 1.95 | 1.95 | 326/325 |
| 36 | 2,2',3',3',4,4',6-Heptachlorobiphenyl | 21.737 | 21.754 | 0.070 | 13.80 | 13.79 | 394/396 |
| 37 | Ben[α]anthracene | 21.805 | 21.806 | 0.033 | 1.18 | 1.23 | 228 |
| 38 | Chrysene-d12 (IS #3) | 21.834 | 21.835 | 0.039 | 0.74 | 0.75 | 240 |
| 39 | 2,2',3,3',4,5,6,6'-Octachlorobiphenyl | 21.903 | 21.904 | 0.036 | 0.94 | 1.04 | 430/428 |
| 40 | Chrysene | 21.928 | 21.928 | 0.034 | 0.80 | 0.80 | 228 |
| 41 | bis(2-Ethylhexyl)phthalate | 22.605 | 22.632 | 0.081 | 9.44 | 10.04 | 149 |
| 42 | cis-Permethrin | 24.692 | 24.720 | 0.077 | 35.23 | 36.04 | 183 |
| 43 | trans-Permethrin | 24.906 | 24.936 | 0.080 | 3.58 | 3.63 | 183 |
| 44 | Benzo[b]fluoranthene | 25.704 | 25.725 | 0.064 | 12.10 | 12.29 | 252 |
| 45 | Benzo[k]fluoranthene | 25.802 | 25.824 | 0.067 | 1.30 | 1.40 | 252 |
| 46 | Benzo[a]pyrene | 26.766 | 26.789 | 0.064 | 12.69 | 13.60 | 252 |
| 47 | Perylene-d12(SS) | 26.966 | 26.991 | 0.067 | 2.60 | 2.74 | 264 |
| 48 | Indeno[1,2,3-c,d]pyrene | 30.434 | 30.470 | 0.077 | 40.75 | 40.99 | 276 |
| 49 | Dibenz[a,h]anthracene | 30.594 | 30.635 | 0.086 | 1.69 | 1.76 | 278 |
| 50 | Benzo[g,h,i]perylene | 31.373 | 31.414 | 0.088 | 7.74 | 8.07 | 276 |

*RSD% was calculated by using the retention times for each compound on the DB-5 ms and DB-5 ms Ultra Inert columns.

Nitro-PAHs

Column: DB-5ms Ultra Inert 15 m x 0.25 mm x 0.25 μm (Agilent part # 122-5512UI)
Carrier: Helium 43.8 cm/sec constant flow
Oven: 70°C (1min) to 310°C at 20°C/min
Inlet: splitless at 250 °C
MSD: 7000A Triple Quadrupole GC/MS, ion source at 300°C and Quadrupole at 150°C



By Frank David
(RIC) and
Matthew Klee
(Agilent)

MRM chromatograms of nitro-PAHs in extract of urban air particulate sample where the concentrations for nitro-naphthalene, nitro-anthracene, nitro-fluranthene, and nitro-pyrene were 21 pg/m³, 10 pg/m³, 77 pg/m³, and 14 pg/m³, respectively.

- Deliver high sensitivity.
- Reliable determination of trace level nitro-PAHs in complex air particulate extract matrices without labor-intensive sample preparation.
- The solutes are selectively detected at pg/ μL level, corresponding to pg/m³ in air.

Column Technology for Reliable Trace Analysis

- Column bleed is only **half the story in trace analysis**.
- Only when a column exhibits **both low bleed** and **low activity** are results reliable.
- Ultra Inert columns set a new industry standard for column inertness QC testing
 - Best columns available for reactive analytes and trace analysis
 - Active analytes adsorbing on the column stops productivity in its tracks.
 - Start with an Ultra Inert Column for guaranteed performance upon installation.

